

Known Defects and Anomalies

Introduction

The Known Defects and Anomalies (KDAs) table consists of system actions that do not meet performance expectations established in VBECS design documents. Some defects and anomalies require user workarounds such as being directed to view information in one report that was expected in another report. Defects and anomalies do not require a workaround if the software performs acceptably. All defects and anomalies are classified as minor (no expectation of injury to the patient, operator, or bystander as a result of software failure) or the workaround sufficiently mitigates the defect to an acceptable risk.

Note: All system errors/shutdowns occur where the user would normally be prohibited from proceeding to process the unit or patient in VBECS; training users to STOP and evaluate the correctness of the action is strongly recommended.

Related Manuals and Materials

- *VistA Blood Establishment Computer Software (VBECS) Installation Guide*
- *VistA Blood Establishment Computer Software (VBECS) Technical Manual-Security Guide*
- *VistA Blood Establishment Computer Software (VBECS) User Guide* (The KDAs are referenced as *Appendix E* but are maintained with a separate revision table).

Printing Requirements

This document is formatted to print on legal (8.5 x 14 inch) paper.

How the Known Defects and Anomalies is Organized

- The table is organized by the option where the issue occurs in VBECS.
- Description of the Issue, Risk Assessment, and Affected Security Role columns provide pertinent information about the defect or anomaly.
- Recommended Workarounds are provided with the Additional Comments column providing more detail as needed.
- Highlighted rows indicate an item that has changed or was added since the last released version of the Known Defects and Anomalies.
- Tentative Schedule column contains the projected release where the anomaly will be corrected.
- Throughout VBECS is the only section where the item may occur in various places in the application not in just one option.

Terms

See the VBECS User Guide Glossary for definitions of other terms and acronyms used in this table.

Security

Six security roles are available in VBECS. Security Levels and User Roles from most to least restrictive are: Blood Bank Technologist (most restrictive), Enhanced Technologist, Lead Technologist, Traditional Supervisor, Enhanced Supervisor, Administrator/Supervisor (least restrictive). In order to simplify this analysis three categories may be used to clarify the AFFECTED USER category: Administrator (Administrative Supervisor and Enhanced Supervisor); Supervisor (Traditional Supervisor and Lead Technologist, and above); All Users (Enhanced Technologist and Technologist, and above). Additional users are a System Administrator (Windows Network Administrator) and a VBECS Administrator (Server Administrator) which are restricted, usually to a primary and a backup person.

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Through-out VBECS	Comments >200 characters in length can be entered in VBECS, but only the first 200 display on the report.	No risk/ Low impact.	All users	The “other” field for comments is free text and accepts more than 200 characters. Users are advised to restrict comments to 200 characters so the entire comment is visible on reports.	None.	CR 2,870 CR 2,871 CR 2,872 CR 2,873	Report Tool Upgrade Project Proposal NSR 20080915
Through-out VBECS	VBECS reports cannot be printed as two sided documents.	No identified risk/ No identified impact.	All users	None available.	The current report generating tool in VBECS does not support two sided documents.	CR 2,992	Report Tool Upgrade Project Proposal NSR 20080915
Through-out VBECS	VBECS map an institution as an orderable location if the VistA institution file is more than 5 digits in length.	Low risk/Low impact.	All users	Set VistA location to 5 digits if this location will be used to map orders to VBECS in the VBECS Administrator.	VBECS can only accept a VistA location 5 digits in length. Sites are advised to limit locations that order blood to 5 digits in length.	CR 3,229 HD 625192	Maintenance TBD
Through-out VBECS	VBECS does not compare the ProcessingID field of an inbound order message to the interface control table. A site could incorrectly configure their VistA test account to send orders to the production VBECS.	Very Low risk/Low impact.	All users	Verify configuration of accounts by processing test orders before releasing the test or production VBECS account for use after system maintenance or mirroring..	None.	CR 3,243 HD 626467	Maintenance TBD
Through-out VBECS	VBECS will display a system error message and will not function if the administrator sets a future date in VBECS settings.	Very Low risk/Low impact.	All users	VBECS Administrators must set the date on the server to a current or past date. Verification of settings is recommended.	This is identified immediately by the first user logged into the system in Test and in Production.	CR 2,153	Maintenance TBD
Through-out VBECS	Outgoing completion messages from VBECS will stay in a status of sending when services are stopped when a message is sending or when an unexpected exception is received from the receiving application.	No identified risk/Low impact.	All users	The messages cannot be resent by a user. Contact product support or create a Remedy ticket. The VBECS team will make the correction.	None.	CR 2,715	Maintenance TBD
Through-out VBECS	Printed reports in VBECS use gray lines to designate the end of patient or unit data. The lines should be darker to make reading the report easier.	No identified risk/ No identified impact.	System Administrator	None required.	The shading on screen and in printed form varies with the monitor and printer settings.	CR 1,514	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Through-out VBECS	A patient name will be truncated in these screens when one of the name fields (Last, First, or Middle) exceed 22 characters: Blood Unit Details tab Patient Specimen Selector control Patient Order Selector control Patient Order Selector for Report control VBECS will only display the first 25 characters of the name.	Low risk/Low impact.	All users	None available.	The full name is displayed elsewhere on the screens identified or it is provided as read only information and is not used to make decisions for the selected patient.	CR 2,658	Maintenance TBD
Through-out VBECS	Patient names are truncated or improperly formatted on these reports when a full name exceeds 25 characters: Outgoing Shipment Invoice Cost Accounting Report Unit History Report Exception Report Inappropriate Transfusion Request Report	Low risk/Low impact.	All users	Users are instructed to correct the names on printed copies of reports to be saved.	None.	CR 2,657	Report Tool Upgrade Project Proposal NSR 20080915
Through-out VBECS	On some reports, VBECS prints the page number but does not include the total number of pages printed (e.g., 1 of 75).	No identified risk/No identified impact.	All users	None required. The total number of report pages display in the print preview and may be used to verify that a complete report has printed.	User may elect to manually update page numbers on impacted reports.	CR 1,967, CR 2,021, CR 2,855	Report Tool Upgrade Project Proposal NSR 20080915
Through-out VBECS	If the user enters a full last name and clicks the ellipsis to search for a patient when entering a unit in Incoming Shipment, VBECS retrieves patients that match the entered name along with others that match the first four letters of the searched name.	Low risk/Low impact.	All users	None required as our feedback from the field sites and other user groups is that the user is trained to enter the full SSN for a patient and not use a pick list to select a patient for any option in VistA or VBECS.	The returned search results are displayed with the best matches at the top of the list. User may have a longer patient list presented to find the required patient.	CR 1,635	Report Tool Upgrade Project Proposal NSR 20080915
Through-out VBECS	When a 10-character last name is entered in the Patient Select Tool, VBECS looks for a specimen UID.	No identified risk/Low impact.	All users	Feedback from the field sites and other user groups is that the user is trained to enter the full SSN for a patient and not use a pick list to select a patient for any option in VistA or VBECS. Search by the initial of the last name and last four digits of the patient ID (SSN), or the full patient ID.	Using standard patient search criteria avoids this potential confusion.	CR 2,004	Maintenance TBD
Through-out VBECS	VBECS records the standard patient name format: the first initial included in the data after the second comma, e.g., last name, first name, middle initial. This conflicts with the VistA standard patient name presentation.	No identified risk/No identified impact.	All users	None required as the patient is correctly identified and name is readable, simply formatted differently from VistA's presentation.	The patient names are not considered a unique identifier. The ID number (DFN, ICN) is used to identify the patient record. The legacy system may record data in the middle-initial field that is not included in the VBECS display of the patient name.	CR 1,903	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Through-out VBECS	When only a few non-specific characters are entered in the patient search field, no matches are found and the option times out without displaying a message to the user.	No identified risk/Low impact.	All users	Enter the full last name, or the full patient ID, last name, last name initial, and last four digits of the patient ID (SSN), or the full patient ID. User will enter standard patient search criteria.	This is not a valid search entry. As this does not result in a patient display, the user must reenter the correct search information.	CR 1,901 CR 2,927	Maintenance TBD
Through-out VBECS	A system error occurs if a user enters non-standard (symbols or punctuation) information into the patient, unit, or product code fields and uses the search.	No Risk/ Low Impact	All Users	Do not enter non-standard information into the patient, unit, or product code fields.	None.	CR 2,768 CR 2,774	Maintenance TBD
Through-out VBECS	When multiple windows are overlaid in the VBECS application, they may appear incompletely drawn.	No risk/ Low impact.	All users	Minimize and maximize the VBECS application or remote desktop connection window to refresh the screen.	The user can easily identify that the data screen is incomplete and cannot use or continue until refreshed.	CR 1,771	Maintenance TBD
Through-out VBECS	Certain combinations of patient data and the printer settings create a report too large to fit on a page leading to a Crystal Report failure (error).	No identified risk/No identified impact	None	Depending on the report: <ul style="list-style-type: none"> Print the report with shorter date ranges. Print the report by section. If you find a section that fails, request the report without that section. Then locate another report that also displays the problem data, view and print. 	Contact your VBECS system administrator to check the printer driver per the Installation Guide instructions. This is typically discovered during local validation and corrected prior to production installation. This has been identified on the Patient History Report and Cost Accounting Report, but may occur in any report.	CR 2,757 HD 372543	Report Tool Upgrade Project Proposal NSR 20080915
ABO/Rh Confirmation	VBECS displays incorrect message to user when they select to perform ABO/Rh confirmation on a unit that has been discard.	No risk/No impact.	All users	None required. The unit is not selectable and only the information message is incorrect.	VBECS displays a message to user stating that the unit is the wrong product type instead of that the unit is in a final status.	CR 2,813 HD 399413	Maintenance TBD
ABO/Rh Confirmation	Confirming more than three units in a batch results in workload multiplication. Increasing the batch confirmed in two unit increments increases the workload count again. (e.g. 3-4 units will have double workload, 4-6 units have triple workload, and so on).	No identified risk/Low impact.	All users	None available.	None.	CR 2,830 HD 400511	Maintenance TBD
ABO/Rh Confirmation	A repeat ABO/Rh Confirmation on a unit that is VBECS does not remove patient assignment from a unit when the unit is quarantined.	No identified risk/No identified impact	All users	None required.	The unit cannot be issued from VBECS while it is quarantined.	CR 2,849	Maintenance TBD
ABO/Rh Confirmation	If units are selected and users change the search criteria, the selected units stay selected even if they do not match the current search criteria.	No risk/No impact.	All users	Deselect units that will not be tested in this option. Units not completely tested are not available for patient selection.	VBECS presents selected units for review prior to performing testing and entering results.	CR 1,578	Testing Grids
ABO/Rh Confirmation	VBECS displays “Abo/Rh” instead of “ABO/Rh” in the ABO/Rh mismatch message.	No risk/No impact.	All users	None required. Format is understandable to users.	None.	CR 1,723	Testing Grids

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
ABO/Rh Confirmation	Patient assignment is not released automatically when a unit is quarantined due to discrepant ABO/Rh retype testing.	Very low risk/No impact.	All users	None required. The unit is quarantined and cannot be selected for a different patient or issued for this patient.	VBECS prevents issuance of a unit with discrepant testing. Quarantined units are clearly marked in VBECS and cannot be issued.	DR 2,982	Testing Grids
ABO/Rh Confirmation	A system error occurs when attempting to save Unit Confirmation tests with an invalid test between valid test results.	No identified risk/Low impact.	All users when attempting to save an invalid confirmation test in the middle of a batch.	Save and exit the option immediately after the invalid unit confirmation test row. Reenter the option and select any remaining units for testing.	User is immediately aware that system data must be reentered.	CR 2,396	Testing Grids
Accept Orders: Accept an Order	A VBECS system error occurs when a user tries to print out an order and the patient location was not included in the order information from CPRS.	No identified risk/Low impact.	All users	Cancel the order without a location. Request a new order using the CPRS GUI version.	The CPRS GUI version requires a patient location for the order to be submitted. The CPRS roll and scroll version allows an order to be place without a patient location. Using the CPRS order dialog roll and scroll version is not recommended.	CR 3,199 HD 540813	Maintenance TBD
Accept Orders: Accept an Order	The VBECS order dialog "Requirements" field that displays the CPRS Modifier is limited to 29 characters.	Low Risk/Low impact.	All users	Limit selectable modifiers in CPRS to 29 characters.	Requirements displayed are CPRS order Modifier entered as part of a blood component order such as Irradiated or Leukoreduced.	DR 3,962	Maintenance TBD
Accept Orders: Accept an Order	Orders that are accessioned in the Laboratory package with a date/time collected of T@U are not accepted by VBECS.	Low Risk/Low impact.	All users	Enter a time and collection date for specimen collection in the Laboratory package to allow VBECS to process the specimen. OR Do not use T@U to accession VBECS orders in the Laboratory package.	None.	CR 3,063	Maintenance TBD
Accept Orders: Accept an Order	VBECS only displays the first 20 characters of the ordering location when the user clicks on the Ordering Details button to see the Ordering Division field.	No identified risk/Low impact.	All users	None available.	Ordering locations in Vista can be shortened or changed if the uniqueness of a location cannot be determined with the first 20 characters.	CR 2,919	Maintenance TBD
Accept Orders: Accept an Order	The display restricted or assigned units button in Accept Orders was disabled for performance.	No identified risk/ No identified impact	All users	The Blood Availability Report displays units assigned and restricted for all patients.	None.	CR 2,903	Maintenance TBD
Accept Orders: Accept an Order	On the Accept Orders screen, clicking the Received and Wanted column headers does not sort on the date and time of the Pending Order List.	No identified risk/ No identified impact.	All users	None available.	All order information is viewable. This is the column sort function which may or may not be used. As sites have stated that they select an order by entering the specimen UID or patient information this is a rarely used search criteria.	CR 2,098	Testing Grids

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Accept Orders: Accept an Order	When the user clicks Print and selects certain date ranges, a shutdown error may occur due to a Crystal Reports bug.	No identified risk/ No identified impact.	All users	Restart VBECS and reprint the report. Select a different date range to include the dates required.	All data is viewable by selecting different ranges. Crystal Reports does not publish specific elements that cause the error.	CR 1,774	Report Tool Upgrade Project Proposal NSR 20080915
Accept Orders: Accept an Order	The emergency order check box is not enabled unless there is a CPRS order and a patient specimen accessioned in VistA.	No identified risk/Low impact.	All users	Accession the component order in VistA without the specimen to use the emergency order function in VBECS. The VistA order can then be changed when the specimen is received.	Alternately, process the emergency issue of units manually, per local policy, until the specimen is received and accessioned. Processing an emergency blood component issue is an uncommon occurrence. Generally, a specimen is collected and may be available at the time though testing may be incomplete at the time of issue depending on the patient population of a facility.	CR 1,604	Testing Grids
Accept Orders: Accept an Order	The Evaluate MSBOS button is disabled for TAS pre-op orders. TAS only orders cannot be marked as inappropriate.	No identified risk/ No identified impact.	All users	None available.	When a surgical procedure has no TAS or component recommendation, VBECS allows the order but does not mark the request for report review. The provider has ordered a TAS when there is no pre-op recommendation.	CR 1,944	Testing Grids
Accept Orders: Accept an Order	A VBECS system error will occur when a user attempts to save a rejection comment that is longer than 255 characters.	No identified risk/ No identified impact.	All users	User must log back into VBECS and enter a comment with a message string that does not exceed 255 characters including the canned comment text.	VBECS configuration of standardized (canned) and free-text comments rarely require this message length.	CR 2,149	Testing Grids
Accept Orders: Accept an Order	A specimen may not be marked unacceptable when Maintain Specimen is accessed during the acceptance of an order.	Low Risk/Low impact.	All users attempting to associate an unacceptable specimen with a VBECS order.	None required. The specimen may be marked unacceptable without relation to the order by using the Maintain Specimen option.	Directly access Maintain Specimen to mark a specimen unacceptable without cancelling the order. This arose from a misunderstanding of the system functionality. The specimen can be marked unacceptable and the order can be cancelled in each one's appropriate option.	DR 2,851	Maintenance TBD
Accept Orders: Pending Order List	The printed version of the Pending Order List differs in appearance from the screen view.	Low risk/Low impact.	All users	None required.	The content of the Pending Order List printout is the same as the screen display.	CR 1,401	Report Tool Upgrade Project Proposal NSR 20080915
Accept Orders: Pending Order List	VBECS requires more than 15 seconds to print a Pending Order List report.	Low risk/Low impact.	All users	None required.	The amount of time required to generate the report is proportional to the number samples processed by the blood bank.	CR 2,622	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Administrative Data Report	The number of transfused units on the Administrative Data Report does not match with the total number of units transfused on other VBECS reports.	Low risk/Low impact.	All users	The total of transfused units includes units based on completion date and time so units may not be included in specific date range reports. The report also includes units that were inactivated.	If the transfusion of the blood unit was invalidated and then the unit was transfused again it will count as 2 transfusions on the report.	CR 2,838 CR 3,079	Report Tool Upgrade Project Proposal NSR 20080915
Administrative Data Report	The Administrative Data Report does not contain all of the plasma product types in the final count of units received via Incoming Shipment.	Low risk/Low impact.	All users	Manual calculations must be used to derive the missing data from the number of units received.	The FP24 and apheresis plasma products are not counted as FFP received.	DR 3,539 HD 338438	Report Tool Upgrade Project Proposal NSR 20080915
Administrative Data Report	The total number of red cell units transfused units on the Administrative Data Report does not match the total transfused on the C:T Ratio Report.	Low risk/Low impact.	All users	Manual calculations must be used to compare the total number of units transfused by using this report and the division C:T Ratio Report.	None.	DR 3,539 HD 338438	Report Tool Upgrade Project Proposal NSR 20080915
Antibodies	The Antigen Negative Compatibility Percentage field allows the entry of a decimal that causes the reversal of the entry (e.g., user entry of "1.5" becomes "51") and may be saved.	Low risk/Low impact.	Administrator	Enter whole numbers; do not enter decimals. Check the accuracy of the entry before saving.	This data is managed by security role and is for information only.	CR 1,842	Maintenance TBD
Antibodies	Anti-A,B is not selectable as a patient antibody.	Low risk/Low impact.	Administrator	Users can enter the Anti-A,B antibody in the patient's Special Instructions.	This antibody is historical from Vista and is not utilized for selection of blood components for transfusion.	CR 1,872	Maintenance TBD
Audit Trail Report	When the user inactivates a unit's ABO/Rh confirmation test results or inactivates a unit record and enters the required comment, VBECS does not print the comment on the Audit Trail Report.	Low risk/Low impact.	Supervisor	Immediately print the Audit Trail Report for this activity and manually complete the comment; save it for review.	The inactivation is recorded and maintained. The unsaved comment does not impact patient testing or transfusion records and is an infrequent occurrence.	CR 1,824 CR 2,353	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	Changes made to a unit's log-in CMV or Sickie Cell status are not displayed on the report.	Low risk/No identified impact	Supervisor when blood unit CMV or Sickie cell changes have been processed.	The initial and updated information in the Unit History Report	The change history is correctly maintained with each unit record. It is standard VBECS behavior to display this type of data change on the Audit Trail Report.	CR 2,318	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	VBECS displays only the ISBT 128 five-digit product code in the sub-headers for the Blood Unit Changes or Blood Unit Financial Changes sections of the report.	No identified risk/No identified impact	Supervisor when a blood unit changes have been processed.	None required. This is the product code without the donation type and divisions which are not required for accurate interpretation of this report entry.	The full eight-digit product code is available in the Unit History Report. Format is understandable to users.	CR 1,824	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Audit Trail Report	Changes to a blood unit's disease marker testing do not appear in the report.	No identified risk/No identified impact .	Supervisor when a blood unit's disease marker status has been changed.	None required.	Changes to the biohazardous marker are displayed in the Edit Unit Information section, testing details are available in the Unit History Report. It is standard VBECS behavior to display this type of data change on the Audit Trail Report.	CR 1,907	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	Special Instructions and Transfusion Requirements comments may not display completely on the report.	Low risk/Low impact.	Supervisor when SI and TR have been changed (removal is security restricted to Supervisor).	None required.	The Patient History Report can be viewed for the entire comment. Although this is an infrequent occurrence, it is non-standard VBECS behavior.	CR 1,968	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	VBECS allows the user to click the report header; it displays as a separate tab.	No identified risk/ No identified impact	All users who look at the print preview of the audit trail report.	None required.	This is a section header. The information is correct.	CR 1,531	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	The report displays changes to the database during MSBOS configuration even though a value was not actually changed by the user. This is due to fields in the database changing from null to NO.	No identified risk/ No identified impact	Supervisor reviewing the Audit Trail after initial configuration.	None required.	This is the day one data update of the database with the standardized default data on the MSBOS. Additional local configuration of the MSBOS is also displayed.	DR 1,312	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	The report does not differentiate reagent types by case size when indicated by a letter only when both have Minimum Reagent Levels have been defined on the same day, e.g. K, k or C, c.	No identified risk/No identified impact.	Supervisor	Update the minimum levels on different days and print the report on each day.	None.	CR 2,730	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	The Audit Trail Report does not display the Date/Time data was originally saved or the user comment entered for the change. Audit Trail Report entries for Units does not display Unit ABO/Rh at log-in, Expiration date/time, Date/Time received, or the user comment entered for the change.	Low risk/Low impact.	Supervisor reviewing the Audit Trail.	Immediately print the Audit Trail Report for this activity and manually complete the comment and other pertinent details; save it for review.	See the Unit History report for the missing data saved in other options.	CR 2,334	Report Tool Upgrade Project Proposal NSR 20080915
Audit Trail Report	The Audit Trail Report does not display the full comments a user can enter when using the Justify ABO/Rh or Inactivate Patient Transfusion Requirement option. The report will only show the first 100 characters.	Low risk/Low impact.	Supervisor reviewing the Audit Trail.	The entire set of comments (max 255 characters) is displayed in the Patient History Report, ABO/Rh Change History Section.	None.	CR 2,766 CR 2,869	Report Tool Upgrade Project Proposal NSR 20080915

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Blood Availability Report	The Blood Availability Report does not include units that went to a final status on the day the user requests the report start and stop on (e.g. start 1/30/11 and end 1/30/11).	Low risk/Low impact.	All users	Expand the date range of the report to include the units in a final status.	None.	CR 3,132	Report Tool Upgrade Project Proposal NSR 20080915
Blood Availability Report	Running the Blood Availability Report for units with no disposition (Available Expired) first then the Daily Inventory (All Units) Report will result in an error on the Daily Inventory (All Units) Report and no values displayed.	No identified risk/ No identified impact.	All Users	Do not run the Blood Availability Report for units with no disposition (Available Expired) before the Daily Inventory (All Units) Report.	None.	CR 3,234	Report Tool Upgrade Project Proposal NSR 20080915
Blood Availability Report	"Units with Final Disposition (Not Transfused)" report, defaults a start date of 1/1/1900 and end date of 12/31/9999. Even though the user manually enters new start/end dates, they are not used when the report data is generated.	Low risk/Low impact.	All users	Use the Custom Report option to generate the report.	The time to return the report will be greater if the user accepts the default dates.	CR 2,463	Report Tool Upgrade Project Proposal NSR 20080915
Blood Availability Report	Blood Availability Report does not work with the report scheduler.	No identified risk/ No identified impact.	All Users	Users must run the report when needed.	This report is printed regularly and reviewed for daily work so system performance is not impacted.	CR 2,793	Report Tool Upgrade Project Proposal NSR 20080915
Blood Products	VBECS allows a user to activate a blood product code for a shipper without an associated HCPCS code or text.	Low risk/Low impact.	All users	Add the HCPCS code after the blood product code is activated.	The user can proceed to activate a unit without the code. This may be preferred functionality when the code is unknown so all blood products can be entered without delay for patient availability.	CR 1,217	Maintenance TBD
Blood Products	A system error occurs when the user attempts to change the shipper FDA Registration Number while a warning message is displayed in the Cost or Return Credit fields.	Low risk/Low impact.	All users	Clear the incorrect data entry prior to changing the shipper information.	This information would be updated very infrequently.	CR 1,921	Maintenance TBD
Canned Comments	System error occurs when an attempt is made to use Maintain Comments if the user who configures a division in VBECS Administrator is not added to the division as a VBECS user.	Low risk/Low impact.	Administrator	The user that configured the division must be added to the users list in VBECS.	None.	CR 2,241	Maintenance TBD
Canned Comments	Inactivation of all the canned comments for a category will prevent a user from completing any process that uses these canned comments.	Low risk/Low impact.	Administrator	Do not delete all of the canned comments.	Comments can be recreated by user with appropriate security.	CR 2,243	Maintenance TBD
Component Classes	VBECS rounds the minutes entered in the maximum transfusion time field to a whole number.	Low risk/Low impact.	Administrator	Check the entry for validity prior to saving.	Explanation for system behavior with no user impact. Transfusion times are reported in minutes.	CR 1,992	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Configure Daily QC	VBECS does not allow a user to save changed values in the Non-Commercial template.	Low risk/Low impact.	Administrator	When editing the Non-Commercial template, the user must first click the Commercial option, then re-select the Non-Commercial option. The OK to save button is then enabled and the user may proceed to change the values and save.	None.	CR 3,066	Maintenance TBD
Configure Daily QC	Screening Cells vial numbers are not displayed on the Audit Trail Report when there is a configuration change.	Low risk/Low impact.	Administrator	Immediately print the Audit Trail Report for this activity and manually complete the comment; save it for review.	Explanation for system behavior with no user impact. Reagents are updated by set, not individual vials.	CR 1,949	Report Tool Upgrade Project Proposal NSR 20080915
Configure Daily QC	The comment details do not display completely on the Audit Trail Report.	Low risk/Low impact.	Administrator	Immediately print the Audit Trail Report for this activity and manually complete the comment; save it for review.	The incompletely displayed comment does not impact patient testing or transfusion records and is an infrequent occurrence.	CR 1,950	Report Tool Upgrade Project Proposal NSR 20080915
Configure Daily QC	In a VBECS configured for multidivision use the lot numbers for testing reagents will not pre-fill day to day when sites are using the same reagent rack designations.	No identified risk/No identified impact	All users performing Daily Reagent QC testing in a multidivision configuration of VBECS.	The sites comprising a multidivision can agree to assign different reagent rack names to prevent overlap. This will allow the reagent lot number to pre-fill with the values from the previous day.	None.	CR 2,690 HD 335126	Daily Reagent Quality Control Upgrade, NSR 20081120
Configure Daily QC	Reverse ABO cell testing results are not displayed on the Testing Worklist Report. Various reagent lot numbers are not displayed including QC kit, Reverse ABO cells, PEG, LISS, or Anti-Human Globulin.	Low risk/Low impact.	All users performing Daily Reagent QC testing. Supervisor who configures and reviews QC output for compliance	Manually record daily reagent QC testing and reagent lot numbers to remain compliant with regulatory requirements.	The incompletely displayed QC lot numbers and results force the blood bank to establish a policy to record results to maintain a complete record for accreditation and regulatory compliance which has been put in place by the sites.	CR 2,385 CR 2,436 HD 249548	Report Tool Upgrade Project Proposal NSR 20080915
Configure Daily QC	The OK button on the Configure Daily QC form is enabled even though the user did not change details on the screen.	Low risk/Low impact.	Administrator	None required.	VBECS evaluates that no changes were made and no changes are saved to the database.	CR 1,309	Daily Reagent Quality Control Upgrade NSR 20081120
Configure Daily QC	VBECS Configure QC option will allow a user to configure a reagent rack to use screening cells 1, 2, and 4, instead of 1, 2, and 3 and save.	Low risk/Low impact.	Administrator	Supervisors must verify the configuration to ensure the correct settings are saved.	Evaluation is performed at initial configuration	DR 1,625	Daily Reagent Quality Control Upgrade, NSR 20081120

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Configure Daily QC	VBECS does not permit configuring daily QC to use either Anti-A or Anti-B test with an O Cell	Low risk/Low impact.	Administrator	None required.	Negative control s are available for both antisera. Explanation for system behavior with no user impact.	DR 2,206	Daily Reagent Quality Control Upgrade, NSR 20081120
Configure Division: Order Alerts	Configure Division/Order Alert setting changes are not included in the Audit Trail Report.	Low risk/Low impact.	Administrator	When changing the order alert setting, document the changes by capturing and saving screen shots.	Access to change settings are restricted by security role.	CR 1,637	Report Tool Upgrade Project Proposal NSR 20080915
Cost Accounting Report	VBECS does not remove cost associated with a unit from the report when a unit is inactivated.	No identified risk/No identified impact	All users can print the report; Supervisor is responsible for its review	None required.	None.	DR 3,792 HD 399931	Report Tool Upgrade Project Proposal NSR 20080915
Cost Accounting Report	The modification time displays in GMT (system time) on the Cost Accounting Report.	No identified risk/No identified impact	All users can print the report; Supervisor is responsible for its review	Calculate the local time, as required.	None.	CR 1,633	Report Tool Upgrade Project Proposal NSR 20080915
Cost Accounting Report	A Reflex ABID test that was entered in error displays on the Cost Accounting Report.	No identified risk/No identified impact	All users can print the report; Supervisor is responsible for its review	None available.	The corrected and entered-in-error entries are included. Extra cost is identifiable on report.	CR 1,994	Report Tool Upgrade Project Proposal NSR 20080915
Cost Accounting Report	If a user places a unit on the outgoing shipment invoice, cancels the invoice, and then places the same unit on another outgoing shipment invoice, VBECS displays a return credit twice on the report.	No identified risk/No identified impact	All users can print the report; Supervisor is responsible for its review	Deselect a unit prior to canceling the invoice to avoid the credit appearing twice on the Cost Accounting Report.	Extra credit is identifiable on report.	CR 1,905	Report Tool Upgrade Project Proposal NSR 20080915
Cost Accounting Report	A unit with an edited return credit in final status will not revert to the original return credit amount when the unit's final status is removed.	No identified risk/No identified impact	All users can print the report; Supervisor is responsible for its review	Select Edit Financial Data and correct the return credit amount.	Extra cost is identifiable on report.	CR 1,871	Report Tool Upgrade Project Proposal NSR 20080915
Cost Accounting Report	The tally of discarded units, waste or credit, may include quarantined units and may not accurately reflect the unit discard as waste vs. credit. When the default credit amount is selected, it appears as \$0.00 on the report.	No identified risk/No identified impact	All users can print the report; Supervisor is responsible for its review.	Units listed on the report were quarantined or discarded for the selected date range. Verify the unit and quarantine status of a unit by checking its Unit History Report.	Extra credit is identifiable on report.	CR 2,448	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
CPRS	Orders that expire in VBECS are updated in CPRS but the laboratory order component is not marked as expired.	No identified risk/ Very Low impact..	All users	None required.	None.	CR 3,184 HD 532136	Maintenance TBD
CPRS	Testing comments appear in duplicate on the CPRS Blood Bank Report.	No identified risk/No identified impact	All users	None required.	The testing comments are saved by specimen and by test in VBECS. Both sets of comments are sent to CPRS resulting in the duplication.	CR 3,143 HD 461287	Maintenance TBD
CPRS	CPRS: Orders are not completed in CPRS and Lab when multiple order completion message identifiers are requested within 1/100 th of a second of each other.	Very Low risk/Very Low impact.	All users	Discontinue the order in CPRS and cancel the order in Lab.	The completed test results display as expected in CPRS. The likelihood of occurrence of this is rare and should not interfere with normal business practices. This would be very infrequent as multiple users would need to complete orders on the same patient simultaneously.	CR 2,397	CPRS v 28 coordinated work.
CPRS	CPRS Order details do not include the identity of the VBECS user who rejected the order.	No identified risk/No identified impact	All users	The user that rejected the order can be retrieved from the Order History Report for the patient in question.	None.	CR 2,776 HD 381217	Maintenance TBD
C:T Ratio Report	When a unit is crossmatched more than once to a patient the C:T Ratio Report will count the transfusion of that unit twice.	Low risk/Low impact.	All users reviewing the report	Disregard the second transfusion instance and recalculate the C:T ratio .	This happens when the unit was crossmatched and released for a patient, then crossmatched with a new specimen to that same patient and transfused.	CR 2,748 HD 421286	Report Tool Upgrade Project Proposal NSR 20080915
C:T Ratio Report	The report will show a C:T Ratio of zero (0) any time units are crossmatch and none are transfused on a patient.	Low risk/Low impact.	All users reviewing the report	Any C:T ratio of zero (0) indicates that the user must take the total number of units crossmatched as the true C:T ratio.	None.	CR 2,780	Maintenance TBD
C:T Ratio Report	Unable to schedule report print time in the future.	No identified risk/No identified impact	All users can print the report, Supervisor is responsible for its review.	None required.	Print report. Administrative report can be printed during the user's session.	CR 2,101	Report Tool Upgrade Project Proposal NSR 20080915
Discard or Quarantine	A VBECS system error occurs when a user tries to discard a unit that was brought in through Incoming Shipment with an inactive shipper.	Low risk/Low impact.	All users	Restart VBECS and print unit history for reference. Inactivate the unit and then bring the unit back through Incoming Shipment with an active blood shipper. Recreate the unit history and then discard.	VBECS does not allow a user to save a unit without an active shipper unless the VistALink connection is lost while the user is trying to process a unit that does not have a shipper defined.	CR 3,048	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Discard or Quarantine	Cannot edit hour: minutes in Discard Date field.	No identified risk/No identified impact	All users	None required.	This option is designed for processing in real time, not retrospective entry. Default is current day/time.	CR 2,074	Maintenance TBD
Discard or Quarantine	Cannot designate a Discard/Quarantine Tech ID though the exception report contains both "Test By" and "Override By" columns.	No identified risk/No identified impact	All users	None required.	Discard/Quarantine option does not accommodate the entry of a second user for retrospective entry.	CR 2,076	Maintenance TBD
Display Order Alerts	The tool tip on the Patient Alert icon will not display after VBECS refreshes and no new updates are detected.	No identified risk/No identified impact	Administrator	None required.	User has no new alerts to view so icon function does not prompt the tool tip.	CR 2,168	Maintenance TBD
Division Workload Report	The Division Workload Report will not print as "Preliminary" when the report contains the current day in the range of data requested.	No identified risk/ No identified impact	All users	None available.	Users are cautioned to pay attention to the date range of this report. Workload is an administrative report that is not required the day the data is created.	CR 2,148	Report Tool Upgrade Project Proposal NSR 20080915
Division Transfusion Report	The report will only print the Division Transfusion Report by physician for the first 50 physician's in the VBECS database.	No identified risk/ No identified impact	All users	Do not run the Division Transfusion Report by all physicians. Selected physician report can be run but should be limited to no more than 50 physicians at a time.	The report search in VBECS is limited to 2000 characters. Requesting a report that exceeds this limit will result in an incomplete report.	CR 3,170	Report Tool Upgrade Project Proposal NSR 20080915
Document ABO Incompatible Transfusions	When transfusion interruption is indicated the OK button enables allowing the user to save prior to changing the transfused amount.	Low risk/Low impact.	Supervisor	Record the amount transfused before saving the transfusion event.	None.	CR 3,062	Maintenance TBD
Edit Financial Data	When a unit is "transferred," VBECS does not enable the Return Credit field in the Edit Unit Financial Data window.	Low risk/Low impact.	All users	Process the unit through Remove Final Status, adjust the return credit amount, and reprocess the outgoing shipment information.	Return credit is set up as a default value so this would only be used for a rare return outside of the vendor contract amount.	CR 1,777	Maintenance TBD
Edit Financial Data	When a Technologist enters a unit ID and product code for a unit that is not in final status, VBECS does not enable the Special Test Cost field.	No identified risk/ No identified impact	Technologist only	An Enhanced Technologist or higher access can change the Special Test Cost field in the Edit Financial Data option.	None.	CR 1,819	Maintenance TBD
Edit Financial Data	Entry of a hyphen (-) or a comma (,) in any of the three editable numeric fields generates an error message: "Value entered is not a valid number."	No identified risk/No identified impact	All users	None required.	None.	CR 1,890	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Edit Unit Information	When Inactivating a unit the incorrect list of canned comments is presented. The list of comments shown corresponds to the canned comment category "Unit Status Removal" instead of "Unit Inactivation".	Low risk/Low impact.	Supervisor	Enter comments selectable during "unit inactivation" in the Canned Comment Category of "Unit Status Removal".	Canned comments are site configurable to include desired comments. "Other" and a free text comment is also available.	CR 2,270	Maintenance TBD
Edit Unit Information	A split unit and associated units may be restricted to and transfused to different patients.	No identified risk/No identified impact	Supervisor reviewing transfusion records	None required.	This may be the desired functionality.	DR 2,218	Maintenance TBD
Edit Unit Information	In the unit volume field, numeric values outside the allowed range can be entered. VBECS changes the entry to an acceptable value without a warning message.	No identified risk/No identified impact	Supervisor	Check the entry for validity prior to confirming the save. Edit as required.	None.	CR 1,917	Maintenance TBD
Edit Unit Information	When two ABO/Rh confirmation tests are present and both require invalidation, VBECS does not allow both tests to be invalidated in the same transaction.	No identified risk/No identified impact	All users	Perform one invalidation action at a time. Invalidate one test, leave the option, and go back to invalidate the second test.	This may be the desired functionality.	CR 1,906	Maintenance TBD
Edit Unit Information	Inactivating the parts of a split units will result in negative workload equal to the number of splits created applied to the total number of units logged in.	No identified risk/Low impact	Supervisor	None available.	Infrequent user action that only affects laboratory workload reports.	CR 2,224	Maintenance TBD
Edit Unit Information	The "Biohazardous?" checkbox can be cleared for Donation Type: For Autologous Use Only, Biohazardous, but the donation type of the unit is not editable. If the unit is re-edited, the checkbox displays as re-checked and disabled. The Unit History Report indicates this field is cleared.	Low risk/Low impact.	All users	When an autologous unit is considered biohazardous because the testing results are incomplete and expected, select the donation type "For Autologous Use Only" and select incomplete disease marker testing which sets the biohazardous indicator. The biohazardous indicator and disease status are editable when testing results are received.	Select the donation type "For Autologous Use Only, Biohazardous" only when the unit testing is not pending	CR 2,262	Maintenance TBD
Edit Unit Information	When attempting to edit a unit in a final status, VBECS displays the wrong message "Blood Unit not found for Unit Id/Product Code" .	No identified risk/ No identified impact.	All users	None required. VBECS does not allow edit of the unit, which is correct.	None.	CR 2,310	Maintenance TBD
Edit Unit Information	The OK button is enabled unnecessarily.	No identified risk/ No identified impact.	All users	None required.	VBECS evaluates that no changes were made and no changes are saved to the database.	CR 2,359	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Enter Daily QC Results	A VBECS system error occurs when a user maximizes the reagent rack selection window while performing QC on the rack.	No identified risk/ Low impact.	All users	Do not maximize the reagent rack window.	User may need to reenter unsaved QC data. Full screen display is not required to view the test grid. The window is set to display in a size that does not lead to this problem. The user must maximize the window purposefully to create this error. Generally this happens to one user during validation and is mitigated by training.	CR 1,911 HD 367084	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	The secondary reagent is Check Cells (CC) for the Daily QC test CC. Changing the AHG, PS lot number does not enable the testing grid cell.	No identified risk/No identified impact.	All users	None required.	None.	CR 3,188	Maintenance TBD
Enter Daily QC Results	Reagents that are not associated with test results when the QC is partially saved are marked as satisfactory on the Testing Worklist Report.	Low risk/Low impact.	All users	Complete all QC testing in one instance saving once. Do not perform partially completed QC testing.	The Testing Worklist Report displays the reagents as unsatisfactory through the indicator appears as satisfactory in QC testing.	CR 3,074	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	Changing the tested with lot number (ex. LISS) does not enable the Enter Daily Reagent QC test grid.	Low risk/Low impact.	All users performing Daily Reagent QC testing.	Change the lot number of the primary reagent as well as the secondary reagent to allow testing of both.	Refer to FAQ: Retesting QC for a listing of primary and secondary reagents.	CR 3,176	Maintenance TBD
Enter Daily QC Results	The user cannot save a partially filled worksheet due to inactivity timeout.	No identified risk/ No identified impact.	All users	The user must maintain activity on their screen/session to prevent the activity timeout, and then the lock on the worksheet will not expire.	Unsaved data must be reentered if the locally configured timeout is exceeded.	CR 2,131	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	Rack Daily QC cannot be saved if one of the QC'd racks was partially QC'd.	No identified risk/ No identified impact.	All users	Perform QC for only one rack at a time or for multiple racks to segregate non-QC'd racks from partially completed racks.	This is a rare issue since normal workflow is to complete all QC testing together or to have each user perform own QC.	CR 2,145	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	If a partially tested rack is designated as not tested for a given day, VBECS displays the reagent rack on both the Partially Tested and Not Tested tabs.	No identified risk/ No identified impact.	All users	None required.	None.	CR 1,761 CR 1,918	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	Twenty-four-hour expiration messages are displayed more frequently than required.	No identified risk/ No identified impact.	All users	Acknowledge the expired reagent-warning message and continue.	Explanation of VBECS behavior. User requirements need further definition.	CR 1,768	Daily Reagent Quality Control Upgrade, NSR 20081120

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Enter Daily QC Results	When a user clicks No on the decision box to not use a reagent/antiserum that is within 24 hours of expiration (23:59 on the expiration date), VBECS moves the cursor to the next reagent lot number field and <i>does not</i> remove the lot number of the previous cell.	No identified risk/ No identified impact.	All users	Change the lot number of the reagent or continue using it until the actual expiration date and time.	Explanation of VBECS behavior. User requirements need further definition.	CR 1,861	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	VBECS antiglobulin testing grids for QC have a title of IAT. The same grid for patient testing is titled AHG.	No identified risk/ No identified impact.	All users	None required.	IAT and AHG are synonyms for the antiglobulin testing and is not a patient safety concern.	DR 1,230	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	PS AHG always appears on the lot number page though it is may not be used.	No identified risk/Low impact.	All users in full service divisions.	Enter the Lot number of the specific AHG used for antibody screen QC. Note in the procedure that the lot number is accurate for the reagent used in QC and testing.	Explanation of VBECS behavior.	DR 2,679	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Daily QC Results	Line items may display in a different order on a display and report after a change to the configuration.	No identified risk/ Low impact.	All users in full service divisions.	None required.	Explanation of VBECS behavior. Re organization of the same line items.	CR 2,437 HD 269544	Daily Reagent Quality Control Upgrade, NSR 20081120
Enter Reflex Test Results	The facilities pick list is not restricted to those marked as testing facilities.	No identified risk/ No identified impact.	All users	None required.	Explanation of VBECS behavior.	DR 1,651	Maintenance TBD
Equipment	The equipment type does not print on the report; column information may be truncated; if the manufacturer city is not entered, the state and zip code do not print.	No identified risk/ No identified impact.	All users when equipment maintenance records are defined in VBECS.	None required.	Report data may be reviewed on screen. This is optional data entry for a site. They may have an existing equipment process that they prefer and do not use the VBECS option to record periodic maintenance.	CR 2,036	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	Thawed plasma product types are displayed as fresh frozen plasma on the Exception Report.	No identified risk/ No identified impact.	All users	Users can view the Unit History Report for the correct product type name.	The Unit History Report and CPRS both show the product code name which is accurately represented.	DR 4,008	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	An Exception Report entry is not created when a unit is resulted as weak D positive or weak D inconclusive and is properly quarantined.	No identified risk/ No identified impact.	All users	None required as the unit is quarantined appropriately and is not available for patient selection.	None.	DR 3,832	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Exception Report	Duplicated exception entries (instead of only one) will appear on the Exception Report when a unit is issued that did not meet the antigen negative requirements of the patient at issue.	No identified risk/ No identified impact.	All users	None required.	The number of duplications that appear correlate to the number of the patient's antigen negative requirements.	CR 2,924	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	"Visual Inspection Information" and the selected processing "User Information" are not included in Exception Report entry.	No identified risk/ No identified impact.	All users	See the Unit History report, Issue/Return section for the "Visual Inspection Information" as an entry detail for Transaction Type: Unit unsatisfactory upon return from issue.	None.	DR 2,479	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	The exception type "Previously recorded results invalidated" displays an incorrect time for the "Date/time results invalidated."	No identified risk/ Low impact.	All users	None required.	The time displayed in the "Date/time results invalidated" field is the time the results were entered.	CR 2,686	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	When invalidating results from Unit Antigen Typing the Exception Report does not show all invalidated values.	No identified risk/ No identified impact.	All users	View the invalidated tests on the Testing Worklist Report and the Unit History Report.	None.	CR 2,022	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	Exception Report does not work with the report scheduler.	No identified risk/ No identified impact.	All Users	Users must run the report when needed.	This report is printed regularly and reviewed for daily work so system performance is not impacted.	CR 2,257	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	The "Expired reagent QC'd" exception type section does not include the rack identification or the phase.	No identified risk/ No identified impact.	All users	The lot numbers and testing phases associated with the DAT testing are included in the Testing Worklist Report in the Miscellaneous testing and QC sections, respectively, of the report. The rack identifier and/or phase may be manually added to the Exception Report prior to or during the review.	None.	CR 1,626 CR 1,636	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	Not all exception types use the date and time of the save as the date and time of the exception.	No identified risk/ No identified impact.	All users	None available.	None.	CR 2,055	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	The previously recorded Results Inactivated exception type is not generated when a crossmatch is inactivated by Invalidate Test Results or by using the red X in the grid.	No identified risk/ No identified impact.	All users	None available.	The Testing Worklist Report contains the details of the invalidated testing and is recommended for daily supervisory review with the Exception Report. Infrequent event that is performed by a user with higher security role if the blood product was issued.	CR 2,035	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Exception Report	The User ID is displayed in the Tested By column (instead of User Name) for exception type: Expired Task Processed.	No identified risk/ No identified impact.	All users	None required.	None.	CR 2,106	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	Modification exceptions do not display full (eight-digit) product codes for ISBT 128 labeled units. The sixth, seventh, and eighth digits are not included.	No identified risk/ No identified impact.	All users	None available.	Obtain the full ISBT 128 product code from the Unit History Report, as necessary. Format is understandable to users.	CR 1,995	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	If patient ABO/Rh results are not entered in the order of performance (current testing is entered before the retrospective data entry), the Exception Report entries for an ABO/Rh interpretation discrepancy are displayed based on the time the data are entered (the OK button is clicked).	No identified risk/ No identified impact.	All users	None available.	The correct blood type system rules are applied. The report columns state Current/Previous, the individual exceptions include the date/time of entry allows the sequencing of testing.	CR 2,048	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	Exception Report and Blood Unit History Report display the testing status code from the database in the "Testing Status" section for autologous units.	No identified risk/ No identified impact.	All users	None required.	The values displayed in "Testing Status" translate to: 1 - NEG - Unit fully tested and negative for all disease markers. 2 - POS - Unit tested and positive for one or more disease markers, Biohazard. 3 - NFT - Unit not fully tested for one or more disease markers, Biohazard.	CR 2,118	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	The user ID of the issuing user appears on the report instead of the user name.	No identified risk/ No identified impact.	All users	None available.	None.	CR 2,201	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	VBECS Exception Report does not display a patient's first name for Unacceptable/Expired Specimen Used exceptions.	No identified risk/ No identified impact	All users	None required.	Patient ID and specimen UID are presented in full with patient last name allowing further investigation of the exception as required.	CR 2,309	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	The Exception Report section "Antigen Testing Phase Change" does not display test interpretations.	No identified risk/ No identified impact.	Supervisor responsible for report review	None required.	The Testing Worklist Report is identified as one of the reports requiring supervisor review. The supervisor reviews the test results as part of their investigation.	CR 2,642	Report Tool Upgrade Project Proposal NSR 20080915
Exception Report	VBECS displays ISBT unit ID and product short name, but not the product code.	No identified risk/ No identified impact.	All users	None available.	See the Unit History report for the missing data saved in other options.	DR 2,323	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Finalize/ Print TRW	Barcode scanning of a specimen UID displays a tilde (~) in front of UID in the pre and post transfusion specimen fields.	No identified risk/ Low impact.	All users	Do not use the barcode scanner to enter a transfusion reaction specimen ID.	None.	CR 3,200	Maintenance TBD
Finalize/ Print TRW	A system error occurs when a user enters serologic testing results for a transfusion reaction and then unchecks the testing box, checks that box again and clicks OK to save.	Low risk/Low impact.	All users	None available.	Once testing is entered there should not be a need to uncheck the testing box and recheck in normal practice.	CR 2,804 HD 381573	Maintenance TBD
Finalize/ Print TRW	Changes to previously selected canned comments associated with an implicated are not saved when a user changes them in a Transfusion Reaction Workup.	Low risk/Low impact.	All users	None available.	None.	CR 2,867	Maintenance TBD
Finalize/ Print TRW	A system error occurs when a user attempts to finalize a TRW with a canned comment and >350 characters of additional text in the comment field.	Low risk/Low impact.	Lead Technologist and above	Limit finalize TRW comments to <350 characters to allow save.	The comment is saved when the system error is generated but VBECS will not be able to print the TRW.	CR 2,852	Maintenance TBD
Finalize/ Print TRW	A system error occurs if a user has a TRW open and then attempts to open a second instance of the same TRW.	Low risk/Low impact.	All users	Restart VBECS session and open one TRW to result and finalize.	None.	CR 2,815	Maintenance TBD
Finalize/ Print TRW	VBECS only displays 350 characters entered into the Transfusion Reaction Details field on the Finalized Transfusion Reaction Report.	No identified risk/ No identified impact.	Supervisor	The details text is available to ~1000 characters on the Transfusion Reaction Count Report (Detailed).	Report is printed for medical director signature and charting. Printed report can be updated manually if required or by using a Vista consult.	CR 2,230	Maintenance TBD
Finalize/ Print TRW	A finalized Transfusion Reaction Workup cannot be corrected once it has been finalized.	Low risk/Low impact.	Supervisor	Request a correction of a finalized Transfusion Reaction Report in the VBECS database by creating a Remedy ticket.	Users are cautioned to double check the workup before finalizing. Report is printed for medical director signature and charting. Printed report can be updated manually if required.	DR 1,633 HD 332336	Maintenance TBD
Finalize/ Print TRW	When finalizing a TRW report and VistaLink is down, a pop up message displays that VistaLink is down and the DSS message cannot be sent. VBECS should resend the DSS message later. However, the message is not resent.	Low risk/Low impact.	Supervisor	None available.	This does not impact patient test reporting. TRW finalization and VistaLink downtime occur infrequently, independently, and rarely concurrently.	CR 2,455	Maintenance TBD
Free Directed Unit For Crossover	VBECS displays the logged on user performing the Free Directed Units for Crossover process even if he selects a different user name in the Removed By field.	Low risk/Low impact.	Enhanced Tech	Do not use Free Directed Units for Crossover during downtime.	When processed in downtime, the selected user is recorded on the downtime form.	CR 1,853	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Free Directed Unit For Cross-over	A duplicate unit record may be created when an ISBT 128 unit is entered with a donation type of "D" and changed to a donation type of "V" during Free Directed Unit option and the unit was also entered during Incoming Shipment with the donation type "V".	Low risk/Low impact.	All users	When a unit is received from the blood supplier with a donation type of "V" and the unit is to be restricted to a patient as a directed donation, the product code must be manually entered with a "D" to allow the restriction in a Full Service Blood Bank. A Transfusion Only facility type will not encounter this problem as the option Free Directed Unit is not enabled.	If the user tried to bring the unit into inventory using both methods, the incorrect unit can be invalidated to prevent the duplicate record. The unit label would reflect the correct unit status per local policy.	CR 2,261	Maintenance TBD
Free Directed Unit For Cross-over	The "Unit is in a final status" VBECS message is missing the "a."	Low risk/Low impact.	Enhanced Tech	None required.	None.	CR 1,448	Maintenance TBD
Incoming Shipment	ISBT expiration date barcodes from the Department of Defense are not scan able due to different embedded prefix characters to indicate date only versus date and time.	Low risk/Low impact.	All users	Manually type in the expiration date and time of the blood unit per the eye readable expiration date information.	This is generally associated with blood units that have a standard 2359 expiration time and the unit may be optionally configured to include time or not.	DR 3,610 HD 360751	Maintenance TBD
Incoming Shipment	VBECS will not accept entry of ISBT units when the first letter of the donor identification number is "ABDIU".	Low risk/Low impact.	All users	Units with these letters cannot be brought into VBECS and must be sent back to your supplier.	The letter generally indicates the unit's country of origin, or a non-collection facility. The U.S. blood collection facilities all use "W".	CR 3,061	Maintenance TBD
Incoming Shipment	VBECS allows a user to change the time in the Date Received field to a future time.	Low risk/Low impact.	All users	VBECS defaults to the current date and time when the incoming shipment screen is opened. Instruct user to not change the time in the Date Received field to a future time.	The invoice date/time is used to determine an entered unit's appropriate maximum storage time.	CR 2,933	Maintenance TBD
Incoming Shipment	A system error occurs when using the restrict to a patient search field and only 5 to 8 characters of the patient's SSN are entered to initiate the search for the patient.	Low risk/Low impact.	All users	Restart VBECS session and log in the unit again and enter the full patient SSN to search for the patient to restrict the unit.	The system error will occur after clicking OK to save the unit. The error is caused by the incomplete patient search conducted to restrict the unit.	CR 2,809	Maintenance TBD
Incoming Shipment	VBECS does not allow the entry of a 10 digit FDA number in Incoming Shipment when the user is prompted to activate the facility.	No identified risk/ Low impact.	All users	Close the Incoming Shipment window and open Tool, Local Facilities to create a new supplier with a 10 digit FDA number.	The entry of a 10 digit FDA number is permitted but must be entered through the Local Facilities option and not directly through Incoming Shipment.	CR 3,131 HD 515950	Maintenance TBD
Incoming Shipment	When a unit is received for the second time from a blood supplier a user can remove the ABO/Rh of the unit and VBECS permits saving without a blood type making the unit unavailable.	Low risk/Low Impact	All Users	Do not remove the blood type when entering the unit back into inventory and review the data entered into VBECS against the unit labels before saving.	If a unit is saved without a blood type inactivate the unit and re-enter it in Incoming Shipment. If a unit is saved without a blood type it cannot be properly confirmed nor selected.	CR 3,191 HD 538535	Maintenance TBD
Incoming Shipment	When the user holds the mouse over the disabled OK button; VBECS displays the tool tip associated with the OK, instead of the tool tip that prevents the activation.	No identified risk/ No identified impact.	All users	Check the tool tips for the active field and OK button. If there is a conflict, comply with the field error message.	None.	CR 1,721	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Incoming Shipment	When a user enters a valid date in the Expiration Date field and tabs out of the field, the save button is enabled. If the user returns to the expiration date, clicks the Delete or Backspace key to delete it, and tabs out of the field, the button remains enabled and the previously entered date and time are saved.	Low risk/Low impact.	All users	Inactivate and reenter the unit to correct the entry.	User is allowed to correct the date before save which would be the usual reason for revisiting the field. The date is presented for review before the user finally accepts.	CR 1,605	Maintenance TBD
Incoming Shipment	The Unit History Report will not show a unit that has been in inventory as CMV negative, shipped out, and then received again through Incoming Shipment as CMV negative.	No identified risk/ Low impact.	All users	None available.	VBECS still holds the CMV negative indicator and treats the unit as CMV negative. User is expected to reenter the CMV test during unit receipt.	CR 1,932	Report Tool Upgrade Project Proposal NSR 20080915
Incoming Shipment	The Unit History report displays results of Antigen Typing on a unit when the additional daily QC rows (POS and NEG) are processed the testing row is displayed three times.	No identified risk/ No identified impact.	All users	Select the one labeled for the Antisera as the result, or look at the Testing Worklist Report.	QC results are also displayed but easily distinguishable from unit testing.	CR 2,181	Report Tool Upgrade Project Proposal NSR 20080915
Incoming Shipment	When entering a response of "YES" to the message that you are attempting to re-enter a unit, a pop up appears "System Error: Unit already exists in the division." Click OK and the expected message about "Are you sure you want to add this unit?" displays and allows the user to proceed normally.	No identified risk/ Low impact	All users where it is not common practice to move blood products in and out of inventory.	None required.	This may occur when a unit is returned to the supplier and then back to the facility for use. User is allowed to proceed with the unit entry. If this is common practice related to inventory maintenance, the user would not have any impact.	CR 3,160	Maintenance TBD
Invalidate Test Results	When a user is prompted to save a comment when invalidating a test result, they can select a comment and then hit the spacebar before saving. Tapping the spacebar removes the comment and still allows the user to save the blank comment. The lack of proper comment text associated with the invalidated test then prevents direct access to the SI and TR option for this patient.	Low risk/Low impact.	All users	Review the comment before saving to confirm the comment is correct.	None.	CR 3,240 HD 720879	Maintenance TBD
Invalidate Test Results	A VBECS system error occurs when a user attempts to invalidate two crossmatch results and the blood unit was assigned, crossmatched, released, reassigned, and crossmatched again on the same specimen.	No identified risk/Low impact.	All users	Invalidate the first crossmatch save and exit the invalidate Patient Results option. Reenter the Invalidate Patient Results option to invalidate the second crossmatch.	Requires security role to complete and would be performed administratively.	CR 1,925	Testing Grids

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Invalidate Test Results	When an antigen typing is invalidated; VBECS selects the Pending Task List check box. If the user does not clear the box, the test is automatically put back on the Pending Task List.	No identified risk/ No identified impact.	All users	Cancel the test on the Pending Task List, as necessary.	The user would be prompted to return to option by the test reappearance.	CR 1,596	Testing Grids
Invalidate Test Results	VBECS does not display the implicated unit product code for a transfusion reaction workup in the Invalidate Patient Test Results option.	No identified risk/ No identified impact.	All users	None required.	The option displays the product short name, which is materially equivalent to the product code.	CR 2,026	Testing Grids
Invalidate Test Results	"The patient has history of justified ABO/Rh change" warning message contains an extra space.	No identified risk/ No identified impact.	All users	None required.	None.	CR 1,597	Testing Grids
Invalidate Test Results	A system error occurs when attempting to invalidate a crossmatch test after the unit has been marked unsatisfactory for issue in Issue Unit.	No identified risk/No identified impact.	All users	If the unit is unsatisfactory for issue because the crossmatch test is incorrect, do not mark the unit unsatisfactory, release the unit from assignment and invalidate crossmatch test.	The steps follow user expectations but the error message should inform the user to use another process. This is an administrative process with security role.	CR 3,190	Testing Grids
Invalidate Test Results	The expected confirmation message does not display to the user informing them that the test was added back to the Pending Task List.	No identified risk/No identified impact.	All users	None required.	None. The inactivated test is added to the Pending Task List as indicated by the user.	CR 2,579	Testing Grids
Issue Blood Components	The area to mark a unit unsatisfactory for issue extends beyond the checkbox. A user that clicks in that area but not the checkbox could inadvertently mark a unit unsatisfactory for issue.	Low risk/Low impact.	All users	A VBECS Confirmation message window appears to verify the user wants to remove patient association and quarantine the unit. The user must click No if they mistakenly indicated the unit was unsatisfactory for issue.	None.	CR 3,208 HD 590527	Maintenance TBD
Issue Blood Components	A system error occurs when a user attempts to issue a unit assigned or crossmatched to a patient with an expired specimen.	No identified risk/Low impact.	All users except Supervisor who is allowed to proceed with override.	Go to Orders, Maintain Specimen, change the specimen expiration date so it is not expired and issue the unit. This may require a temporary change of the Maximum Specimen Expiration Date in the Configure Division option.	The steps follow user expectations but the error message should inform the user to use another process. This is an administrative process with security role.	CR 2,092	Select Unit Project Proposal
Issue Blood Components	When the user clicks Cancel, VBECS closes without displaying a confirmation message that data are lost.	No identified risk/Low impact.	All users	If canceled in error, recreate the data.	The steps follow user expectations to reduce the number of user responses.	CR 1,598	Select Unit Project Proposal
Issue Blood Components	Plasma can be issued frozen without thawing when associated with the order type is "OTHER."	No identified risk/Low impact.	All users	"OTHER" products require special consideration by local policy.	Most Plasma product codes are processed in an FFP order not in the OTHER orderable. "OTHER" products are handled with local policies as described in the User Guide section.	CR 2,254	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Issue Blood Components	Missing a space between the sentences: "This unit had a crossmatch result of incompatible or the crossmatch was not completed. This unit cannot be issued."	No identified risk/ No identified impact.	None	None required.	Format is understandable to users.	CR 2,085	Maintenance TBD
Issue Blood Components	Due to an unjustified discrepant ABO, the assigned O negative units appear on the Emergency Issue tab, not the assigned tab, without an information message to explain why.	No identified risk/ No identified impact.	None	None required.	The unit may be issued. The unit may meet crossmatch requirements but VBECS still invokes Emergency Issues rules since other testing criteria have not been met.	CR 1,794	Select Unit Project Proposal
Issue Blood Components	When retrospective data is entered into VBECS after a more recent specimen test result, the blood type of the retrospective entry is compared to the most recent test result.	Low risk/Low impact.	All users when entering retrospective testing out of chronologic order of execution.	None required.	Enter retrospective testing in chronological order so the current specimen is the last test entry for the patient. Any testing discrepancy requires supervisor review and action.	CR 2,011	Select Unit Project Proposal
Issue Blood Components	VBECS may display an inactivated unit in the unit search screen without any indication of its inactivated record status.	No identified risk/ No identified impact.	All users when inactivated unit records are available in the division.	None required.	The best practice for issue is to scan the unit information. The unit relocation cannot proceed if an inactivated unit is selected off the list. Users are directed by User Guide and local policy to use the unit for selection.	CR 1,877	Select Unit Project Proposal
Issue Blood Components	When VBECS opens the Issue Blood Components window the focus (cursor) is not in the unit ID field.	No identified risk/ No identified impact.	All users	A user must first click in the Unit ID field to bring the focus to that field before scanning or entering a unit number.	No errors are expected to result.	CR 2,129	Select Unit Project Proposal
Issue Blood Components	When VBECS warns that a unit is assigned to another patient; the user is prompted to process an override and enter a comment to proceed.	Low risk/Low impact. The comment can be viewed with the historical record.	All users	Complete the override to proceed.	The entered comment is not recorded or displayed on the Exception Report.	CR 1,798	Select Unit Project Proposal
Issue Blood Components	The message displayed reads "Original and repeat ABO/Rh interpretation do not match." The 2nd sentence in the designed message "You must resolve the discrepancy before units can be issued." is not displayed.	Low risk/Low impact.	All users	None required.	VBECS will not allow the user to proceed with selecting a unit until the discrepancy is resolved thus mitigating the hazard.	CR 2,203	Select Unit Project Proposal
Issue Blood Components	The workload recorded for Issue Units is doubled on the VBECS Workload Report.	Low risk/Low impact.	All users	None required as the VistA workload reported is accurate and correct.	Workload validation clarification.	CR 2,242	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Issue Blood Components	Surgical initiative status is sent as assigned not compatible or issued.	Low risk/Low impact.	All users	None required.	The unit is available for bedside verification to the surgical initiative only after the unit has been physically issued from the blood bank regardless of the unit status in VBECS. Unit is correctly handled between VBECS and Surgery so only units that are issued by VBECS are transfused.	CR 2,333	Maintenance TBD
Issue Blood Components	Blood is issued to the VistA Hospital location associated with the VBECS division, not a VistA hospital location associated with a mapped associated VistA Institution hospital location in VBECS Administrator.	Low risk/Low impact.	All users	None required.	Optionally, use the remote storage location to add details of exactly where the products for transfusion were delivered to at the remote site. Clarification of design to allow more details as to transfusion locations.	DR 2,881	Select Unit Project Proposal
Issued/Returned Units Report	The report returns more information than is requested. The Issued/Returned Report will include the transactions from the day prior to the start date requested by the user.	No identified risk/ No identified impact.	All users	None required.	The report provides an additional day of information.	CR 2,736	Report Tool Upgrade Project Proposal NSR 20080915
Issued/Returned Units Report	The date and time of unit return is not included on this report if it is different from the time the returned unit is selected (retrospective entry) at the time of data entry. The report does not include issue and return comments, as expected.	Low risk/Low impact.	All users	None available.	The report does not include the selected date/ time or processing date/ time fields. This information is displayed on the Unit History Report in the Issue Information section. Comment fields are not available at the time of issue or return, so none are available for the report.	CR 2,044	Report Tool Upgrade Project Proposal NSR 20080915
Issued/Returned Units Report	If a unit is entered into Incoming Shipment as one product code and then issued, and after issue modified to a new product code, the Issued/Returned Units Report will not show the original product code. Only the new modified product code displays.	Low risk/Low impact.	All users	A user can determine the original product by viewing the Unit History Report.	This is performed by a user with proper security role.	CR 2,185	Report Tool Upgrade Project Proposal NSR 20080915
Issued/Returned Units Report	If a unit is assigned to a patient, released from assignment, and then assigned to any patient, duplicate issue records will display for both patients on the Issued/Returned Units Report.	No identified risk/ No identified impact.	All users	None required.	The report displays the correct issue information with a time and date stamp. A user can easily determine the last issue event.	CR 2,596 HD 443058	Report Tool Upgrade Project Proposal NSR 20080915
Issued/Returned Units Report	Report title does not print when requested using the report scheduler.	No identified risk/ No identified impact.	All users	Do not use the report scheduler to print the Issued/Returned Report.	The report title prints when the report is printed on demand rather than queued.	CR 2,512	Report Tool Upgrade Project Proposal NSR 20080915
Issued/Returned Units Report	The Issue/Returned Report does not display the processing user and time of entry when results are entered retrospectively.	Low risk/Low impact.	All users	None required.	Retrospective data entry has supporting downtime documentation to capture the user information.	DR 2,072	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Justify ABO/Rh Change	When a database conversion blood type is one of the two blood type results involved in a justification, the database conversion blood type does not display on the Audit Trail Report entry of the justification.	Low risk/Low impact.	Administrator	None required.	The information is available on the Patient Record Report. Details are readily available to the reviewer.	CR 1,913	Report Tool Upgrade Project Proposal NSR 20080915
Justify ABO/Rh Change	When initially displayed, the VBECS patient name field in the Justify ABO/Rh Change option may not display the expanded field.	Low risk/Low impact.	Administrator	Maximize the window to display the expanded (full name) field for long patient names.	Details are readily available to the user.	CR 1,914	Maintenance TBD
Log Into VBECS and VistA	A cancelled VistA Logon - Authorization window at initial sign-on reopens automatically several times before staying closed.	Low risk/Low impact.	All users	None available.	None.	CR 1,837	Maintenance TBD
Log Into VBECS and VistA	System error occurs when a user tries to log in and their role has been inactivated in that division instead of presenting the message "Your role within the division <Name> was inactivated. Please contact your system administrator."	Low risk/Low impact.	All users	None required.	The user is not able to access the division.	CR 2,361	Maintenance TBD
Maintain Minimum Levels	The minimum stock level for the reagents field allows the entry of a decimal that causes the reversal of the entry (e.g., user entry of "1.5" becomes "51") and may be saved.	Low risk/Low impact.	Administrator	Enter whole numbers; do not enter decimals. Check the accuracy of the entry before saving.	There is no location to allow decimal entries but may occur through misuse of the option. Entries are not tied to patient testing and are used for reagent maintenance.	CR 1,686	Daily Reagent Quality Control Upgrade, NSR 20081120
Maintain Minimum Levels	Multiple shipments of the same lot number of a reagent type can cause a discrepancy in the total number available displayed on the Reagent Inventory Report and what is in inventory.	No identified risk/ No identified impact	All users	Refer to the most recent entry of a reagent for the number of vials available when the shipment was received.	VBECS can only display the data entered and does not manage inventory.	CR 2,599 HD 301562	Daily Reagent Quality Control Upgrade, NSR 20081120
Maintain Minimum Levels	Update reagents shows different information for a lot number than the Reagent Inventory report. The Reagent Inventory report displays each individual entry but the update reagent lot number view displays only the last entry for the lot number.	Low risk/Low impact.	Administrator	Enter with * at the beginning and end of the lot number. User may also re-enter lot number if changing the invoice number. See the Reagent Inventory Report.	The report displays all changes in full.	CR 2,209	Daily Reagent Quality Control Upgrade, NSR 20081120
Maintain Minimum Levels	In a multidivisional environment, the reagent report for a division will print one line item for each division that has set minimum levels for that reagent type. For example, if you have minimum levels for Reagent A, and two other divisions have minimum levels for Reagent A, that reagent information will print three times on your report.	No identified risk/ No identified impact	Administrator	There is no crossover of data, and no safety issue. If minimum levels are not set in a division, then it will not factor into the display. In a single-division environment, there is no such problem.	The report accurately displays the inventory at each division.	CR 2,233	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Maintain Specimen	When extending a specimen expiration date the override message that the patient has been transfused in the past 3 months appears even though the patient was transfused more than 3 months prior to the specimen expiration date.	No Risk/ Low Impact	All Users when extending a specimen expiration date.	None required. Complete the override as requested to proceed per local policy and procedure.	None.	CR 2,802 HD 394421	Maintenance TBD
Maintain Specimen	A VBECS system error occurs when a user allows two consecutive system timeouts.	Low risk/Low impact.	All users	Once in the option, perform specimen processing immediately.	This occurs when the user repeatedly enters the option and does not complete the window.	CR 875	Testing Grids Project Proposal
Maintain Specimen	For a component order that does not require a specimen, VBECS calculates the order's expiration date and time to the "minute" from the collection time, not the appropriate day with an expiration time of 23:59.	Low risk/Low impact.	All users	Request a new component order, as needed.	VBECS expires the order at the time of day 10 days after the order was accepted. The order is simply expired hours earlier on the 10 th day.	CR 1,930	Testing Grids Project Proposal
Maintain Specimen	When a specimen expiration is more than 72 hours in the future and a blood product is issued, the specimen expiration is recalculated to 72 hours from the time of issue.	Low risk/Low impact.	All users	None required.	The expiration date is recalculated at blood product issue to expire hours earlier on the 3 rd day when it would normally expire.	CR 2,485 HD 382525	Testing Grids Project Proposal
Medication Profile	A Crystal Reports Forms Viewer Error may appear: "Query Engine Error."	Low risk/Low impact.	All users	None available. Reenter the report search.	Use the patient tool bar link to the medication profile report to avoid this error.	CR 1,560	Report Tool Upgrade Project Proposal NSR 20080915
Medication Profile	VBECS does not display the "No medications were found for the date range" message when the report does not return data.	Low risk/Low impact.	All users	None required.	A blank report is displayed showing the date range and "End of Report." Business processes and patient safety are not impacted; display of the blank report is inconsistent with report messaging. This is a presentation inconsistency.	CR 1,577	Report Tool Upgrade Project Proposal NSR 20080915
Modify Units (not Pool or Split)	The Volume Reduce Unit form contains fields for discarded plasma volume and original unit volume. When either field contains data entry errors, the OK button is enabled, which allows the user to save the form.	Low risk/Low impact.	All users	None available. Check the accuracy of the entry before saving.	The volumes are saved as entered, The volume total may not match the original default volume.	CR 1,924	Maintenance TBD
Modify Units (not Pool or Split)	When a user enters a value outside the allowable range in the Target Unit Volume field; VBECS accepts the entry but disables the OK button.	Low risk/Low impact.	All users	Correct the entry within the acceptable volume range to enable the OK button and save the modification.	Tool tip on the OK button directs the user to correct their entry.	CR 1,920	Maintenance TBD
Modify Units (not Pool or Split)	VBECS truncates the product name.	No identified risk/No identified impact.	All users	View the Unit History Report header.	The full product name is not normally required for modification and is not required by a user already in the modification process.	CR 1,698	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Modify Units (not Pool or Split)	When an inappropriate product type is selected for THAW, the user is allowed to proceed but cannot fully complete the modification as there is no target product.	No identified risk/ No identified impact.	All users	None required.	Re-enter the unit with the proper modification method. Code should have prevented the selection of the unit in this modification.	CR 2,286	Maintenance TBD
Modify Units: Pool Units	User cannot add units to a pooled unit created in VBECS with a Sterile Connecting Device (SCD).	Low risk/Low impact.	All users	Recreate the pool.	<ol style="list-style-type: none"> 1. Release the patient assignment , if any. 2. Inactivate the pooled unit (Blood Unit, Edit Unit Information). 3. Remove the final status of each original unit in the pool (Supervisor, Remove Final Status). 4. Create the new pool including all selected units. 	CR 2,875	Maintenance TBD
Modify Units: Pool Units	VBECS permits a biohazardous unit to be added to a pool and does not automatically designate the pool as biohazardous. By system design, a biohazardous unit is not automatically marked quarantined to allow processing of a unit that has to be given regardless of this biohazardous state.	Very Low risk/Very Low impact.	All users receiving biohazardous blood unit. (extremely rare).	The user has the option to make the unit quarantined and biohazardous. This includes AUTOLOGOUS.	Unit Short Name displays "QUAR" and the long name includes it as well. Blood centers do not ship biohazardous products unless there is a documented medical necessity. When a biohazardous unit must be added to inventory the blood bank has specific storage and handling policies to segregate this unit from the main blood inventory.	DR 1,984	Maintenance TBD
Modify Units: Pool Units	During the Pool Unit function the assigned to patient information is not displayed even though the assigned to patient has been selected for the pooled unit.	No identified risk/No identified impact	All users in a division where pooling modification has been enabled.	None required.	Before the pool is saved and proceeding to label verification, VBECS requires confirmation of the patient assignment. This occurs after the patient has been selected by the user during pool modification.	CR 1,629	Maintenance TBD
Modify Units: Pool Units	When a user chooses to add or remove units from a pooled unit , the patient assignment is duplicated on the pool. Leading to additional complications: <ul style="list-style-type: none"> • Users are also prevented from selecting additional units for a patient that has a duplicate assignment until the duplicate is removed. • A VBECS system error occurs when a user attempts to issue this unit with duplicate assignment. 	Low risk/Low impact.	All users in a division where pooling modification has been enabled.	Remove one of the assignments by using the Release Unit(s) from Patient Assignment option to allow issue of the pooled unit.	None.	CR 2,650	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Modify Units: Pool Units	VBECS allows a user to add a unit to a pool and transfer assignment to the pooled unit when both units are already assigned to the same patient. This creates a unit with dual assignment and the unit cannot be issued.	Low risk/Low impact. Platelets are not generally assigned to patients prior to pooling.	All users	Verify that the pooled unit being edited is not assigned to the same patient or do not move assignments when editing a pool.	Should a unit become assigned in duplicate to a patient and therefore unavailable for issue, use the Release Units From Patient Assignment option to remove the second assignment.	CR 2,387	Maintenance TBD
Modify Units: Split a Unit	The system does not allow split modification for product types E065 = Washed Thawed Apheresis PLATELETS and E067 = Reconstituted FRESH FROZEN PLASMA	Low risk/Low impact.	All users	None.	None.	CR 3,248	Maintenance TBD
Modify Units: Split a Unit	A system error occurs when a user attempts to Split, Discard, or Quarantine a unit and the unit has multiple antigen typing records of the same antigen (e.g., from Incoming Shipment and from testing).	Low risk/Low impact.	All users	In order to perform the process, Split, Discard, or Quarantine, the user must remove the typing record that occurred in Incoming Shipment or Edit Unit Information by editing the unit, and then perform the split, discard, or quarantine process.	This is associated with a unit record and not a patient.	CR 2,096	Maintenance TBD
Modify Units: Split a Unit	When red cell products are ABO/Rh confirmed and then split, the confirmation is not inherited by the target units. The split units are not available for selection until a confirmation test is entered for each unit.	Low risk/Low impact.	All users When split modification is enabled at the facility.	Selected for the patient prior to split modification, the split units are issuable to that patient without workaround.	User must enter ABO/Rh confirmation test results on the split units to select them. Unit is available for issue if selected for the patient prior to modification.	CR 2,295	Maintenance TBD
Modify Units: Split a Unit	When red cell products are antigen typed and subsequently divided, the antigen typing information is not inherited by the target units. When selected in Select Units, the message that the units are not antigen negative displays.	Low risk/Low impact.	All users When split modification is enabled at the facility.	Prepare the original unit for all antigen negative requirements prior to selection and modification. The split units are issuable without override. Ensure that the split units are properly labeled and that the BTRF is correctly printed.	User must enter the antigen typing information on each of the split units before selecting them. Normal business process has the split as the last step prior to issue to maximize the time available for transfusion, unless a sterile connecting device is used.	CR 2,354	Maintenance TBD
Modify Units: Split a Unit	When red cell products are divided the CMV negative and SC negative status are not inherited by the target units.	Low risk/Low impact.	All users when a unit with CMV and SC negative units are split. When Split modification is enabled at the facility.	VBECS alerts the user to the missing requirement. The user may proceed, completing the required override documentation.	The user should ensure that the unit is labeled CMV negative per local policy by viewing the Incoming Shipment section of the unit's Unit History Report. The unit record cannot be changed without removing the assignment which then leads to CR 2,295 and CR 2,354 (Split units in available status, but unavailable for selection. Split units do not have antigen negative status.). VBECS notifies the user of the lack of computer documentation and allows continued processing for patient use.	CR 2,363	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Modify Units: Split a Unit	An incorrectly formatted date may be entered in the expiration date field; the user is not warned; VBECS saves its calculated expiration date/time.	Low risk/Low impact.	All users when they opt to change the system calculated expiration date of a modified product.	Enter a date/time in the correct format, mm/dd/yyyy hh:mm.	The user may refer to the unit history report for the saved expiration date/time. There is no reason to change the calculated date/time in modification.	CR 2,494	Maintenance TBD
Order History Report	The Order History Report only displays the first 20 characters of an ordering location.	No identified risk/Low impact.	All users	None available.	Ordering locations in Vista can be shortened or changed if the uniqueness of a location cannot be determined with the first 20 characters.	CR 2,914	Report Tool Upgrade Project Proposal NSR 20080915
Order History Report	A VBECS error loop occurs when a user attempts to print the Order History Report and selects an end date before the start date.	No identified risk/Low impact.	All users	Select a date range for the report that has an end date after the start date of the report.	Close VBECS and restart session if the error loop occurs as a result of the problem described in the description.	CR 2,989 HD 465511	Report Tool Upgrade Project Proposal NSR 20080915
Order History Report	Canceled orders are not designated on the Order History Report (summary report).	Low risk/Low impact.	Supervisors (When compiling local statistics)	The Single Order History Report (detailed) includes the canceled order information.	Enhancement request to the summary report format. When the order is also cancelled in Vista, which would be the norm, there is a lab report that can be referenced.	CR 1,945	Report Tool Upgrade Project Proposal NSR 20080915
Order History Report	The Single Order History Report details the events of each order, rather than all orders, placed for the specimen.	Low risk/Low impact.	Supervisors (When compiling local statistics)	Create an inclusive report for all tests performed on a specimen by viewing or printing each order's history report.	Enhancement request to the single order history repeat to organize by specimen rather than ordered test or component.	CR 1,956	Report Tool Upgrade Project Proposal NSR 20080915
Order History Report	VBECS does not provide a message that it is compiling the report.	No identified risk/Low impact.	All users	None available.	The report takes a long time to process (A few minutes). It is not instantaneously generated. This appears to be a Crystal Reports bug. The report information is available through other reports that do not have this delay.	CR 1,846	Report Tool Upgrade Project Proposal NSR 20080915
Order History Report	Order History Report does not work with the report scheduler.	No identified risk/ No identified impact.	All users	Users must run the report when needed.	The report is available for on demand printing.	CR 2,584	Report Tool Upgrade Project Proposal NSR 20080915
Order History Report	Single Order History Report contains duplicated workload information for TAS and DAT orders.	No identified risk/ No identified impact.	All users	None required.	Correct information is duplicated.	CR 2,159	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Order History Report	Specimen History section of the Order History Report displays historic specimen acceptability incorrectly, in that an acceptable specimen is displayed as unacceptable and vice versa.	No identified risk/Low impact.	All users when a specimen has multiple updates regarding its acceptability.	None required as the current specimen status is correct.	View the record of specimen changes in the Patient History Report, Patient Specimen section where all activities are audited. The current specimen status is correct. The historic information is reversed.	CR 2,337	Report Tool Upgrade Project Proposal NSR 20080915
Order History Report	The Single Order History Report displays the patient's middle name as the provider's middle name too.	No identified risk/ No identified impact.	All users	None available.	None.	CR 2,778 HD 383473	Report Tool Upgrade Project Proposal NSR 20080915
Order Reflex Tests	VBECS does not allow a user to reflex to a weak D test from an ABO/Rh test.	Low risk/Low impact.	All users	A user can order an ABID then the weak D test and cancel the ABID to tie the weak D test to the ABO/Rh test.	None.	DR 3,643	Maintenance TBD
Order Reflex Tests	The Orderable Reflex test for DAT displays as a "DAT" instead of a "Repeat DAT."	No identified risk/ No identified impact.	All users	None required.	By definition Reflex tests are secondary tests only orderable within VBECS.	CR 1,158	Testing Grids
Order Reflex Tests	VBECS displays the Xg(a) blood group antigen incorrectly as "XgA."	No identified risk/No identified impact.	All users	None required.	Format is understandable to users.	CR 1,881	Testing Grids
Outgoing Shipment	The outgoing shipping invoice is missing required information regarding a blood unit's CMV Negative status, biohazardous, and autologous unit testing status. When assigned or restricted, the patient first name and ID are not printed on the outgoing shipment document to maintain patient privacy.	No identified risk/No identified impact.	All Users	Handwrite pertinent information on the VBECS invoice, when used to ship blood products.	None.	CR 2,734	Report Tool Upgrade Project Proposal NSR 20080915
Outgoing Shipment	VBECS allows the user to hide columns in the Selected Units section.	No identified risk/No identified impact.	All users	Exit and return to the outgoing shipment invoice; the columns are resized to default width.	If the user wishes to redisplay the column, they simple close and open the option.	CR 116	Maintenance TBD
Outgoing Shipment	The Expiration Date column in the list of units included in the shipment includes only the expiration date.	No identified risk/No identified impact.	All users	None required.	The unit expiration time is displayed correctly on the mock unit label during invoice creation. This does not impede ship out process.	CR 1,980	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Outgoing Shipment	When the return credit field is opened; VBECS enables the OK button although no change is made.	No identified risk/No identified impact.	All users	None required.	VBECS evaluates that no changes were made and no changes are saved to the database. OK button is enabled early which is an inconsistency.	CR 1,981	Maintenance TBD
Patient History Report	An expired order override exception report entry does not appear on the Patient History Report when the Exception Report is requested.	Low risk/Low impact.	All users	Users can view the expired specimen override information on the Exception Report.	None.	CR 3,204	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	An exception item that includes the product name displays only the first four characters of the product name.	Low risk/Low impact.	All users	Refer to the Unit History report for the full product name of the unit if that information is required.	None.	CR 2,886	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	A VBECS automatic update of the Decision Support System (DSS) appears as an update to the transfusion record on the Patient History and Unit History Reports.	No identified risk/No identified impact.	All users	None required.	None.	CR 3,007 HD 479169	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	A target blood unit volume may display as zero during batch processing to irradiate or thaw units. If it does, the target unit volume displays as zero on the Unit History Report.	No identified risk/No identified impact.	All users	Perform single unit modifications when irradiating or thawing units. Alternately, edit the target's volume in Edit Unit Information.	VBECS volume is a default and does not indicate a transfused volume. This is displayed only within VBECS and is updatable.	CR 1,999	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	The 2nd digit of the minutes under the Date Processed column of the Demographics section is difficult to read when a 30 character patient name is displayed.	No identified risk/No identified impact.	All users	None required.	None.	CR 2,610	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	When units are restricted for a patient in Incoming Shipments, the restricted units do not display on the Patient History Report.	Low risk/Low impact.	All users	Create a custom look up using the Blood Availability report searching for restricted units with minimal criteria selected, or as applicable to view restricted units in inventory.	Restricted units must be assigned to the patient to display on this report. VBECS forces the user to select this unit with a corresponding component order.	DR 1,643	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	VBECS displays a transfused volume of "0" on the Patient History Report when the defaults are unchanged.	No identified risk/No identified impact.	All users	When a transfusion was marked as completed in Enter Post-transfusion Information, the volume of the original or modified blood product is available from the Unit History Report.	VBECS volume is a default and does not indicate a transfused volume. This is displayed only within VBECS and is updatable.	CR 1,991	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	The VBECS domain user name (exp. VHATESTVBECSCLUSTER) will appear in the Processed By Field as having added patients to the database on the Patient History Report.	No identified risk/No identified impact.	All users	None required.	This is accurate but inconsistent as the user's name is usually displayed.	DR 2,580	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient History Report	The Patient History Report displays the NT logon of the user in the Processed By space for database conversion data.	No identified risk/No identified impact.	All users	None required.	The user is identifiable. This is accurate but inconsistent as name is usually displayed.	CR 1,908	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	"The" is misspelled in the "A patient must be selected for the report" tool tip.	No identified risk/No identified impact .	All users	None required.	Typographical error.	CR 1,703	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	The section separator between Patient Association & Unit Testing is missing.	No identified risk/No identified impact .	All users	None required.	A line is missing between the labeled sections.	CR 2,123	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	The Exception Type: Expired Task Processed does not display in the Patient History Report.	No identified risk/Low impact.	All users when using the Patient History is used to look up a processed exception rather than the Exception Report.	None required.	The Exception Report is recommended for daily supervisor review and save.	CR 2,133	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	The exception for antigen typing testing phase changes does not appear on the Patient History Report.	Low risk/Low impact.	All users	Users are directed to review the Exception Report for this information.	The Exception report is identified as one of the reports requiring supervisory review. The impact to the antigen typing result is investigated and addressed at that time not as a result from an entry on the unit or history report.	CR 2,684	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	When an incoming unit is associated to a restricted patient, a duplicate line displays in the Demographics section of the report.	No identified risk/No identified impact .	All users	None required.	Correct information is duplicated.	CR 1,639	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient History Report	In a multidivisional setting, if a patient had a historical ABO/Rh coming from a Transfusion Only (TO) division and the subsequent ABO/Rh testing on that patient is performed in the full-service division, the ABO/Rh result from the TO division is treated as a historical ABO/Rh. When displaying "ABO Discrepant unit issue" or "Unit issued with unsatisfied Transfusion Requirement" exception types on the Exception Report, VBECS calculates the historical system ABO/Rh from the test results entered not the interpretation entered. However, since a TO records only interpretations, these results do not exist and the report displays incomplete information instead of "N/A."	Low risk/Low impact when at least ONE TO division exists in a consolidated database and the most recent blood typing was performed at the TO facility. No identified risk/No identified impact for other division configurations.	All users in a full service division when at least ONE TO division exists in a consolidated database and the most recent blood typing was performed at the TO facility.	None required.	Refer to the Patient History Report or a Patient Selection window to view the full blood type. VBECS uses the TO interpretation correctly as a historical type, the Exception Report does not include the TO typing result. The user will have to refer to the patient's history report to determine if the typing was performed at the other facility.	CR 2,202	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	The report displays the (sub) section header not the heading (sub) section on following pages when a (sub) section prints over multiple pages.	No identified risk/No identified impact	All users	None required.	The header is displayed on the first page where the subsection begins.	CR 2,225	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	Post-transfusion information is printed twice on the Patient History Report "Transfusions" section.	No identified risk/No identified impact	All users	None required.	Correct information is duplicated.	CR 2,227	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	Implicated Units display same information 5 times on report.	No identified risk/No identified impact	All users	None required.	Correct information is duplicated.	CR 2,226	Report Tool Upgrade Project Proposal NSR 20080915
Patient History Report	Patient History Report: Antigen Positive/Untested Unit issued exception is listed twice.	No identified risk/No identified impact	All users	None required.	Correct information is duplicated.	CR 2,408	Report Tool Upgrade Project Proposal NSR 20080915
Patient Information Toolbar	VBECS does not display transfusion reactions that are not finalized.	Low risk/Low impact.	All users when a TRW is pending completion and this option is used to look up the information.	Check the PTL for non-finalized transfusion reaction workups.	This test is incomplete without the finalization and is not ready for evaluation using this option.	CR 1,658	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Information Toolbar	The Recent Orders option does not display pending orders (not accepted) or completed transfusion reaction workups.	Low risk/Low impact.	All users when a TRW is completed and this option is used to look up the information rather than the appropriate use of reports that contain completed test results.	Access transfusion reaction workup information from the patient's Transfusion Reaction History Report or the Finalize/Print TRW option.	This may result in a repeat order of the TRW, which would be cancelled when received and the patient history checked during processing.	CR 1,927	Maintenance TBD
Patient Merge	VBECS displays duplicate merge events to user for merge.	Very Low risk/Low impact.	All users	Users can process one of the merge events manually and clear both merge alerts.	The Master Veteran Index team is working on a patch to discontinue sending a duplicate message.	CR 3,099 MVI CR 232	Master Veterans Index Patch
Patient Merge	A merge alert displays but has no discernable differences between the patient information in VBECS.	Very Low risk/Low impact.	All users	None available.	If the merge from matches the merge to file a Remedy ticket for Tier 3 support. A merge of DFN to another DFN is not rare but the patient information associated to those DFN values being identical is uncommon. VBECS does not display a change in the DFN identifier when two records are merged. It is possible that a site would merge two identical patient records with only the DFN being different.	CR 3,060	Maintenance TBD
Patient Testing	VBECS does not correctly re-display AHG Only crossmatch testing grid as originally selected for the component order.	Very Low risk/Low impact.	All users	Place a new CPRS RBC order and select the AHG only or All Phases grid format desired for the additional units' crossmatch tests.	The grid appears enabled for all phases but only the AHG and CC fields allow data entry. The AHG only grid pattern does not include enabled CC fields. VBECS does not allow "not tested" entry for the CC.	CR 2,812 HD 395766	Maintenance TBD
Patient Testing	A user is permitted to select two partially completed tests with different test phases (e.g. AHG only and All Phases) along with a not started test. The testing phase for the not started test will present as AHG only even though results for all phases are required by VBECS. Entering only AHG test results for the specimen that was not started and saving will result in a partially completed test that cannot be completed or invalidated.	Very Low risk/Low impact.	All users	Select partially completed orders that have been previously saved with the same grid test pattern. Optionally, a not started test may be selected with these but must be tested with the same grid pattern to avoid the problem.	Specimens that cannot be completed or invalidated require the attention of Tier 3 support. File a Remedy ticket for correction.	CR 3,118	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Testing	When the user clicks the red X button in the upper right corner and clears the test results; VBECS clears the grid but does not disable the OK button when there are results in a second grid.	No identified risk/No identified impact	All users	None required.	The grid is cleared to allow for reentry of results. If the OK button is clicked only the results of the second grid are saved. The test remains pending for reentry.	CR 1,617	Testing Grids
Patient Testing	VBECS does not clear the test interpretation cell when the user clicks the red X in a Transfusion Only division and is in any of the testing grids in Patient Testing.	No identified risk/Low impact.	All users at TO facilities.	None required.	Invalidation is accepted only in Invalidate Patient Testing. The intent of the red x is to clear the grid for re entry of reaction results. The user may clear the interpretation cell and re enter the interpretation. VBECS will perform row validation between the entered reaction results and the interpretation enforcing various system rules.	CR 1,644 CR 1,645	Testing Grids
Patient Testing	A VBECS system error occurs when a user tries to save a result of "H" (hemolysis) in the Patient Antigen Typing testing grids.	No identified risk/ No identified impact.	All users	Users are instructed not to enter "H" in the Patient Antigen Typing test grids. Any attempt to save a result of "H" will cause a system error and clear the testing grid of that result when the user re-enters VBECS.	The key for acceptable entries in the testing grid does not include "H" but it is permitted to enter until saving.	CR 2,685	Testing Grids
Patient Testing	The pattern IS-N 37-X AHG-P CC-N is considered valid for antibody screen testing.	Low risk/Low impact.	All users	None required.	This is technically acceptable for a prewarmed crossmatch or sites using PeG enhancement. This is allowed with the processing of an override for Nonstandard test methods.	DR 2,140	Testing Grids
Patient Testing	A system error may result when the user tries to enter unacceptable characters or tab where not allowed.	No identified risk/ No identified impact.	All users	None required.	Difficult to reproduce requiring significant manipulation to trigger. Do not tab around needlessly while in the patient testing grids.	CR 2,367	Testing Grids
Patient Testing	A user is allowed to enter more than 50 characters in the Off-Site Location field when entering antibody ID results without a message from VBECS that only the first 50 characters will be saved.	Very Low risk/Very Low impact.	All users	Do not exceed 50 characters for the name of the off-site testing location.	None.	CR 2,626	Testing Grids
Patient Testing	When a site is defined as "full service" and daily QC was not performed the Exception Report exception type: "QC not performed on rack used for testing" entry is not saved when the associated patient has a middle initial(MI).	Very Low risk/Very Low impact.	All users when Daily QC was not performed and the patient with an MI is tested.	Standard practice is not to override the QC not performed warning message. Verify the QC was performed daily by reviewing the Testing Worklist Report.	The user is warned that QC has not been performed in all instances. This exception is accurately collected and displayed for unit testing and patients without a middle initial. The user is warned at each testing episode when QC has not been performed and the override is collected in all but a rare situation. There is some flexibility in the 24 hour clock regarding QC testing.	CR 2,402	Testing Grids

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Testing	Comments entered in the DAT POS and NEG control cell rows are not displayed on any report in VBECS or CPRS.	Low risk/Low impact.	All users	Enter comments only in the patient test row.	Display should not have multiple comment rows. Comments related to the patient testing are entered in the patient row. This is a training issue. Comments related to the patient test should only be entered in the patient test comment cell.	CR 2,421	Testing Grids
Patient Testing	Users cannot order a reflex weak D test for an ABO/Rh test as ordered from CPRS.	No identified risk /Very Low impact.	All users	Order a reflex antibody ID (ABID) and the weak D test at the same time. Then cancel the ABID.	None.	DR 3,643	Testing Grids
Patient Testing: Pending Task List	A system error occurs when a user cancels an order adding canned and free text comments that exceed 255 characters total.	No identified risk/ No identified impact.	All users	Users must not exceed a total of 255 characters of combined canned and free text when cancelling an order.	None.	CR 2,874	Testing Grids
Patient Testing: Pending Task List	When a user enters random alphanumeric characters (nonsense); VBECS may or may not return a list of patients.	No identified risk/ No identified impact.	All users	Enter viable search criteria: full name or first initial of last name and last four digits of the patient ID.	Test sites state that they also use the full SSN for patient identification.	CR 1,891	Testing Grids
Patient Testing: Pending Task List	VBECS does not display the patient's middle initial on the printed Pending Task List.	No identified risk/ No identified impact.	All users	The patient's middle initial is displayed in the Pending Task List user interface.	The online display is normally used.	CR 1,775	Testing Grids
Patient Testing: Pending Task List	VBECS does not include a date range search field.	No identified risk/ No identified impact.	All users	None required.	This is an enhancement request as this is not a normal search usage per the field test sites.	DR 2,586	Testing Grids
Patient Testing: Pending Task List	A system error occurs when the search option "Results Corrected" is used.	No identified risk/ No identified impact.	All users	None required.	There are other fields available for use that will present the orders associated for a single patient. Field sites report that they require entry of the patient full SSN or specimen UID to select an order for processing.	CR 2,490 HD 264904	Testing Grids
Patient Testing: Record Patient ABO/Rh	Test entries are not being cleared when canceling out of the Invalid Results message.	Low risk/Low impact.	All users	Correct the testing entries or click the red X to return to the PTL and retest.	None.	CR 2,042	Testing Grids

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Testing: Record a Cross-match	When the user attempts to save a negative or not tested check cell result, VBECS displays the "Invalid result. Check cells must have a positive result. Repeat test," not "Interpretation does not match your results. Please correct" message.	Low risk/Low impact.	All users	None required.	The message displayed to the user is misleading; VBECS requires entry correction to proceed. The user must enter a valid result.	CR 1,862	Testing Grids
Patient Testing: Record a Cross-match	A currently crossmatched unit is selected for a different patient. The unit status displays as "Crossmatched."	No identified risk/ No identified impact.	All users when the local practice is to crossmatch a single blood unit for multiple patients.	None required.	VBECS displays the unit status in a hierarchy when a crossmatch is present on any patient or unit and the unit is available for selection, assignment, or crossmatch to another patient.	CR 1,601	Testing Grids
Patient Testing: Record a Cross-match	When returning to partially completed crossmatch tests or testing additional units, the originally selected crossmatch grid configuration is applied.	No identified risk/ No identified impact.	Any user that has saved incomplete test results and returns to enter the remaining results.	None required.	The test is partially completed with a test pattern setting and must be completed in that format. This is a training issue.	CR 1,972	Testing Grids
Patient Testing: Record a Cross-match	A user can bypass the entry of test results in the IS field when performing a crossmatch.	No identified risk/ No identified impact.	All users	None required.	VBECS saves the entered results but does not allow the entry of a crossmatch interpretation while the testing fields are incomplete, as required to accommodate saving partially completed tests. The user may enter their results in any order but may not save the entry until completed.	CR 1,620	Testing Grids
Patient Testing: Record a Direct Antiglob-ulin Test	DAT grid does not properly calculate the QC status of the PS AHG reagent when multiple lot numbers are used on the same day.	Low risk/Low impact.	All users when multiple lots of the same reagent type are in use on that day.	View the QC data for the day for the lot number in question from the Testing Worklist Report and re-enter the results for the test grid (or repeat the testing).	The user may repeat the reagent QC results or look them up and reenter them. The requirement stated that the reagent should not require repeat on the same day. Some sites require retesting of the reagent at each use.	CR 2,179	Testing Grids
Patient Testing: Record a Direct Antiglob-ulin Test	The Anti-Human Globulin reagent lot number entered is not saved when QC is recorded with the patient test. The Patient Testing Worklist Report displays the positive control lot number.	Low risk/Low impact.	All users	Record the reagent lot number for the Anti-Human Globulin in the test comment field.	When the primary AHG reagent is tested with daily QC, there is no issue. This is an issue with subsequent AHG testing associated with an antibody identification where the user must manually add the reagent lot number in addition to the displayed entry.	CR 2,020	Testing Grids

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Testing: Record a Direct Antiglob-ulin Test	VBECS displays “IgG” instead of “Anti-IgG” in the Tested With column of the Patient Testing Worklist Report.	No identified risk/ No identified impact .	All users	None required.	Reagent name is clear as this is the only reagent with this name.	CR 1,638	Testing Grids
Patient Testing: Record a Patient Antibody Screen	VBECS displays a system error message when the user clicks a tab for antigen testing, does not enter any test results, and clicks Cancel and Yes to close the window.	No identified risk/ No identified impact.	All users	Do not cancel out of testing without entering test results.	The user has indicated to exit this option. The user must log into VBECS to continue.	CR 1,797	Testing Grids
Patient Testing: Record a Transfusion Reaction Workup (TRW)	When a second implicated unit is selected prior to the completion of the first unit, VBECS does not record the associated data entry for the second unit.	Low risk/Low impact.	All users	Add the first unit completely, and then add subsequent units.	Details affected are the Bag Returned information, hemolysis, Checks OK, and Comments. This does not occur when the user selects units for entry from left to right. The user may invalidate the TRW testing prior to finalization and enter the data.	CR 2,053	Testing Grids
Patient Testing: Record a Transfusion Reaction Workup (TRW)	A system error occurs when a user attempts to access a transfusion reaction workup that is already in progress and locked by another user.	No identified risk/ No identified impact .	All users	None available.	Simultaneous data entry on the same order is not allowed to avoid data corruption. They must log into VBECS to continue with other orders.	CR 2,215	Testing Grids
Patient Testing: Record a Transfusion Reaction Workup (TRW)	Date Reaction Noted and Date Reaction Investigated is updated and reset to the current division date/time each time the TRW is opened for editing.	No identified risk/ No identified impact .	All users	None required.	The Date Reaction Noted is recorded correctly and displayed as first entered on the Patient History Report. This is a display issue only as the original Date Reaction Noted is correctly recorded in the database and on the report.	CR 2,218	Testing Grids

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Testing Worklist and Testing Worklist Reports	The Rack QC and Reagent QC sections display only the first 100 characters of comments entered in the POS or NEG control comment fields as saved in a Patient or Unit Antigen Typing test.	No identified risk/ No identified impact .	All users	Limit comments to less than 100 characters to allow the full comment to display on the report.	None.	CR 2,868	Maintenance TBD
Patient Testing Worklist and Testing Worklist Reports	Testing Worklist Report is missing its Title and cannot be scheduled.	No identified risk/ No identified impact.	All users attempting to print the report; Supervisor is responsible for reviewing the report.	Print the report and immediately write the title of the report on the pages.	The report may be printed on demand. The title may be written on the report manually, if required. This is the only report that contains the reaction results and is clearly the Testing Worklist Report.	CR 2,070	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	The Rack QC Testing Worklist Report section of the Testing Worklist Report only displays the testing tech and not the identity of the logged in tech that entered the results.	Low risk/Low impact.	Supervisor responsible for report review	Maintain the original testing records for the QC that was performed offline with the testing tech information. Do not change the testing tech name when entering results in VBECS so that the logged in tech that enters the results is included on the Testing Worklist report.	Local policy dictates the storage of downtime records. Storage of the original work in addition to the computer entry is customary.	DR 3,556	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	Crossmatch tests that were invalidated and retested are not marked on the testing report.	No identified risk/ No identified impact .	Supervisor responsible for report review	The sequence of test performance by date and time on the report identifies the most recent result.	This is a training issue. The tests are displayed in the order of testing.	CR 2,034	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	The Testing Worklist Report does not display the "Entered in Error" comment when a Unit ABO/Rh Confirmation Test is invalidated.	No identified risk/ No identified impact .	Supervisor responsible for report review	None required.	The unit has a status of "Limited" and is not available for patient selection. The unit is properly restated with the test invalidation. The invalidated result is included on the report for review. The supervisory review would require minimal additional time.	CR 2,083	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	The Testing Worklist Report only displays the first 218 characters entered in the ABID results comment field.	No identified risk/ No identified impact .	Supervisor responsible for report review	The full comment is viewable in the Patient History Report.	Users are limited to entering 255 characters in the ABID results comment field.	CR 2,627	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Testing Worklist and Testing Worklist Reports	"Automated Instrument" does not display on the report.	Low risk/Low impact.	Supervisor responsible for report review	Add instrument information from the instrument output during the report's review.	The results from an automated instrument must be entered by a VBECS user whose name displays. This is not incorrect but does not include the secondary information that the test was performed with a specific instrument.	CR 2,046	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	The generation of testing reports may be prolonged.	No identified risk/ No identified impact .	Supervisor responsible for report review	The report can be scheduled to print.	This report is printed for review. This is dependent on the amount of data being retrieved.	CR 1,867	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	Rack IDs are not consistently displayed in the Patient Testing Worklist Report.	No identified risk/ No identified impact .	Supervisor responsible for report review	None available.	1) The one letter rack IDs are supposed to appear in the first line of the testing entry. In many cases, it is appearing in seemingly random line number positions, sometimes appearing multiple times in one testing entry. 2) When an XM test is invalidated, the Rack ID is not displayed. The missing rack IDs are an artifact when the same rack is used for the series of tests.	CR 2,063	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	The weak D report format includes line items for IS and RT although these phases are disabled in the testing grid and may never have entered results.	Low risk/Low impact.	Supervisor responsible for report review	None required.	The report displays these phases with no results when none have been entered.	CR 2,040	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	The Patient Testing Worklist Report displays "Inc" for test interpretations that have different meaning. In the Interpretation column for crossmatch tests, "Inc" means incompatible; for other diagnostic tests, "Inc" means "inconclusive" interpretation.	No identified risk/ No identified impact.	Supervisor responsible for report review	None available.	Crossmatch test does not include an "inconclusive" interpretation. Other diagnostic tests do not have an "incompatible" interpretation. The definition is associated with the particular test and is not confusing to the blood bank staff.	CR 2,032	Report Tool Upgrade Project Proposal NSR 20080915

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Testing Worklist and Testing Worklist Reports	VBECS displays an extra comma after the last specimen UID when a patient has multiple associated specimen UIDs.	No identified risk/ No identified impact .	All users	None required.	Extra comma does not interfere with the displayed data.	CR 1,579	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	Testing Worklist Report: Unit and Patient Testing Details window reflects the reagent rack's testing status on the current date not on previous dates (retrospective data entry) and may require an override to proceed.	No identified risk/ No identified impact .	Supervisor responsible for report review	None Available.	Verify that QC was performed on the selected date(s) using the Testing Worklist Report. This is a report pulled for review and part of that review is to review that QC was properly performed on each date selected for the report.	CR 2,384	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	The Testing Worklist Report does not display the last name of the tester if their user name includes a site abbreviation [e.g. First Name, Middle Initial, (Division Abbreviation Last Name)].	No identified risk/ No identified impact .	All users	There is enough of the user name present on the report to identify the individual that performed the testing.	The VBECS report generating tool is displaying the first, middle, and division abbreviation, for the users that have a division abbreviation included in their user name.	CR 2,694	Report Tool Upgrade Project Proposal NSR 20080915
Patient Testing Worklist and Testing Worklist Reports	Testing Worklist Report does not provide a comprehensive list of reagent lot numbers entered for daily reagent rack QC.	Low risk/Moderate impact.	All users	Continue to record all reagent lot numbers on a rack with existing process (hard copy or spreadsheet). Retain with Testing Worklist review records.	The user must record some lot numbers manually as they are not displayed on the report. Field sites have opted to continue current manual recording of all lot numbers to minimize confusion.	CR 2,385	Report Tool Upgrade Project Proposal NSR 20080915 or Daily Reagent Quality Control Upgrade NSR 20081120
Patient Updates	VBECS performance will slow when a patient has ten or more active ordered components because it must calculate the active orders for each patient checked.	Low risk/Low impact.	All users	None Available.	It is unlikely that this will occur unless the patient has multiple duplicate orders for the same blood component.	CR 2,177	Maintenance TBD
Patient Updates	The last update date time is updated on all Patient Updates displayed each time a new update is viewed.	Low risk/Low impact.	All users	None required.	None.	CR 2,711	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Patient Updates	Patient updates display in VBECS for active patients when the data change is unrelated to VBECS data.	No identified risk/ No identified impact.	All users	None Available.	The changing demographics do not affect VBECS or the patient identity so the update will not reflect any changes. However, the patient update event is still logged and appears on the Patient History Report and the Transfusion Requirements Report.	CR 2,580 CR 3,249	Maintenance TBD
Patient Updates	When a Vista patient update removes a patient's middle name VBECS displays the Patient Name and the Previous name in the Patient Name column. The Previous Name column displays the previous name correctly.	Low risk/Low impact.	All users	Close and re-open the Patient Update window. The subsequent display of the patient information is correct.	This is a display issue only. The VBECS database, reports and any future transactions have the correct information.	CR 2,709	Maintenance TBD
Post – Transfusion Information	VBECS allows the user to save a transfused volume greater than the default average volume associated with the product code in Blood Products with a transfusion record.	No identified risk/ No identified impact.	All users	None required.	This may actually be a benefit as an accurate transfused volume may be entered. The default product volume is an average, not the maximum.	CR 1,655 CR 2,607	Maintenance TBD
Post – Transfusion Information	The Presumed Transfused indicator does not clear when post-transfusion information is updated.	No identified risk/ No identified impact .	All users	None required.	The transfusion end date is maintained properly in the database; the Presumed Transfused indicator is retained. The updated transfusion data are displayed in the Unit History Report. This may actually be a benefit as the record is clearly marked and indicates that the update occurred over 48 hours after administration and may require investigation as to why the information was not earlier made available.	CR 1,936	Maintenance TBD
Post – Transfusion Information	The OK button is enabled on the Post Transfusion Information window when a Traditional Supervisor or above selects a unit that was marked transfused.	No identified risk/Low impact.	Traditional Supervisor and above.	None required.	The OK button should only be enabled after a change has been made to saved information.	CR 2,550	Maintenance TBD
Post – Transfusion Information	No warning displays for missing workload process when entering post-transfusion information.	No identified risk/Low impact.	Supervisor responsible for workload reporting.	None required.	Associate workload process with Enter Post-transfusion Data. The application does behave correctly when a workload process is defined.	CR 2,263	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Post – Transfusion Information	A system error occurs when entering post-transfusion information if the tabs are selected out-of-order.	Low risk/Low impact If VBECS error occurs no data is saved and the information must be re-entered.	All users	Enter the tabs in order presented to avoid causing VBECS error.	If a VBECS error occurs no data is saved and the information must be re-entered.	CR 2,208, CR 2,879, HD 426258	Maintenance TBD
Print Unit Caution Tag & Transfusion Record Form	A VBECS system error occurs when a user enters a crossmatch for a patient and attempts to print caution tags or a blood transfusion record form while another user is performing testing on the same patient.	Low risk/Low impact.	All users	Reprint the caution tags from the Print Unit Caution Tags & Transfusion Record Form option.	None.	CR 2,922	Maintenance TBD
Print Unit Caution Tag & Transfusion Record Form	VBECS displays a message that the caution tag(s) were successfully printed, but the printer has failed and printing was not successful.	No identified risk/ No identified impact.	All users	Physically verify that the requested tag is printed. The "Success" message is a message that the print has queued to the print queue, and does not verify that the printing has actually occurred.	If the tag didn't print, the user would reprint the tag. The tag usually has printed successfully without incident.	CR 2,059	Report Tool Upgrade Project Proposal NSR 20080915
Print Unit Caution Tag & Transfusion Record Form	A VBECS system error occurs when a Caution Tag is requested and the printer is not correctly configured.	No identified risk/ No identified impact.	All users	Configure the printer as detailed in the VBECS installation guide with subsequent testing prior to use.	This will be discovered during local validation and corrected prior to production installation.	CR 1,660	Report Tool Upgrade Project Proposal NSR 20080915
Print Unit Caution Tag & Transfusion Record Form	VBECS can only print ten blank caution tags per request. A user request for more than ten blank caution tags will result in only ten tags printed.	No identified risk/ No identified impact.	All users	Perform several requests for ten blank caution tags to fulfill the quota if more than ten caution tags are required.	None.	CR 2,624	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Print Unit Caution Tag & Transfusion Record Form	The printed Caution Tag displays "Not Applicable" in the crossmatch field when a unit does not require a crossmatch. The Blood Transfusion Record Form (BTRF) displays "Not Required" in the same situation.	No identified risk/ No identified impact.	All users	None required.	The messages are materially equivalent. If the user feels that this is unacceptable they would manually change the text on one of the documents to match letter for letter with the other.	CR 2,007	Report Tool Upgrade Project Proposal NSR 20080915
Print Unit Caution Tag & Transfusion Record Form	When an antigen negative requirement was entered in Special Instructions and Transfusion Requirements, the message and tool tip includes the antigen negative requirement in a sentence structured for an antibody history insertion.	No identified risk/ No identified impact.	All users	None required.	The statement is grammatically incorrect but conveys that the unit is untested or positive for the required antigen to a Blood Bank user.	CR 1,814	Maintenance TBD
Print Unit Caution Tag & Transfusion Record Form	VBECS cannot print a BTRF documenting retrospective crossmatch compatibility.	No identified risk/ No identified impact	Supervisor reviewing retrospective data entry related to emergency issue of blood products	None required.	The information is available in the various reports associated with testing documentation. The user may want to write up a document "for the re cords" based on paper record driven policy for documentation for retrospective test entry related to the issuance of blood products with incomplete pretransfusion testing. This form is not required as one has already been written and used to document the blood product administration.	DR 1,956	Report Tool Upgrade Project Proposal NSR 20080915
Print Unit Caution Tag & Transfusion Record Form	Remove the ">" from the message: >This order must be emergency issued at this time. This unit does not qualify for emergency issue. Tag(s) cannot be printed.	No identified risk/ No identified impact	All users	None required.	The statement is grammatically incorrect but is clear to a Blood Bank user.	CR 2,086	Maintenance TBD
Print Unit Caution Tag & Transfusion Record Form	When blank Caution Tags are printed in preparation for down time and the print job is canceled, the number printed is not correct in the message.	No identified risk/ No identified impact	All users	None available.	The user may count the number of tags printed to determine how many blank tags should be printed.	CR 2,003	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Product Modifications (Division Configuration)	When a user enters an amount less than 0.00 or greater than 9999.99 in a Cost field and presses the Enter key, VBECS changes the amount to 0.00 and 9999.99, respectively, and enables the OK button.	Low risk /Low impact.	Administrator	Check the cost field entries for validity prior to saving.	This is a configuration option where an upper level user is defining the cost of the process. VBECS will save the cost at the maximum allowed number.	CR 1,732	Maintenance TBD
Prolonged Transfusion Time Report	Units issued to a remote storage location will appear on the report when the transfusion start time is greater than 30 minutes after issue time.	No identified risk/ No identified impact	All users can print the report, but a Supervisor would be responsible for review	None available.	Units issued to a remote storage location should not appear to have a prolonged issue time as this setting indicates that the unit is issued to a monitored location maintained at the proper temperature for the product type.	CR 2,803	Report Tool Upgrade Project Proposal NSR 20080915
Prolonged Transfusion Time Report	The Prolonged Transfusion Time Report closes after the user prints it.	No identified risk/ No identified impact	All users	Open report again if on screen review is still desired.	None.	CR 3,100	Report Tool Upgrade Project Proposal NSR 20080915
Prolonged Transfusion Time Report	When the user clicks OK after selecting the printer, the print preview view of the Prolonged Transfusion Report closes.	No identified risk/ No identified impact	All users can print the report, but a Supervisor would be responsible for review	None available.	The user has requested the hard copy and no longer requires this view. This is not consistent with other reports in the application.	CR 2,231	Report Tool Upgrade Project Proposal NSR 20080915
Prolonged Transfusion Time Report	The delayed start time presented is the total time from issue to start time; the prolonged transfusion presented is the time from issue to the transfusion end time.	No identified risk/ No identified impact	All users can print the report, but a Supervisor would be responsible for review	None available.	The times are not inaccurate but are not in the preferred format where a delayed start would calculate only 30 minutes after issue to the start for its calculation (example, 34 minutes, not 4 minutes). The prolonged transfusion would be from the start to end time that exceeds the maximum transfusion time (MTT) set for the component class, for example, 500 minutes, not 30 minutes when maximum MTT is 470. This is a report for a transfusion committee. The data is accurate, but does not present as the VBECS business rule is written.	CR 2,499	Report Tool Upgrade Project Proposal NSR 20080915
Reagents	A system error occurs when a user enters reagents, clicks Ok to save and then NO to not commit the changes.	Low risk /Low impact.	All users	No save occurs with the system error event. A user only need re-enter reagents if they truly meant to save instead of cancelling the save.	None.	CR 2,810	Daily Reagent Quality Control Upgrade, NSR 20081120

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Reagents	Checking all of the reagents types to view on the Reagent Inventory Report will only display reagents with at least one vial in inventory.	Low risk/Low impact.	All users	Select specific reagents individually to view and print. Print and save a copy of the report prior to zeroing out any expired product lot numbers.	Reagents that have zero vials remaining will not display.	CR 3,133 HD 514507	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	Comments entered for unsatisfactory reagents do not appear on any VBECS reports.	Low risk /Low impact.	All users	The comments are saved in VBECS but are not viewable. No workaround is available.	If the comments must be verified or viewed the user can file a request for this information with the National Help Desk.	DR 4,346 HD 541145	Report Tool Upgrade Project Proposal NSR 20080915
Reagents	Once tripped, the message warning of low reagent inventory continues to alert the user in Update Reagents even when the inventory is above the minimum.	No identified risk/ No identified impact	All users	Set the "Minimum Stock Vial Level" in Maintain Reagent Minimum Inventory to zero to disable the warning message.	This is an annoyance that is corrected by resetting the configured number to zero for the specific reagent lot number.	CR 1,505	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	When reagents are received, updated and marked unsatisfactory in the same transaction; comments and details entered are concatenated and are not displayed on the Reagent report.	Low risk /Low impact.	All users can print the report, but a Supervisor would be responsible for review	Document the reason a reagent was unsatisfactory upon receipt on the manufacturer's invoices for future reference.	None.	CR 1,875	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	VBECS does not maintain a reagent lot number's history when the quantity of an expired reagent is set to zero.	Low risk /Low impact.	All users can print the report, but a Supervisor would be responsible for review	Do not set the reagent lot number quantity to zero when the record is to be maintained.	This is an expired reagent and would not be used routinely. The lot number continues to be available for use should that be the case.	CR 1,958 CR 1,931 HD 359714	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	A system error occurs when sorting reagents by invoice if the invoice contains any non-numeric characters	No identified risk/ No identified impact .	All Users	Do not enter non-numeric characters for the invoice number.	None.	CR 2,528	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	A warning icon appears when the user selects "No" to the entry confirmation message.	No identified risk/ No identified impact .	All users can print the report, but a Supervisor would be responsible for review	Click the Clear button to continue.	The user must click clear to close the window. This is a nuisance in a non-patient care option in the application.	CR 1,952	Daily Reagent Quality Control Upgrade, NSR 20081120

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Reagents	The Vials Received per Lot Number field in Log In Reagents allows the entry of a decimal that causes the reversal of the entry (e.g., user entry of 1.5 becomes 51).	No identified risk/ No identified impact .	All users can print the report but a Supervisor would be responsible for review	Enter whole numbers; do not enter decimals. Check the accuracy of the entry before saving.	This is a nuisance in a non-patient care option in the application.	CR 2,067	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	Previously entered manufacturers' names may not persist in the manufacturers' drop-down list in Log In Reagents.	No identified risk/ No identified impact .	All users can print the report but a Supervisor would be responsible for review	Enter the manufacturer's name.	This is a nuisance in a non-patient care option in the application.	CR 1,501	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	VBECS allows the entry of only free-text comments after revising fields when the user selects "Other" from the Reason for Change drop-down list.	No identified risk/ No identified impact.	All Users	None required.	This is consistent with the selection of the comment "OTHER" in the rest of the application.	CR 1,486	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	The warning message displayed when the user exceeds the upper limit of the Vials Received field in Log In Reagents states the range as 0–999 when it is actually 1–999.	No identified risk/No identified impact.	All users can print the report but a Supervisor would be responsible for review	None required.	The lower range is incorrect in the documentation. The upper limit is the same in both cases.	CR 1,910	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	When the inclusive dates entered for Expiration Date Before/After include the current date, reagents that expire on the current date may not be included.	No identified risk/No identified impact.	All users can print the report but a Supervisor would be responsible for review	None required.	The reagent is not expired on the date the report is requested. The reagent is available.	CR 1,966	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	Anti Pk, Anti-PP1, Anti-I, and Anti-I(int) are available as reagent types with no corresponding active test.	No identified risk/No identified impact	All users can print the report, but a Supervisor would be responsible for review	None required.	There are no associated test for these antigens in VBECS.	CR 1,882	Daily Reagent Quality Control Upgrade, NSR 20081120

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Reagents	VBECS truncates text in the print preview and printed versions of the Reagent Inventory Report. When VBECS does not find a match for selections entered in the Reagent Search window, it creates a report with "No Reagents Found" instead of presenting a message that no match was found.	No identified risk/No identified impact	All users can print the report, but a Supervisor would be responsible for review	None required.	This is an inconsistent response.	CR 1,515	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	The time defaults to the time when the user enters the Reagent, Update Inventory window, not the time when a lot number is selected.	No identified risk/No identified impact	All users can print the report, but a Supervisor would be responsible for review	None required.	As the report is generated after the user enters the option, this would not be a significant time difference that would impact the information in the report.	CR 2,068	Daily Reagent Quality Control Upgrade, NSR 20081120
Reagents	In a multi-divisional environment, the reagent report for a division will print one line item for each division that has set minimum levels for that reagent type. For example, if you have minimum levels for Reagent A, and 2 other divisions on your VBECS have minimum levels for Reagent A, your division's reagent information will print three times on your report.	No identified risk/Low impact in a consolidated division No identified risk/No identified impact in a single division database And in consolidated database where not all facilities set minimum reagent levels.	All users can print the report, but a Supervisor would be responsible for review	None required.	There is no crossover of data and no safety issue. If minimum levels are not set in a division, then it will not factor into the display. In a single-division environment, there is no problem. Correct information is duplicated. This occurs only in consolidated divisions.	CR 2,233	Daily Reagent Quality Control Upgrade, NSR 20081120
Recent Orders	Cancelled patient orders do not appear in the Recent Orders menu option.	No identified risk/No identified impact	All users	Users can view information about cancelled orders on a patient by using the Order History Report.	Cancelled orders have never appeared in the Recent Order view.	CR 3,241	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Release Units From Patient Assignment	The scanner icon is currently displayed but is not working.	No identified risk/No identified impact	All users	None required.	This is an inconvenience. The user may opt to release units using the patient route rather than the unit driven route.	CR 2,200	Maintenance TBD
Release Units From Patient Assignment	In the Release Units menu option, a system error occurs if the user selects a comment, then deletes it by pressing the space bar and clicks the OK button to save.	No identified risk /No identified impact.	All users	None required.	The release of units requires a comment to save.	CR 2,769	Maintenance TBD
Return Issued Units To Blood Bank	When all assigned units are released, the component order status is returned to "Not Started" on the Pending Task List.	No identified risk /No identified impact.	All users	None required.	The order remains available for use and expires normally based on division settings. This is viewed as a training issue. The order status is calculated based on the units currently associated with the order.	CR 2,051	Select Unit Project Proposal
Return Issued Units To Blood Bank	VBECS may display an inactivated unit in the unit search screen without any indication of its inactivated record status.	No identified risk /Low impact.	All users	None required.	All unit information is correct.	CR 1,877	Select Unit Project Proposal
Return Issued Units To Blood Bank	Visual Inspection Information is not included in Exception Report.	No identified risk /No identified impact.	All users	View the Unit History Report for the visual inspection associated with an exception recorded for a blood unit relocation.	The display of the visual inspection response is included in the unit's Unit History report which may be opened during the investigation of the exception by the supervisor.	CR 2,119	Select Unit Project Proposal
Return Issued Units To Blood Bank	The selected unit marked unsatisfactory for return is removed from the Select Units to Return list view, but the information in the Selected Unit pane for the unsatisfactory unit remains unchanged.	No identified risk /No identified impact.	All users	None required.	The display changes when another unit is selected and updates the view.	CR 2,089	Select Unit Project Proposal
Select Units	A system error occurs when a user enters a unit product code and presses "Enter" rather than "Tab" key.	Low risk /Low impact.	All users	When a user is entering unit ID and product code via the keyboard, press the "Tab" key or use the mouse to move to the next field.	Recommended business practice is to use the hand held barcode scanner not type in the barcoded information. The error occurs because the user presses "Enter" which clicks OK without submitting the unit information causing the crash.	CR 3,227	Select Unit Project Proposal

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Select Units	A system error occurs when a user selects a unit and then filters the available units search screen for a product different than the unit selected and uncheck the unit selected.	Very Low risk / Very Low impact.	All users	Uncheck all selected units before filtering available units.	The steps to get to the system error are convoluted in normal practice and may never be encountered.	CR 2,821	Select Unit Project Proposal
Select Units	In a multidivisional database, restricted units residing in another division will be displayed but are not selectable. The tool tip (mouse over) message wording incorrectly states, "Unit is not in division ###."	No identified risk /No identified impact.	All users in consolidated VistA with multiple VBECS Divisions	Read this message without the word NOT in the sentence. The division shown is the division where the unit is physically located.	The unit is not selectable as it is physically in a different location. This occurs only in consolidated divisions which are limited in number and with restricted units which are not a high volume product.	CR 2,107	Select Unit Project Proposal
Select Units	VBECS displays ISBT unit ID and product short name, but not the product code. A split ISBT unit's uniqueness is determined by the seventh and eighth digit of the product code; therefore, VBECS cannot identify the difference.	Very Low risk / Very Low impact. No identified risk/No identified impact when SPLIT modification not enabled at the facility.	All users in divisions where split unit modification is enabled.	None available.	The user scans/enters the product information and does not use the pick list.	DR 1,729	Select Unit Project Proposal
Select Units	VBECS blood product table does not include the ICCBBA classes Thawed POOLED PLASMA (E055), Thawed POOLED FRESH FROZEN PLASMA (E056), WASHED GRANULOCYTES (E057) and Liquid POOLED PLASMA (E058).	No identified risk/No identified impact.	None.	None required.	The product codes associated with these blood product classes are available for use and the class changes are for mapping tables that are not visible to the user.	DR 1,703	Select Unit Project Proposal
Select Units	The message used to state that the patient has a history of Anti-K and the unit is positive for the antigen K, references the antigen type as anti-K and the antibody as K as well as all other antibody/antigen pairs.	No identified risk/No identified impact.	All users when the patient has an identified irregular antibody and selected an untested or positive unit for the patient.	None required.	The terms are flipped but the content of the message is clear to a blood bank user.	CR 2,093	Select Unit Project Proposal
Select Units	Blood product 00950 E001 WB CPD is listed as WB ACD-B in the short name and 27223 E002 RBC CPDA-1 is listed as RBC CP2D in the short name.	No identified risk/No identified impact.	All users	None available. The core conditions and expiration dates are correct for these products.	The short names include the incorrect anticoagulant listing. VBECS rules are not dependent on this information and enforce correct handling.	DR 2,703	Select Unit Project Proposal

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Select Units	Autologous, directed, or designated or dedicated blood units located in one division of a multidivisional database do not display in another division if the status is assigned or crossmatched.	Low risk /Low impact.	Limited to all users in consolidated Vista with multiple VBECS Divisions	None required.	Release the unit from assignment; the patient is no longer in the facility. Transfer any autologous or directed and assigned units with the patient. This assumes that the patient has been moved and physically in the other site while the restricted unit is actively selected for the patient. This requires a current specimen and its testing. This is a very unusual circumstance as autologous and directed units are generally associated with elective surgical procedures where a patient transfer is not common.	CR 2,050	Select Unit Project Proposal
Select Units	When there is a patient blood typing discrepancy and the unit does not meet requirements for selection, VBECS displays a warning with text that is not worded exactly as specified. System message displayed "VBECS- Not selectable unit [unit number]. This order must be emergency issued. This unit doesn't qualify for emergency issue and may not be selected."	No identified risk/No identified impact.	All users when an ABO discrepancy is found during the course of testing a selected patient.	None required.	VBECS is working correctly and enforcing the emergency issue compatibility rules and issuing universally compatible unit(s). ABO discrepancy is an uncommon occurrence in a transfusion service regardless of the root cause. The patient ABO/Rh discrepancy information is available in the Patient History or Exception Report. The message should read: "This patient had a previous ABO/Rh discrepancy and entry of a justified blood type with the following comment: <insert type of discrepancy from record>. Perform ABO/Rh on the current specimen to continue using normal rules and policies or follow emergency issue rules and policy as related to patient ABO/Rh retrieval and blood compatibility."	CR 1,548	Select Unit Project Proposal
Select Units	The Unit Expiration field is blank prior to the unit selection confirmation.	No identified risk/No identified impact.	All users	None required.	The unit expiration displays when a unit is selected. The expiration date is always displayed on the mock face label. The information is accurately displayed on the screen, just not in this specific field.	CR 1,899	Select Unit Project Proposal
Select Units	When a user selects an antigen positive unit for a patient with an antigen negative requirement; VBECS displays the product name instead of product type in the warning message.	No identified risk/No identified impact.	All users	None required.	The long product name contains more information than the product type for the unit.	CR 2,010	Select Unit Project Proposal
Select Units	VBECS displays a Codabar product code with an appended donation type code in the tree view of Select Units for Patients.	No identified risk/No identified impact.	All users	None required.	Codabar product codes do not require the appending of the donation type. This appended donation type information is correct for the unit.	CR 2,030	Select Unit Project Proposal

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Select Units	When the patient's blood type is unknown or inconclusive and the user uses the unit search option, VBECS displays available cryoprecipitate units for selection without grouping the units by blood type.	No identified risk/No identified impact.	All users	None required.	Available units are displayed correctly. Click the ABO/RH column header of the available unit list to reorganize the unit list by blood type, as desired.	CR 1,727	Select Unit Project Proposal
Select Units	When associating a specimen with a patient in the "Select Unit for a Patient" window; the "Expires" field within the "Associate with Specimen" panel, displays the time in AM/PM format.	No identified risk/No identified impact.	All users	VBECS normally displays date/time fields in the "military" time format. Workaround is to look in Maintain Specimen and check the expiration time military time format.	This is an inconsistent presentation which does not impede the user.	CR 2,120	Select Unit Project Proposal
Select Units	The display of Unit ID/Product Code for Codabar labeled units may display an ISBT 128 type code at the end of the number.	No identified risk/No identified impact.	All users	Select the tree view of the unit information where the full screen displays correctly.	None	CR 2,030	Select Unit Project Proposal
Select Units	The message displayed reads "Original and repeat ABO/Rh interpretation do not match." The 2nd sentence in the designed message "You must resolve the discrepancy before units can be issued." is not displayed.	No identified risk/No identified impact.	All users	None required.	VBECS will not allow the user to proceed with specimen processing until the discrepancy is resolved.	CR 2,203	Select Unit Project Proposal
Select Units	The Lab Order number does not display when the clipboard icon tool is selected.	No identified risk/No identified impact.	All users	User can view the Lab Order number using the Pending Order List or the Order History Report.	The lab order is available on those screens where it is needed. This is an additional display of the information that is not essential to the user.	CR 2,190	Select Unit Project Proposal
Select Units	When the patient's historic blood type from database conversion or previous testing is Rh Positive and the current specimen is unknown VBECS is not generating the message "Patient is Rh negative and the selected unit is Rh positive."	No identified risk/No identified impact.	All users when selecting an RH positive unit for a currently untested specimen where the database conversion blood type is Rh positive.	None required.	When the current testing is performed and the patient is Rh negative in conflict with their historic Rh information, VBECS does document the problem and alert the user. Issue of Rh Positive units to a patient where incomplete testing exists regardless of the historic typing does require override processing by the user at issue. This is a preliminary warning to the user, the mitigation is at Issue Blood Product.	CR 2,219	Select Unit Project Proposal

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Select Units	When a Rh positive unit is selected for a Rh neg patient with anti-D, the Enhanced Technologist does not receive the warning that unit is antigen positive.	No identified risk/No identified impact.	Technologists when selecting Rh positive blood products for an Rh negative patient.	None available. The warning message for selection of an Rh positive unit for an Rh negative patient does appear.	A Lead Technologist is required to issue this unit to the patient and an exception is collected. Issue Unit is handling this situation correctly as designed.	CR 2,406 DR 3,020	Select Unit Project Proposal
Select Units	A system error occurs when attempting to associate a specimen in Select Unit with an order when no unit is selected.	No identified risk/No identified impact.	All users, retrospective data entry when specimen has not been associated during Accept Order.	Ensure that a unit is selected when associating a specimen UID with units processed without a specimen for emergency issue.	Alternately, associate the specimen UID in Accept Order when possible. The user is attempting to associate the specimen to nothing. A unit must be selected. This is in place to allow users to associate a specimen with blood products for retrospective testing after the unit has been processed in emergency circumstances.	CR 2,411	Select Unit Project Proposal
Select Units	A system error occurs when the user clicks very quickly on the multiple messages in select unit regarding eXM eligibility and printing tags/forms.	No identified risk/No identified impact.	All users who are clicking faster than the messages can be read and responded to.	View, read, and respond to the presented messages.	This has been identified is a training issue and possible enhancement to consolidate redundant messages for the selected units. The user is "counting" and clicking ahead of VBECS message presentation. The original information is saved, but VBECS crashes due to the database constraint violation. The user is clicking faster than the messages can be read and responded to.	CR 2,501 HD 354125	Select Unit Project Proposal
Server System Administrator	Microsoft Windows limits the number of non-interactive processes on the server. On rare occasions when the limit is exceeded the VBECS patch process will fail. User will be notified by message: "Error: occurred during VBECS application lock processing on the server." Repeat the patch process.	No identified risk/ No identified impact.	System Administrator	Reboot prior to the report patch process.	This may occur during the patch process when users are currently off line due to the patch installation. The administrator would reboot the system. A slightly longer downtime may occur.	CR 2,234	Maintenance TBD
Server System Administrator	Locking: a system error occurred when the user opened the Invalidate Patient Test Results window.	No identified risk/ No identified impact.	All users	None available.	This is a random Microsoft bug with form loading. Inconsistent and difficult to reproduce.	CR 1,670	Maintenance TBD
Server System Administrator	A VBECS system error will occur when a database timeout occurs. A Database timeout error should be handled by the software.	Low risk /Low impact.	All users	The user should restart VBECS and continue where they left off.	If the database receives numerous queries or save requests simultaneously or there is a high-volume of network traffic, there is the possibility for a database timeout. This means that the database cannot respond to all of the requests in the allotted response time and VBECS error message is generated.	CR 2,113	Maintenance TBD

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Special Instructions & Transfusion Requirements	When opening the Special Instructions and Transfusion Requirements (SI and TR) option, there is no audible alert for patients with existing entries.	No identified risk/ No identified impact.	All users	None required.	The user is in the option proper, this was a misunderstanding of the requirement for an audible alert associated with the existence of SI and TR.	DR 2,152	Maintenance TBD
Special Instructions & Transfusion Requirements	After a patient has been selected and the SI/TR entry screen appears, the OK button remains enabled.	No identified risk/ No identified impact.	All users	None required.	None	CR 2,078	Maintenance TBD
Supplies: Log in Supplies	VBECS does not allow a user to select an expired supply item with an override during modification.	No identified risk/ No identified impact.	All users, when selecting an expired supply item to associate with a unit modification.	Do not use an expired supply item.	Standard of practice is to use only in-date supplies. This would be an unusual situation. The lot number can be entered with a different expiration date as a work around if in-date supplies must be used. Supplementary documentation related to the use of an expired supply are required, including a local risk assessment.	CR 1,807	Maintenance TBD
Supplies: Log in Supplies	Log-In Supplies Inventory list view does not sort when column headers are clicked.	No identified risk/ No identified impact.	All users when attempting to sort the list using the column headers.	None required.	None	CR 2,100	Maintenance TBD
Transfusion Complications Report	VBECS cannot be queried for data in a date range occurring prior to VBECS installation.	No identified risk/ No identified impact.	All users can print the report, but a Supervisor would be responsible for review.	The user must include a Vista report request for dates prior to VBECS installation.	The existing Vista report is used for data retrieval prior to the production use of VBECS. Vista blood bank reports remain available for retrieval.	CR 2,255	Report Tool Upgrade Project Proposal NSR 20080915

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Transfusion Effectiveness Report	Lab Test Name displays the VistA short name for the test.	No identified risk/ No identified impact.	All users can print the report, but a Supervisor would be responsible for review.	None required.	The user should be familiar with this test name. If they are not, they would query VistA.	CR 2,135	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Effectiveness Report	When requesting the Transfusion Effectiveness Report when no data is found the message "No information is available, per entered search criteria" displays.	No identified risk/ No identified impact.	All users can print the report, but a Supervisor would be responsible for review.	None required.	This is consistent with design.	DR 2,539	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Effectiveness Report	The report scheduler does not schedule prints on the same date.	No identified risk/ No identified impact.	All users can print the report, but a Supervisor would be responsible for review.	Select to print on a future day. For example, a user must queue a report to print data from the previous day.	The report is available for on demand printing. This prevents the user from creating a partial report for the date requested.	CR 2,136	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Reaction Count Report	The report only displays text entered in the "details" field of the Finalize TRW screen. The canned comments entered are suppressed on the report.	Low risk /Low impact.	All users can print the report, but a Supervisor would be responsible for review.	Enter all comments in the "details" field for them to print on the report.	None.	CR 2,860 HD 415291	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Reaction Count Report	Transfusion Reaction Count Rpt: Transfusion Information section: Patient. The patient prints twice in the Transfusion Information section on the report. Once with the correct treating specialty, (e.g., HEMATOLOGY/ONCOLOGY) and once with a blank treating specialty.	No identified risk/ No identified impact.	All users	None required.	Correct information is duplicated.	CR 2,282	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Reaction Count Report	Ordered and pending transfusion reaction workups are not included in the report.	No identified risk/ No identified impact.	All users can print the report but a Supervisor would be responsible for review.	None required.	These are incomplete tests. Only finalized TRW are included in the report.	CR 1,974	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Transfusion Reaction Count Report	Transfusion Reaction Count Report will include one or all Divisions for the report.	No identified risk/ No identified impact.	All users can print the report but a Supervisor would be responsible for review.	None required.	The report includes a single division or all of the divisions in the consolidated database. The option to select a subset of the consolidated divisions is not available.	DR 2,602	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Reaction Count Report	Hardcopy of Summary Report does not print "Page x of y" at bottom right page.	No identified risk/ No identified impact.	All users can print the report, but a Supervisor would be responsible for review.	View the Report on screen. "Page x of y" appears at the bottom right of page.	The page numeration is not centered	CR 2,216	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Reaction Count Report	Date Range uses the date the Transfusion Reaction Work Up (TRW) was entered, rather than the Date Reaction Noted as set in Patient Testing: Enter Transfusion Reaction Workup.	No identified risk/ Low impact.	All users can print the report, but a Supervisor would be responsible for review.	The search terms for the Transfusion Reaction Count Report use the last date the TRW was updated, and the "Date Reported" is found in the details of the report.	The TRW may or may not be included in the count if it was updated after the date reported and that is used for the report selection.	CR 2,220	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Reaction Count Report	The first page of the report is blank except for the header, title and section header.	No identified risk/ No identified impact.	All users can print the report but a Supervisor would be responsible for review.	None required.	The report information begins on the second page and is accurately displayed.	CR 2,221	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Reaction Count Report	Hovering a mouse over a report produces tool tips that are meaningless.	Low risk /Low impact.	All users can print the report, but a Supervisor would be responsible for review.	None required.	No incorrect information is displayed. The tool tips are related to Crystal Reports section names of the report. This does not interfere with the contents of the report.	CR 2,223	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Requirements Report	Transfusion Requirements Report does not work with the report scheduler.	No identified risk/ No identified impact.	All Users can print the report, but a Supervisor would be responsible for review.	Users must run the report when needed.	The report is available for on demand printing.	CR 2,257	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Transfusion Requirements Report	The partial report displays patient blood types and Transfusion Requirements. The cumulative report does not display ABO/Rh when there is no other requirement.	Low risk /Low impact.	All users can print the report but a Supervisor would be responsible for review.	To create a report with all patients' blood types and Transfusion Requirements, print a partial report with the date range from VBECS implementation through the current date.	The user must request the report using the workaround. If a cumulative report format is selected only patients with transfusion requirements and antibodies are displayed, not all patients on file.	CR 1,971	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Requirements Report	A report printed with an end date of today does not print "Preliminary" in the header.	No identified risk/ No identified impact.	All users can print the report but a Supervisor would be responsible for review.	Check the subsequent pages of the report and write "PRELIMINARY" on the front page of the report prior to filing.	The information included in the report is accurate to the time printed. There is a possibility of a gap of information should the user not overlap the print requests. Generally, this report is printed for the day prior or time period that does not include the current day.	CR 1,969 DR 2,584	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Requirements Report	The Transfusion Requirements report prints all divisions; no option to print selected divisions.	No identified risk/ No identified impact.	All users can print the report but a Supervisor would be responsible for review.	None required.	The patient files are available at this level to all divisions in the consolidated database. This is a disconnect between the written rule and the desired information handling.	DR 2,605	Report Tool Upgrade Project Proposal NSR 20080915
Transfusion Requirements Report	The Transfusion Requirements report will not print if the user selects the date VistA records were loaded into the VBECS database.	No identified risk/ No identified impact.	All users can print the report but a Supervisor would be responsible for review.	Access the VistA Blood Bank Reports option and use Patient Antibody Report (short list) to retrieve the list of transfusion requirements for the patients that were imported during conversion.	Aside from antigen negative requirements, other transfusion requirements are not set in VBECS during the database conversion.	CR 2,762	Report Tool Upgrade Project Proposal NSR 20080915
Transmit Workload Data	User will get the "No workload code defined. Workload credit cannot be applied to this transaction." warning message (if appropriate) when you cancel an order.	No identified risk/ No identified impact.	None	None required.	Setup requires that workload codes be assigned. Canceling an order has no specific workload associated with it. This does not impact billable tests.	CR 2,236	Maintenance TBD
Unit Antigen Typing	When a user is completing antigen typing results for a unit and the control results were previously saved a user can navigate to the QC comments field but is not allowed to save new comments.	No identified risk/ No identified impact.	All users	None required.	None.	CR 2,888	Testing Grids
Unit Antigen Typing	VBECS inserts the negative control vial number in the positive control result field if the user uses the down arrow key.	Low risk /Low impact.	All users	Navigate using the mouse to prevent this problem in the antigen typing test grid.	None.	CR 3,014	Testing Grids

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Unit Antigen Typing	If the database connection is lost in the milliseconds between testing result save and testing worklist update, the worklist may not be updated as Completed. If the worklist has been Completed but not updated, attempting to access it from the worklist listing will cause a VBECS system error, and the user will need to restart VBECS, and should not invalidate this worklist from the Unit Antigen typing function.	Very Low risk /Low impact.	All users	When invalidating a worklist, it is suggested that the user access the worklist to confirm that it is actually Incomplete before invalidating it.	Rare occurrence. Invalidations should be made in the Edit Unit Information function.	CR 2,109	Testing Grids
Unit Antigen Typing	When a user cancels out of an antigen testing worklist without performing any testing and then creates a new worklist for the same antiserum lot number; VBECS does not display the control cells on this worklist.	No identified risk/ No identified impact.	All users	When the intent is to recreate the worklist for the same antiserum lot number, perform the control test when it is displayed whether or not the selected units are tested.	The user invalidates this worklist and begins again.	CR 1,585	Testing Grids
Unit Antigen Typing	A system error occurs when a user enters reagent lot numbers, uses the Backspace key to erase one of multiple entries, and clicks OK.	Low risk /Low impact.	All users	Reenter the reagent lot numbers.	The user must log back into VBECS. This has been identified as a training issue. When the user highlights and retypes the lot number this does not occur.	CR 1,572	Testing Grids
Unit Antigen Typing	VBECS displays "Unit was issued and results indicate a possible patient conflict. Notify the physician of this potential problem IMMEDIATELY, per hospital policy!" when a unit was emergency issued prior to the completion of testing and an antigen negative requirement is subsequently identified. The expected message is "Units have been emergency issued. Notify the physician of this potential problem IMMEDIATELY, per hospital policy!"	No identified risk/ No identified impact.	All users	None required.	The intent of the message is the same. The message is grammatically different from the written business rule but is materially equivalent.	CR 1,575	Testing Grids
Unit Antigen Typing	VBECS does not display testing comments entered for the control cells when partially completed testing is recalled.	No identified risk/ No identified impact.	All users	Testing comments are properly displayed in the Testing Worklist Report.	None	CR 1,744	Testing Grids
Unit Antigen Typing	When a user applies a sort order when selecting units for antigen typing that order is not inherited by the testing grid.	No identified risk/ No identified impact.	All users	Organize physical tests materials according to the order of units presented on the testing grid.	The units sort order is based on the last update date/time which is not displayed to the user.	CR 2,646	Testing Grids

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Unit Antigen Typing	When a duplicate antigen test is performed on a unit that does not match the original, VBECS displays a message that does not exactly match the expected message.	No identified risk/ No identified impact.	All users	None required.	VBECS correctly displays the message for antigen typing discrepancies when the discrepancy is due to an update as well as login or testing. The message is grammatically different from the written business rule but is materially equivalent.	CR 1,830	Testing Grids
Unit Antigen Typing	Selecting a large number of units (~38) for antigen typing takes a long time for VBECS to process (~3 minutes). During that time the screen appears distorted.	Low risk /Low impact.	All users	None required. However, selecting fewer units for a batch avoids this problem.	The order of units does change from the selection order and requires verification.	CR 2,645	Testing Grids
Unit History Report	When a user retroactively updates a unit status through the Discard or Quarantine option, the Unit History Report displays the updated information, but does not display the date the change was made.	No identified risk/ No identified impact.	All users	None available.	None	DR 2,020	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	Unit Tests section of the Unit History Report incorrectly displays the QC results for unit antigen typing. The control cells display the unit test result not their own results.	No identified risk/ No identified impact.	All users	View or print the Testing Worklist Report, Miscellaneous Reagent QC section where the correct antigen typing QC results are displayed.	None	CR 2,639	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The Antigen Testing Phase Change Exception does not appear on the Unit History Report.	No identified risk/ No identified impact.	All users	The Antigen Testing Phase Change Exception appears on the Exception Report which should be reviewed daily by the supervisor.	None	CR 2,643	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The last digit of the minute printed for the Date Processed field is difficult to read.	No identified risk/ No identified impact.	All users	None available.	None	CR 2,545	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The Unit History report Costs column displays the incorrect cost if units are pooled.	No identified risk/ No identified impact.	All users	See the Cost Accounting Report for cost of Splits/Divides and Pools.	None	CR 2,158	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	If a unit is entered in inventory, shipped via Outgoing Shipment, and reentered in inventory via Incoming Shipment; the Unit History Report includes incomplete Outgoing Shipment information on the Incoming Shipment page.	No identified risk/ No identified impact.	All users	The correct information appears on the Unit History Report but the duplicated display of outgoing shipment information may require investigation.	None	CR 1,607	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	When a unit is added, the Unit History Report for a pooled unit displays duplicate targets for that unit.	No identified risk/ No identified impact.	All users	None required.	Although duplicated, the correct modification information is included in the report.	CR 1,880	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Unit History Report	The Unit History report does not indicate the CMV status correctly if a unit was re-entered into VBECS without a CMV negative status. This only occurs if the unit was originally entered with the CMV negative status.	Very Low risk / Very Low impact.	All users when re-entering a CMV negative unit.	None required.	VBECS correctly retrieves the unit's original information and considers the unit CMV negative for selection and issue. If the unit CMV testing status changed and is no longer considered negative at the time of reentry, clear the CMV negative special testing check box in Edit Unit Information.	CR 1,932	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	VBECS displays the five-digit blood product code in the Report Criteria section of the report.	No identified risk/ No identified impact.	All users	None required.	The full product code is displayed properly in the header of the report.	CR 1,989	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The Crossmatch Tests section of the Unit History Report will display testing without interpretations if the unit selected was tested in a batch of units and results were not entered.	No identified risk/ No identified impact.	All users	Enter crossmatch testing results for all units selected for crossmatch.	The Unit History Report correctly displays crossmatch results when the testing is complete or if the testing fields for the unit are left blank and not edited during testing.	CR 2,716	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	VBECS displays the user ID rather than the user name in the Transfusion Information Processed By column.	No identified risk/ No identified impact.	All users	None required.	The user is identifiable though this is not the usual presentation.	CR 1,943	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The exception detail captured for the exception type Previously Recorded Results Invalidated do not appear on the Unit History Report.	Low risk/Low impact.	All users	View the Exception Details on the Exception Report for this information for the date the exception was generated.	The Exception report is identified as one of the reports requiring supervisory review. The impact to the antigen typing result is investigated and addressed at that time not as a result from an entry on the Unit History Report.	CR 2,682	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The exception details captured for Antigen Testing Phase Change do not appear on the Unit History Report.	Low risk/Low impact.	All users	View the Exception Details on the Exception Report for this information for the date the exception was generated.	The Exception report is identified as one of the reports requiring supervisory review. The impact to the antigen typing result is investigated and addressed at that time not as a result from an entry on the Unit History Report.	CR 2,689	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The column header reads "Crossmatched By" instead of "Released By."	No identified risk/ No identified impact.	All users	Information displayed in this column is clearly unit release information when taken in the context of the entire report.	This is misleading as a header but follows the issue information so is logically located for the user.	CR 2,008	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	Unit History Report section of the report has the word "antigen" spelled incorrectly.	No identified risk/ No identified impact.	All users	None required.	There is a typographical error.	CR 2,114	Report Tool Upgrade Project Proposal NSR 20080915

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Unit History Report	In the Unit History Report, Patient Association section, the Crossmatched to Patient Name and Assigned to Patient Name may not have a space between the first and middle name depending on the length of the patient's names.	No identified risk/ No identified impact.	All users	None required.	There is a typographical error.	CR 2,116	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	When there are inactivated blood units in the VBECS system, an additional redundant Select Units window appears when attempting to create a Unit History Report.	No identified risk/ No identified impact.	All users	None required.	The second window must be closed to proceed.	CR 2,084	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	Tested By column is blank for Unit Antigen Typing.	No identified risk/ No identified impact.	All users	See the Testing Worklist Report for the tested by user name.	This is a secondary report for this information.	CR 2,180	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	Antigen typing comments for QC do not appear on the Unit History Report. The testing comments entered for a unit appear for the positive and negative controls instead.	No identified risk/ No identified impact.	All users	Testing comments for the antigen typing QC can be viewed as entered on the Testing Worklist Report, Miscellaneous Reagent QC section.	None.	CR 2,641	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	Antigen typing interpretations from Incoming Shipment, Edit Unit Information, or testing are presented as the antigen status only on pick lists throughout VBECS.	No identified risk/ No identified impact.	All users	None required.	Unit testing information is found on the Unit History Report and the Testing Worklist Report. This relates to pick lists not the blood availability report which would be used to search for type specific units. Assigning units to match patient requirements is not impacted.	DR 1,640	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	The Unit History Report does not display when a unit has been Returned to Shipper or Shipped to Another Facility.	No identified risk/ No identified impact.	All users	The Cost Accounting Report contains this information.	This is a secondary report for this information.	DR 1,641	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	Unit History Report: Transfusion Information section: Patient. The patient prints twice in the Transfusion Information section. Once with the correct treating specialty, (e.g., HEMATOLOGY/ONCOLOGY) and once with a blank treating specialty.	No identified risk/ No identified impact.	All users	None required.	Correct information is duplicated.	CR 2,281	Report Tool Upgrade Project Proposal NSR 20080915
Unit History Report	Pooled ISBT units of differing ABO/Rh types appear on the Unit History Report as "Mx" rather than "Pooled" designation for the ABO/Rh.	Low risk/Low impact. .	All users	None required.	This is a display only; there is no system-generated blood product label. The nomenclature is described in user documentation.	CR 1,628	Report Tool Upgrade Project Proposal NSR 20080915

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
Update Reagent Inventory	The list of selectable comments in Reason for Change field is not associated with a canned comment category.	No identified risk/No identified impact	None	Select a comment or OTHER and enter a free text comment to enter additional information.	These appear to be part of the Reagent and Supply comment category but are not part of that list.	CR 2,752	Maintenance TBD
VBECS Administrator	Active Directory name changes cannot be imported into VBECS.	Low risk/Low impact.	All users	There is no way for a system administrator to update the database record with the new name. A name update requires a manual change to the database by development.	None.	DR 3,808 HD 415139	Maintenance TBD
VBECS Administrator	Users are prompted with a VBECS Confirmation box when exiting the VBECS Administrator that asks if the user wants to abandon changes and continue when no change was made.	No identified risk/ No identified impact	System Administrator	Click Yes on the message to exit and continue. VBECS does not alter the records that were not changed.	Clicking No on the message will remove the message and leave the user in the VBECS Administrator.	CR 2,841	Maintenance TBD
VBECS Administrator	When the user attempts to exit the VBECS Administrator by clicking the X button on the window it will not close. The user then clicks Cancel on the Vista Logon window and a system error occurs.	No identified risk/ No identified impact	System Administrator	Click on File then Exit to close the application.	None.	CR 2,862	Maintenance TBD
VBECS Administrator	Facility Name Field does not allow a full view of the registration number and full name of the selected facility.	No identified risk/ No identified impact	System Administrator	None required.	The user can read the FDA registration number which is the unique identifier for the facility. The FDA registration number is not editable where the Name may be edited. The FDA registration number is the unique identifier for a facility.	CR 2,340	Maintenance TBD
VBECS Administrator	System Administrator can save the Patient Update and Patient Merge tabs for the VBECS Interface without entering a value in the Facility ID field. Failure to enter a facility ID will prevent patient updates and merge alerts from being processed.	Low risk/Low impact.	System Administrator	Enter the facility ID in the Facility ID field on the Patient Update and Patient Merge tabs of VBECS Administrator.	VBECS System Administrator is instructed to enter the facility ID in the Facility ID field in the VBECS Technical Manual-Security Guide when VBECS is configured for use.	CR 2,573	Maintenance TBD
VBECS Administrator	Clicking the Clear button after saving any changes still displays a warning about discarding changes that were made to the division details.	Low risk /Low impact.	System Administrator	None required.	The user must reject this errant message and save their changes. The user has the option to reject the message and save the details.	CR 2,344	Maintenance TBD
VBECS Administrator	A system error occurs in the VBECS Administrator when the system administrator edits a user of VBECS and minimizes the VBECS Administrator window.	No identified risk/No identified impact	System Administrator	Restart VBECS Administrator.	None.	CR 2,616	Maintenance TBD
VBECS Administrator	The tool tip message for the disabled Clear button is incorrect. It states "Click to erase changes and reset screen." Expected to state "The data on the screen were not changed".	No identified risk/No identified impact	System Administrator	None required.	The message is incorrect but is rarely viewed by a user. The button is disabled and cannot be selected in either case.	CR 2,348	Maintenance TBD

*CR (Code Request Tracking System Number), DR (Document Request Tracking System Number), HD (Remedy Help Desk Ticket Number), CQ (CPRS ClearQuest Tracking System Number)

Option Where the Issue Occurs	Description	Risk Assessment /Impact to Patient Care	Security Role Mitigations (Affected User)	Recommended Workaround	Additional Comments	Tracking System Number*	Tentative Schedule
VBECS Administrator	A user can configure a VistA Institution as a division and an associated institution. Orders placed at that CBOC to be rejected at acknowledgment. The user is alerted to the error in configuration by this message: "Unable to find valid Associated Institutions information. Please check configuration."	Low risk /No identified impact.	System Administrator	None required.	This is corrected during configuration verification during validation testing.	CR 2,376	Maintenance TBD
Workload Codes (Division Configuration)	Invalidating split units will result in negative workload equal to the number of splits created applied to the total number of units logged in.	No identified risk/ No identified impact.	All users When split modification is enabled at the facility.	None available.	This occurs when split units are incorrectly processed in VBECS.	CR 2,224	Maintenance TBD
Workload Codes (Division Configuration)	Workload totals for VBECS processes may not match the VistA workload report totals. Repeat orders are known only to VBECS so if workload is generated as a result of a repeat order test the VBECS workload report will include that workload in its totals but the VistA workload report will not.	No identified risk /Low impact.	All users	None required.	Compare the workload code number totals that are equivalent.	CR 2,240	Maintenance TBD
Workload Codes (Division Configuration)	LMIP/NLT associated CPT codes are not changed by checking or unchecking the boxes.	No identified risk/ No identified impact.	Administrator	Adjust the LMIP/NLT code to CPT code associations in VistA Lab.	This has been identified as a training issue.	CR 2,438	Maintenance TBD

Revision History

Date	Revision	Description	Author
7-21-06	1.0	Initial Version	BBM team
6-22-07	2.0	<p>Deleted CRs 1,051, 1,373, 1,441, 1,625, 1,691, 1,709, 1,736, 1,756, 1,791, 1,801, 1,804, 1,826, 1,847, 1,855, 1,859, 1,922, 1,939, 1,957, 1,982, 1,983, 1,986, 1,987, 1,988, 2,031, 2,057, 2,060, and 2,061.</p> <p>Deleted CR 457 and CR 1,512, included in error in the 1.0 version of this document.</p> <p>Added CRs 1,858, 2,022, 2,030, 2,070, 2,074, 2,076, 2,078, 2,083, 2,084, 2,085, 2,086, 2,089, 2,092, 2,093, 2,096, 2,098, 2,100, 2,101, 2,103, 2,105, 2,106, 2,107, 2,109, 2,110, 2,111, 2,114, 2,115, 2,116, 2,119, 2,120, 2,123, 2,129, 2,130, 2,131, 2,133, 2,134, 2,135, 2,136, 2,145, 2,148, 2,149, 2,153, 2,158, 2,159, 2,168, 2,177, 2,179, 2,180, 2,181, 2,185, 2,186, 2,188, 2,200, 2,201, 2,202, 2,203, 2,209, 2,215, 2,216, 2,217, 2,218, 2,219, 2,220, 2,221, 2,223, 2,224, 2,225, 2,226, 2,227, 2,230, 2,231, 2,233, 2,236, 2,237, 2,240, 2,241, 2,242, 2,243, 2,244, 2,247, 2,248, 2,251, 2,252, 2,253, and 2,254, and 2,255.</p> <p>Added DRs 1,230, 1,312, 1,322, 1,625, 1,632, 1,633, 1,640, 1,641, 1,643, 1,651, 1,661, 1,703, 1,704, 1,707, 1,729, 1,942, 1,956, 1,967, 1,984, 2,020, 2,072, 2,133, 2,135, 2,140, 2,152, 2,206, 2,218, 2,255, 2,323, 2,479, 2,536, 2,539, 2,257, 2,580, 2,584, 2,586, 2,602, 2,605, 2,616, 2,618, 2,681, and 2,703.</p>	BBM team
2/25/08	3.0	<p>Updated for the release of VBECS 1.2.0.0:</p> <p>Combined Known Defects and Anomalies section into one and added a column to identify anomaly vs. defect.</p> <p>Added column for CR/DR numbers. Separated CR and DR numbers from the description.</p> <p>Added a column to separate workaround and comments.</p> <p>Changed all No workaround available to None available as the column only contains workarounds.</p> <p>Consolidated Configure Initial Division Parameters, Configure Interfaces, Configure System Administrators and Define Interface Control Parameters under VBECS Administrator section.</p> <p>Consolidated CR 2,115 with throughout VBECS audible alert.</p> <p>Deleted CRs 1,352, 1,805, 1,844, 1,874, 1,904, 1,947, 1,977, 1,979, 2,009, 2,091, 2,110, 2,134, 2,155, 2,251, 2,186, 2,188, 2,248, 2,253, 1,857, 2,356, 1,533, 1,883, 1,109, 1,856, 1,870, 1,161, 1,602, 1,651, 1,858.</p> <p>Deleted DRs 2,135, 2,255, 2,257.</p> <p>Added CRs 1,828, 2,065, 2,095, 2,113, 2,118, 2,190, 2,193, 2,195, 2,199, 2,208, 2,233, 2,234, 2,257, 2,259, 2,261, 2,262, 2,263, 2,270, 2,273, 2,281, 2,282, 2,286, 2,295, 2,309, 2,310, 2,318, 2,330, 2,337, 2,340, 2,344, 2,347, 2,348, 2,353, 2,354, 2,356, 2,359, 2,361, 2,363, 2,366, 2,367, 2,368, 2,370, 2,372, 2,376, 2,385, 2,387.</p> <p>Added DRs 1,707, 2,672, 2,679, 2,690, 2,735, 2,821, 2,851, 2,881, 2,883, 2,891, 2,937, and 2,982.</p> <p>Made multiple grammatic and format changes throughout the document.</p> <p>Edited workaround text in CR 1,824, CR 2,334, CR 2,353, CR 1,949, CR 2,185, CR 1,976, CR 2,286, CR 2,337, CR 2,387, CR 2,363, CR 1,902 and CR 1,950.</p> <p>Combined CR 2,353 and CR1,824</p> <p>Revised the description of CR 2, 334, CR 2,310, CR 2,131, CR 1,897, DR 2,536, CR 2,333, DR 2,881, CR 2,361, CR 1,797, CR 2,053, CR 2,070, CR 1,515, CR 1,548, CR 2,120, CR 2,041, CR 2,083, CR 1,961, CR 1,963, CR 2,003, CR 2,051, CR 2,220</p> <p>Revised all DR 2,891, DR 2,681, DR 2,133 CR 1,853, DR 1,643, CR 2,032,</p> <p>Moved CR 2130 to the section "CPRS".</p> <p>Moved CR 2042 to the section "Record Patient ABO/Rh". Changed the heading of the section from Record Patient ABO/Rh to "Patient Testing: Record Patient ABO/Rh". Moved the section to Patient Testing section.</p>	BBM team

Date	Revision	Description	Author
		Moved DR 1,322 to the Incoming Shipment section. Separated the comment from the workaround for CR 2148, DR 1,230 Combined CR 2,240 and CR 2,244 Combined Special Instructions/Transfusion Requirements sections. Removed comment from CR 1,890, CR 2,085 Added comment to CR 2,359, CR 2,254, CR 1,846 “Not Related to the User Interface” section changed to “Server System Administrator”.	
4/01/08	4.0	Updated for the release of VBECS 1.3.0.0. Added CR 2,384, CR 2,387, CR 2,396, CR 2,397, CR 2,401, CR 2,402, CR 2,406, CR 2,408, CR 2,410, CR 2,411.	BBM Team
4/22/08	5.0	Updated for the release of VBECS 1.3.0.0. Added CR 2,419, CR 2,422, DR 3,016, DR 3,041, DR 3,046 Edited DR 2,891 workaround and comments.	BBM Team
6/17/08	6.0	Updated for the release of VBECS 1.4.0.0. Added CR 2,417, CR 2,433, CR 2,435, CR 2,436, CR 2,437, CR 2,438, CR 2,443, CR 2,448 and CR 2,454. Edited CR 2,422. Removed CR 1,904, CR 2,419, DR 3,016, CR 2,130, CR 2,237	BBM Team
7/11/08	7.0	Updated for the release of VBECS 1.4.0.0. Added DR 3,124.	BBM Team
7/17/08	8.0	Updated for the release of VBECS 1.4.0.0: Added CR 2,462.	BBM Team
9/19/08	9.0	Updated for the full install release: Added CR 2,421, CR 2,451, CR 2,455, CR 2,456, CR 2,459, CR 2,461, CR 2,463, CR 2,466, CR 2,467, CR 2,483, CR 2,484, CR 2,485, CR 2,490, CR 2,494, CR 2,500, CR 2,501 CR 2,385, CR 2,489, CR 2,499, CR 2,502, CR 2,425, CR 2,105. Added DR 3,098. Removed DR 2,618, DR 3,124, DR 3,041 Added the hyphen to Post-Transfusion for consistency	BBM Team
1/13/09	10.0	Added additional introduction. Per feedback from the Patient Safety Office review: added columns for Risk Assessment/Impact to Patient Care, Security Role Mitigations (Affected User), and Tentative Schedule for correction. Moved explanatory text from the workaround column to the comments column and added additional comments. Options put in alphabetical order. CR 1,514, moved from Server System Administrator to Throughout VBECS section. CR 2,262, Recommended Workaround, removed space before “When”. CR 1,212, Additional Comments, changed “The use cannot” to “The user cannot” CR 2,092, Security Role Mitigations, Reworded. CR 1,628, Additional Comments, added a period. CR 1,925, Additional Comments, removed space before “Requires” DR 2,020, Description, changed “of” to “or”. CR 1,501, and CR 1,514, Additional Comments, added the hyphen to “non-patient”. Removed DR 2,347, and 2,672.	BBM Team
9/1/09	11.0	Updated for the release of VBECS 1.5.0.0: Added by section: Throughout VBECS: CR 2,536, CR 2,657, CR 2,658, CR 2,708, CR 2,715 ABO/Rh Confirmation: CR 2,707	BBM Team

Date	Revision	Description	Author
		<p>Accept Orders: Pending Order List: CR 2,622 Accept Orders: Accept an Order: CR 2,712 and CQ#18690 Administrative Data Report : DR 3,539 Audit Trail Report : CR 2,614, CR 2,714 and HD 345870 Configure Daily QC: CR 2,690, HD 335126, CR 2,696, added HD 249548 to CR 2,385 CPRS: DR 3,510 CPRS CQ#18494, CR 2,713 and CQ#18690 Discard or Quarantine: CR 2,706 Enter Daily QC Results: added HD 269544 to CR 2,437 Exception Report: CR 2,642, CR 2,686 Finalize / Print TRW: added HD 332336 and DR 2,084 to DR 1,633 Incoming Shipment: DR 3,542 HD 339722 Invalidate Test Results: CR 2,579, CR 2,674 Issued/Returned Units Report: CR 2,512, CR 2,596 Maintain Minimum Levels: CR 2,599 HD 301562 Modify Units: Pool Units: CR 2,650, CR 2,691 Modify Units (not Pool or Split): CR 2,612, CR 2,656, CR 2,678 Order History Report: CR 2,584 Patient History Report: CR 2,610, CR 2,684 Patient Testing: CR 2,576 HD 288163, CR 2,626, CR 2,685 Patient Testing: Pending Task List: CR 2,490 HD 264904 Patient Testing Worklist and Testing Worklist Reports: CR 2,627, CR 2,694 HD 336190, DR 3,556 Patient Updates: CR 2,580, CR 2,711 Post-Transfusion Information: CR 2,550, CR 2,589, CR 2,607 Print Unit Caution Tag & Transfusion Record Form: CR 2,493, CR 2,624, edited CR 2,454 Reagents: CR 2,528 Transfusion Complications Report: added HD 337715 to CR 2,252 Transfusion Effectiveness Report: CR 2,695, HD 335307 Transfusion Requirements Report: CR 2,581 Unit Antigen Typing: CR 2,646 Unit History Report: CR 2,545, CR 2,639, CR 2,641, CR 2,643, CR 2,682, CR 2,689, CR 2,716 VBECS Administrator: CR 2,573, CR 2,616 Workload Codes (Division Configuration): CR 2,438, CR 2,700 HD 339277</p> <p>Removed by section: Throughout VBECS: CR 2,065, CR 2,443, CR 2,467 Accept Orders: Pending Order List: CR 1,828, CR 2,095 Accept Orders: Accept an Order: CR 1,912 Blood Availability Report: CR 2,006 Discard or Quarantine: CR 2,417 Document ABO Incompatible Transfusions: CR 2,456 Edit Unit Information: CR 2,330 Incoming Shipment: CR 2,466, CR 2,459, DR 1,322, DR 3,098</p>	

Date	Revision	Description	Author
		<p>Modify Units: Pool Units: CR 1,976 Patient History Report: CR 2,433 Patient Testing Worklist and Testing Worklist Reports: CR 2,025 Patient Testing: Antibody Identification: CR 2,489 Print Unit Caution Tag & Transfusion Record Form: CR 2,483 Post-Transfusion Information: CR 1,653, CR 1,654, CR 1,961, CR 1,963, CR 2,461, CR 2,500 Select Units: CR 2,502 Unit Antigen Typing: CR 1,984 VBECS Administrator: CR 1,902, CR 2,366, CR 2,451, CR 2,259, CR 2,462 Workload Codes (Division Configuration): CR 2,425</p> <p>Added Printing Requirements section to address DR 3,390</p> <p>In addition, the Tentative Schedule column was updated to reflect current project proposals.</p>	
4/7/10	12.0	<p>Updated for the release of VBECS 1.5.0.0 after field testing:</p> <ul style="list-style-type: none"> The Known Defects and Anomalies Addendum for 10.0, Version 1.0 included the updates made in Known Defects and Anomalies, Version 11.0 for VBECS 1.5.0.0. The Known Defects and Anomalies, Version 12.0 incorporates defects and anomalies previously released in the Known Defects and Anomalies Addendum for 10.0 Version 1.0 along with other newly identified defects and anomalies. <p>How the Known Defects and Anomalies is Organized, fourth bullet was updated to explain that the highlighting now indicates changes from the last released version of the Known Defects and Anomalies.</p> <p>Added by section:</p> <p>Throughout VBECS: CR 2,757, CR 2,768 Accept Orders: Pending Order List: CR 2,784 Audit Trail Report : CR 2,730, CR 2,766 Blood Availability Report: CR 2,793 CPRS: DR 3,616 CQ 18842 HD 364290, DR 3,617 CQ 18841, CR 2742, CR 2,776 HD 381217, CR 2,782 HD 382825 HD 383708 CQ 19139, DR 2,883 C:T Ratio Report: CR 2,748, CR 2,780 Incoming Shipment: DR 3,605, DR 3,607, DR 3,610 HD 360751, DR 3,683, CR 2,759 Issue Blood Components: CR 2,758 Issued/Returned Report: CR 2,736 Modify Units (not Pool or Split): CR 2,587, CR 2,735 Modify Units: Split a Unit: CR 2,756 Order History Report: CR 2,778 HD 383473 Outgoing Shipment: CR 2,722, CR 2,734 Patient Testing: DR 3,643, CR 2,777 Patient Updates: CR 2,709 Post Transfusion Information: CR 2,761 Recent Orders: CR 2,785 HD 388392 Release Units From Patient Assignment: CR 2,769</p>	BBM Team

Date	Revision	Description	Author
		<p>Transfusion Requirements Report: CR 2,762 Unit Antigen Typing: CR 2,645, CR 2,724 HD 353065, CR 2,754, CR 2,764 HD 377080 Update Reagent Inventory: Added CR 2,752 Workload Report to Austin: DR 3,601, CR 2,737</p> <p>Removed by section: Accept Orders: Pending Order List: DR 1,967 Accept Orders, Accept an Order: CQ 18690 CPRS: CR 2,713, CQ 18494 Patient Testing: Pending Task List: DR 1,967 Unit Antigen Typing: CR 2,484</p> <p>Moved by section: Post Transfusion Information: Removed CR 2,589 from this section and added it to the CPRS section. Reworded the Additional Comments. Configure Daily QC: Removed CR 2,696 from this section and added it to the Enter Daily QC Results section.</p> <p>Edits by section: Updated revision history of version 11.0 of this document by adding reference to changes made for DR 3,390 (DR 3,670) CPRS: Added 1.6.0.0 to the tentative schedule for DR 2,821 Modify Units: Pool Units: Changed tentative schedule for DR 2,133 to 1.6.0.0 Removed '#' from CQ numbers Server System Administrator: Changed tentative schedule for CR 2,273 to 1.5.1.0 Revised DR 3,539, CR 2,614, CR 2,385, CR 2,642, DR 1,633, CR 2,691</p>	
12/1/11	13.0	<p>Updated for the release of VBECS 1.6.0.</p> <p>Added by section: Throughout VBECS: CR 2,870, CR 2,871, CR 2,872, CR 2,873, CR 2,855, CR 2,992 ABO/Rh Confirmation: CR 2,813, HD 399413, CR 2,830, HD 400511. CR 2,849 Accept Orders: Accept an Order: CR 2,903, CR 2,919, CR 3,063, DR 3,962 Administrative Data Report: CR 2,838, CR 3,079 Blood Availability Report: CR 3,132 Configure Daily QC: CR 3,066 Cost Accounting Report: DR 3,792, HD 399931 CPRS: CR 3,143, HD 461287, CR 3,184, HD 532136 Discard or Quarantine: CR 3,048 Division Transfusion Report: CR 3,170 Document ABO Incompatible Transfusion: CR 3,062 Enter Daily QC Results: CR 3,074, CR 3,188 Exception Report: CR 2,924, DR 3,832, DR 4,008 Finalize/Print TRW: CR 2,852, CR 2,867, CR 2,804, HD 381573, CR 2,815 Incoming Shipment: CR 2,809, CR 2,933, CR 3,061, DR 3,610, CR 3,131, HD 515950, CR 3,191, HD 538535 Maintain Specimen: CR 2,802, HD 394421</p>	BBM Team

Date	Revision	Description	Author
12/1/11	13.0	<p>Exception Report: CR 2,924, DR 3,832, DR 4,008 Finalize/Print TRW: CR 2,852, CR 2,867, CR 2,804, HD 381573, CR 2,815 Incoming Shipment: CR 2,809, CR 2,933, CR 3,061, DR 3,610, CR 3,131, HD 515950, CR 3,191, HD 538535 Maintain Specimen: CR 2,802, HD 394421 Modify Units: Pool Units: CR 2,875, CR 3,135 Order History Report: CR 2,914, CR 2,989, HD 465511 Order Reflex Tests: DR 3,643 Patient History Report: CR 2,886, CR 3,007, HD 479169 Patient Merge: CR 3,060, CR 3,099, MVI CR 232 Patient Testing: CR 2,812, HD 395766, CR 3,118 Patient Testing: Pending Task List: CR 2,874 Patient Testing Worklist and Testing Worklist Reports: CR 2,868 Post Transfusion Information: CR 2,879, HD 426258 Print Unit Caution Tag & Transfusion Record Form: CR 2,932, CR 2,922 Prolonged Transfusion Time Report: CR 2,803, CR 3,100 Reagents: CR 2,810, CR 3,133, HD 514507, DR 4,346, HD 531145 Select Units: CR 2,821 Special Instructions & Transfusion Requirements: CR 2,998 Transfusion Reaction Count Report: CR 2,860, HD 415291 Unit Antigen Typing: CR 2,888, CR 3,014 Unit History Report: CR 3,192, HD 538535 VBECS Administrator: CR 2,841, DR 3,808, HD 415139, CR 2,862</p> <p>Removed by section: Throughout VBECS: CR 1,174, CR 2,536, CR 2,033, CR 1,893, CR 2,045, CR 2,038, CR 2,014, CR 2,017, CR 2,019, CR 2,115, CR 2,193, CR 2,195, CR 2,199, CR 2,708 ABO/Rh Confirmation: CR 2,707 Accept Orders: Accept an Order: CR 2,712 Accept Orders: Pending Order List: CR 2,784 Administrative Data Report: CR 2,103, CR 2,105 Audit Trail Report: CR 1,671, CR 2,614, CR 2,714, HD 345870 CPRS: CPRS CQ 18841, DR 3,617, CR 2,782, HD 382825, HD 383708, CPRS CQ 19139, CPRS CQ 18493, DR 3,616, CPRS CQ 18842, HD 364290, DR 2,821, CR 2,368, CR 2,742, CR 2,435, CR 2,589, DR 2,883 Discard or Quarantine: CR 2,706 Enter Daily QC Results: CR 2,696 Exception Report: CR 1,897 Finalize/Print TRW: DR 2,084, DR 2,681 Free Directed Unit For Crossover: DR 2,735 Incoming Shipment: DR 3,683, CR 1212, CR 1,536, CR 1,591, CR 1,610, CR 1,719, CR 1,895, CR 1,935, DR 2,536, DR 3,542, HD 339722, DR 3,605, DR 3,607, CR 2,759 Invalidate Test Results: CR 2,674 Issue Blood Components: CR 2,758</p>	BBM Team

Date	Revision	Description	Author
12/1/11	13.0	<p>Modify Units (not Pool or Split): CR 2,656, DR 1,632, DR 1,942, DR 2,616, CR 2,587, CR 2,678, CR 2,612, CR 2,735 Modify Units: Pool Units: DR 2,133, CR 1,594, CR 2,691 Modify Units: Split a Unit: CR 2,756 Order Reflex Tests: CR 1,843 Outgoing Shipment: CR 2,722 Patient Testing: CR 2,576, HD 288163, CR 2,777 Patient Testing: Pending Task List: CR 2,056 Patient Testing Worklist and Testing Worklist Reports: HD 336190 Post-Transfusion Information: CR 1,656, DR 2,891, CR 2,761, CR 2,372 Print Unit Caution Tag & Transfusion Record Form: CR 2,454. CR 2,493 Select Units: DR 1,661, DR 1,707 Server System Administrator: CR 2,273, CR 1,926 Transfusion Complications Report: CR 2,247, CR 2,252, HD 337715 Transfusion Effectiveness Report: CR 2,695, HD 335307 Transfusion Reaction Count Report: DR 2,690 Transfusion Requirements Report: CR 2,581 Unit Antigen Typing: CR 1,574, CR 2,764, HD 377080, CR 2,724, HD 353065, CR 2,754 VBECS Administrator: CR 2,370, DR 2,937, CR 1,990, CR 2,410 Workload Codes (Division Configuration): CR 2,111, CR 2,244, CR 2,700, HD 339277 Workload Report to Austin: CR 2,737, DR 3,601</p> <p>Moved by section: Modify Units: Pool Units: Moved CR 1,628 to Unit History Report section.</p> <p>Edits by section: Throughout VBECS: Added CR 2,927 to CR 1,901, Added HD 372543 to CR 2,757 Audit Trail Report: Added CR 2,869 to CR 2,766 and reworded description Blood Availability Report: Changed tentative schedule for CR 2,463 to Report Tool Upgrade Project Proposal Configure Daily QC: Changed tentative schedule for CR 2,385, CR 2,436 to Report Tool Upgrade Project Proposal C:T Ration Report: Added HD 421286 to CR 2,748 Enter Daily QC Results: Added HD 367084 to CR 1,911 Finalize/Print TRW: Changed tentative schedule for CR 2,230 to Maintenance TBD, Changed tentative schedule for CR 2,455 to Maintenance TBD Incoming Shipment: Added HD 360751 to DR 3,610, Renumbered CR 2,401 to CR 3,160. Invalidate Test Results: Renamed DR 3,046 to CR 3,190 Issue Blood Components: Changed tentative schedule for CR 2,242 to Maintenance TBD Issued/Returned Units Report: Reworded description of CR 2,596 (DR 3953) and added HD 443058 Maintain Minimum Levels: Changed tentative schedule for CR 2,233 to Maintenance TBD Maintain Specimen: Added HD 382525 to CR 2,485 Modify Units: Pool Units: Reworded Description and Additional Comments for CR 2,650 Patient Information Toolbar: Reworded CR 1,927 (DR 3,952) Post Transfusion Information: Changed tentative schedule for CR 1,655, CR 2,607, CR 1,936, CR 2,263, CR 2,208, CR 2,879 to Maintenance TBD Reagents: Added HD 359714 to CR 1,958 and CR 1,931 Select Units: Added HD 354125 to CR 2,501</p>	BBM Team

Date	Revision	Description	Author
12/1/11	13.0	Transfusion Reaction Count Report: Added CR 2,282 from Unit History Report (Task 1286) Unit History Report: Reworded CR 2,281 and removed CR 2,282	BBM Team
6/27/12	14.0	<p>Updated for the release of VBECS 1.6.0.</p> <p>Added by section: Throughout VBECS: CR 3,229, HD 625192 Accept Orders: Accept an Order: CR 3,199, HD 540813 Blood Availability Report: CR 3,234 Finalize/Print TRW: CR 3,200 Invalidate Test Results: CR 3,240, HD 720879 Issue Blood Components: CR 3,208, HD 590527 Patient History Report: CR 3,204 Recent Orders: CR 3,241 Select Units: CR 3,227</p> <p>Removed by section: Print Unit Caution Tag & Transfusion Record Form: DR 1,704, CR 2,932 Recent Orders: CR 2,785, HD 388392 Outgoing Shipment: CR 2,041 & CR 3,192</p> <p>Edits by section: Incoming Shipment CR 3,191: Updated Additional Comments Outgoing Shipment CR 2,734: Updated Description</p>	BBM Team
		<p>Updated for the release of VBECS 1.6.0. Moved the revision history to the end of the document.</p> <p>Added by section: Throughout VBECS: CR 3,243, HD 626467 Enter Daily QC Results: CR 3,176 Modify Units: Split a Unit: CR 3,248</p> <p>Removed by section: Special Instructions & Transfusion Requirements: CR 2,998 Modify Units: Pool Units: CR 3,135</p> <p>Edits by section: Patient Updates: Added CR 3,249 to CR 2,580</p> <p>Changed tentative schedule from Type III Blood Product Table Update to Maintenance TBD for these CRs and DRs: CR 2,100, CR 2,254, CR 1,924, CR 1,920, CR 1,698, CR 2,286, CR 1,629, CR 2,650, CR 2,387, CR 3,248, CR 2,096, CR 2,295, CR 2,354, CR 2,363, CR 2,494, CR 1,732, CR 1,807, DR 1,984</p>	

Date	Revision	Description	Author
8/20/12	15.0	Updated the tentative schedule for Daily Reagent Quality Control Upgrade to include NSR 20081120 Updated the tentative schedule for Report Tool Upgrade Project Proposal to include NSR 20080915	BBM Team