



**Update \_2\_0\_16**

**Clinical Reminders**

**VA-WH MAMMOGRAM SCREENING**

**Install Guide**

**October 2017**

Product Development  
Department of Veterans Affairs

# Contents

---

<b>Introduction .....</b>	<b>3</b>
<b>Pre-Installation .....</b>	<b>6</b>
<b>Install Details .....</b>	<b>6</b>
<b>Install Example.....</b>	<b>7</b>
<b>Post-Install Set-up Instructions .....</b>	<b>17</b>

## Introduction

The mammogram screening reminder was originally released in 2005 by the Office of Women's Health. In May 2017, the VHA adopted the American Cancer Society (ACS) breast cancer screening guidelines that give women a choice to begin screening at age 40. The guidelines also recommend starting yearly mammograms by age 45 and then every other year from age 55. The patch releases these edits to the existing national mammogram screening reminder.

### Revisions:

1. Previously two reminder definitions and reminder dialogs were deployed; one for patients aged 50-74 and another aged 40-49. The 2017 revision will replace the two definitions and dialogs with one reminder definition and one reminder dialog. The reminder dialog will use branching logic to display three separate dialogs to the user based on the age of the patient; one for patients 40-44, one for patients 45-74, and one for patients >74.
2. The reminder definition baseline frequencies have been revised to the following:
  - a. 0Y for patients 39 and younger
  - b. 99Y for patients 40-44
  - c. 1Y for patients 45-54
  - d. 2Y for patients 55 and older
3. New health factors are deployed and appropriate terms and logic changes have been added to the reminder definition.
  - a. VA-WH BR CA 40-44 BEGIN AGE 45
    - i. Indicates a patient aged 40-44yo has chosen to begin annual mammogram screening at age 45. This health factor has a condition that will allow it to be evaluated as true only until the patient reaches age 45.
  - b. VA-WH BR CA 40-44 WANTS SCREENING
    - i. Indicates a patient has chosen to begin annual mammogram screening between ages 40-44. This health factor has a condition that will allow it to be evaluated as true only until the patient reaches age 45.
4. The condition within the currently deployed health factor WH BR CA 40-49 WANTS SCREENING, so that it will only be evaluated as true for patients until the patient reaches age 45. At age 45 the reminder will become applicable annually.
5. The 4-month frequency option has been removed from the reminder dialogs. The frequency finding was not removed from the reminder definition logic, so the reminder will still become applicable at 4 months for the appropriate patients. Providers can then choose a different frequency.
6. Logic added to the reminder definition to make the reminder applicable to patients aged 40-44, so the provider and patient can decide to begin screening at age 45 or between 40-44yo.
7. Logic added to have reminder to continue to become applicable after age 74, so provider and patient can decide to continue mammogram screening if patient is in good health and has life expectancy greater than 10 years.
8. Logic added to have reminder applicable annually for patients 45-54 and every two years for patients > 54yo. Providers will have opportunity to change frequency to annually for patients 55-74.

9. The reminder definition VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49 will be inactivated. The print name will be changed to “ZZDiscuss Breast CA (Discontinued)” to indicate it should not be used.
10. The reminder dialog VA-WH DISCUSS BREAST CA WOMAN 40-49 will be disabled.

**Two reminder Definitions are being distributed with this update:**

VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49  
VA-WH MAMMOGRAM SCREENING

UPDATE\_2\_0\_16 contains 1 Reminder Exchange entry:

UPDATE\_2\_0\_16 VA-WH MAMMOGRAM SCREENING

**The exchange file contains the following components:**

**TIU TEMPLATE FIELD**

BLANK TEXT FIELD FOR DIALOGS  
VA-WH MAMMOGRAM GUIDANCE

**HEALTH FACTORS**

VA-REMINDER UPDATES  
VA-UPDATE\_2\_0\_16  
WH MAMMOGRAM  
WH BR CA SCREEN N/A 5 YRS-COMORBIDITIES  
WH BR CA SCREEN N/A 5 YRS-LE<5YRS  
WH MAMMOGRAM SCREEN NOT INDICATED  
VA LIFE EXPECTANCY  
VA LIMITED LIFE EXPECTANCY  
WH BILATERAL MASTECTOMY  
VA-WH BR CA DEFER 1Y  
VA-WH BR CA DEFER 6M  
WH ORDER MAMMOGRAM SCREEN HF  
WH OUTSIDE INCOMPLETE MAMMOGRAM  
WH OUTSIDE ABNL MAMMOGRAM  
WH OUTSIDE NORMAL MAMMOGRAM  
VA-WH BR CA 40-44 WANTS SCREENING  
VA-WH BR CA 40-44 BEGIN AGE 45  
WH MAMMOGRAM SCREEN FREQ - 2Y  
WH MAMMOGRAM SCREEN FREQ - 1Y  
WH MAMMOGRAM SCREEN FREQ - 6M  
WH MAMMOGRAM DEFERRED  
WH MAMMOGRAM DECLINED  
WH BR CA 40-49 WANTS SCREENING  
WH MAMMOGRAM SCREEN FREQ - 4M  
REMINDER FACTORS  
INACTIVATE BREAST CANCER SCREEN

WH ORDER MAMMOGRAM UNILAT HF  
WH ORDER MAMMOGRAM BILAT HF  
WH MAMMOGRAM OUTSIDE  
WH BR CA 40-49 DEFER 5Y  
WH BR CA 40-49 DEFER 1Y  
WH BR CA 40-49 DEFER 6M  
WH BR CA 40-49 BEGIN AGE 50

**REMINDER SPONSOR**

Women Veterans Health Program  
Women Veterans Health Program and National Center for Health Promotion and Disease Prevention (NCP)  
Mental Health Services and Women Veterans Health Program  
Office of Patient Care Services  
VHA National Center for Health Promotion and Disease Prevention (NCP) and VHA Office of Women's Health  
National Center for Health Promotion and Disease Prevention (NCP)

**REMINDER TAXONOMY**

VA-BREAST TUMOR  
VA-MAMMOGRAM/SCREEN  
VA-TERMINAL CANCER PATIENTS  
VA-WH BILATERAL MASTECTOMY  
VA-MASTECTOMY

**REMINDER TERM**

VA-UPDATE\_2\_0\_16  
VA-AGE >39&<45  
VA-AGE >74  
VA-WH MAMMOGRAM AGE<55  
VA-WH BR CA 40-44 WANTS SCREEN TERM  
VA-WH BR CA 40-44 BEGIN AGE 45  
VA-WH MAMMOGRAM ORDER  
VA-WH HX BREAST CANCER/ABNORMAL MAM  
VA-WH MAMMOGRAM SCREEN FREQ - 2Y  
VA-WH MAMMOGRAM SCREEN FREQ - 1Y  
VA-WH MAMMOGRAM SCREEN FREQ - 6M  
VA-WH MAMMOGRAM SCREEN FREQ - 4M  
VA-WH MAMMOGRAM UNSATISFACTORY IN RAD/WH PKG  
VA-WH MAMMOGRAM SCREEN DEFER  
VA-WH MAMMOGRAM SCREEN NOT INDICATED  
VA-WH BREAST CARE ORDER HEALTH FACTOR  
VA-WH MAMMOGRAM SCREEN DONE  
VA-WH MAMMOGRAM SCREEN IN RAD PKG  
VA-WH MAMMOGRAM SCREEN IN WH PKG  
VA-TERMINAL CANCER PATIENT

VA-WH BILATERAL MASTECTOMY  
**REMINDER DEFINITION**  
VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49  
VA-WH MAMMOGRAM SCREENING

**REMINDER DIALOG**  
VA-WH DISCUSS BREAST CA WOMAN 40-49  
VA-WH MAMMOGRAM SCREENING

## Pre-Installation

---

1. Create a copy of the reminder inquiry output for the current VA-WH MAMMOGRAM SCREENING and VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49 installed in VistA to use as a reference for mapping to the new reminder definition terms.
2. Create a copy of the detailed display of the VA-WH MAMMOGRAM SCREENING and VA-WH DISCUSS BREAST CA WOMAN 40-49 reminder dialogs to use as a reference for mapping findings to the new reminder dialog groups and elements.
3. In reminder Dialog manager, Dialog manager, Group Vies, Inquire on Group VA-WH GP ORDERS-MAMMOGRAM SCREEN and note the order used as a finding. You will need to map this order in the post install section.

## Install Details

---

This update is being distributed as a web host file. The address for the host file is:  
**[https://vaww.va.gov/reminders/docs/UPDATE\\_2\\_0\\_16.PRD](https://vaww.va.gov/reminders/docs/UPDATE_2_0_16.PRD)**

The file will be installed using Reminder Exchange, programmer access is not required.

*This patch can be installed with users on the system, but it should be done during non-peak hours. Estimated Installation Time: 5-10 minutes*

## Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: http://vista.med.va.gov/reminders/UPDATE_2_0_16.PRD
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

[https://vaww.va.gov/reminders/docs/UPDATE\\_2\\_0\\_16.PRD](https://vaww.va.gov/reminders/docs/UPDATE_2_0_16.PRD)

```
http://vista.med.va.gov/reminders/UPDATE_2_0_16.PRD successfully loaded.
```

You should see a message at the top of your screen that the file successfully loaded.

Search and locate an entry titled UPDATE\_2\_0\_16 VA-WH MAMMOGRAM SCREENING in reminder exchange

```
+Item  Entry                                     Source                                     Date Packed
61  UPDATE_2_0_16 VA-WH MAMMOGRAM SCREENING  GRIFFITH@SALT LAKE CI  09/15/2017@05:59
62  UPDATE_2_0_17 VA-WH PRECONCEPT ION CARE  GRIFFITH@SALT LAKE CI  11/14/2016@09:36
63  UPDATE_2_0_18 MISC UPDATE 2             GRIFFITH@SALT LAKE CI  10/20/2016@12:42
64  UPDATE_2_0_19 VA-CSC INITIAL CAREGIVER ASSESSMENT  GRIFFITH@SALT LAKE CI  01/31/2017@12:14
65  UPDATE_2_0_2 CAREGIVER ELIG             GRIFFITH@SALT LAKE CI  02/18/2016@13:10
66  UPDATE_2_0_21 VA-BATTLEFIELD ACUPUNCTURE  GRIFFITH@SALT LAKE CI  04/18/2017@05:39
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History
Select Action: Next Screen// IFE  Install Exchange File Entry
Enter a list or range of numbers (1-220): 61
```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry titled UPDATE\_2\_0\_16 VA-WH MAMMOGRAM SCREENING (in this example it is entry 61, it will vary by site)

Component	Category	Exists
Source: GRIFFITH,ELIZABETH A at SALT LAKE CITY		
Date Packed: 09/15/2017@05:59:09		
Package Version: 2.0P35		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER TERM		
VA-UPDATE_2_0_16		
REMINDER DEFINITION		
VA-WH MAMMOGRAM SCREENING		
VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49		
Keywords:		
Components:		
+ Enter ?? for more actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// IA		

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.

You will see several prompts, for all new entries you will choose **I to Install**. For existing entries you will **Overwrite** or **MERGE** if it is a reminder term. Here are the prompts you may see:

```
Select Action: Next Screen// IA Install all Components
TIU TEMPLATE FIELD entry named BLANK TEXT FIELD FOR DIALOGS already
exists and the packed component is identical, skipping.
TIU TEMPLATE FIELD entry named VA-WH MAMMOGRAM GUIDANCE already exists
but the packed component is different, what do you want to do?
  Select one of the following:
    C      Create a new entry by copying to a new name
    O      Overwrite the current entry
    U      Update
    Q      Quit the install
    S      Skip, do not install this entry
Enter response: O// overwrite the current entry
Are you sure you want to overwrite? N// y YES
HEALTH FACTORS entry named VA-REMINDER UPDATES already exists and the
packed component is identical, skipping.
HEALTH FACTORS entry VA-UPDATE_2_0_16 is NEW,
what do you want to do?
  Select one of the following:
    C      Create a new entry by copying to a new name
    I      Install
    Q      Quit the install
    S      Skip, do not install this entry
Enter response: I// nstall
HEALTH FACTORS entry named WH MAMMOGRAM already exists and the packed
component is identical, skipping.
HEALTH FACTORS entry named WH MAMMOGRAM SCREEN FREQ - 1Y already
exists and the packed component is identical, skipping.
```



HEALTH FACTORS entry named WH MAMMOGRAM SCREEN FREQ - 6M already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BR CA SCREEN N/A 5 YRS-COMORBIDITIES already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BR CA SCREEN N/A 5 YRS-LE<5YRS already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH MAMMOGRAM SCREEN NOT INDICATED already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named VA LIFE EXPECTANCY already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named VA LIMITED LIFE EXPECTANCY already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BILATERAL MASTECTOMY already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry VA-WH BR CA DEFER 1Y is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry VA-WH BR CA DEFER 6M is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry named WH ORDER MAMMOGRAM SCREEN HF already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH OUTSIDE INCOMPLETE MAMMOGRAM already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH OUTSIDE ABNL MAMMOGRAM already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH OUTSIDE NORMAL MAMMOGRAM already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry VA-WH BR CA 40-44 WANTS SCREENING is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry VA-WH BR CA 40-44 BEGIN AGE 45 is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall

HEALTH FACTORS entry named WH MAMMOGRAM DECLINED already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH MAMMOGRAM SCREEN FREQ - 2Y already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH MAMMOGRAM DEFERRED already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BR CA 40-49 WANTS SCREENING already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH MAMMOGRAM SCREEN FREQ - 4M already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named REMINDER FACTORS already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named INACTIVATE BREAST CANCER SCREEN already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH ORDER MAMMOGRAM UNILAT HF already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH ORDER MAMMOGRAM BILAT HF already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH MAMMOGRAM OUTSIDE already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BR CA 40-49 DEFER 5Y already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BR CA 40-49 DEFER 1Y already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BR CA 40-49 DEFER 6M already exists and the packed component is identical, skipping.  
HEALTH FACTORS entry named WH BR CA 40-49 BEGIN AGE 50 already exists and the packed component is identical, skipping.  
REMINDER SPONSOR entry named Women Veterans Health Program already exists and the packed component is identical, skipping.  
REMINDER SPONSOR entry named Women Veterans Health Program and National Center for Health Promotion and Disease Prevention (NCP) already exists and the packed component is identical, skipping.  
REMINDER SPONSOR entry named Mental Health Services and Women Veterans Health Program already exists and the packed component is identical, skipping.  
REMINDER SPONSOR entry named Office of Patient Care Services already exists and the packed component is identical, skipping.  
REMINDER SPONSOR entry named VHA National Center for Health Promotion and Disease Prevention (NCP) and VHA Office of Women's Health already exists and the packed component is identical, skipping.  
REMINDER SPONSOR entry named National Center for Health Promotion and Disease Prevention (NCP) already exists and the packed component is identical, skipping.  
REMINDER TAXONOMY entry named VA-MASTECTOMY already exists but the packed component is different, what do you want to do?  
Select one of the following:  
C Create a new entry by copying to a new name  
O Overwrite the current entry  
U Update  
Q Quit the install  
S Skip, do not install this entry

Enter response: **O// verwrite the current entry**  
Are you sure you want to overwrite? **N// YES**  
REMINDER TAXONOMY entry named VA-BREAST TUMOR already exists but the packed component is different, what do you want to do?  
Select one of the following:  
C Create a new entry by copying to a new name  
O Overwrite the current entry  
U Update  
Q Quit the install  
S Skip, do not install this entry

Enter response: **O// verwrite the current entry**  
Are you sure you want to overwrite? **N// YES**  
REMINDER TAXONOMY entry named VA-TERMINAL CANCER PATIENTS already exists but the packed component is different, what do you want to do?  
Select one of the following:  
C Create a new entry by copying to a new name  
O Overwrite the current entry  
U Update  
Q Quit the install  
S Skip, do not install this entry

Enter response: **O// verwrite the current entry**  
Are you sure you want to overwrite? **N// YES**  
REMINDER TAXONOMY entry named VA-MAMMOGRAM/SCREEN already exists but the packed component is different, what do you want to do?  
Select one of the following:  
C Create a new entry by copying to a new name  
O Overwrite the current entry  
U Update  
Q Quit the install  
S Skip, do not install this entry

Enter response: **O// verwrite the current entry**  
Are you sure you want to overwrite? **N// YES**  
REMINDER TAXONOMY entry named VA-WH BILATERAL MASTECTOMY already exists but the packed component is different, what do you want to do?  
Select one of the following:  
C Create a new entry by copying to a new name  
O Overwrite the current entry  
U Update  
Q Quit the install  
S Skip, do not install this entry

Enter response: **O// verwrite the current entry**  
Are you sure you want to overwrite? **N// y YES**  
REMINDER TERM entry VA-UPDATE\_2\_0\_16 is NEW,  
what do you want to do?  
Select one of the following:  
C Create a new entry by copying to a new name  
I Install  
Q Quit the install  
S Skip, do not install this entry

Enter response: **I// nstall**  
REMINDER TERM entry VA-AGE >39<45 is NEW,  
what do you want to do?  
Select one of the following:

C Create a new entry by copying to a new name  
I Install  
Q Quit the install  
S Skip, do not install this entry

Enter response: **I// nstall**

REMINDER TERM entry VA-AGE >74 is NEW,  
what do you want to do?

Select one of the following:

C Create a new entry by copying to a new name  
I Install  
Q Quit the install  
S Skip, do not install this entry

Enter response: **I// nstall**

REMINDER TERM entry VA-WH MAMMOGRAM AGE<55 is NEW,  
what do you want to do?

Select one of the following:

C Create a new entry by copying to a new name  
I Install  
Q Quit the install  
S Skip, do not install this entry

Enter response: **I// nstall**

REMINDER TERM entry VA-WH BR CA 40-44 WANTS SCREEN TERM is NEW,  
what do you want to do?

Select one of the following:

C Create a new entry by copying to a new name  
I Install  
Q Quit the install  
S Skip, do not install this entry

Enter response: **I// nstall**

REMINDER TERM entry VA-WH BR CA 40-44 BEGIN AGE 45 is NEW,  
what do you want to do?

Select one of the following:

C Create a new entry by copying to a new name  
I Install  
Q Quit the install  
S Skip, do not install this entry

Enter response: **I// nstall**

REMINDER TERM entry named VA-WH MAMMOGRAM ORDER already exists but  
the packed component is different, what do you want to do?

Select one of the following:

C Create a new entry by copying to a new name  
M Merge findings  
O Overwrite the current entry  
U Update  
Q Quit the install  
S Skip, do not install this entry

Enter response: **O// Merge findings**

REMINDER TERM entry named VA-WH HX BREAST CANCER/ABNORMAL MAM already  
exists and the packed component is identical, skipping.

REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN FREQ - 2Y already  
exists and the packed component is identical, skipping.

REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN FREQ - 1Y already  
exists and the packed component is identical, skipping.

REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN FREQ - 6M already exists and the packed component is identical, skipping.  
REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN FREQ - 4M already exists and the packed component is identical, skipping.  
REMINDER TERM entry named VA-WH MAMMOGRAM UNSATISFACTORY IN RAD/WH PKG already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- M Merge findings
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: **O// Merge findings**

REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN DEFER already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- M Merge findings
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: **O// Merge findings**

REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN NOT INDICATED already exists and the packed component is identical, skipping.  
REMINDER TERM entry named VA-WH BREAST CARE ORDER HEALTH FACTOR already exists and the packed component is identical, skipping.  
REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN DONE already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- M Merge findings
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: **O// Merge findings**

REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN IN RAD PKG already exists and the packed component is identical, skipping.  
REMINDER TERM entry named VA-WH MAMMOGRAM SCREEN IN WH PKG already exists and the packed component is identical, skipping.  
REMINDER TERM entry named VA-TERMINAL CANCER PATIENT already exists and the packed component is identical, skipping.  
REMINDER TERM entry named VA-WH BILATERAL MASTECTOMY already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- M Merge findings
- O Overwrite the current entry
- U Update

Q Quit the install  
 S Skip, do not install this entry  
 Enter response: **O// Merge findings**  
 REMINDER TERM entry named VA-WH BR CA 40-49 WANTS SCREEN TERM  
 already exists and the packed component is identical, skipping.  
 REMINDER DEFINITION entry named VA-WH MAMMOGRAM SCREENING already  
 exists but the packed component is different, what do you want to do?  
 Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: **O// overwrite the current entry**  
 Are you sure you want to overwrite? **N// YES**  
 REMINDER DEFINITION entry named VA-WH DISCUSS BREAST CA SCREEN WOMAN  
 40-49 already exists but the packed component is different, what do  
 you want to do?  
 Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: **O// overwrite the current entry**  
 Are you sure you want to overwrite? **N// YES**

After the components have installed,

**You will be prompted to install the first reminder dialog component:**

Packed reminder dialog: VA-WH MAMMOGRAM SCREENING

Item	Seq.	Dialog Findings	Type	Exists
1		VA-WH MAMMOGRAM SCREENING	dialog	X
2	5	VA-WH MAMMOGRAM SCREENING HEADER* Finding: *NONE*	group	
3	5.5	VA-WH GP MAM SCREEN HEADER* Finding: *NONE*	group	X
4	5.5.5	VA-WH MAM HEADER1 - SCREEN Finding: *NONE*	element	X
5	5.5.10	VA-WH GP MAM SCREEN HEADER-EXPANDED Finding: *NONE*	group	X
6	5.5.10.5	VA-WH MAM HEADER2 - SCREEN Finding: *NONE*	element	X
7	5.5.15	VA-WH DOTTED LINE Finding: *NONE*	element	X
8	5.5.20	VA-WH GP MAMMOGRAM SCREENING	group	X
<b>+ + Next Screen - Prev Screen ?? More Actions</b>				
DD	Dialog Details	DT	Dialog Text	IS Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU Quit
DS	Dialog Summary	IA	Install All	

Select Action: Next Screen// IA

At the **Select Action** prompt, type **IA** to install the dialog – **VA-WH MAMMOGRAM SCREENING**

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// YES

```
Packed reminder dialog: VA-WH MAMMOGRAM SCREENING
VA-WH MAMMOGRAM SCREENING (reminder dialog) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-WH MAMMOGRAM SCREENING	dialog	X
2	5	VA-WH MAMMOGRAM SCREENING HEADER*	group	X
		Finding: *NONE*		
3	5.5	VA-WH GP MAM SCREEN HEADER*	group	X
		Finding: *NONE*		
4	5.5.5	VA-WH MAM HEADER1 - SCREEN	element	X
		Finding: *NONE*		
5	5.5.10	VA-WH GP MAM SCREEN HEADER-EXPANDED	group	X
		Finding: *NONE*		
6	5.5.10.5	VA-WH MAM HEADER2 - SCREEN	element	X
		Finding: *NONE*		
7	5.5.15	VA-WH DOTTED LINE	element	X
		Finding: *NONE*		
8	5.5.20	VA-WH GP MAMMOGRAM SCREENING	group	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

**You will then be prompted to install the next Dialog component:**

```
Packed reminder dialog: VA-WH DISCUSS BREAST CA WOMAN 40-49
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-WH DISCUSS BREAST CA WOMAN 40-49	dialog	X
2	5	VA-BR CA 40-44 DISCUSS INTRO	group	
		Finding: *NONE*		
3	5.5	VA-WH GP MAM SCREEN HEADER-EXPANDED	group	X
		Finding: *NONE*		
4	5.5.5	VA-WH MAM HEADER2 - SCREEN	element	X
		Finding: *NONE*		
5	6	VA-WH DOTTED LINE	element	X
		Finding: *NONE*		
6	10	VA-BR CA 40-44 DISCUSS DOCUMENTING RESULTS	group	
		Finding: *NONE*		
7	10.5	VA-BR CA 40-44 DISCUSS BEGIN AGE 45	element	

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-WH DISCUSS BREAST CA WOMAN 40-49**

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// YES

```

Packed reminder dialog: VA-WH DISCUSS BREAST CA WOMAN 40-49
VA-WH DISCUSS BREAST CA WOMAN 40-49 (reminder dialog) installed from exchange fi
Item  Seq.  Dialog Findings                                     Type  Exists
  1    VA-WH DISCUSS BREAST CA WOMAN 40-49                dialog  X
  2    5 VA-BR CA 40-44 DISCUSS INTRO                      group   X
      Finding: *NONE*
  3    5.5 VA-WH GP MAM SCREEN HEADER-EXPANDED            group   X
      Finding: *NONE*
  4    5.5.5 VA-WH MAM HEADER2 - SCREEN                   element  X
      Finding: *NONE*
  5    6 VA-WH DOTTED LINE                                 element  X
      Finding: *NONE*
  6    10 VA-BR CA 40-44 DISCUSS DOCUMENTING RESULTS      group   X
      Finding: *NONE*
  7    10.5 VA-BR CA 40-44 DISCUSS BEGIN AGE 45           element  X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```

Component                                     Category  Exists
Source:      GRIFFITH,ELIZABETH A at SALT LAKE CITY
Date Packed: 09/15/2017@05:59:09
Package Version: 2.0P35

Description:
The following Clinical Reminder items were selected for packing:
REMINDER TERM
  VA-UPDATE_2_0_16

REMINDER DEFINITION
  VA-WH MAMMOGRAM SCREENING
  VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49

Keywords:

Components:
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components      IS  Install Selected Component
Select Action: Next Screen// Q

```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install Completed



# Post-Install Set-up Instructions

---

## 1. Map local health factors to new terms.

### A. Term Name: **VA-WH BR CA 40-44 BEGIN AGE 45**

- a. Map any local findings that represent a patient has decided to begin annual mammogram screening at age 45.

### B. Term name: **VA-WH BR CA 40-44 WANTS SCREEN TERM**

- a. Map any local finding(s) that represent a patient has decided to begin annual screening at age 40 and continue to age 44.

## 2. Confirm following dialog groups and elements have appropriate local orders embedded.

### A. Group Name: **VA-WH GP ORDERS - MAMMOGRAM SCREEN**

- a. Map local mammogram screening order to the FINDING ITEM field. This group was not revised in this update, but sites should confirm the appropriate order is still embedded.

NAME: VA-WH GP ORDERS - MAMMOGRAM SCREEN  
DISABLE:  
CLASS: LOCAL//  
SPONSOR: Women Veterans Health Program//  
REVIEW DATE:  
RESOLUTION TYPE:  
ORDERABLE ITEM:  
FINDING ITEM: **Order should be placed here**

## 3. Confirm the reminder definition **VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49** has been removed from the CPRS cover sheet list and that it is inactive.

- A. The following should display when accessing the reminder definition from the reminder definition manager menu.

Select Reminder Definition: **VA-WH DISCUSS BREAST CA SCREEN WOMAN 40-49  
NATIONAL (INACTIVE)**

***NOTE:** If you require further technical assistance, please notify your local IT support to log a national CA Service Desk Manager (SDM) ticket (previously a Remedy™ ticket) or contact the VA Service Desk at 1-888-596-4357 and have them submit a national CA ticket to the Incident Area: **NTL.APP.VISTA.CLINICAL REMINDERS 2\_0** and someone will contact you.*