



# **PROSTHETICS NPPD DETAIL DISPLAY PATCH RMPR\*3\*71 INSTALLATION GUIDE**

Version 3.0

**November 2002**

Department of Veterans Affairs  
VISTA System Design and Development

# Patch RMPR\*3\*71 Installation Guide

## Extract the file

---

**Introduction** With Patch RMPR\*3\*71, Prosthetics offers the **NPPD Detail Display** feature, which is the first GUI windows-based functionality available to users. This ***Installation Guide*** will help IRM and users to download this patch.

---

**Anonymous directory** The file for this patch can be obtained from the ANONYMOUS.SOFTWARE directory at one of the OI Field Offices. The preferred method is to FTP the file from DOWNLOAD.VISTA.MED.VA.GOV, which will transmit the file from the first available server.

Alternatively, you may elect to retrieve the file from a specific OI Field Office. The documentation can be read via the Adobe Acrobat Reader browser program. The .EXE file is binary file and must be transferred using binary file transfer tools.

---

**Steps** To begin the install process and unzip the files, follow these steps:

---

Step	Action
1	From the Anonymous Directory, download the file <b>RMPR_30_71.EXE</b> .
2	Double click the <b>RMPR_30_71.exe</b> icon.

---


*Continued on next page*

## Extract the file, Continued

**Unzipping files** The file RMPR\_30\_71.EXE is a zip file that contains the “setup” files needed for the automated installation. The files will unzip to the C:\PROS folder for setup purposes. As an alternative to extracting these files to every client computer, IRM staff can unzip this file to a shared folder (directory) on a network server and then run the installation “SETUP.EXE” remotely.

**Recommendation:** Do not change the file folder, as future patches will need to overwrite these files.

### WinZip dialog box



### Steps (continued)

To continue to unzip the files, follow these steps:

Step	Action
3	In the <b>WinZip Self-Extractor</b> dialog box, click the <b>Unzip</b> button (as shown above).
4	A successful message box displays when the unzip process has completed. (See below.)
5	Click the <b>OK</b> button on the <b>WinZip Self-Extractor</b> confirmation box.
6	The <b>WinZip Self-Extractor</b> dialog box continues to display (as shown above).
7	Click the <b>Close</b> button.

### Successful message box



# Install the NPPD File

## Automatic Wizard


Once you have extracted the files onto your c-drive, you can perform an installation using an automatic wizard that is provided for you. This feature has combined several steps to make installation easy.

## Steps

To begin the installation, follow these steps:

Step	Action
1	Open the <b>My Computer</b> icon on your desktop by double clicking it.
2	Double click the C-drive in your <b>Windows Explorer</b> view.
3	Double click the folder that contains the files: <b>PROS</b> .

## Windows Explorer



*Continued on next page*


## Install the NPPD File, Continued

### Steps (continued)

To run the installation program, follow these steps:

Step	Action
4	Double click the <b>Setup.exe</b> file in the C:\PROS folder shown on your <b>Windows Explorer</b> view.

### Windows Explorer



*Continued on next page*


## Install the NPPD File, Continued

### Steps (continued)

To run the wizard for the **NPPD Detail Display** installation, follow these steps:

Step	Action
5	Click the <b>Next</b> button on the first window of the <b>Install Shield Wizard</b> for the <b>NPPD Detail Display</b> .

### NPPD Detail Display Setup



*Continued on next page*


## Install the NPPD File, Continued

### Steps (continued)

To continue the installation, follow these steps:

Step	Action
6	Click the <b>Yes</b> button on the License Agreement if you agree to the terms of the license.  <b>Note:</b> If you choose <b>No</b> , the setup will automatically close.

### License Agreement



*Continued on next page*

## Install the NPPD File, Continued

---

**Recommendation** The destination folder “\Program Files\Vista\PROS” has been established by VA policy so as not to conflict with other VA programs. It is recommended that you do NOT change this portion of the path.

---

**Steps  
(continued)**


To continue the installation, follow these steps:

---


Step	Action
7	Click the <b>Next</b> button on the <b>Destination Location</b> dialog box.
8	Click <b>OK</b> on the confirmation box. Notice that it states that the PROS file folder was created.

---

**File Destination  
Location**



**Confirmation  
message**



*Continued on next page*




## Install the NPPD File, Continued

### Steps (continued)

To continue the installation, follow these steps:

Step	Action
9	Click the <b>Next</b> button on the <b>Start Copy Files</b> dialog box.

### Start Copying Files



*Continued on next page*


## Install the NPPD File, Continued

### Steps (continued)

To continue to install the Prosthetics feature onto your desktop, follow these steps:

Step	Action
10	The <b>Install Shield Wizard</b> will be complete with this step. Click the <b>Finish</b> button.

### Install Shield Wizard



### Installation Complete

The application is now installed on C:\Program Files\Vista\PROS. An icon has been placed on your desktop to run the application NPPD Detail Display. The application can also be run from **Start → Programs → NPPD Detail Display**.

### Uninstall procedures

To uninstall the application, go to the **Control Panel → Add/Remove Programs** option and select the **NPPD Detail Display** item. Click the **REMOVE** button and the program will automatically uninstall. To complete the full removal, the setup files will need to be deleted from the file location you selected.