

Update \_2\_0\_119

## **Clinical Reminders**

# **MENTAL HEALTH REMINDER UPDATES 7.0**

### **Install Guide**

October 2020

Product Development Office of Information Technology Department of Veterans Affairs

Table of Contents	
INTRODUCTION	
PRE-INSTALL	2
INSTALL DETAILS	6
INSTALL EXAMPLE	6
POST INSTALLATION	

# Introduction

Starting January 1, 2021, performance measures will be updated to include the following changes:

- 1. Depression screening will use PHQ-2 instead of PHQ-2+I9
- 2. PTSD screening will use PC PTSD instead of PC-PTSD-5+I9
- 3. Primary suicide screening will be a new stand-alone clinical reminder/dialog pair that will use Mental Health instrument C-SSRS.
- 1. Alcohol screening will continue to use the AUDC as the screening tool, but the instrument will now "flex" based on the patient sex, that is, different questions will be asked. Formerly sites had the choice to link the screening and follow-up combined dialog to the alcohol screening reminder. The screening reminder will have its own dialog and the follow up for positive alcohol will have its own dialog, there will no longer be a combined dialog.
- 2. Follow up for positive alcohol reminder has been simplified to remove distinctions for Alcohol Use Dependency (AUD) and Substance Use Disorder (SUD) clinic visits as part of the branching logic for the dialog. This reminder will be due for any patient who scores 5 or greater on the AUDC tool

Other changes included in this Update:

Reminder Taxonomy VA-DEMENTIA DIAGNOSIS was updated to include A81.00, A81.01, A81.09, A81.2, A81.89, F01.50, F01.51, F02.80, F02.81, F03.90, F03.91, F10.27 F19.97, G23.1, G30.0, G30.8, G31.01, G31.09, G31.83, G90.3 When installing the update and are prompted, choose OVERWRITE for VA-DEMENTIA DIAGNOSIS taxonomy.

Reminder Location List VA-MH NO SHOW APPT CLINICS LL had the following telephone stop codes added: 527,528,530,542,545,546,579 and 584

### Six reminder definitions and their associated dialog are being distributed with this update

VA-ALCOHOL USE SCREEN (AUDIT-C) VA-ALCOHOL AUDIT-C POSITIVE F/U EVAL VA-DEPRESSION SCREENING VA-PTSD SCREENING VA-PRIMARY SUICIDE SCREEN VA-FOLLOW-UP POS PTSD/DEP COMBINATION

UPDATE\_2\_0\_119 contains 1 Reminder Exchange entry: UPDATE\_2\_0\_119 MH UPDATES 7.0 The exchange file contains the following components:

#### **REMINDER DIALOG LINK TYPE**

VA-ALCOHOL POSITIVE WARM HANDOFF CHECK VA-CSSRS RN WARM HANDOFF CHECK VA-PTSD POSITIVE WARM HANDOFF CHECK VA-DEPRESSION POSITIVE WARM HANDOFF CHECK

#### TIU TEMPLATE FIELD

VA-SUD CPG BLANK TEXT FIELD FOR DIALOGS VA-URL MANAGE AUD TX IN PRIMARY CARE VA-PTSD CPG VA-PTSD PCL INFO BLANK SPACE1 VA-URL PHQ9 VA-DEP CPG

#### **EDUCATION TOPICS**

MEDICAL PROBLEMS OF ALCOHOL (SCREENING) ALCOHOL USE AND MEDICAL PROBLEMS

#### **HEALTH FACTORS**

VA-REMINDER UPDATES VA-UPDATE 2 0 119 MENTAL HEALTH VA-WARM HANDOFF ALCOHOL POSITIVE VA-NURSE/LPN-PATIENT COG IMPAIRMENT ALCOHOL USE **REFUSED ALCOHOL USE SCREENING** SEVERE CHRONIC COGNITIVE IMPAIRMENT **UNABLE TO SCREEN - ACUTE ILLNESS REMINDER FACTORS** LIFE EXPECTANCY < 6 MONTHS ALCOHOL - ADVISE ON SAFE LIMITS **REFER FOR ALCOHOL TREATMENT** ALCOHOL - ADVISE TO ABSTAIN ALCOHOL USE WITHIN SAFE LIMITS VA-ALCOHOL ABSTAIN OTHER ADVICE **REFUSED REFERRAL TO SUBSTANCE ABUSE** VA-ALCOHOL OTHER VA-ALCOHOL REFER TO SUD PROGRAM VA-ALCOHOL REFER TO PCMHI VA-ALCOHOL TREAT IN PRIMARY CARE VA-ALCOHOL SELF HELP/AA MEETINGS NON-DRINKER (NO ALCOHOL FOR >1 YR)

VA-ALCOHOL USE AND MEDICAL PROBLEMS ALCOHOL - PRIOR TREATMENT ALCOHOL COUNSELING - PT. RESPONSE ALCOHOL - RESPONSE TO COUNSELING ALCOHOL - ADDL ASSESSMENT OF USE ALCOHOL - ASSESS USE ALCOHOL - NO PRIOR TREATMENT VA-WARM HANDOFF PTSD/DEP POSITIVE **REFUSED PTSD SCREEN** DEP/PTSD/SI ASSESSMENT COMPLETE PTSD ASSESSMENT COMPLETE DEPRESSION ASSESS COMPLETE VA-SUICIDE RISK EVALUATION TYPE [C] VA-SUICIDE RISK UPDATE EVALUATION VA-SUICIDE RISK NEW EVALUATION VA-CSSRS RN COGNITIVE IMPAIRMENT VA-CSSRS DECLINED **REFUSED DEPRESSION SCREENING** 

#### **REMINDER SPONSOR**

Mental Health Services Office of Public Health and Environmental Hazards National Clinical Practice Guideline Council Office of Patient Care Services Mental Health and Behavioral Science Strategic Group Office of Quality & Performance NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION (NCP) Office of Mental Health Services

#### **REMINDER TAXONOMY**

VA-TERMINAL CANCER PATIENTS VA-AUD DIAGNOSIS VA-DEMENTIA DIAGNOSIS VA-PTSD DX OUTPT VISIT VA-DEPRESSION DX OUTPT VISIT

#### **REMINDER LOCATION LIST**

VA-MH NO SHOW APPT CLINICS LL VA-SUD CLINICS

#### **REMINDER TERM**

VA-REMINDER UPDATE\_2\_0\_119 VA-USER CLASS FOR COGNITIVE IMPAIRMENT VA-BL NURSE/LPN COG IMPAIRMENT VA-MH ACUTE ILLNESS

VA-POS ALCOHOL SCREEN VA-COGNITIVE IMPAIRMENT VA-REFUSED ALCOHOL SCREENING VA-LIFE EXPECTANCY < 6 MONTHS VA-ALCOHOL USE SCREEN VA-ALC USE INTERVENTIONS VA-BL AUDIT-C >7 VA-SEX VA-AUD DIAGNOSIS VA-SUD CLINIC VISIT VA-ALCOHOL NONE PAST 1YR VA-ALC MEDICAL PROBLEMS VA-ALCOHOL USE INFORMATION VA-PC-PTSD NEGATIVE VA-DEMENTIA DIAGNOSIS VA-VETERAN PATIENT VA-PTSD DIAGNOSIS VA-SUICIDE SCREEN POSITIVE PTSD VA-SUICIDE SCREEN NEGATIVE PTSD VA-PTSD/DEPRESSION ASSESSMENT COMPLETE VA-REFUSED PTSD SCREEN VA-PTSD SCREEN POSITIVE VA-PTSD SCREEN NEGATIVE VA-PTSD SCREEN VA-BL C-SSRS DONE TODAY VA-C-SSRS DONE TODAY POSITIVE QUESTION 8 VA-C-SSRS DONE TODAY POSITIVE QUESTION 5 VA-C-SSRS DONE TODAY POSITIVE QUESTION 4 VA-C-SSRS DONE TODAY POSITIVE QUESTION 3 VA-BL C-SSRS DONE TODAY POSITIVE VA-USER CLASS NURSE VA-CSSRS DECLINED VA-CSSRS ACUTE ILLNESS **VA-PHO-2 NEGATIVE** VA-SUICIDE SCREEN POSITIVE DEPRESSION VA-SUICIDE SCREEN NEGATIVE DEPRESSION VA-REFUSED DEPRESSION SCREENING VA-DEPRESSION SCREEN POSITIVE VA-DEPRESSION SCREEN NEGATIVE VA-DEPRESSION DIAGNOSIS

#### **REMINDER DEFINITION**

VA-BL NURSE/LPN COG IMPAIRMENT VA-BL C-SSRS DONE TODAY POSITIVE VA-ALCOHOL USE SCREEN (AUDIT-C) VA-ALCOHOL AUDIT-C POSITIVE F/U EVAL VA-PTSD SCREENING VA-PRIMARY SUICIDE SCREEN VA-DEPRESSION SCREENING VA-FOLLOW-UP POS PTSD/DEP COMBINATION

#### HEALTH SUMMARY COMPONENT

MHA Score PCE HEALTH FACTORS SELECTED

#### HEALTH SUMMARY TYPE

OB AUDC VA-CSRE DONE TODAY PC-PTSD PCL+I9 PHQ2+I9

#### **HEALTH SUMMARY OBJECTS**

AUDC (TIU) VA-CSRE DONE TODAY (TIU) PC-PTSD PCL+I9 (TIU) PHQ2+I9 (TIU)

#### **TIU DOCUMENT DEFINITION**

AUDIT-C VA-CSRE DONE TODAY PC-PTSD PCL+I9 PHQ2+I9

#### **REMINDER DIALOG**

VA-ALCOHOL USE SCREENING (AUDIT-C) VA-AUD-C POS BRIEF ALCOHOL COUNSELING VA-PRIMARY SUICIDE SCREEN C-SSRS VA-FOLLOW-UP POS PTSD/DEP COMBINATION VA-PTSD SCREENING VA-DEPRESSION SCREEN

### **Pre-Install**

1. Sites may want to check for quick order(s) in the VA-FOLLOW-UP POS PTSD/DEP COMBINATION dialog and have those available when installing this dialog. Even though the name of the elements leads one to think this is only for PTSD, this should be a generic consult to MH for either a depression and/or a PTSD issue. Below are the elements you can check for Quick Orders.

VA-PTSD SCREEN POS MH EMERGENT EVAL VA-PTSD SCREEN POS MH NON-EMERGENT EVAL

- 2. Sites may want to check for quick order(s) in the VA-AUD-C POS BRIEF ALCOHOL COUNSELING reminder dialog and have available when installing this dialog. Below is the element you can check for Quick Orders.
  - VA-ALC REFER TO SPECIALTY SUD PROG
- 3. Reminder update 37 (UPDATE\_2\_0\_37), update 53 (UPDATE\_2\_0\_53), update 63 (UPDATE\_2\_0\_63) should be considered a prerequisite to Update 119.

## **Install Details**

This update is being distributed as a web host file. The address for the host file is: <u>https://REDACTED/UPDATE\_2\_0\_119.PRD</u>

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

This update can be loaded with users on the system. Installation will take less than 20 minutes.

### **Install Example**

To Load the exchange file, use LWH.	The URL is <u>https://REDACTED</u>
/UPDATE 2 0 119.PRD	

÷	+ Next Screen - Pre	v Screen	??	Nore	Actions	>>>
CHF	Create Host File		LHH	Load	MailMan Message	
CHH	Create MailMan Message		LR	List	Reminder Definiti	ions
DFE	Delete Exchange File Entry		LHH	Load	Heb Host File	
IFE	Install Exchange File Entr	y	RI	Remir	nder Definition In	nquiry •
IH	Installation History		RP	Repa	:k	
Sele	ct Action: Next Screen// LI	H Load	<b>Heb</b>	Host I	ile	
Inpu	t the URL for the .prd file	: https:	11			'UPDATE_2_0_119.PRD
Added Reminder Exchange entry UPDATE_2_0_119 MH UPDATES 7.0						

If you use LWH and it works, you will see something similar to this message. Make sure the file was loaded successfully.

https://

UPDATE\_2\_0\_119.PRD successfully loaded.

Most components will be the same and will be skipped. Here are general rules for prompts to which you will have to respond.

If a component exists, but is different, use OVERWRITE action

If a component is new, use **INSTALL** action For Reminder Terms below, use **OVERWRITE** action

> VA-DEPRESSION SCREEN NEGATIVE VA-DEPRESSION SCREEN POSITIVE VA-PTSD SCREEN VA-PTSD SCREEN NEGATIVE VA-PTSD SCREEN POSITIVE VA-SUICIDE SCREEN NEGATIVE PTSD VA-SUICIDE SCREEN POSITIVE PTSD VA-PC-PTSD NEGATIVE VA-SUICIDE SCREEN NEGATIVE DEPRESSION VA-SUICIDE SCREEN POSITIVE DEPRESSION

Search and locate an entry titled UPDATE\_2\_0\_119 MH UPDATES 7.0 in reminder exchange. If you were a test site, you may have more than 1 entry with this name. Be sure and choose the entry with date of 11/11/2020



At the <u>Select Action</u> prompt, enter <u>IFE</u> for Install Exchange File Entry Enter the number that corresponds with your entry titled UPDATE\_2\_0\_119 MH UPDATES 7.0 *(in this example it is entry 76, it will vary by site)* that is dated 11/11/2020.

At the <u>Select Action</u> prompt, type <u>IA</u> for Install all Components and hit enter.

You will see several prompts:

For all new entries you will choose I to Install

For reminder definitions, Taxonomies, TIU Template fields, Health Summary Type, Reminder Location List entries that exist, but the packed component is different, use OVERWRITE action. For Reminder Terms, use MERGE action except for the ones below: For these, use OVERWRITE action.

VA-DEPRESSION SCREEN NEGATIVE

VA-DEPRESSION SCREEN POSITIVE VA-PTSD SCREEN VA-PTSD SCREEN NEGATIVE VA-PTSD SCREEN POSITIVE VA-SUICIDE SCREEN NEGATIVE PTSD VA-SUICIDE SCREEN POSITIVE PTSD VA-PC-PTSD NEGATIVE VA-SUICIDE SCREEN NEGATIVE DEPRESSION VA-SUICIDE SCREEN POSITIVE DEPRESSION

Dialog installation may not necessarily appear in the order listed in the install guide



You will be promped to install the first reminder dialog component:

Packed	reminder dialog: VA-ALCOHOL USE SCREENING (AUDIT-C)		
Item	Seq. Dialog Findings	Туре	Exists
1	VA-ALCOHOL USE SCREENING (AUDIT-C)	dialog	Х
2	1 VA-ALCOHOL SCREENING VERSION Finding: *NONE*	group	x
3	1.5 VA-TEXT ALCOHOL MISUSE SCREEN INSTRUCTIONS	element	Х
	Finding: *NONE*		
4	1.10 VAL-GP ALC SCREEN AUDC	group	Х
5	Finding: *NONE* 1.10.5 VA-GP MH AUDC	group	х
6	1.10.5.5 VA-MH AUDC	element	X
7	PYRM AUDO RESULT GROUP	rs aroun	Ŷ
8	PXRM AUDC RESULT ELEMENT 1	rs.element	x
9	PXRM AUDC RESULT ELEMENT 3	rs.element	X
+	+ Next Screen - Prev Screen ?? More Actions		
DD D	ialog Details DT Dialog Text IS 1	Install Selecte	d
DF D	ialog Findings DU Dialog Usage QU (	Quit	
DS D	ialog Summary IA <mark>I</mark> nstall All		
Select	Action: Next Screen// IA		

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – VA-ALCOHOL USE SCREEN (AUDIT-C)

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// YES

Pac	ke	ed remi	nder di	alog: \	/A-ALCO	HOL	USE SCI	REEN	ING (AU	DIT-C	;)			
VA-	AL	COHOL	USE SCR	EENING	(AUDIT	-C)	(remin	der (	dialog)	inst	alled f	⁼rom	exchar	nge fil
It	en	ı Seq.	Dialog	Findir	ngs							Тур	be E	xists
	1	VA-7	ALCOHOL	. USE SC	REENIN	G (A	UDIT-C	)				dia	Log	×
	2	1 VA	-ALCOHC Findin	DL SCREE	NING V	ERSI	ON					gro	oup	×
	З	1.5	VA-TEXT	ALCOHO	DL MISU	SE S	CREEN	INST	RUCTION	S	e	leme	ent	×
	4	1.10	Findin VAL-GP Findin	IG: *NON ALC SC	NE* CREEN A	UDC						gro	pup	×
	5	1.10	.5 VA-G	P MH AL	JDC							gro	oup	×
	6	1.10	Findin .5.5 VA Findin	IG: *NON A-MH AUE	IE* )С ; (МН Т	ESTS	AND SI	URVE	Y)		e	eleme	ent	× ×
	7	PXRM	AUDC F	ESULT C	ROUP						rs	s.gro	bup	×
	8	PXRM	AUDC F	ESULT E	ELEMENT	1					rs.e	leme	ent	×
	9	PXRM	AUDC F	ESULT E	ELEMENT	3					rs.e	eleme	ent	×
+		+	Next S	creen	- Pre	v Sc	reen	??	More Ac	tions	;			
DD		Dialog	Detail	S	DT	Dial	.og Tex	t		IS	Install	Se]	Lected	
DF		Dialog	Findin	igs	DU		og Usa	ge		QU	Quit			
DS	-		Summar	y t Canad		Inst	ALL AL	1						
Sel	ec	T ACTI	on: Nex	t Scree	en// Q									

When the dialog has completed installation, you will then be returned to this screen. At the <u>Select Action</u> prompt, type <u>Q</u>.

### You will be promped to install the next reminder dialog component:

Packe	d reminder dialog: VA-AUD-C POS BRIEF ALC	COHOL COUNSELING	[NATIONAL	DIALOG]
Item	Seq. Dialog Findings		Туре	Exists
1	VA-AUD-C POS BRIEF ALCOHOL COUNSELING		dialog	Х
2	1 VA-POS ALCOHOL FOLLOW UP VERSION Finding: *NONE*		element	Х
3	25 TEXT AUDC PRIOR RESULTS Finding: *NONE*		element	Х
4	30 VA-AUDC SCORE GREATHER THAN/EQUAL TO Finding: *NONE*	0 5	element	
5	35 VA-ALC SCREEN POS MEDICAL PROBLEMS Finding: *NONE*		element	x
6	VA-ALCOHOL BAC FEEDBACK OPTIONS		prompt	X
+	+ Next Screen - Prev Screen ?? M	Nore Actions		
DD	Dialog Details DT Dialog Text	IS Insta	all Selecte	ed
DF	Dialog Findings DU Dialog Usage	QU Quit		
DS	Dialog Summary IA Install All			
Selec	t Action: Next Screen// IA			

# At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – VA-AUD-C POS BRIEF ALCOHOL COUNSELING

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// YES

Pac	ke	d reminder dialog: V	A-AUD-C POS E	BRIEF AL	COHOL COU	NSELING	[NATIONAL	DIALOG]
VA-	AU	D-C POS BRIEF ALCOHO	L COUNSELING	(remind	er dialog	) instal	led from	exchange
It	em	<u>Seq. Dialog Findin</u>	gs				Туре	Exists
	1	VA-AUD-C POS BRIE	F ALCOHOL COU	JNSELING			dialog	Х
	2	1 VA-POS ALCOHOL F Finding: *NON	OLLOW UP VER: E*	SION			element	х
	3	25 TEXT AUDC PRIOR Finding: *NON	RESULTS E*				element	x
	4	30 VA-AUDC SCORE G Finding: *NON	REATHER THAN, E*	/EQUAL T	05		element	x
	5	35 VA-ALC SCREEN P Finding: *NON	OS MEDICAL PF E*	ROBLEMS			element	X
	6	VA-ALCOHOL BAC FEE	DBACK OPTION	S			prompt	X
+		+ Next Screen	- Prev Scree	en ??	More Acti	ons		
DD DF DS Sel	.ec	Dialog Details Dialog Findings Dialog Summary t Action: Next Scree	DT Dialog DU Dialog IA Instal: n// Q	Text Usage 1 All	IS QU	Īnsta Quit	ll Select	ed

When the dialog has completed installation, you will then be returned to this screen. At the <u>Select Action</u> prompt, type Q.

#### You will be promped to install the next reminder dialog component:

Packed	reminder dialog: VA-PRIMARY SUICIDE SCREEN C-SSRS	[NATIONAL DIALOO	3]
Item	Seq. Dialog Findings	Туре	Exists
1	VA-PRIMARY SUICIDE SCREEN C-SSRS	dialog	
2	3 VA-PRIMARY SUICIDE SCREEN C-SSRS VERSION NUMBER Finding: *NONE*	element	
3	5 VA-PRIMARY SUICIDE SCREEN C-SSRS HEADER Finding: *NONE*	group	
4	10 VA-GP PRIMARY SUICIDE SCREEN C-SSRS* Finding: *NONE*	group	
5	10.5 VA-SUICIDE RISK ASSESSMENT C-SSRS* Finding: C-SSRS (MH TESTS AND SURVEY)	element	X X
6	PXRM C-SSRS RESULT GROUP 3	rs.group	Х
7	PXRM C-SSRS RESULT ELEMENT 3	rs.element	Х
8	PXRM C-SSRS RESULT GROUP 4	rs.group	X
+	+ Next Screen - Prev Screen ?? More Actions	3	
DD D:	ialog Details DT Dialog Text IS	Install Selected	k
DF D:	ialog Findings DU Dialog Usage QU	Quit	
DS D:	ialog Summary IA Install All		
Select	Action: Next Screen// IA Install All		

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – VA-PRIMARY SUICIDE SCREEN C-SSRS

Select Action: Next Screen // IA Install All

Install reminder dialog and all components with no further changes: Y// YES

You will see these prompts:

Reminder Dialog VA-PRIMA

RY SUICIDE SCREEN C-SSRS is not linked to a reminder.

Select Reminder to Link: VA-PRIMARY SUICIDE SCREEN// HIT ENTER TO LINK REMINDER

Packe	d reminder dialog: VA-PRIMARY SUICIDE SCREEN C-SSRS	NATIONAL DIALO	3]
VA-PF	IMARY SUICIDE SCREEN C-SSRS (reminder dialog) install	Led from exchang	ge file.
Item	n Seq. Dialog Findings	Туре	Exists
1	VA-PRIMARY SUICIDE SCREEN C-SSRS	dialog	Х
2	3 VA-PRIMARY SUICIDE SCREEN C-SSRS VERSION NUMBER Finding: *NONE*	element	Х
3	5 VA-PRIMARY SUICIDE SCREEN C-SSRS HEADER Finding: *NONE*	group	х
4	10 VA-GP PRIMARY SUICIDE SCREEN C-SSRS* Finding: *NONE*	group	X
5	10.5 VA-SUICIDE RISK ASSESSMENT C-SSRS* Finding: C-SSRS (MH TESTS AND SURVEY)	element	X X
6	PXRM C-SSRS RESULT GROUP 3	rs.aroup	X
7	PXRM C-SSRS RESULT ELEMENT 3	rs.element	X
8	PXRM C-SSRS RESULT GROUP 4	rs.group	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DF DS Selec	Dialog Findings DU Dialog Usage QU C Dialog Summary IA Install All et Action: Next Screen// Q Quit	Duit	

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

#### You will be promped to install the next reminder dialog component: Packed reminder dialog: VA-FOLLOW-UP POS PTSD/DEP COMBINATION [NATIONAL

Tacke	a reminder dialog. VA-rollow-or roo rrob/ber conbination		DIALOGI
Item	Seq. Dialog Findings	Туре	Exists
1	VA-FOLLOW-UP POS PTSD/DEP COMBINATION	dialog	Х
2	2 VA-DEPRESSION/PTSD FOLLOW UP VERSION Finding: *NONE*	element	х
З	5 VA-PTSD/DEPRESSION COMBINATION MASTER GP	group	Х
	Finding: DEP/PTSD/SI ASSESSMENT COMPLETE (HEALTH	FACTOR)	X
4	5.5 VA-GP POS PTSD/DEP HEADER	group	X
	Finding: *NONE*		
5	5.5.5 VA-GP POS PTSD/DEP OBJECTS	group	Х
	Finding: *NONE*		
6	5.5.5.5 VA-POS PTSD/DEP OBJECTS	element	Х
	Finding: *NONE*		
7	5.10 VA-PTSD/DEP SCREEN POS DISPOSTIONS GP	group	Х
	Finding: *NONE*		
+	+ Next Screen - Prev Screen ?? More Actions		
DD	Dialog Details DT Dialog Text IS Inst	tall <mark>S</mark> electe	d
DF	Dialog Findings DU Dialog Usage QU Quit	t	
DS	Dialog Summary IA Install All		
Selec	t Action: Next Screen// IA		

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – VA-FOLLOW-UP POS PTSD/DEP COMBINATION

#### Select Action: Next Screen // IA Install All

Install reminder dialog and all components with no further changes: Y// YES

-		-						
Packe	ed reminder dialog: VA-FOLLOW-UP POS PTSD/DEP COMBINATIO	N [NATIONAL	_ DIALOG]					
VA-FO	DLLOW-UP POS PTSD/DEP COMBINATION (reminder dialog) inst	alled from	exchange					
Iter	n Seq. Dialog Findings	Туре	Exists					
1	VA-FOLLOW-UP POS PTSD/DEP COMBINATION	dialog	Х					
2	2 VA-DEPRESSION/PTSD FOLLOW UP VERSION Finding: *NONE*	element	х					
З	5 VA-PTSD/DEPRESSION COMBINATION MASTER GP	group	x					
	Finding: DEP/PTSD/SI ASSESSMENT COMPLETE (HEALTH	FACTOR)	Х					
4	5.5 VA-GP POS PTSD/DEP HEADER	group	Х					
	Finding: *NONE*	· · ·						
5	5.5.5 VA-GP POS PTSD/DEP OBJECTS	group	Х					
	Finding: *NONE*							
6	5.5.5.5 VA-POS PTSD/DEP OBJECTS	element	Х					
	Finding: *NONE*							
7	5.10 VA-PTSD/DEP SCREEN POS DISPOSTIONS GP	group	Х					
	Finding: *NONE*							
+	+ Next Screen - Prev Screen ?? More Actions							
DD	Dialog Details DT Dialog Text IS Ins	tall <mark>S</mark> elect	ed					
DF	Dialog Findings DU Dialog Usage QU Qui	t						
DS	Dialog Summary IA Install All							
Seled	Select Action: Next Screen// Q							

When the dialog has completed installation, you will then be returned to this screen. At the <u>Select Action</u> prompt, type Q.

DTAL OG L

You will be promped to install the next reminder dialog component: Packed reminder dialog: VA-PTSD SCREENING [NATIONAL DIALOG]

Iten	1 Seq. Dialog Findings Type E	xists
1	VA-PTSD SCREENING dialog	Х
2	1 VA-PTSD SCREENING VERSION group	Х
	Finding: *NONE*	
3	1.5 VA-PTSD SCREEN group	Х
	Finding: *NONE*	
4	1.5.10 VA-MH PC-PTSD+I9 SCREEN element	Х
	Finding: PC-PTSD-5 (MH TESTS AND SURVEY)	Х
5	PXRM PTSD RESULT GROUP rs.group	Х
6	PXRM PTSD RESULT ELEMENT 1 rs.element	Х
7	PXRM PTSD RESULT ELEMENT 2 rs.element	Х
8	1.5.20 VA-GP PTSD UNABLE TO SCREEN group	Х
	Finding: *NONE*	
9	1.5.20.5 VA-HF ACUTE ILLNESS element	Х
	Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	Х
+	+ Next Screen - Prev Screen ?? More Actions	
DD	Dialog Details DT Dialog Text IS Install Selected	
DF	Dialog Findings DU Dialog Usage QU Quit	
DS	Dialog Summary IA Install All	
Selec	ct Action: Next Screen// IA	

At the Select Action prompt, type IA to install the dialog - VA-PTSD SCREENING Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// YES

Packe	d reminder dialog: VA-PTSD SCREENING [NATIONAL DIALOG]					
VA-P1	SD SCREENING (reminder dialog) installed from exchange	file.				
Iten	Seq. Dialog Findings	Туре	Exists			
1	VA-PTSD SCREENING	dialog	Х			
2	1 VA-PTSD SCREENING VERSION Finding: *NONE*	group	x			
3	1.5 VA-PTSD SCREEN	group	X			
	Finding: *NONE*					
4	1.5.10 VA-MH PC-PTSD+I9 SCREEN	element	X			
	Finding: PC-PTSD-5 (MH TESTS AND SURVEY)		X			
5	PXRM PTSD RESULT GROUP	rs.group	Х			
6	PXRM PTSD RESULT ELEMENT 1	rs.element	X			
7	PXRM PTSD RESULT ELEMENT 2	rs.element	X			
8	1.5.20 VA-GP PTSD UNABLE TO SCREEN	group	X			
	Finding: *NONE*					
9	1.5.20.5 VA-HF ACUTE ILLNESS	element	Х			
	Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEAL	TH FACTOR)	X			
+	+ Next Screen - Prev Screen ?? More Actions					
DD	Dialog Details DT Dialog Text IS Ins	stall Selecte	d			
DF	Dialog Findings DU Dialog Usage QU Qu	it				
DS	Dialog Summary IA Install All					
Selec	Select Action: Next Screen// Q					

When the dialog has completed installation, you will then be returned to this screen. At the Select Action prompt, type Q.

#### You will be promped to install the last reminder dialog component: Packed reminder dialog: VA-DEPRESSION SCREEN [NATIONAL DIALOG]

Iten	n Seq. Dialog Findings Type E	xists			
1	VA-DEPRESSION SCREEN dialog	Х			
2	1 VA-DEPRESSION SCREENING VERSION group Finding: *NONE*	X			
3	1.5 VA-DEP PHQ2 group Finding: *NONE*	X			
4	1.5.5 VA-PHQ-2 element				
	Finding: PHQ-2 (MH TESTS AND SURVEY)	Х			
5	PXRM PHQ2 RESULT GROUP rs.group	Х			
6	PXRM PHQ2 RESULT ELEMENT 1 rs.element	Х			
7	PXRM PHQ2 RESULT ELEMENT 2 rs.element	X			
8	1.5.20 VA-GP DEP UNABLE TO SCREEN group	Х			
	Finding: *NONE*				
9	1.5.20.5 VA-HF ACUTE ILLNESS element	X			
	Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	X			
+	+ Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details DT Dialog Text IS Install Selected				
DF	Dialog Findings DU Dialog Usage QU Quit				
DS	Dialog Summary IA Install All				
Select Action: Next Screen// IA					

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – VA-DEPRESSION SCREENING Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// YES

<b>D</b>			-			
Раск	ed reminder dialog: VA-DEPRESSION SCREEN [NAT.	IONAL L	DIALOG			
VA-D	EPRESSION SCREEN (reminder dialog) installed 1	from ex	change file.			
Ite	m Seq. Dialog Findings		Туре	Exists		
1	VA-DEPRESSION SCREEN		dialog	Х		
			J. J			
2	1 VA-DEPRESSION SCREENING VERSION		group	х		
	Finding: *NONE*		0 1			
3	1.5 VA-DEP PHQ2		group	х		
	Finding: *NONE*		<b>3</b> 1			
4	1.5.5 VA-PHQ-2		element	х		
	Finding: PHQ-2 (MH TESTS AND SURVEY)			X		
5	PXRM PHQ2 RESULT GROUP		rs.group	x		
6	PXRM PHQ2 RESULT ELEMENT 1		rs.element	х		
7	PXRM PHQ2 RESULT ELEMENT 2		rs.element	X		
8	1.5.20 VA-GP DEP UNABLE TO SCREEN		group	x		
	Finding: *NONE*		0.			
9	1.5.20.5 VA-HF ACUTE ILLNESS		element	X		
	Finding: UNABLE TO SCREEN - ACUTE ILL	NESS (H	HEALTH FACTOR)	X		
+	+ Next Screen - Prev Screen ?? More	Action	าร			
DD	Dialog Details DT Dialog Text	IS	Install Selecte	d		
DF	Dialog Findings DU Dialog Usage	QU	Quit			
DS	Dialog Summary IA Install All					
Sele	Select Action: Next Screen// Q					

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

This will complete the Dialog Install section.

Component	Category Exists
Source: Date Packed: 11/11/2020@07:14:27 Package Version: 2.0P46	
Description: The following Clinical Reminder items were selected REMINDER DIALOG VA-DEPRESSION SCREEN VA-PTSD SCREENING VA-FOLLOW-UP POS PTSD/DEP COMBINATION VA-FOLLOW-UP POS PTSD/DEP COMBINATION VA-PRIMARY SUICIDE SCREEN C-SSRS VA-AUD-C POS BRIEF ALCOHOL COUNSELING VA-ALCOHOL USE SCREENING (AUDIT-C)	for packing:
REMINDER LOCATION LIST VA-MH NO SHOW APPT CLINICS LL	
+ Enter ?? for more actions IA Install all Components IS Install Select Action: Next Screen// Q	>>> L Selected Component

Type Q to exit the installation

### **Post Installation**

- 1. Reminder Term VA-USER CLASS NURSE has two mapped items, both are CF.VA-ASU USER CLASS. Make sure you have your exact RN and LPN user class names in the CF parameter field. If you have more than 2 user classes that represent RN and LPN, you will have to add a new instance of the computed finding VA-ASU USER CLASS, this is done by using quotes: "CF.VA-ASU USER CLASS". This term is used in branching logic to display the warm handoff option in the primary suicide screen reminder dialog, if you have additional user classes that may complete the primary suicide screen and should have this option in the dialog, you can add them to this term.
- 2. Ensure the new VA-PRIMARY SUICIDE SCREEN reminder (print name is "Suicide Screen") is added to your coversheet for users.
- 3. Add your local Quick Orders back into the following elements from Pre-Install
  - a. VA-PTSD SCREEN POS MH EMERGENT EVAL
    - b. VA-PTSD SCREEN POS MH NON-EMERGENT EVAL
    - c. VA-ALC REFER TO SPECIALTY SUD PROG

Install complete.