



**Update \_2\_0\_173M**

**Clinical Reminders**

**VA-COVID-19 IMMUNIZATION REMINDER UPDATE  
Install Guide**

**October 2021**

Product Development  
Office of Information Technology  
Department of Veterans Affairs

Table of Contents

<b>INTRODUCTION</b> .....	<b>1</b>
<b>INSTALL DETAILS</b> .....	<b>2</b>
<b>INSTALL EXAMPLE</b> .....	<b>2</b>

# Introduction

---

## DESCRIPTION:

This is an update to taxonomies only. No reminder dialog or reminder definition changes.

These taxonomies are updated:

VA-SARS-COV-2 IMM MODERNA SERIES BOOSTER	0064A & 91306
VA-SARS-COV-2 IMM PFIZER SERIES BOOSTER	0004A & 91300
VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 1	0051A & 91305
VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 2	0052A & 91305
VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 3	0053A & 91305
VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES BOOSTER	0054A & 91305
VA-SARS-COV-2 IMM NOVAVAX SERIES 1	0041A
VA-SARS-COV-2 IMM NOVAVAX SERIES 2	0042A

Please note – the recent Lexicon patch did not contain the new code for the Janssen 2<sup>nd</sup> /booster dose so that code could not be added to this update. For questions about that code, please contact your HIMS Department.

Thank you to the following test sites who installed and tested this reminder update:

James J Peters VA Medical Center (Bronx)  
Chillicothe VA Medical Center  
VA Nebraska Western Iowa Health Care System

UPDATE\_2\_0\_173M contains 1 Reminder Exchange entry:

## **UPDATE\_2\_0\_173M VA-COVID-19 IMMUNIZATION REMINDER**

### **HEALTH FACTORS**

REMINDER FACTORS [C]  
VA-UPDATE\_2\_0\_173M

### **REMINDER SPONSOR**

NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION  
(NCP)

### **REMINDER TAXONOMY**

VA-SARS-COV-2 IMM NOVAVAX SERIES 2  
VA-SARS-COV-2 IMM NOVAVAX SERIES 1  
VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES BOOSTER  
VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 3  
VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 2

VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 1  
VA-SARS-COV-2 IMM PFIZER SERIES BOOSTER  
VA-SARS-COV-2 IMM MODERNA SERIES BOOSTER

## REMINDER TERM

VA-REMINDER UPDATE\_2\_0\_173M

## Pre-Install

You should have installed update 173L prior to installing this update.


## Install Details

This update is being distributed as a web host file. The address for the host file is:  
**https://Redacted/UPDATE\_2\_0\_173M.PRD**

The file will be installed using Reminder Exchange, programmer access is not required. This update can be loaded with users on the system. Installation will take less than 5 minutes

## Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista



```
+          + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry          LHF  Load Host File
CHF  Create Host File                    LMM  Load MailMan Message
CMM  Create MailMan Message              LWH  Load Web Host File
DFE  Delete Exchange File Entry          RI   Reminder Definition Inquiry
IFE  Install Exchange File Entry         RP   Repack
IH   Installation History

Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: https://[redacted] UPDATE_2_0_173M.PRD
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

[https://Redacted/UPDATE\\_2\\_0\\_173M.PRD](https://Redacted/UPDATE_2_0_173M.PRD)

You should see a message that the file was successfully loaded.



```
https://[redacted] UPDATE_2_0_173M.PRD successfully loaded.
```

+Item	Entry	Source	Date Packed
423	UPDATE_2_0_173M VA-COVID-19 IMMUNIZATION REMINDER		10/28/2021@10:55
424	UPDATE_2_0_174 VAAES TEMPLATE UPDATES		11/12/2020@11:26
425	UPDATE_2_0_175 VA-CLINICAL APPEAL DECISION		11/20/2020@10:46
426	UPDATE_2_0_176 VA-ONC LUNG AND PROSTATE MOLECULAR TESTING		02/05/2021@06:36
427	UPDATE_2_0_177 VA-HOMELESSNESS FOOD		02/09/2021@13:33
+    + Next Screen    - Prev Screen    ?? More Actions    >>>			
CFE	Create Exchange File Entry	LHF	Load Host File
CHF	Create Host File	LMM	Load MailMan Message
CMM	Create MailMan Message	LWH	Load Web Host File
DFE	Delete Exchange File Entry	RI	Reminder Definition Inquiry
IFE	Install Exchange File Entry	RP	Repack
IH	Installation History		
Select Action: Next Screen// IFE    Install Exchange File Entry			
Enter a list or range of numbers (1-835): 423			

Search and locate an entry titled **UPDATE\_2\_0\_173M VA-COVID-19 IMMUNIZATION REMINDER**.

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry  
Enter the number that corresponds with your entry **UPDATE\_2\_0\_173M VA-COVID-19 IMMUNIZATION REMINDER**. (in this example it is entry 423 it will vary by site). The date of the exchange file should be 10/28/2021.

```

Component                                     Category   Exists
Source:
Date Packed: 10/28/2021@10:55:04
Package Version: 2.0P42

Description:
The following Clinical Reminder items were selected for packing:
REMINDER TAXONOMY
  VA-SARS-COV-2 IMM NOVAVAX SERIES 1
  VA-SARS-COV-2 IMM NOVAVAX SERIES 2
  VA-SARS-COV-2 IMM MODERNA SERIES BOOSTER
  VA-SARS-COV-2 IMM PFIZER SERIES BOOSTER
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 1
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 2
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 3
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES BOOSTER

+      Enter ?? for more actions                >>>
IA  Install all Components                      IS  Install Selected Component
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for **Install all Components** and hit enter. Any items that are NEW in any of the below categories, you should **INSTALL SKIP** any HEALTH FACTORS if you are prompted to overwrite **OVERWRITE** any REMINDER TAXONOMIES that you are prompted to overwrite

```

Component                                     Category   Exists
Source:
Date Packed: 10/28/2021@10:55:04
Package Version: 2.0P42

Description:
The following Clinical Reminder items were selected for packing:
REMINDER TAXONOMY
  VA-SARS-COV-2 IMM NOVAVAX SERIES 1
  VA-SARS-COV-2 IMM NOVAVAX SERIES 2
  VA-SARS-COV-2 IMM MODERNA SERIES BOOSTER
  VA-SARS-COV-2 IMM PFIZER SERIES BOOSTER
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 1
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 2
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES 3
  VA-SARS-COV-2 IMM PFIZER TRIS/SUC SERIES BOOSTER

+      + Next Screen  - Prev Screen  ?? More Actions                >>>
IA  Install all Components                      IS  Install Selected Component
Select Action: Next Screen// Q

```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.