

Update _2_0_182

Clinical Reminders

VA-OSP SUICIDE SAFETY PLAN UPDATE Install Guide

August 2021

Product Development Office of Information Technology Department of Veterans Affairs

Table of Contents	
INTRODUCTION	
INSTALL DETAILS	4
Pre-Install	4
INSTALL EXAMPLE	5
POST INSTALLATION	9

Introduction

DESCRIPTION:

The Suicide Prevention Safety Plan and Safety Plan Review/Decline notes have been updated with the goal of streamlining and clarifying language and access to prior Safety Plans as well as to capturing Veteran contact data, specifically:

- No major additions in content
- Added 31B "button" functionality to streamline the presentation of 'Tips' geared toward providers embedded in the template
- Added 31B "button" functionality to display the most recent Suicide Prevention Safety Plan completed (if present) to quickly and easily display prior safety Plan details to providers *this is display only and does not appear in the note*
- A few clarifying language changes: to 'Tips' for providers; to step 6.
- Update of the lethal means firearm section in step 6.
- Additions to document current phone number and physical address
- There are several new health factors added to existing choices that did not previously have a health factor, these were items where the NO response had a health factor and the Yes or NA responses did not, health factors were added to the Yes and NA responses. The naming convention of the health factors was kept consistent with the No choices.
- The hyperlinks throughout should all be updated to the correct site.

The Suicide Behavior and Overdose Report template has one txml template field that is being updated to a new website address for the printable version. The dialog component will not need to be installed.

UPDATE_2_0_182 contains 1 Reminder Exchange entry: UPDATE_2_0_182 VA-OSP SUICIDE SAFETY PLAN UPDATE

The exchange file contains the following components:

REMINDER GENERAL FINDINGS

VIEW PROGRESS NOTE TEXT

TIU TEMPLATE FIELD

VA-OSP DATE EG VA-OSP BLANK SBOR VA-OSP DATE(R) VA-OSP EDIT 40 R IND VA-OSP CITY ST VA-OSP STREET ADD VA-OSP SRM CONSULT EMAIL VA-OPIOID USE SPANISH VA-OPIOID USE ENGLISH VA-OPIOID SAFETY SPANISH

VA-OPIOID SAFETY ENGLISH VA-OSP MEANS SAFETY VA-OSP FIREARM SAFETY VA-OSP MEANS SAFETY VET VA SP EB30 VA SP 3 LINES VA-OSP LMS POCKET CARD SAFETY PLAN NUMBER VA SP EB25 VA SP EB30 IND VA SP CITY ST VA SP EB40 VA SP CRISIS CHAT VA-OSP VHB URL VA-OSP SAFETY PLAN BLANK CLINICIAN VA-OSP SAFETY PLAN BLANK

HEALTH FACTORS

VA-REMINDER UPDATES [C] VA-UPDATE 2 0 182 VA-SUICIDE BEHAVIOR AND OVERDOSE REPORT [C] VA-SBOR EVENT DATE VA-SUICIDE SAFETY PLAN [C] VA-SP SAFETY PLAN REVIEWED VA-SP EMERG CONTACT NEEDS UPDATE VA-SP EMERG CONTACT ARE UP-TO-DATE VA-SP REFUSED SAFETY PLANNING VA-SP CG PARTICIPATE IN SAFETY PLAN N/A VA-SP CG PARTICIPATE IN SAFETY PLAN NO VA-SP CG PARTICIPATE IN SAFETY PLAN YES VA-SP CG GIVEN COPY NOT APPLICABLE VA-SP CG GIVEN COPY NO VA-SP CG GIVEN COPY YES VA-SP VET GIVEN COPY NO VA-SP VET GIVEN COPY YES VA-NALOXONE [C] VA-NALOXONE NOT NEEDED VA-PATIENT DECLINED NALOXONE VA-NALOXONE RX REQUEST PROVIDER NOTIFIED VA-NALOXONE RX ORDERED VA-SP OPIOID DISCUSSION NO VA-SP OPIOID DISCUSSION YES VA-SP OPIOID ACCESS YES VA-SP OPIOID ACCESS NO VA-SP FIREARM ACCESS OTHER ENGAGE OWNER VA-SP FIREARM STORAGE OTHER

VA-SP FIREARM STORAGE OUTSIDE OF HOME VA-SP FIREARM STORAGE DISASSEMBLED VA-SP FIREARM STORAGE UNLOADED WITH AMMO VA-SP FIREARM STORAGE UNLOADED SEPARATE FROM AMMO VA-SP FIREARM STORAGE LOADED W AMMO VA-SP FIREARM STORAGE UNLOCKED VA-SP FIREARM STORAGE GUN SAFE OR LOCKBOX VA-SP FIREARM STORAGE EXTERNAL LOCK VA-CSP VET ASSESSMENT [C] VA-SP GUNLOCK PROVIDER DID NOT OFFER VA-SP GUNLOCK OFFERED N/A SAFE STORAGE IN PLACE VA-SP GUNLOCK VETERAN REQUESTED NONE AVAILABLE VA-SP GUNLOCK OFFERED VETERAN DECLINED VA-SP GUNLOCK OFFERED YES GIVEN VA-SP FIREARM DISCUSSION NO VA-SP FIREARM DISCUSSION YES VA-SP FIREARM ACCESS YES HOME OR ELSEWHERE VA-SP FIREARM ACCESS NO HOME OR ELSEWHERE VA-SP DISTRESS HELP NOT DISCLOSED VA-SP DISTRESS HELP NO VA-SP DISTRESS HELP YES VA-SP SOCIAL CONTACTS NO VA-SP SOCIAL CONTACTS YES

REMINDER SPONSOR

MENTAL HEALTH SERVICES

REMINDER TERM

VA-REMINDER UPDATE_2_0_182

HEALTH SUMMARY TYPE

VA-OSP SAFETY PLAN NOTE

HEALTH SUMMARY OBJECTS

VA-OSP SAFETY PLAN NOTE (TIU)

TIU DOCUMENT DEFINITION

VA-OSP SAFETY PLAN NOTE

REMINDER DIALOG

VA-OSP SUICIDE SAFETY PLAN REVIEW/DECLINE VA-OSP SUICIDE SAFETY PLAN VA-GP SBOR EVENT DATE

Install Details

This update is being distributed as a web host file. The address for the host file is: https://Redacted /UPDATE_2_0_182.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

This update can be loaded with users on the system. Installation will take less than 5 minutes.

Pre-Install

Inquire and save the following reminder dialog element for review post install, use the menu REMINDER DIALOG MANAGER, REMINDER DIALOGS and change View to E for element:

1. VA-EL SAFETY PLAN STEP 6 OPIOID KIT YES RX – you will add the order back post install.

Install Example

T	o Lo	ad the W	leb Host File.	Navig	ate to F	Remir	nder (excha	ange in V	/ista	
	+	+ N	lext Screen	- Prev	Screen	??	More	Acti	ons		>>>
	CFE	Create E	Exchange File	Entry		LHF	Load	Host	File		
	CHF	Create H	lost File			LMM	Load	Mail	Man Messa	ge	
	СНИ	Create M	lailMan Messag	e		LR	List	Remi	nder Defi	nitions	
	DFE	Delete B	Exchange File	Entry		LWH	Load	Web I	Host File		
	IFE	Install	Exchange File	Entry		RI	Remi⊓	nder I	Definitio	n Inquiry	
	IH	Installa	tion History			RP	Repa	ck 🧹			
	Seled	ct Actior	n: Next Screer	I// LWH	Load	Web I	Host I	File			
	Inpu	t the URL	for the .prd	file:	https:/	11				/UPDATE	_2_0_182.PRD

At the **Select Action:** prompt, enter **LWH** for Load Web Host File At the **Input the url for the .prd file:** prompt, type the following web address: https://Redacted/UPDATE 2 0 182.PRD

You should see a message at the top of your screen that the file loaded successfully

https://	<pre>/UPDATE_2_0_182.PRD successfully loaded.</pre>	
Search and locate an entry title	ed UPDATE_2_0_182 VA-OSP SUICIDE SAFETY PLAI	N

+Item	Entry	Source	Date Packed
168	UPDATE_2_0_182 VA-OSP SUICIDE		08/24/2021@10:07
	SAFETY PLAN UPDATE	the second se	
169	UPDATE_2_0_184 VA-COVID-19		01/14/2021@06:21
	TEMPLATE UPDATES		
170	UPDATE_2_0_186 VA-COVID-19		02/09/2021009:57
474	UPDATE 0 A 10 VA COC INITIAL		01/01/0017010.14
171	CAREGIVER ASSESSMENT		01/31/2017@12:14
172	UPDATE 2 0 190 VA-MATERNITY		07/28/2021806:34
	CARE COORDINATION UPDATE		0172072021000101
+	+ Next Screen - Prev Scree	n ?? More Actions	>>>
CFE	Create Exchange File Entry	LHF Load Host File	
CHF	Create Host File	LMM Load MailMan Me	ssage
CMM	Create MailMan Message	LR List Reminder D	efinitions
DFE	Delete Exchange File Entry	LWH Load Web Host F.	ile
IFE	Install Exchange File Entry	RI Reminder Defini	tion Inquiry
IH	Installation History	RP Repack	
Selec	t Action: Next Screen// IFE Inst	tall Exchange File Ent	гу
Enter	a list or range of numbers (1-43)): 168	

At the <u>Select Action</u> prompt, enter <u>IFE</u> for Install Exchange File Entry Enter the number that corresponds with your entry **UPDATE_2_0_182 VA-OSP SUICIDE SAFETY PLAN UPDATE.** *(in this example it is entry 168 it will vary by site).* The date of the exchange file should be 08/24/2021.

Component		Category	Exists
Source: Date Packed: 08/24/2021@10:07:47 Package Version: 2.0P42 Description: The following Clinical Reminder items we REMINDER DIALOG VA-OSP SUICIDE SAFETY PLAN VA-OSP SUICIDE SAFETY PLAN REVIEW/DECL VA-GP SBOR EVENT DATE	re selected INE	for packing:	
REMINDER TERM VA-REMINDER UPDATE_2_0_182 Keywords:			
+ Enter ?? for more actions IA Install all Components Select Action: Next Screen// IA	IS Instal	l Selected Compone	>>> ent

At the <u>Select Action</u> prompt, type <u>IA</u> for Install all Components and hit enter. Select Action: Next Screen// <u>IA</u> Install all Components You will see several prompts, for all new entries you will choose <u>I to Install</u> If you are prompted to overwrite any items, choose O to Overwrite.

You will be promped to install the FIRST reminder dialog component, YOU DO NOT NEED TO INSTALL THIS DIALOG, hit Q to quit:

Dialo	g Components	Jan 27, 2021@14:26:11	Page: 1	of 1
Packe	d reminder dialog: V	A-GP SBOR EVENT DATE [NATIONAL	DIALOG]	
Item	<u>Seq. Dialog Findin</u>	igs	Туре	Exists
1	VA-GP SBOR EVENT	DATE	group	×
	Finding: VA-S	BOR EVENT DATE (HEALTH FACTOR)		X
2	PXRMZ BLANK		prompt	X
_				
3	10 VA-EL SBOR EVEN	IT DATE APPROX	element	X
	Finding: *NON	IE*		
	+ Noxt Sopoop	Brov Sonoon 22 Mono Action		
DD	Dialog Dataila	DT Dialog Taxt	Thetall Salaata	vd.
DE	Dialog Eindings	DI Dialog Usago Ol	Oui+	a
	Dialog Summany	TA Install All	Guite	
	t Astissi Ouit// O			
Selec	t Action: Quit// Q			

At the Selection Action prompt, type Q to quit.

You will be promped to install the NEXT reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

Packe	d reminder dialog: VA-OSP SUICIDE SAFETY PLAN REVIEW/DECL	INE [NATIO	NAL DIAL
Item	Seq. Dialog Findings	Туре	Exists
1	VA-OSP SUICIDE SAFETY PLAN REVIEW/DECLINE	dialog	Х
2	1 VA-SUICIDE SAFETY PLAN REVIEW NOTE VERSION Finding: *NONE*	element	
3	5 SAFETY PLAN VIEW PRIOR REVIEW	element	
	Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL	FINDING)	Х
4	PXRM GF VIEW BUTTON	prompt	Х
5	10 VA-GP SAFETY PLAN DECLINE/REVIEW CHOICE Finding: *NONE*	group	x
6	10.10 VA-GP SAFETY PLAN REFUSES OFFER	group	Х
	Finding: VA-SP REFUSED SAFETY PLANNING (HEALTH FAC	TOR)	Х
7	PXRM COMMENT	prompt	Х
8	10.10.20 VA-GP DECLINE INFORMATION	group	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DD	Dialog Details DT Dialog Text IS Insta	11 Selecte	d
DF	Dialog Findings DU Dialog Usage QU Quit		
DS	Dialog Summary IA <mark>I</mark> nstall All		
Selec	t Action: Next Screen// IA		

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – <u>VA-OSP SUICIDE SAFETY</u> <u>PLAN REVIEW/DECLINE</u>

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// Yes

Packe	d reminder dialog: VA-OSP SUICIDE SAFETY PLAN REVIEW/DE	CLINE [NAT	IONAL DIAL
VA-0S	P SUICIDE SAFETY PLAN REVIEW/DECLINE (reminder dialog)	installed ⁻	from excha
Item	Seq. Dialog Findings	Туре	Exists
1	VA-OSP SUICIDE SAFETY PLAN REVIEW/DECLINE	dialog	Х
2	1 VA-SUICIDE SAFETY PLAN REVIEW NOTE VERSION Finding: *NONE*	element	x
3	5 SAFETY PLAN VIEW PRIOR REVIEW	element	x
	Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENER	AL FINDING) X
4	PXRM GF VIEW BUTTON	prompt	Х
5	10 VA-GP SAFETY PLAN DECLINE/REVIEW CHOICE Finding: *NONE*	group	х
6	10.10 VA-GP SAFETY PLAN REFUSES OFFER	group	Х
	Finding: VA-SP REFUSED SAFETY PLANNING (HEALTH F	ACTOR)	Х
7	PXRM COMMENT	prompt	Х
8	10.10.20 VA-GP DECLINE INFORMATION	group	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DD	Dialog Details DT Dialog Text IS Ins	tall Selec	ted
DF	Dialog Findings DU Dialog Usage QU Qui	.t	
DS	Dialog Summary IA Install All		
Selec	t Action: Next Screen// Q		

After completing this dialog install, you will type Q.

You will be promped to install the next reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

Packe	d reminder dialog: VA-OSP SUICIDE SAFETY PLAN [NATIONAL D	IALOG]	
Item	Seq. Dialog Findings	Туре	Exists
1	VA-OSP SUICIDE SAFETY PLAN	dialog	Х
2	1 VA-SUICIDE SAFETY PLAN NOTE VERSION Finding: *NONE*	element	Х
3	5 SAFETY PLAN VIEW PRIOR	element	
	Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL	FINDING)	Х
4	PXRM GF VIEW BUTTON	prompt	Х
5	10 VA-GP SAFETY PLAN INFORMATION Finding: *NONE*	group	Х
6	10.10 VA-EL SAFETY PLAN DECLINE/REVIEW INFO Finding: *NONE*	element	Х
7	10.18 SAFETY PLAN INTRODUCTION TIP	element	
	Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL	FINDING)	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DD I	Dialog Details DT Dialog Text IS Insta	ll Selecte	d
DF I	Dialog Findings DU Dialog Usage QU Quit		
DS I	Dialog Summary IA Install All		
Select	t Action: Next Screen// IA		

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – <u>VA-OSP SUICIDE SAFETY</u> <u>PLAN</u>

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// Yes

Packed	reminder dialog: VA-OSP SUICIDE SAFETY PLAN [NATIONAL	DIALOG]	
VA-OSP	SUICIDE SAFETY PLAN (reminder dialog) installed from e	exchange fil	e.
Item	Seq. Dialog Findings	Туре	Exists
1	VA-OSP SUICIDE SAFETY PLAN	dialog	X
2	1 VA-SUICIDE SAFETY PLAN NOTE VERSION Finding: *NONE*	element	х
3	5 SAFETY PLAN VIEW PRIOR	element	х
	Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERA	L FINDING)	Х
4	PXRM GF VIEW BUTTON	prompt	Х
5	10 VA-GP SAFETY PLAN INFORMATION Finding: *NONE*	group	X
6	10.10 VA-EL SAFETY PLAN DECLINE/REVIEW INFO Finding: *NONE*	element	x
7	10.18 SAFETY PLAN INTRODUCTION TIP	element	х
	Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERA	L FINDING)	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DD D	ialog Details DT Dialog Text IS Inst	all Selecte	d
DF D	ialog Findings DU Dialog Usage QU Quit		
DS D	ialog Summary IA Install All		
Select	Action: Next Screen// Q		

After completing this dialog install, you will type Q.

```
Category
      Component
                                                                          Exists
Source:
Date Packed: 08/24/2021@10:07:47
Package Version: 2.0P42
Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-OSP SUICIDE SAFETY PLAN
 VA-OSP SUICIDE SAFETY PLAN REVIEW/DECLINE
 VA-GP SBOR EVENT DATE
REMINDER TERM
 VA-REMINDER UPDATE_2_0_182
Keywords:
          + Next Screen - Prev Screen ?? More Actions
                                                                              >>>
IA
     Install all Components
                                         IS
                                              Install Selected Component
Select Action: Next Screen// Q
You will then be returned to this screen. At the Select Action prompt, type Q.
```

Install complete.

Post Installation

1. Format Prior Safety Plan health summary object display text.

To format the display text for easier reading, edit the Alternate Progress Note Text field of the elements SAFETY PLAN VIEW PRIOR <u>AND</u> SAFETY PLAN VIEW PRIOR REVIEW. Replace the current text |VA-OSP SAFETY PLAN NOTE| with the following text: 'FMT{|VA-OSP SAFETY PLAN NOTE|}FMT.

2. Add the button prompt to the element SAFETY PLAN VIEW PRIOR REVIEW. Edit the reminder dialog element SAFETY PLAN VIEW PRIOR REVIEW and add the following prompt so that the button will display properly:

```
SELECT SEQUENCE: 5

ADDITIONAL PROMPT/FORCED VALUE: PXRM GF VIEW BUTTON PROMPT

NATIONAL

...OK? YES// (YES)

OVERRIDE PROMPT CAPTION: SAFETY PLAN

START NEW LINE:

EXCLUDE FROM PN TEXT:

REQUIRED:

DISABLE PROMPT:
```

3. Update Health Summary Type VA-OSP SAFETY PLAN NOTE and add your local note title used for SUICIDE PREVENTION SAFETY PLAN NOTE.

SELECT BUILD HEALTH SUMMARY TYPE MENU <TEST ACCOUNT> OPTION: 1 CREATE/MODIFY HEALTH SUMMARY TYPE SELECT HEALTH SUMMARY TYPE: VA-OSP SAFETY PLAN NOTE VA-OSP SAFETY PLAN NOTE OK? YES// VA-OSP SAFETY PLAN NOTE WARNING: YOU ARE ABOUT TO EDIT A HEALTH SUMMARY TYPE THAT IS BEING USED BY A HEALTH SUMMARY OBJECT. CHANGING THE STRUCTURE OF THIS HEALTH SUMMARY TYPE WILL ALTER HOW THE OBJECT WILL DISPLAY. DO WANT TO CONTINUE? NO// YES

NAME: VA-OSP SAFETY PLAN NOTE REPLACE TITLE: SUPPRESS PRINT OF COMPONENTS WITHOUT DATA: SUPPRESS SENSITIVE PRINT DATA: LOCK: OWNER: GRIFFITH,ELIZABETH A//

DO YOU WISH TO COPY COMPONENTS FROM AN EXISTING HEALTH SUMMARY TYPE? YES// N NO SELECT COMPONENT: SPN PROGRESS NOTES SELECTED SPN SUMMARY ORDER: 5// 5 OCCURRENCE LIMIT: 1 TIME LIMIT: HEADER NAME: SELECTED PROG NOTES// RECENT SAFETY PLAN

NO SELECTION ITEMS CHOSEN.

SELECT NEW ITEMS ONE AT A TIME IN THE SEQUENCE YOU WANT THEM DISPLAYED. YOU MAY SELECT ANY NUMBER OF ITEMS.

SELECT SELECTION ITEM: ADD YOUR NOTE TITLE HERE

YOU SHOULD CLOSE AND RE-OPEN CPRS TO DISPLAY THE OBJECT PROPERLY.

4. Review template field VA-OSP SP MATERIAL Confirm the Address is: REDACTED

If that is not the address displaying, update what is in the address to field and choose Apply: REDACTED

5. Review from pre-install, add your site's quick order (or ordering menu) for ordering Naloxone as an additional finding to the following element: VA-EL SAFETY PLAN STEP 6 OPIOID KIT YES RX