



Update _2_0_192

Clinical Reminders

MENTAL HEALTH REMINDER UPDATES 8.0

Install Guide

June 2021

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION	1
PRE-INSTALL	1
INSTALL DETAILS	2
INSTALL EXAMPLE	2
POST INSTALLATION	19

Introduction

This reminder update basically is serving as a Mental Health miscellaneous update. There are 10 items that were updated to better server our Veterans.

Six reminder definitions and their associated dialog are being distributed with this update

VA-ALCOHOL USE SCREEN (AUDIT-C)
VA-ALCOHOL AUDIT-C POSITIVE F/U EVAL
VA-DEPRESSION SCREENING
VA-PTSD SCREENING
VA-PRIMARY SUICIDE SCREEN
VA-FOLLOW-UP POS PTSD/DEP COMBINATION

UPDATE_2_0_192 contains 1 Reminder Exchange entry:

UPDATE_2_0_192 MH UPDATES 8.0

The exchange file contains the following components:

Pre-Install

1. Sites may want to check for quick order(s) in the VA-FOLLOW-UP POS PTSD/DEP COMBINATION dialog and have those available when installing this dialog. Even though the name of the elements leads one to think this is only for PTSD, this should be a generic consult to MH for either a depression and/or a PTSD issue. Below are the elements you can check for Quick Orders.
VA-PTSD SCREEN POS MH EMERGENT EVAL
VA-PTSD SCREEN POS MH NON-EMERGENT EVAL
2. Sites may want to check for quick order(s) in the VA-AUD-C POS BRIEF ALCOHOL COUNSELING reminder dialog and have available when installing this dialog. Below is the element you can check for Quick Orders.
VA-ALC REFER TO SPECIALTY SUD PROG
3. Reminder update 37 (UPDATE_2_0_37), update 53 (UPDATE_2_0_53), update 63 (UPDATE_2_0_63), update 119 (UPDATE_2_0_119) should be considered a prerequisite to Update 192.

Install Details

This update is being distributed as a web host file. The address for the host file is:
https://REDACTED/UPDATE_2_0_192_1.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 20 minutes.

Install Example

To Load the exchange file, use LWH. The URL is https://Redacted/UPDATE_2_0_192_1.PRD

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CHF  Create Host File          LMM  Load MailMan Message
CMM  Create MailMan Message    LR   List Reminder Definitions
DFE  Delete Exchange File Entry LWH  Load Web Host File
IFE  Install Exchange File Entry RI   Reminder Definition Inquiry
IH   Installation History      RP   Repack
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: https://[REDACTED]/UPDATE_2_0_192_1.PRD
Added Reminder Exchange entry UPDATE_2_0_192 MH UPDATES 8.0
```

If you use LWH and it works, you will see something similar to this message. Make sure the file was loaded successfully.

```
https://[REDACTED]/UPDATE_2_0_192_1.PRD successfully loaded.
```

Most components will be the same and will be skipped. Here are general rules for prompts to which you will have to respond.

For Health Factors that exist, but the component is different, you should use the **SKIP** action. This is due to patch PX*1.0*211 being installed.

For Health Summary Types that exist, but the component is different, you should use the **OVERWRITE** action.

For Taxonomies that exist, but the component is different, you should use the **OVERWRITE** action

For all Terms, *except* VA-USER CLASS FOR COGNITIVE IMPAIRMENT, you should use the **SKIP** action. For VA-USER CLASS FOR COGNITIVE IMPAIRMENT, you should use the **MERGE** action.

For all Education Topics that exist, but the component is different, you should use the **SKIP** action

If a component is new, use **INSTALL** action

Search and locate an entry titled UPDATE_2_0_192 MH UPDATES 8.0 in reminder exchange. If you were a test site, you may have more than 1 entry with this name. Be sure and choose the entry with date of 05/04/2021

Item	Entry	Source	Date Packed
164	UPDATE_2_0_192 MH UPDATES 8.0		05/04/2021@12:38
165	UPDATE_2_0_193 VA-OSP TH INTAKE ASSESSMENT		03/19/2021@09:25
166	UPDATE_2_0_194 VA-COVID-19 CPRS STATUS VERSION 6		04/20/2021@05:54
167	UPDATE_2_0_195 VAAES SKIN INSPECTION/ASSESSMENT UPDATE		04/06/2021@13:49
168	UPDATE_2_0_197 VA-AES NON-VENTILATOR ORAL CARE UPDATE		04/09/2021@08:30

Action	Description
CFE	Create Exchange File Entry
CHF	Create Host File
CMM	Create MailMan Message
DFE	Delete Exchange File Entry
IFE	Install Exchange File Entry
IH	Installation History
LHF	Load Host File
LHM	Load MailMan Message
LR	List Reminder Definitions
LWH	Load Web Host File
RI	Reminder Definition Inquiry
RP	Repack

Select Action: Next Screen// IFE Install Exchange File Entry
Enter a list or range of numbers (1-403): 164

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry. Enter the number that corresponds with your entry titled UPDATE_2_0_192 MH UPDATES 8.0 (in this example it is entry 164, it will vary by site) that is dated 05/04/2021.

Component	Category	Exists
Source:		
Date Packed:	05/04/2021@12:38:12	
Package Version:	2.0P46	
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-GP ED TRIAGE SUICIDE IDEATION CSSRS		
REMINDER TERM		
VA-REMINDER UPDATE_2_0_192		
REMINDER DEFINITION		
VA-FOLLOW-UP POS PTSD/DEP COMBINATION		
VA-DEPRESSION SCREENING		
VA-PRIMARY SUICIDE SCREEN		
Enter ?? for more actions		
IA	Install all Components	
IS	Install Selected Component	

Select Action: Next Screen// IA

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.

You will see several prompts:

For Health Factors that exist, but the component is different, you should use the **SKIP** action.

This is due to patch PX*1.0*211 being installed.

For Health Summary Types that exist, but the component is different, you should use the **OVERWRITE** action.

For Taxonomies that exist, but the component is different, you should use the **OVERWRITE** action

For all Terms, *except* VA-USER CLASS FOR COGNITIVE IMPAIRMENT, you should use the **SKIP** action. For VA-USER CLASS FOR COGNITIVE IMPAIRMENT, you should use the **MERGE** action.

For all Education Topics that exist, but the component is different, you should use the **SKIP** action

If a component is new, use **INSTALL** action

Dialog installation may not necessarily appear in the order listed in the install guide

You will be prompted to install the first reminder dialog component:

Packed reminder dialog: VA-GP ED TRIAGE SUICIDE IDEATION CSSRS

Item	Seq.	Dialog Findings	Type	Exists
1		VA-GP ED TRIAGE SUICIDE IDEATION CSSRS Finding: *NONE*	group	X
2	5	VA-GP ED TRIAGE PERFORM C-SSRS Finding: *NONE*	group	X
3	5.5	VA-SUICIDE RISK ASSESSMENT C-SSRS ED NEW Finding: C-SSRS (MH TESTS AND SURVEY)	element	X
4		PXRM C-SSRS RESULT GROUP 3 ED	rs.group	
5		PXRM C-SSRS RESULT ELEMENT 3 ED	rs.element	
6		PXRM C-SSRS RESULT GROUP 4 ED	rs.group	
7		PXRM C-SSRS RESULT ELEMENT 4 ED	rs.element	
8		PXRM C-SSRS RESULT GROUP 5 ED	rs.group	
9		PXRM C-SSRS RESULT ELEMENT 5 ED	rs.element	
10		PXRM C-SSRS RESULT GROUP 8 ED	rs.group	
11		PXRM C-SSRS RESULT ELEMENT 8 ED	rs.element	

+ + Next Screen - Prev Screen ?? More Actions

DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		

Select Action: Next Screen// IA Install All

At the **Select Action** prompt, type **IA** to install the dialog group– VA-GP ED TRIAGE SUICIDE IDEATION CSSRS

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: **Y// YES**

Packed reminder dialog: VA-GP ED TRIAGE SUICIDE IDEATION CSSRS

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-GP ED TRIAGE SUICIDE IDEATION CSSRS Finding: *NONE*	group	X		
2	5	VA-GP ED TRIAGE PERFORM C-SSRS Finding: *NONE*	group	X		
3	5.5	VA-SUICIDE RISK ASSESSMENT C-SSRS ED NEW Finding: C-SSRS (MH TESTS AND SURVEY)	element	X		
4		PXRM C-SSRS RESULT GROUP 3 ED	rs.group	X		
5		PXRM C-SSRS RESULT ELEMENT 3 ED	rs.element	X		
6		PXRM C-SSRS RESULT GROUP 4 ED	rs.group	X		
7		PXRM C-SSRS RESULT ELEMENT 4 ED	rs.element	X		
8		PXRM C-SSRS RESULT GROUP 5 ED	rs.group	X		
9		PXRM C-SSRS RESULT ELEMENT 5 ED	rs.element	X		
10		PXRM C-SSRS RESULT GROUP 8 ED	rs.group	X		
11		PXRM C-SSRS RESULT ELEMENT 8 ED	rs.element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// Q						

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the next reminder dialog component:

Packed reminder dialog: VA-FOLLOW-UP POS PTSD/DEP COMBINATION

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-FOLLOW-UP POS PTSD/DEP COMBINATION	dialog	X		
2	2	VA-DEPRESSION/PTSD FOLLOW UP VERSION Finding: *NONE*	element	X		
3	5	VA-PTSD/DEPRESSION COMBINATION MASTER GP Finding: DEP/PTSD/SI ASSESSMENT COMPLETE (HEALTH FACTOR)	group	X		
4	5.5	VA-GP POS PTSD/DEP HEADER Finding: *NONE*	group	X		
5	5.5.5	VA-GP POS PTSD/DEP OBJECTS Finding: *NONE*	group	X		
6	5.5.5.5	VA-POS PTSD/DEP OBJECTS Finding: *NONE*	element	X		
7	5.10	VA-PTSD/DEP SCREEN POS DISPOSTIONS GP Finding: *NONE*	group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA Install All						

At the **Select Action** prompt, type **IA** to install the dialog – **VA-FOLLOW-UP PTSD/DEP COMBINATION**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// YES**

Packed reminder dialog: VA-FOLLOW-UP POS PTSD/DEP COMBINATION
VA-FOLLOW-UP POS PTSD/DEP COMBINATION (reminder dialog) installed from exchange

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-FOLLOW-UP POS PTSD/DEP COMBINATION	dialog	X		
2	2	VA-DEPRESSION/PTSD FOLLOW UP VERSION Finding: *NONE*	element	X		
3	5	VA-PTSD/DEPRESSION COMBINATION MASTER GP Finding: DEP/PTSD/SI ASSESSMENT COMPLETE (HEALTH FACTOR)	group	X		
4	5.5	VA-GP POS PTSD/DEP HEADER Finding: *NONE*	group	X		
5	5.5.5	VA-GP POS PTSD/DEP OBJECTS Finding: *NONE*	group	X		
6	5.5.5.5	VA-POS PTSD/DEP OBJECTS Finding: *NONE*	element	X		
7	5.10	VA-PTSD/DEP SCREEN POS DISPOSTIONS GP Finding: *NONE*	group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// Q						

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the next reminder dialog component:

Packed reminder dialog: VA-PRIMARY SUICIDE SCREEN C-SSRS [NATIONAL DIALOG]

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-PRIMARY SUICIDE SCREEN C-SSRS	dialog	X		
2	3	VA-PRIMARY SUICIDE SCREEN C-SSRS VERSION NUMBER Finding: *NONE*	element	X		
3	5	VA-PRIMARY SUICIDE SCREEN C-SSRS HEADER Finding: *NONE*	group	X		
4	10	VA-GP PRIMARY SUICIDE SCREEN C-SSRS* Finding: *NONE*	group	X		
5	10.5	VA-SUICIDE RISK ASSESSMENT C-SSRS* Finding: C-SSRS (MH TESTS AND SURVEY)	element	X		
6		PXRM C-SSRS RESULT GROUP 3	rs.group	X		
7		PXRM C-SSRS RESULT ELEMENT 3	rs.element	X		
8		PXRM C-SSRS RESULT GROUP 4	rs.group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// IA Install All						

At the **Select Action** prompt, type **IA** to install the dialog – **VA-PRIMARY SUICIDE SCREEN C-SSRS**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// YES**

VA-PRIMARY SUICIDE SCREEN C-SSRS (reminder dialog) installed from exchange file.

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-PRIMARY SUICIDE SCREEN C-SSRS	dialog	X		
2	3	VA-PRIMARY SUICIDE SCREEN C-SSRS VERSION NUMBER Finding: *NONE*	element	X		
3	5	VA-PRIMARY SUICIDE SCREEN C-SSRS HEADER Finding: *NONE*	group	X		
4	10	VA-GP PRIMARY SUICIDE SCREEN C-SSRS* Finding: *NONE*	group	X		
5	10.5	VA-SUICIDE RISK ASSESSMENT C-SSRS* Finding: C-SSRS (MH TESTS AND SURVEY)	element	X		
6		PXRM C-SSRS RESULT GROUP 3	rs.group	X		
7		PXRM C-SSRS RESULT ELEMENT 3	rs.element	X		
8		PXRM C-SSRS RESULT GROUP 4	rs.group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// Q						

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the next reminder dialog component:

Packed reminder dialog: VA-PTSD SCREENING [NATIONAL DIALOG]

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-PTSD SCREENING	dialog	X		
2	1	VA-PTSD SCREENING VERSION Finding: *NONE*	group	X		
3	1.5	VA-PTSD SCREEN Finding: *NONE*	group	X		
4	1.5.10	VA-MH PC-PTSD+I9 SCREEN Finding: PC-PTSD-5 (MH TESTS AND SURVEY)	element	X		
5		PXRM PC-PTSD-5 RESULT GROUP	rs.group			
6		PXRM PC-PTSD-5 RESULT ELEMENT 1	rs.element			
7		PXRM PC-PTSD-5 RESULT ELEMENT 2	rs.element			
8	1.5.20	VA-GP PTSD UNABLE TO SCREEN Finding: *NONE*	group	X		
9	1.5.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA						

At the **Select Action** prompt, type **IA** to install the dialog – **VA-PTSD SCREENING**

Select Action: Next Screen// IA **Install All**

Install reminder dialog and all components with no further changes: **Y// YES**

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

VA-PTSD SCREENING (reminder dialog) installed from exchange file.						
Item	Seq.	Dialog Findings	Type	Exists		
1		VA-PTSD SCREENING	dialog	X		
2	1	VA-PTSD SCREENING VERSION Finding: *NONE*	group	X		
3	1.5	VA-PTSD SCREEN Finding: *NONE*	group	X		
4	1.5.10	VA-MH PC-PTSD+I9 SCREEN Finding: PC-PTSD-5 (MH TESTS AND SURVEY)	element	X		
5		PXRM PC-PTSD-5 RESULT GROUP	rs.group	X		
6		PXRM PC-PTSD-5 RESULT ELEMENT 1	rs.element	X		
7		PXRM PC-PTSD-5 RESULT ELEMENT 2	rs.element	X		
8	1.5.20	VA-GP PTSD UNABLE TO SCREEN Finding: *NONE*	group	X		
9	1.5.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// Q						

You will be prompted to install the next reminder dialog component:

```
Packed reminder dialog: VA-PTSD SCREENING [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-PTSD SCREENING	dialog	X
2	1	VA-PTSD SCREENING VERSION Finding: *NONE*	group	X
3	1.5	VA-PTSD SCREEN Finding: *NONE*	group	X
4	1.5.10	VA-MH PC-PTSD+I9 SCREEN Finding: PC-PTSD-5 (MH TESTS AND SURVEY)	element	X
5		PXRM PTSD RESULT GROUP	rs.group	X
6		PXRM PTSD RESULT ELEMENT 1	rs.element	X
7		PXRM PTSD RESULT ELEMENT 2	rs.element	X
8	1.5.20	VA-GP PTSD UNABLE TO SCREEN Finding: *NONE*	group	X
9	1.5.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details      DT Dialog Text      IS Install Selected
DF Dialog Findings    DU Dialog Usage    QU Quit
DS Dialog Summary     IA Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – VA-PTSD SCREENING

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// YES**

```
Packed reminder dialog: VA-PTSD SCREENING [NATIONAL DIALOG]
VA-PTSD SCREENING (reminder dialog) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-PTSD SCREENING	dialog	X
2	1	VA-PTSD SCREENING VERSION Finding: *NONE*	group	X
3	1.5	VA-PTSD SCREEN Finding: *NONE*	group	X
4	1.5.10	VA-MH PC-PTSD+I9 SCREEN Finding: PC-PTSD-5 (MH TESTS AND SURVEY)	element	X
5		PXRM PTSD RESULT GROUP	rs.group	X
6		PXRM PTSD RESULT ELEMENT 1	rs.element	X
7		PXRM PTSD RESULT ELEMENT 2	rs.element	X
8	1.5.20	VA-GP PTSD UNABLE TO SCREEN Finding: *NONE*	group	X
9	1.5.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details      DT Dialog Text      IS Install Selected
DF Dialog Findings    DU Dialog Usage    QU Quit
DS Dialog Summary     IA Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the last reminder dialog component:

```
Packed reminder dialog: VA-DEPRESSION SCREEN [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-DEPRESSION SCREEN	dialog	X		
2	1	VA-DEPRESSION SCREENING VERSION Finding: *NONE*	group	X		
3	1.5	VA-DEP PHQ2 Finding: *NONE*	group	X		
4	1.5.5	VA-PHQ-2 Finding: PHQ-2 (MH TESTS AND SURVEY)	element	X		
5		PXRM PHQ2 RESULT GROUP	rs.group	X		
6		PXRM PHQ2 RESULT ELEMENT 1	rs.element	X		
7		PXRM PHQ2 RESULT ELEMENT 2	rs.element	X		
8	1.5.20	VA-GP DEP UNABLE TO SCREEN Finding: *NONE*	group	X		
9	1.5.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// IA						

At the **Select Action** prompt, type **IA** to install the dialog – **VA-DEPRESSION SCREENING**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// YES**

Screen shot of the terminal window showing a list of elements, groups and dialogs with the text with IA for "Install All" entered at the "Select Action:" Prompt. The correct entry is

described in the text. **THIS is for the VA-DEPRESSION SCREEN.**

VA-DEPRESSION SCREEN (reminder dialog) installed from exchange file.						
Item	Seq.	Dialog Findings	Type	Exists		
1		VA-DEPRESSION SCREEN	dialog	X		
2	1	VA-DEPRESSION SCREENING VERSION Finding: *NONE*	group	X		
3	1.5	VA-DEP PHQ2 Finding: *NONE*	group	X		
4	1.5.5	VA-PHQ-2 Finding: PHQ-2 (MH TESTS AND SURVEY)	element	X		
5		PXRM PHQ2 RESULT GROUP	rs.group	X		
6		PXRM PHQ2 RESULT ELEMENT 1	rs.element	X		
7		PXRM PHQ2 RESULT ELEMENT 2	rs.element	X		
8	1.5.20	VA-GP DEP UNABLE TO SCREEN Finding: *NONE*	group	X		
9	1.5.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// Q						

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the last reminder dialog component:

Packed reminder dialog: VA-ALCOHOL USE SCREENING (AUDIT-C) [NATIONAL DIALOG]

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-ALCOHOL USE SCREENING (AUDIT-C)	dialog	X		
2	1	VA-ALCOHOL SCREENING VERSION Finding: *NONE*	group	X		
3	1.5	VA-TEXT ALCOHOL MISUSE SCREEN INSTRUCTIONS Finding: *NONE*	element	X		
4	1.10	VAL-GP ALC SCREEN AUDC Finding: *NONE*	group	X		
5	1.10.5	VA-GP MH AUDC Finding: *NONE*	group	X		
6	1.10.5.5	VA-MH AUDC Finding: AUDC (MH TESTS AND SURVEY)	element	X		
7		PXRM AUDC RESULT GROUP	rs.group	X		
8		PXRM AUDC RESULT ELEMENT 1	rs.element	X		
9		PXRM AUDC RESULT ELEMENT 3	rs.element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// IA						

At the **Select Action** prompt, type **IA** to install the dialog – VA-ALCOHOL USE SCREENING (AUDIT-C)

Select Action: Next Screen// IA **Install All**

Install reminder dialog and all components with no further changes: **Y// YES**

VA-ALCOHOL USE SCREENING (AUDIT-C) (reminder dialog) installed from exchange fil

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-ALCOHOL USE SCREENING (AUDIT-C)	dialog	X		
2	1	VA-ALCOHOL SCREENING VERSION Finding: *NONE*	group	X		
3	1.5	VA-TEXT ALCOHOL MISUSE SCREEN INSTRUCTIONS Finding: *NONE*	element	X		
4	1.10	VAL-GP ALC SCREEN AUDC Finding: *NONE*	group	X		
5	1.10.5	VA-GP MH AUDC Finding: *NONE*	group	X		
6	1.10.5.5	VA-MH AUDC Finding: AUDC (MH TESTS AND SURVEY)	element	X X		
7		PXRM AUDC RESULT GROUP	rs.group	X		
8		PXRM AUDC RESULT ELEMENT 1	rs.element	X		
9		PXRM AUDC RESULT ELEMENT 3	rs.element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// Q						

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the last reminder dialog component:

```
Packed reminder dialog: VA-AUD-C POS BRIEF ALCOHOL COUNSELING [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-AUD-C POS BRIEF ALCOHOL COUNSELING	dialog	X
2	1	VA-POS ALCOHOL FOLLOW UP VERSION Finding: *NONE*	element	X
3	25	TEXT AUDC PRIOR RESULTS Finding: *NONE*	element	X
4	30	VA-AUDC SCORE GREATER THAN/EQUAL TO 5 Finding: *NONE*	element	X
5	35	VA-ALC SCREEN POS MEDICAL PROBLEMS Finding: *NONE*	element	X
6		VA-ALCOHOL BAC FEEDBACK OPTIONS	prompt	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-AUD-C POS BRIEF ALCOHOL COUNSELING**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// YES**

VA-AUD-C POS BRIEF ALCOHOL COUNSELING (reminder dialog) installed from exchange			
Item	Seq.	Dialog Findings	Type Exists
1		VA-AUD-C POS BRIEF ALCOHOL COUNSELING	dialog X
2	1	VA-POS ALCOHOL FOLLOW UP VERSION Finding: *NONE*	element X
3	25	TEXT AUDC PRIOR RESULTS Finding: *NONE*	element X
4	30	VA-AUDC SCORE GREATER THAN/EQUAL TO 5 Finding: *NONE*	element X
5	35	VA-ALC SCREEN POS MEDICAL PROBLEMS Finding: *NONE*	element X
6		VA-ALCOHOL BAC FEEDBACK OPTIONS	prompt X
+ + Next Screen - Prev Screen ?? More Actions			
DD	Dialog Details		DT Dialog Text IS Install Selected
DF	Dialog Findings		DU Dialog Usage QU Quit
DS	Dialog Summary		IA Install All
Select Action: Next Screen// Q			

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

This will complete the dialog installation section

Component	Category	Exists	
Source: [REDACTED]			
Date Packed: 05/04/2021@12:38:12			
Package Version: 2.0P46			
Description:			
The following Clinical Reminder items were selected for packing:			
REMINDER DIALOG			
VA-GP ED TRIAGE SUICIDE IDEATION CSSRS			
REMINDER TERM			
VA-REMINDER UPDATE_2_0_192			
REMINDER DEFINITION			
VA-FOLLOW-UP POS PTSD/DEP COMBINATION			
VA-DEPRESSION SCREENING			
VA-PRIMARY SUICIDE SCREEN			
+ + Next Screen - Prev Screen ?? More Actions >>>			
IA	Install all Components		IS Install Selected Component
Select Action: Next Screen// Q			

Type **Q** to exit the installation

Post Installation

1. Add your local Quick Orders back into the following elements from Pre-Install
 - a. VA-PTSD SCREEN POS MH EMERGENT EVAL
 - b. VA-PTSD SCREEN POS MH NON-EMERGENT EVAL
 - c. VA-ALC REFER TO SPECIALTY SUD PROG

Install complete.