



Update_2_0_200

Clinical Reminders

**VA-COVID-19 SCREEN UPDATE
Install Guide**

April 2021

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION 1
PRE-INSTALL..... 3
INSTALL DETAILS 3
INSTALL EXAMPLE 4
POST INSTALLATION 7

Introduction

DESCRIPTION:

The national Infectious Disease and Infection Control experts, along with input from other program offices and in review of feedback from the field, have updated the COVID-19 Screening template. The symptoms reflect the questions posed by the CDC more precisely. These changes are also being created in a tool for Cerner sites, which will be coming through the Cerner channels.

Immunization is being removed from the COVID-19 Screening template to provide consistency with the Cerner tool. It is being added in Update 194 as additional information in the banner.

The screening template will be updated in each of the locations where it is currently in use (Screen, ED Triage, CLC, SCI/D) as well as the symptoms being updated in additional templates where they currently exist (Telephone follow-up and Antibody testing).

The CLC and SCI/D have an unable to screen option for use with patients that are unable to answer the screening questions.

Reminder Dialog Template: SCREENING COVID

Coronavirus Disease 2019 (COVID-19) Screen version 7.1

Record the patient's responses to the following questions: (check all that apply, or None at the end)

In the last 14 days have you had new onset of any of the following symptoms?

- Chills Comment:
- Cough Comment:
- Diarrhea Comment:
- Fatigue Comment:
- Fever Comment:
- Headache Comment:
- Loss of taste or smell Comment:
- Muscle pain (myalgias) Comment:
- Nausea Comment:
- Runny nose (Rhinorrhea) Comment:
- Shortness of Breath (Dyspnea) Comment:
- Sore throat Comment:
- Vomiting Comment:
- No symptoms

Within the last 14 days, have you had:

- Close exposure (within 6 feet for more than 15 minutes) to someone with a febrile/respiratory illness
Comment:
- Close exposure (within 6 feet for more than 15 minutes) to someone with known or suspected case of COVID-19
Comment:
- No known exposure

Any symptom or exposure equal to positive screen

- Patient has a POSITIVE symptom or exposure and requires further evaluation Nurse/Provider/Other notified
- Screen is negative

- Patient is waiting on COVID-19 test results. Comment:
- Patient reports prior COVID-19 Diagnosis. Comment:

UPDATE_2_0_200 contains 1 Reminder Exchange entry:
UPDATE_2_0_200 VA-COVID-19 SCREEN UPDATE

The exchange file contains the following components:

TIU TEMPLATE FIELD

VA-COVID-COVID-19 AB INFORMATION SHEET

HEALTH FACTORS

VA-TRAVEL SCREEN [C]
VA-EXPOSURE TO COVID NO
VA-EXPOSURE TO COVID YES
VA-COVID-19
VA-EXPOSURE TO FEVER/RESP ILL 14D
VA-COVID-19 NO SYMPTOMS
VA-SYMPTOM VOMITING
VA-SYMPTOM SORE THROAT
VA-SYMPTOM SHORTNESS OF BREATH
VA-SYMPTOM RUNNY NOSE
VA-SYMPTOM NAUSEA
VA-SYMPTOM MUSCLE PAIN
VA-SYMPTOM LOSS OF TASTE OR SMELL
VA-SYMPTOM HEADACHE
VA-SYMPTOM FEVER
VA-SYMPTOM FATIGUE
VA-SYMPTOM DIARRHEA
VA-SYMPTOM COUGH
VA-SYMPTOM CHILLS
VA-COVID-19 OUTSIDE PCR SPEC COLLECTED
VA-COVID-19 NON PCR AG LAB OUTSIDE NEG
VA-COVID-19 UNK OUTSIDE LAB TEST NEG
VA-COVID-19 PCR LAB OUTSIDE NEGATIVE
VA-COVID-19 NON PCR AG LAB OUTSIDE POS
VA-COVID-19 UNK OUTSIDE LAB TEST POS
VA-COVID-19 PCR LAB OUTSIDE POSITIVE
VA-COVID-19 PRIOR DIAGNOSIS
VA-COVID WAITING TEST RESULTS
VA-COVID SCREEN NEGATIVE
VA-COVID SCREEN POSITIVE
VA-TRAVEL 21 DAYS UNABLE TO ANSWER
VA-NO HEADACHE LOSS OF TASTE OR SMELL
VA-SYMPTOM LOSS OF TASTE
VA-SYMPTOM LOSS OF SMELL

REMINDER SPONSOR

Office of Geriatric Extended Care

REMINDER TERM

VA-REMINDER UPDATE_2_0_200
VA-COVID-19 SCREEN COMPLETE TODAY

REMINDER DIALOG

VA-COVID-19 CLC SCREEN
VA-GP COVID-19 SCREEN UPDATE 7

Pre-Install

Using the Vista menu Reminder Dialog Manager, Reminder Dialog, change your view to GROUP and choose inquire/print for the reminder dialog group, VAL-COVID-19 CLC SCREEN LOCAL GROUP.

Save this info for the post install, it will make adding information back post -install easier if your site has made local edits. Not all sites will have made local edits to this dialog groups.

Install Details

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 5 minutes. You will need to allow additional time for the pre and post install steps.

This update is being distributed as a web host file. The address for the host file is:

[https://REDACTED /docs/UPDATE_2_0_200.PRD](https://REDACTED/docs/UPDATE_2_0_200.PRD)

The file will be installed using Reminder Exchange, programmer access is not required.

Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry       RI   Reminder Definition Inquiry
IH   Installation History              RP   Repack
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: https://[REDACTED]/docs/UPDATE_2_0_200.PRD
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:
https://REDACTED /UPDATE_2_0_200.PRD

You should see a message at the top of your screen that the file loaded successfully

```
https://[REDACTED]/UPDATE_2_0_200.PRD successfully loaded.
```

Search and locate an entry titled **UPDATE_2_0_200 VA-COVID-19 SCREEN UPDATE** in reminder exchange.

```
+Item  Entry                                     Source                                     Date Packed
272    UPDATE_2_0_200 VA-COVID-19
      SCREEN UPDATE
273    UPDATE_2_0_21 VA-BATTLEFIELD
      ACUPUNCTURE
274    UPDATE_2_0_22 VA-GATEWAY TO
      HEALTHY LIVING
275    UPDATE_2_0_23 VA-NALOXONE USE
276    UPDATE_2_0_24 MISC UPDATE 3
277    UPDATE_2_0_26 VA-HEAD AND
      NECK CANCER DIALOGS
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry       RI   Reminder Definition Inquiry
IH   Installation History              RP   Repack
Select Action: Next Screen// ife  Install Exchange File Entry
Enter a list or range of numbers (1-511): 272
```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry **UPDATE_2_0_200 VA-COVID-19 SCREEN UPDATE** (in this example it is entry 272 it will vary by site). The date of the

exchange file should be 04/01/2021.

```

Component                                     Category      Exists
Source:
Date Packed: 04/01/2021@04:33:02
Package Version: 2.0P46

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-COVID-19 CLC SCREEN
  VA-GP COVID-19 SCREEN UPDATE 7

REMINDER TERM
  VA-REMINDER UPDATE_2_0_200

Keywords:

Components:
+      Enter ?? for more actions                >>>
IA  Install all Components                      IS  Install Selected Component
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.
 Select Action: Next Screen// **IA Install all Components**
 You will see several prompts, for all new entries you will choose **I to Install**
*For components that already exists but the packed component is different, you will choose **S to Skip**.*

You will be prompted to install the FIRST reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

```

Packed reminder dialog: VA-GP COVID-19 SCREEN UPDATE 7 [NATIONAL DIALOG]

```

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-GP COVID-19 SCREEN UPDATE 7 Finding: *NONE*	group			
2	5	VA-GP COVID SYMPTOMS AND TRAVEL Finding: *NONE*	group	X		
3	5.5	VA-GP COVID UPDATED SCREENING QUESTIONS Finding: *NONE*	group			
4	5.5.5	VA-COVID SYMPTOM CHILLS Finding: VA-SYMPTOM CHILLS (HEALTH FACTOR)	element	X		
5	5.5.10	VA-COVID SYMPTOM COUGH Finding: VA-SYMPTOM COUGH (HEALTH FACTOR)	element	X		
6	5.5.15	VA-COVID SYMPTOM DIARRHEA Finding: VA-SYMPTOM DIARRHEA (HEALTH FACTOR)	element	X		
7	5.5.20	VA-COVID SYMPTOM FATIGUE Finding: VA-SYMPTOM FATIGUE (HEALTH FACTOR)	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		

```

Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-GP COVID-19 SCREEN UPDATE 7**

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// Yes

```
Packed reminder dialog: VA-GP COVID-19 SCREEN UPDATE 7 [NATIONAL DIALOG]
VA-GP COVID-19 SCREEN UPDATE 7 (group) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-GP COVID-19 SCREEN UPDATE 7 Finding: *NONE*	group	X
2	5	VA-GP COVID SYMPTOMS AND TRAVEL Finding: *NONE*	group	X
3	5.5	VA-GP COVID UPDATED SCREENING QUESTIONS Finding: *NONE*	group	X
4	5.5.5	VA-COVID SYMPTOM CHILLS Finding: VA-SYMPTOM CHILLS (HEALTH FACTOR)	element	X
5	5.5.10	VA-COVID SYMPTOM COUGH Finding: VA-SYMPTOM COUGH (HEALTH FACTOR)	element	X
6	5.5.15	VA-COVID SYMPTOM DIARRHEA Finding: VA-SYMPTOM DIARRHEA (HEALTH FACTOR)	element	X
7	5.5.20	VA-COVID SYMPTOM FATIGUE Finding: VA-SYMPTOM FATIGUE (HEALTH FACTOR)	element	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

After completing this dialog install, you will type Q.

You will be prompted to install the next reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

```
Packed reminder dialog: VA-COVID-19 CLC SCREEN [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-COVID-19 CLC SCREEN	dialog	X
2	5	VA-COVID-19 CLC SCREEN VERSION Finding: *NONE*	element	X
3	10	VA-COVID-19 CLC COVID SCREEN HEADER Finding: *NONE*	group	X
4	10.3	VAL-COVID-19 CLC TEMP OBJECT Finding: *NONE*	element	
5	10.5	VA-COVID-19 CLC TEMPERATURE Finding: TEMPERATURE (GMRV VITAL TYPE)	group	X
6	10.10	VA-GP COVID-19 SCREEN CLC Finding: *NONE*	group	X
7	10.10.10	VA-GP COVID-19 SCREEN COMPLETE CLC Finding: *NONE*	group	

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-COVID-19 CLC SCREEN**

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// Yes


```

Packed reminder dialog: VA-COVID-19 CLC SCREEN [NATIONAL DIALOG]
VA-COVID-19 CLC SCREEN (reminder dialog) installed from exchange file.
Item  Seq.  Dialog Findings                                     Type  Exists
  1      VA-COVID-19 CLC SCREEN                               dialog  X
  2      5 VA-COVID-19 CLC SCREEN VERSION                       element  X
        Finding: *NONE*
  3      10 VA-COVID-19 CLC COVID SCREEN HEADER                 group   X
        Finding: *NONE*
  4      10.3 VAL-COVID-19 CLC TEMP OBJECT                     element  X
        Finding: *NONE*
  5      10.5 VA-COVID-19 CLC TEMPERATURE                       group   X
        Finding: TEMPERATURE (GMRV VITAL TYPE)
  6      10.10 VA-GP COVID-19 SCREEN CLC                       group   X
        Finding: *NONE*
  7      10.10.10 VA-GP COVID-19 SCREEN COMPLETE CLC          group   X
        Finding: *NONE*
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// Q

```

After completing this dialog install, you will type Q.

```

Component                                     Category  Exists
Source:
Date Packed: 04/01/2021@04:33:02
Package Version: 2.0P46

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-COVID-19 CLC SCREEN
  VA-GP COVID-19 SCREEN UPDATE 7

REMINDER TERM
  VA-REMINDER UPDATE_2_0_200

Keywords:

Components:
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components      IS  Install Selected Component
Select Action: Next Screen// Q

```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.

Post Installation

1. Open/review the COVID-19 Screen Template (VA-SCREEN COVID) and confirm the version is 7.1
2. Open/review the CLC COVID-19 SCREEN (VA-COVID-19 CLC SCREEN) and confirm the version is 4.1
3. Open/review the COVID-19 Telephone Follow-Up (VA-COVID-19 TELEPHONE FOLLOW-UP) and confirm the version is 6.1.

4. **Open/review the SCI/D COVID-19 SCREEN (VA-COVID-19 SCI/D SCREEN) and confirm the version is 4.1**
5. **Open/review the COVID-19 ANTIBODY/PCR PANEL TESTING template (VA-COVID-19 ANTIBODY/PCR PANEL TESTING) and confirm the version is 2.1.**
6. **Open/review the EMERGENCY ROOM TRIAGE template and confirm the COVID-19 screen embedded in the template is version 7.1**
7. **Add local content to reminder dialog groups and elements if needed.**
(Optional)

Review the pre-install for the reminder dialog group, VAL-COVID-19 CLC SCREEN LOCAL GROUP and make any changes that are appropriate for your site.