

Update\_2\_0\_216

## **Clinical Reminders**

## VA-COLORECTAL CANCER SCREENING UPDATE

## **Install Guide**

July 2021

Product Development Office of Information Technology Department of Veterans Affairs

# **Table of Contents**

INTRODUCTION	1
PRE-INSTALL	5
INSTALL DETAILS	5
INSTALL EXAMPLE	6
POST INSTALLATION	17

### Introduction

The four colorectal screening reminders were released in 2017 as Clinical reminder Update\_2\_0\_9 and updated with Update\_2\_0\_94. This reminder update includes updates/additions to the VA-CRC Average Risk reminder definition.

#### The following changes are in this update:

1. The VA-CRC AVERAGE RISK baseline frequency changed to 45yo to meet new guidelines to begin CRC screening at age 45.

neminder Inquiry: Avg Risk Colorectal Cancer Screen					
Baseline Frequency:					
Do In Advance Time Frame: Sex Specific:	Do if DUE within 2 months				
Ignore on N/A: Frequency for Age Range:	l year for ages 45 to 75				
No Match Text:					

2. All reminder dialog options with branching logic based on age have been updated to the new guidelines. The options below will now display for patients 35yo and older, instead of beginning at age 40. These options are in the VA-CRC COLONOSCOPY AVG RISK, VA-CRC WAITING FOR RECOMMEND GAP, VA-CRC COLONOSCOPY FOLLOW UP REMINDER, VA-CRC COLON SCREENING DOCUMENTATION reminder dialogs.

Note: The age range has changed in the following terms, however, the reminder term name <u>has not</u> changed. Changing the name to match the new age range would have installed new terms, which would have required remapping throughout the dialogs. The decision was made to keep the existing names to reduce the amount of remapping at local sites.

VA-AGE <50 VA-AGE <40!>74 VA-AGE <40!>65 VA-AGE <40!>66 VA-AGE <40!>67 VA-AGE <40!>68 VA-AGE <40!>69 VA-AGE <40!>70 VA-AGE <40!>71 VA-AGE <40!>72 VA-AGE <40!>73 VA-AGE <40

₱ Establish a new due date for average risk colorectal cancer screening or Colonoscopy reminders (justification required)
**Caution** Do not use this option for refusals or deferrals. Use only if this reminder is not set correctly for this patient.
Justification:
•
Set average risk screening reminder (Prompts a provider to offer screening options such as FIT or Colonoscopy)
**Critical Note** Reminders will be set from TODAY, JUN 01, 2021 not from the date of the procedure/test.
C Average Risk screening due now.
C 6 months from TODAY, JUN 01, 2021
🖸 1 year from TODAY, JUN 01, 2021
🖸 1.5 years from TODAY, JUN 01, 2021
🖸 2 years from TODAY, JUN 01, 2021
C 2.5 years from TODAY, JUN 01, 2021
🖸 3 years from TODAY, JUN 01, 2021
C 4 years from TODAY, JUN 01, 2021
🖸 5 years from TODAY, JUN 01, 2021
C 6 years from TODAY, JUN 01, 2021
C 7 years from TODAY, JUN 01, 2021
🖸 8 years from TODAY, JUN 01, 2021
□ 9 years from TODAY, JUN 01, 2021
□ 10 years from TODAY, JUN 01, 2021

3. In this reminder dialog option to set the average risk screening reminder, the option has been changed from "currently 50 yo" to "currently 45 yo." <u>This option will only display</u> <u>for patients less than 35 yo.</u> This option is in the VA-CRC COLONOSCOPY AVG RISK, VA-CRC WAITING FOR RECOMMEND GAP, VA-CRC COLONOSCOPY FOLLOW UP REMINDER, VA-CRC COLON SCREENING DOCUMENTATION.



4. In this reminder dialog option to set the colonoscopy reminder, the option has been changed from "currently 50 yo" to "currently 45 yo." <u>This option will only display for</u> <u>patients less than 35 yo.</u> This option is in the VA-CRC COLONOSCOPY AVG RISK, VA-CRC WAITING FOR RECOMMEND GAP, VA-CRC COLONOSCOPY FOLLOW UP REMINDER, VA-CRC COLON SCREENING DOCUMENTATION.

```
♥ Establish a new due date for average risk colorectal cancer screening or Colonoscopy reminders
 (justification required)
**Caution** Do not use this option for refusals or deferrals. Use only if this reminder is not set
 correctly for this patient.
Justification:
     *
  Set average risk screening reminder (Prompts a provider to offer screening options such as FIT or
    Colonoscopy)
  🗈 Set Colonoscopy reminder (Specifically prompts a Colonoscopy)
    **Critical Note**
    Reminders will be set from TODAY, JUN 01, 2021 not from the date of the procedure/test.
    C 3 months or less from TODAY, JUN 01, 2021
    □ 6 months from TODAY, JUN 01, 2021
    □ 1 year from TODAY, JUN 01, 2021
    □ 1.5 years from TODAY, JUN 01, 2021
    □ 2 years from TODAY, JUN 01, 2021
    C 2.5 years from TODAY, JUN 01, 2021
    C 3 years from TODAY, JUN 01, 2021
    C 4 years from TODAY, JUN 01, 2021
    □ 5 years from TODAY, JUN 01, 2021
    C 6 years from TODAY, JUN 01, 2021
    C 7 years from TODAY, JUN 01, 2021
    C 8 years from TODAY, JUN 01, 2021
    C 9 years from TODAY, JUN 01, 2021
    □ 10 years from TODAY, JUN 01, 2021
    C Repeat colonoscopy at the age meeting USPSTF guidelines, currently 45yo.
```

5. In the reminder dialog option below, the option in the patient letter portion of the VA-CRC COLON SCREENING DOCUMENTATION has been changed from "age 50" to "age 45." <u>This option will only display for patients less than 45 yo.</u>

```
Patient letter: After review of pathology
Dear C Mr. C Mrs. BEAVIS, MISTER:
Thank you for completing your *
                                                                               on *
                                                                                               - at the
**Edit the VA-VAMC NAME TIU template field** VA Medical Center. I am writing to give you the pathology
 results for the tissue removed during the procedure.
🖻 Click here if communicating polyp biopsy results.
  \Box I am happy to report that the tissue we removed is not cancer. This is good news.
    ***If CANCER then you should call the patient and document this communication in a CPRS telephone note***
  The tissue removed included:
    □ Small (Less that 10 millimeters) hyperplastic polyp(s).
    One or more tubular adenoma.
    One or more Advanced Adenoma.
    □ One or more sessile serrated adenoma/polyp.
    Only normal tissue
    C Other: *
    \square --Click here to add free text comment
  After reviewing your results I recommend:
   C You discuss with your provider options for completing colon cancer screening at age 45.
     O You complete a follow-up colonoscopy in approximately:
    🖸 I do not recommend any future colon cancer screening because of your age and/or other medical
       conditions.
```

# It will be very important that you follow the install guide carefully during the install.

The exchange file contains the following components:

#### **REMINDER DIALOGS**

VA-CRC AVG RISK UPDATE JUNE 2021

**REMINDER DEFINITIONS** VA-CRC AVERAGE RISK

#### **REMINDER TERMS**

VA-REMINDER UPDATE\_2\_0\_216 VA-AGE <50 VA-AGE <40!>74 VA-AGE <40!>65 VA-AGE <40!>66 VA-AGE <40!>67 VA-AGE <40!>68 VA-AGE <40!>69 VA-AGE <40!>70 VA-AGE <40!>71 VA-AGE <40!>72 VA-AGE <40!>73 VA-AGE <40

#### **Pre-Install**

- 1. Using the Vista menu, Reminder Managers Menu, inquire about Reminder Definition, run an inquire/print on the following reminder. **Remapping should not be required if the appropriate items are selected for installation.** However, if inappropriate terms are accidentally overwritten, a copy of the definition will be helpful.
  - a. VA-CRC AVERAGE RISK
- Using the Vista menu, Reminder Dialog Management, Reminder Dialogs, change your view to dialog view and run an inquire/print on the following reminder dialogs.
   Remapping should not be required if the appropriate items are selected for installation. However, if inappropriate groups or elements are accidentally overwritten, a copy of the dialog will be helpful.
  - a. VA-CRC COLONOSCOPY AVG RISK

### **Install Details**

This update is being distributed as a web host file. The address for the host file is: <u>https://REDACTED/UPDATE\_2\_0\_216.PRD</u>

This update can be loaded with users on the system. Installation should take less than 5 minutes.

# Follow the instructions below carefully. Only a few items in the exchange file will be installed. DO NOT INSTALL ALL.

At the <u>Select Action</u>: prompt, enter <u>LWH</u> for Load Web Host File At the <u>Input the url for the .prd file</u>: prompt, type the following web address:

#### https://REDACTED/UPDATE 2 0 216.PRD

Search and locate an entry titled UPDATE\_2\_0\_216 VA-CRC AVERAGE RISK UPDATE

+Item	Entry	Source	e Date Packed
200	UPDATE_2_0_216 VA-CRC AVERAG	E	06/01/2021@08:27
	RISK UPDATE		
201	UPDATE_2_0_22 VA-GATEWAY TO		05/12/2017@09:36
	HEALTHY LIVING		
202	UPDATE_2_0_23 VA-NALOXONE US	E	08/18/2017005:23
203	UPDATE_2_0_23 VA-NALOXONE US	E	08/25/2017@06:21
204	UPDATE_2_0_23 VA-NALOXONE US	E	10/13/2017@12:09
205	UPDATE_2_0_24 MISC UPDATE 3		05/08/2017@04:44
206	UPDATE_2_0_26 VA-HEAD AND		08/22/2017013:17
	NECK CANCER DIALOGS		
+	+ Next Screen - Prev Sc	reen ??	More Actions >>>
CFE	Create Exchange File Entry	LHF	Load Host File
CHF	Create Host File	LMM	Load MailMan Message
CMM	Create MailMan Message	LR	List Reminder Definitions
DFE	Delete Exchange File Entry	LWH	Load Web Host File
IFE	Install Exchange File Entry	RI	Reminder Definition Inquiry
IH	Installation History	RP	Repack
Selec	t Action: Next Screen// IFE	Install Ex	xchange File Entry
Enter	a list or range of numbers (1	-450): 200	0

At the <u>Select Action</u> prompt, enter <u>IFE</u> for Install Exchange File Entry Enter the number that corresponds with your entry UPDATE\_2\_0\_216 VA-CRC AVERAGE **RISK UPDATE.** *(in this example it is entry 200 it will vary by site).* The date of the exchange file should be 06/01/2021.

Follow the screenshots below to install selected items. **Do not select Install All.** 

Follow the instructions below carefully. Only a few items in the exchange file will be installed. DO NOT INSTALL ALL.

Component	Category E	Exists
Source:		
Date Packed: 06/01/2021@08:27:11		
Package Version: 2.0P42		
Description.		
The following Clinical Reminder items were se	elected for packing:	
REMINDER DIALOG		
VA-CRC WAITING FOR RECOMMEND GAP		
VA-CRC COLONOSCOPY AVG RISK		
VA-CRC COLONOSCOPY DIAGNOSTIC REMINDER		
VA-CRC AVG RISK UPDATE JUNE 2021		
REMINDER TERM		
VA-AGE <40		
VA-AGE <40!>74		
VA-AGE <40!>73		
+ Enter ?? for more actions		>>>
IA Install all Components	Tretall Colocted Component	
Select Action: Next Screen// Calast Enter Day	mot Tractal1 A 11	
Select Effer. Do I	not ilistali Ali.	

+ Co	mponent	Category	Exists
VA-AGE	<40!>72		1.1.1
VA-AGE	<40!>71		
VA - AGE	<40!>70		
VA-AGE	<40!>69		
VA-AGE	<40!>68		
VA-AGE	<40!>67		
VA-AGE	<40!>66		
VA - AGE	<40!>65		
VA - AGE	<50		
VA-REM	INDER UPDATE_2_0_216		
REMINDER	DEFINITION		
VA-CRC	AVERAGE RISK		
Non-exch	angeable TIU object(s):		
TIU Obje	ct: TODAY'S DATE		
+	+ Next Screen - Prev Screen ?? More Actions		>>>
IA Ins	tall all Components IS Install Selec	ted Component	
Select A	ction: Next Screen// Select Enter Do not Install A1		
	Sciect Enter. Do not instant Al.		



Excha	nge File Components Jun 01, 2021@12:2	0:23 Page:	4 of 18
+	Component	Category	Exists
1	VA-CRC NCP HYPERLINK		X
2	VA-CRC RESULT TYPES		X
3	VA-CRC RESULT TYPES C		X
4	TEXT 40 CHAR (REQ)		X
5	VA-OTHER COMMENT (EXC FROM NOTE)		X
EDUCA	TION TOPICS VA-CRC DIAGNOSTIC BENEFITS		x
7	VA-CRC SCREEN BENEFITS/RISK DECLINING	Install item 9 from this s	creen.
HEALT	H FACTORS	Ensure to select the corre	ect item
8	VA-REMINDER UPDATES [C]	much as fee the VA	
9	VA-UPDATE_2_0_216	number for the VA-	
10	VA-CHC HF CATEGORY [C] VA-CHC AVG RISK SET FUTURE USPTF AGE	UPDATE_2_0_216 heal	th factor.
12	VA-CRC PROMPT COLONOSCOPY USPTF AGE		X
+	+ Next Screen - Prev Screen ?? Mo	re Actions	>>>
IA Selec	Install all Components IS In t Action: Next Screen//	stall Selected Compone	nt

At the HEALTH FACTORS section, select **IS**, then choose the item number 9 from the screenshot above.

Excha	inge File Components Jun 01, 2021012:21:17	Page:	5 of	18
+	Component	Category	Exi	sts
13	VA-CRC DIAGNOSTIC WAITING		X	
14	VA-CRC DIAGNOSTIC COLON NOT INDICATED		×	
15	VA-CRC DIAGNOSTIC DEFER 1 M		×	
16	VA-CRC DIAGNOSTIC DECLINED		X	
17	VA-CRC DIAGNOSTIC OUTSIDE CARE		×	
18	VA-CRC PROVIDE OUTSIDE RESULT		×	
19	VA-CRC DIAGNOSTIC SCHEDULED SCRN		×	
20	VA-CRC COLON SCREEN NOT INDICATED		×	
21	VA-CRC DEFER OTHER		×	
22	VA-CRC MALIGNANCY CONFIRMED BY PATH		X	
23	VA-CRC ORDER DCBE PREP KIT		×	
24	VA-CRC GAP ORDER DCBE		×	
25	VA-CRC ORDER CTC PREP KIT		×	
26	VA-CRC GAP ORDER CTC		×	
27	VA-CRC TISSUE SENT TO PATHO		×	
28	VA-CRC PROMPT COLONOSCOPY 2.5YR		×	
+	+ Next Screen - Prev Screen ?? More Actions			>>>
IA	Install all Components IO Install Colors	tod Compone	in th	
Selec	t Action: Next Screen// Select Enter. Do not install any ite	ems displaye	d on this	scree
i i				

+	Compone	nt						
+	Compone	nt						
	MA-CBC					Category	Exi	.sts
29	WA-Ono	PROMPT	COLONOSCOPY	1.5YR			>	
30	VA-CRC	PROMPT	COLONOSCOPY	′ 10YR			>	(
31	VA-CRC	PROMPT	COLONOSCOPY	9YR			>	(
32	VA-CRC	PROMPT	COLONOSCOPY	' 8YR			>	(
33	VA-CRC	PROMPT	COLONOSCOPY	7YR			>	(
34	VA-CRC	PROMPT	COLONOSCOPY	6YR			>	(
35	VA-CRC	PROMPT	COLONOSCOPY	′ 5YR			>	(
36	VA-CRC	PROMPT	COLONOSCOPY	′ 4YR			>	(
37	VA-CRC	PROMPT	COLONOSCOPY	' 3YR			>	( ) ( ) ( ) ( ) ( ) ( ) ( ) ( ) ( ) ( )
38	VA-CRC	PROMPT	COLONOSCOPY	2YR			>	t i
39	VA-CRC	FOLLOW	UP ORDER CO	LONOSCOPY			>	ί.
40	VA-CRC	PROMPT	COLONOSCOPY	′ 1YR			>	( ) ( )
41	VA-CRC	FOLLOW	UP DECLINED	)			>	( ) ( )
42	VA-CRC	PROMPT	COLONOSCOPY	′ 6M			>	( ) ( ) ( ) ( ) ( ) ( ) ( ) ( ) ( ) ( )
43	VA-CRC	PROMPT	COLONOSCOPY	′ 3M			>	( ) ( ) ( ) ( ) ( ) ( ) ( ) ( ) ( ) ( )
44	VA-CRC	DIAGNOS	STIC RECORD	PRIOR			>	(
+	+ N	lext Scr	reen - Pre	v Screen	?? More Actio	ns		>>>
IA 1	Install	all Con	nponents		IS Install Se	lected Compone	nt	
Select	t Action	1: Next	Screen// S	elect Enter.	. <u>Do not</u> install an	y items displayed	d on this	scree

Excha	nge File ComponentsJun 01, 2021@12:22:05	Page:	7 of	18
+	Component	Category	Exis	sts
45	VA-CRC RECORD COLONOSCOPY		X	
46	VA-CRC RECORD SIGMOID		x	
47	VA-CRC RECORD CT COLONOGRAPHY		x	
48	VA-CRC AVG RISK DUE NOW		x	
49	VA-CRC AVG RISK 6 M		x	
50	VA-CRC AVG RISK 2.5 YR		×	
51	VA-CRC AVG RISK 1.5 YR		X	
52	VA-CRC SET FUTURE AVG RISK STOP SCREEN		×	
53	VA-CRC AVG RISK 10 YR		X	
54	VA-CRC AVG RISK 9 YR		x	
55	VA-CRC AVG RISK 8 YR		X	
56	VA-CRC AVG RISK 7 YR		X	
57	VA-CRC AVG RISK 6 YR		×	
58	VA-CRC AVG RISK 5 YR		X	
59	VA-CRC AVG RISK 4 YR		X	
60	VA-CRC AVG RISK 3 YR		X	
+	+ Next Screen - Prev Screen ?? More Actions			>>>
IA	Install all Components	tod Compone		
Selec	t Action: Next Screen// Select Enter. Do not install any it	ems displaye	d on this	scree

Excha	ange File Components Jun 01, 2021012:22:29	Page:	8 of 18	
+	Component	Category	Exists	
61	VA-CRC AVG RISK 2 YR		х	
62	VA-CRC AVG RISK 1 YR		X	
63	VA-CRC DIAGNOSTIC ORDER COLONOSCOPY		X	
64	VA-CRC REQ ECONSULT		X	
65	VA-CRC STOP SCREEN TOLERATE		X	
66	VA-CRC STOP SCREEN OTHER		×	
67	VA-CRC <5 YR LIFE EXPECT		X	
68	REMINDER FACTORS [C]	×	×	
69	LIFE EXPECTANCY < 6 MONTHS		×	
70	VA-CRC TOTAL COLECTOMY		X	
71	VA-CRC AVG RISK DEFER 3M		X	
72	VA-CRC AVG RISK WAITING		X	
73	VA-CRC AVG RISK SCRN DECLINED		X	
74	VA-CRC CTC ORDER REASON: NOT COMPLETED		X	
75	VA-CRC CTC ORDER REASON: RESEARCH		X	
76	VA-CRC CTC ORDER REASON: PT DECLINED		X	
+	+ Next Screen - Prev Screen ?? More Actions		>>>	
IA	Install all Components	od Compone	nt	
Seled	ct Action: Next Screen// Select Enter, Do not install any ite	ms displayed	l on this scre	en

Exch	ange File Components Jun 01, 2021012:22:52	Page:	9 of	18
+	Component	Category	EXI	ists
77	VA-CRC CTC ORDER REASON: CONTRAINDICATED		>	<
78	VA-CRC ORDER CT COLONOGRAPHY		>	<
79	VA-CRC AVG RISK MT-SDNA POSITIVE		)	<
80	VA-CRC AVG RISK MT-SDNA NEGATIVE		)	<
81	VA-CRC AVG RISK FOBT/FIT ABNORMAL		)	<
82	VA-CRC AVG RISK FOBT/FIT NORMAL		)	ζ
83	VA-CRC ORDER SIGMOID		)	< C
84	VA-CRC AVG RISK ORDER FOBT/FIT		)	<
85	VA-CRC OUTSIDE CARE		)	<
86	VA-CRC SCHEDULED SCREENING		)	<
87	VA-CRC AVG RISK ORDER COLONOSCOPY		)	<
88	VA-CRC GAP DEFER 1M		)	<
89	VA-CRC GAP DEFER 3M		)	<
90	VA-CRC ADDL TESTS RECOMMENDED		)	¢
91	VA-CRC RECORD OTHER		)	< .
92	VA-CRC GAP ADDL TEST DECLINED		)	<
+	+ Next Screen - Prev Screen ?? More Actions			>>>
IA	Install all Components	+	- +	
Sele	ct Action: Next Screen// Select Enter Do not install any i	tems displaye	d on thi	s scree
	Sciect Enter. Do not install ally i	como anopiaje	o on m	5 50100

Exchai	ge File Components Jun 01, 2021012:23:16	Page:	10 01	18
+	Component	Category	E	xists
93	VA-CRC OUTSIDE CARE COMMENT			X
94	VA-CRC RECORD BE			×
95	VA-CRC GAP ORDER FOBT/FIT			×
96	VA-CRC FOLLOW UP SCHEDULED SCRN			x
97	VA-CRC FOLLOW UP OUTSIDE CARE			×
98	VA-CRC FOLLOW UP DEFER 4 M			×
99	VA-CRC FOLLOW UP WAITING			×
100	VA-CRC AVG RISK SET FUTURE 50 YEARS OLD			×
101	VA-CRC PROMPT COLONOSCOPY 1M			×
102	VA-CRC AVG RISK DEFER 1M			×
REMIN	DER SPONSOR			
103	VHA GI FIELD ADVISORY COMMITTEE			×
104	OFFICE OF NURSING SERVICE			х
105	VA NATIONAL CENTER FOR HEALTH PROMOTION AND			х
	DISEASE PREVENTION (NCP)			
+	+ Next Screen - Prev Screen ?? More Actions			>>>
IA	nstall all Components			
Selec	Action: Next Screen// Select Enter. Do not install any ite	ms displaye	d on th	is scree

Exchange File Components Jun 01, 2021012:23:40 Page	: 11 of 18
+ Component Catego	ory Exists
106 OFFICE OF PATIENT CARE SERVICES	×
107 NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE	x
PREVENTION (NCP)	
REMINDER COMPUTED FINDINGS	
VA - AGE	x
VA-REMINDER DEFINITION	x
REMINDER TAXONOMY	
108 VA-CRC DIAGNOSTIC COLONOSCOPY	×
109 VA-CRC CT COLONOGRAPHY	×
110 VA-CRC COLONOSCOPY AND SIGMOIDOSCOPY	×
111 VA-CRC ADDITIONAL TESTS	×
112 VA-CRC BARIUM ENEMA	x
113 VA-TERMINAL CANCER PATIENTS	x
+ + Next Screen - Prev Screen ?? More Actions	>>>
IA Install all Components Te Tratall Salasted Com	aanant
Select Action: Next Screen// Select Enter. Do not install any items disp	layed on this scree

Excha	nge File ComponentsJun	01, 2021@12:24	4:43	Page:	12 01	18
+	Component			Category	E	Exists
REMIN	DER TERM	1				
114	VA-REMINDER UPDATE_2_0_216					X
115	VA-AGE <50					X
116	VA-AGE <40!>65					X
117	VA-AGE <401>65					×.
110	VA-AGE <401>69					÷
120	VA-AGE <401>69					Ŷ
121	VA-AGE <401>70					
122	VA-AGE <40!>71		Install items 1	14 - 126 -	from t	his
123	VA-AGE <40!>72		moun nemo n		uom u	
124	VA-AGE <40!>73		screen. Ensur	e to select	the co	orrect
125	VA-AGE <40!>74		item numbers	for these	terms.	
126	VA-AGE <40					
127	VA-CRC VETERANS CHOICE SCRE	ËN				X
128	VA-CRC GAP TRIGGERS					Х
+	+ Next Screen - Prev :	Screen ?? Moi	re Actions	d Company		>>>
IA	Install all Components	IS Ins	stall Selecte	a compone	ant	
Selec	t Component(s): (1-208): 11	A-126	ted component			
20100	(1-208). (1-208). II	4-120				

At the REMINDER TERM section, select **IS**, then choose the item numbers 114-126 from the screenshot above.

## Overwrite the terms.

Excha	ange File Components Jun 01, 2021012:29:23	Page:	13 of	18
+	Component	Category	Exi	sts
129	VA-CRC BL TESTING WI PAST YEAR		>	(
130	VA-AGE >64		>	(
131	VA-CRC ADDITIONAL TESTS		>	(
132	VA-CRC GAP ORDER CTC		>	ζ
133	VA-CRC GAP ORDER DCBE		>	ζ
134	VA-CRC MALIGNANCY CONFIRMED BY PATH		>	ζ
135	VA-CRC BL MALIGNANCY		>	ζ
136	VA-CRC COLONOSCOPY COMPLETED		>	٢
137	VA-CRC BL ADDITIONAL TEST COMPLETE		>	٢
138	VA-CRC BL ADDITIONAL TEST ORDERED		>	ζ
139	VA-CRC TISSUE SENT TO PATHO		>	(
140	VA-CRC BL TISSUE SENT TO PATHOLOGY		>	ζ
141	VA-CRC BL GAP A		>	(
142	VA-CRC MT-SDNA LAB TEST POSITIVE		>	(
143	VA-CRC AVG RISK DUE NOW		>	ζ
144	VA-CRC FOLLOW-UP APPLICABLE/NOT RESOLVED		>	(
+	+ Next Screen - Prev Screen ?? More Actions			>>>
IA	Install all Components			
Sele	ct Action: Next Screen// Select Enter. Do not install any ite	ms displaye	d on this	screet

Excha	Exchange File Components Jun 01, 2021@12:30:27 Page: 14 of 18							
+	Component	Category	Exi	sts				
145	VA-CRC AVG RISK 2.5 YR		×	:				
146	VA-CRC AVG RISK 1.5 YR		×					
147	VA-CRC AVG RISK SET FUTURE USPTF AGE		×					
148	VA-CRC DIAGNOSTIC COLONOSCOPY PENDING		×					
149	VA-CRC CONFIRM OUT OF STUDY		×					
150	VA-CRC CONFIRM IN STUDY		×	5				
151	VA-CRC OUTSIDE CARE		×	(				
152	VA-CRC SCHEDULED SCREENING	Terrer 11 Sec. 156 - 4154	. d 177					
153	VA-CRC GAP DEFER 1M	Install item 150, which	is the VP	A-				
154	VA-CRC GAP ORDER COL/SIG/FIT/FOBT	CRC AVG RISK RESE	T from f	his				
155	VA-CRC BARTUM ENEMA	CITC ITTO INDIVIDUO						
156	VA-CRC AVG RISK RESET	screen. Ensure to select	the corre	ect				
157	VA-CRC POSITIVE MI-SONA HEALTH FACTORS	item number for this re-	minder					
158	VA-CRC POSITIVE FOBT/FIT HEALTH FACTORS	field humber for this fer	minuer					
159	VA-CRC DIAGNOSTIC COLONOSCOPY OUTCOMES	term.						
160	VA-CRC POSITIVE FOBT/FIT LAB TESTS							
TA	Testall all Components	stall Palastad Company	ant					
Roloo	t Action: Next Senson// IS In	tod Component	enc					
80100	t Component(s): (1 208): 156	ted component						
Selec	(1-208): 156							

At the REMINDER TERM section, select **IS**, then choose the item number 156 from the screenshot above.

# <mark>Merge</mark> the term.

REMINDER TERM entry named VA-CRC AVG RISK RESET already exists but the packed component is different, what do you want to do?					
Select one of the following:					
С	Create a new entry by copying to a new name				
М	Merge findings				
Ο	Overwrite the current entry				
U	Update				
Q	Quit the install				
S	Skip, do not install this entry				
Enter response: O//	Merge findings				

Exch	ange File Components Jun 01, 2021@12:44:07	Page:	15 of 18
+	Component	Category	Exists
161	VA-CRC GAP DECLINED		X
162	VA-CRC GAP DEFER 3M		x
163	VA-CRC DEFER OTHER		X
164	VA-CRC SCREENING NOT INDICATED		x
165	VA-CRC COLONOSCOPY FOLLOW UP		x
166	VA-CRC GAP REMINDER DUE		x
167	VA-CRC AVG RISK DEFER 3M		x
168	VA-CRC AVG RISK WAITING		X
169	VA-CRC AVG RISK SCRN DECLINED		x
170	VA-CRC AVG RISK 6 M		X
171	VA-CRC GAP REMINDER APPLICABLE/NOT RESOLVED		x
172	VA-LIFE EXPECTANCY < 6 MONTHS		x
173	VA-CRC AVG RISK 10 YR		X
174	VA-CRC AVG RISK 9 YR		x
175	VA-CRC AVG RISK 8 YR		X
176	VA-CRC AVG RISK 7 YR		x
+	+ Next Screen - Prev Screen ?? More Actions		>>>
IA	Install all Components IS Install Select	ed Compone	ant
Sele	ct Action: Next Screen// Select Enter. Do not install any iter	ms displaye	ed on this scre

Excha	inge File Components	Jun 01, 202	21012:44:31	Page:	16 of	18
+	Component			Category	Exi	sts
177	VA-CRC FOBT/FIT				>	
178	VA-CRC AVG RISK 5 YR				>	4
179	VA-CRC AVG RISK 4 YR				>	¢
180	VA-CRC AVG RISK 3 YR				>	(
181	VA-CRC AVG RISK 2 YR				>	(
182	VA-CRC AVG RISK 1 YR				>	(
183	VA-CRC AVG RISK ORDER	FOBT/FIT			>	t in the second s
184	VA-CRC AVG RISK 6 YR				>	(
185	VA-CRC AVG RISK ORDER	COL/SIG/CTC			>	
REMIN	DER DEFINITION					
186	VA-CRC BL TESTING WIT	HIN LAST YEAR			>	(
187	VA-CRC BL MALIGNANCY				>	(
188	VA-CRC BL ADDITIONAL	TEST COMPLETE			>	(
189	VA-CRC OBJ GAP TRIGGE	R			>	
190	VA-CRC BL ADDITIONAL	TEST ORDERED			>	
+	+ Next Screen -	Prev Screen	?? More Actions			>>>
IA	Install all Components		IS Install Select	ed Compone	ent	
Selec	t Action: Next Screen/	/ Select Enter.	Do not install any iter	ns displaye	d on this	scree



Select IS, then the item number for the VA-CRC AVERAGE RISK reminder.

## Overwrite the reminder definition.

Excha	nge File Components Jun 01, 2021@12	:46:44 Page:	18 of 18	
+	Component	Category	/ Exists	
199	VAL-CRC FOBT;3;2Y (TIU)		X	
200	VAL-CRC FIT;1;2Y (TIU)		X	
201	VA-CRC GAP TRIGGER DETAIL (TIU)		×	
TIU D	OCUMENT DEFINITION			
	TODAY'S DATE		x	
202	VA-CRC FOBT;3;3Y	Trantall items 205 from thi		
203	VA-CRC FIT;3;3Y	Install items 200 from the	s screen.	
204	VA-CRC GAP TRIGGER DETAIL	Ensure to select the corre	ct item	
REMIN	DER DIALOG	number for the VA-CRC	AVG	
205	VA-CRC AVG RISK UPDATE JUNE 2021	RISK LIDDATE ILINE 20	021	
200	VA-CRC COLONOSCOPY DIAGNOSTIC REMINDER	KISK OFDATE JUNE 20	521	
207	VA-CRC COLONOSCOPY AVG RISK	reminder dialog.		
208	VA-CRC WAITING FOR RECOMMEND GAP	0		
	Nevt Concer Drev Concer 00			
TA	Testall all Components	More Actions	>>>	2
Soloo	t Action: Quit// IS Install Soloctod C	omponent	ienc	
Selec	t Component(s): (1-208): 205	omponent		
Serec	(1-208). 205			

Select IS, then the item number for the VA-CRC AVG RISK UPDATE JUNE 2021 reminder dialog.

Dialo	g Components	Jun 01, 2021@12:49:55	Page: 1	of 1
Packe	d reminder dialog: V	A-CRC AVG RISK UPDATE JUNE 2021		
Item	<u>Seq. Dialog Findin</u>	gs	Туре	Exists
1	VA-CRC AVG RISK U	PDATE JUNE 2021	dialog	X
2	5 VA-CRC RECOMMEND	FUTURE RISK AT USPTF GUIDELINES	element	X
	Finding: VA-C	RC AVG RISK SET FUTURE USPTF AGE	(HEALTH	x
3	PXRM VISIT DATE FO	RCED TODAY	forced	х
4	10 VA-CRC RECOMMEN	D FUTURE COLO AT USPTF GUIDELINES	element	x
	Finding: VA-C	RC PROMPT COLONOSCOPY USPTF AGE (F	HEALTH FACTOR)	X
5	PXRM VISIT DATE FO	RCED TODAY	forced	X
6	15 VA-CRC RESULT L	ETTER REC CRC 50	group	X
	Finding: VA-C	RC AVG RISK SET FUTURE USPTF AGE (	(HEALTH	X
	FACTOR)			
	+ Next Screen	- Prev Screen ?? More Actions		
DD	Dialog Details	DT Dialog Text IS In	nstall Selecte	d
DF	Dialog Findings	DU Dialog Usage QU QU	uit	
DS	Dialog Summary	IA Install All		
Selec	t Action: Quit// IA	Install All		

The VA-CRC AVG RISK UPDATE JUNE 2021 dialog will display.

Select IA for Install All. At the Install reminder dialog and all components with no further changes, enter Yes.

Packed	d remi	nder dialog: V/	A-CRC	AVG RISK	UPDATE JUNE	2021			
VA-CRC	C AVG I	RISK UPDATE JU	NE 202	1 (remind	er dialog)	instal	lled fr	om exchang	ge file.
Item	Seq.	Dialog Finding	js					Туре	Exists
1	VA-0	CRC AVG RISK UP	PDATE	JUNE 2021				dialog	x
2	5 VA	- CRC RECOMMEND	FUTUR	RE RISK AT	USPTF GUI	DELINE	3	element	X
		Finding: VA-CF FACTOR)	RC AVG	RISK SET	FUTURE USP	PTF AGE	E (HEAL	тн	x
3	PXRM	VISIT DATE FOR	RCED T	ODAY				forced	x
4	10 V	A-CRC RECOMMEN	Ο Γυτι	IRE COLO A	T USPTF GU		s	element	x
		Finding: VA-CF	RC PRC	MPT COLON	OSCOPY USP	FF AGE	(HEALT	H FACTOR)	X
5	PXRM	VISIT DATE FOR	RCED T	ODAY				forced	×
6	15 V/	A-CRC RESULT LE	ETTER	REC CRC 5	0			group	×
		Finding: VA-CF	RC AVO	A RISK SET	FUTURE USP	PTF AGE	E (HEAL	TH	x
		TAOTON)							
	+	Next Screen	- Pre	ev Screen	?? More A	Actions	6		
DD D	Dialog	Details	DT	Dialog Te	xt	IS	Instal	1 Selected	d
DF D	Dialog	Findings	DU	Dialog Us	age	QU	Quit		
DS D	Dialog	Summary	IA	Install A	11				
Select	Select Action: Quit// QUIT								

After the dialog installs, you will then be returned to this screen. At the <u>Select Action</u> prompt, type  $\underline{O}$ .

+	Component	Category	Exists
тти			
110	TODAY'S DATE		×
~ ~ ~ ~	TUDAY S DATE		X
202	VA-CHC GAP TRIGGER DETAIL		X
203	VA-CRC FOBT;3;3Y		X
204	VA-CRC FIT;3;3Y		Х
DENT			
REMI	NDER DIALOG		
205	VA-CRC AVG RISK UPDATE JUNE 2021		Х
206	VA-CRC COLONOSCOPY AVG RISK		Х
207	VA-CRC WAITING FOR RECOMMEND GAP		Х
208	VA-CRC COLONOSCOPY DIAGNOSTIC REMINDER		Х
	+ Next Screen - Prev Screen ?? More Actions		>>>
TA	Install all Commonents IS Install Select	ed Commonent	
So 1 o	ct Action: Quit// Quit	ou componente	
OCIC	CC ACTON, WUIL// WUIL		

After the dialog installs, you will then be returned to this screen. At the <u>Select Action</u> prompt, type  $\underline{\mathbf{O}}$ .

Install Complete

## **Post Install**

- 1. Confirm the changes in the reminder dialogs and the VA-CRC AVERAGE RISK reminder definition as described in the installation description.
- 2. Confirm the appropriate text changed in the VA-CRC DOCUMENTATION dialog if local modifications were made to this dialog.