



Update _2_0_229

Clinical Reminders

MST SCREENING UPDATE

Install Guide

August 2021

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION	1
PRE-INSTALL	1
INSTALL DETAILS	2
INSTALL EXAMPLE	3
POST INSTALLATION	6

Introduction

This installation guide is for changes being made to the CPRS Military Sexual Trauma (MST) clinical reminder that is used to screen patients for experiences of MST. These changes are being made to align the clinical reminder with the definition of MST and MST eligibility requirements in recently passed legislation, to streamline the screening process, and to promote consistency across CPRS and the Cerner Electronic Health Record (EHR). Most of the revisions in this update are minor text changes to update language in the screening or changes to the logic of when and for whom the clinical reminder will be activated. The core aspects of the MST screening and referral process as included in the prior version of the MST clinical reminder largely remain the same.

The specific changes that have been made in this version of the MST Clinical Reminder are:

- The clinical reminder will be automatically set as due for certain other former Service members (FSMs), along with continuing to be set as due for Veterans.
- Patients with a current positive screen will no longer have the clinical reminder automatically set as due after a new period of military service.
- The text of the screen has been modified to include a revised screening question 1 (updating the definition of MST per new legislation), to be inclusive of both Veterans and FSMs, and to include additional information about the availability of care for MST-related physical health conditions.
- The text that is automatically pasted into the CPRS progress note stating the results of the screen has been simplified and shortened.

UPDATE_2_0_229 contains 1 Reminder Exchange entry:
UPDATE_2_0_229 VA-MST SCREENING UPDATE

The exchange file contains the following components:

TIU TEMPLATE FIELD

VA-MST FACTSHEET

HEALTH FACTORS

VA-REMINDER UPDATES [C]
VA-UPDATE_2_0_229
MST CATEGORY [C]
MST DECLINES TO ANSWER
VA-MST DECLINES MH REFERRAL
VA-MST CURRENTLY ENROLLED IN MH
VA-MST REQUESTS MH REFERRAL

MST YES REPORTS
MST NO DOES NOT REPORT

REMINDER SPONSOR

MENTAL HEALTH SERVICES
OFFICE OF QUALITY & PERFORMANCE
MENTAL HEALTH SERVICES AND WOMEN VETERANS HEALTH PROGRAM

REMINDER TERM

VA-REMINDER UPDATE_2_0_229
VA-MST OTH ELIGIBILITY
VA-IRAQ/AFGHAN PERIOD OF SERVICE
VA-MST DECLINES REPORT
VA-MST NEGATIVE REPORT
VA-MST POSITIVE REPORT

REMINDER DEFINITION

VA-MST SCREENING

REMINDER DIALOG

VA-MST SCREENING

Pre-Install

1. Find MST Referral order currently in use.

- Using the Vista Menu, Reminder Dialog Manager, DI- Reminder Dialog, change view to Group view and Inq/Print the reminder dialog group VA-MST REQUESTS MH REFERRAL. You will need to identify the local order used as an additional finding and add this order back post install.

Install Details

This update is being distributed as a web host file. The address for the host file is:
https://Redacted/UPDATE_2_0_229.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 30 minutes.

Install Example

To Load the exchange file, use LWH. The URL is
`https://Redacted/UPDATE_2_0_229.PRD`

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LWH  Load Web Host File
DFE  Delete Exchange File Entry       RI   Reminder Definition Inquiry
IFE  Install Exchange File Entry       RP   Repack
IH   Installation History
Select Action: Next Screen// LWH    Load Web Host File
Input the URL for the .prd file: https://[REDACTED]/UPDATE_2_0_229.PRD
```

If you use LWH and it works, you will see something like this message. Make sure the file was loaded successfully.

```
https://[REDACTED]/UPDATE_2_0_229.PRD successfully loaded.
```

Search and locate an entry titled UPDATE_2_0_229 VA-MST SCREENING UPDATE in reminder exchange.

Item	Entry	Source	Date Packed
206	UPDATE_2_0_229 VA-MST SCREENING UPDATE		08/19/2021@11:35
207	UPDATE_2_0_23 VA-NALOXONE USE		10/13/2017@12:09
208	UPDATE_2_0_24 MISC UPDATE 3		04/18/2017@14:58
209	UPDATE_2_0_27 VA-CHOICE DIALOGS		06/14/2017@11:33
210	UPDATE_2_0_28 VA-MENINGOCOCCAL IMMUNIZATIONS		06/26/2017@05:56
211	UPDATE_2_0_29 VA-HOMELESSNESS FOOD INSECURITY		07/18/2017@07:56

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry       RI   Reminder Definition Inquiry
IH   Installation History              RP   Repack
Select Action: Next Screen// IFE    Install Exchange File Entry
Enter a list or range of numbers (1-430): 206
```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry
Enter the number that corresponds with your entry titled UPDATE_2_0_229 VA-MST
SCREENING UPDATE (*in this example it is entry 206, it will vary by site*) that is dated
08/19/2021.

```
Component Category Exists
Source:
Date Packed: 08/19/2021@11:35:56
Package Version: 2.0P42

Description:
The following Clinical Reminder items were selected for packing:
REMINDER TERM
  VA-REMINDER UPDATE_2_0_229

REMINDER DEFINITION
  VA-MST SCREENING

Keywords:

Components:

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA   Install all Components          IS   Install Selected Component
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.
For all new entries in any category you will choose **I to Install**
Skip any Health Factors that you are prompted to Overwrite.
Merge any Reminder Terms that exist but are different
OVERWRITE the reminder definition

You will then be prompted to install the first dialog component. You will choose IA to install all even though some items will already have an X to indicate they exist.

```
Packed reminder dialog: VA-MST SCREENING
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-MST SCREENING	dialog	X
2	5	VA-GP MST MASTER Finding: *NONE*	group	
3	5.2	VA-MST VERSION NUMBER Finding: *NONE*	element	
4	5.5	VA-TEXT MST STATUS QUESTIONS Finding: *NONE*	group	X
5	5.5.5	VA-HF MST YES Finding: MST YES REPORTS (HEALTH FACTOR)	group	X
6	5.5.5.5	VA-MST FACTSHEET Finding: *NONE*	group	X
7	5.5.5.5.5	VA-MST REQUESTS MH REFERRAL Finding: VA-MST REQUESTS MH REFERRAL (HEALTH FACTOR)	group	X
8	5.5.5.5.10	VA-MST DECLINES MH REFERRAL	element	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-MST SCREENING**

```
Packed reminder dialog: VA-MST SCREENING
VA-MST SCREENING (reminder dialog) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-MST SCREENING	dialog	X
2	5	VA-GP MST MASTER Finding: *NONE*	group	X
3	5.2	VA-MST VERSION NUMBER Finding: *NONE*	element	X
4	5.5	VA-TEXT MST STATUS QUESTIONS Finding: *NONE*	group	X
5	5.5.5	VA-HF MST YES Finding: MST YES REPORTS (HEALTH FACTOR)	group	X
6	5.5.5.5	VA-MST FACTSHEET Finding: *NONE*	group	X
7	5.5.5.5.5	VA-MST REQUESTS MH REFERRAL Finding: VA-MST REQUESTS MH REFERRAL (HEALTH FACTOR)	group	X
8	5.5.5.5.10	VA-MST DECLINES MH REFERRAL	element	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Component	Category	Exists
Source: [REDACTED]		
Date Packed: 08/19/2021@11:35:56		
Package Version: 2.0P42		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER TERM		
VA-REMINDER UPDATE_2_0_229		
REMINDER DEFINITION		
VA-MST SCREENING		
Keywords:		
Components:		
+ + Next Screen - Prev Screen ?? More Actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// Q		

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.

Post Installation

1. **Add MST Referral identified in pre-install to reminder dialog.** Using the Vista Menu, Reminder Dialog Manager, DI- Reminder Dialog, change view to Group view edit the reminder dialog group VA-MST REQUESTS MH REFERRAL, add the quick order identified in the pre-install as an additional finding.
2. **Open the MST reminder** on a test patient and confirm the version number displayed at the top of the screen is 2.1.
3. **Inquire into the reminder definition VA-MST SCREENING** and confirm the cohort logic matches what is shown below:

Customized PATIENT COHORT LOGIC to see if the Reminder applies to patient:

(FI(4)!FI(9))!((FI(4)!FI(9))&FF(2))