

Update _2_0_235

Clinical Reminders

VA-DEMENTIA WARNING SIGNS

March 2022

Product Development Office of Information Technology Department of Veterans Affairs

Table of Contents	
INTRODUCTION	1
PRE-INSTALL	2
INSTALL DETAILS	2
INSTALL EXAMPLE	3
POST-INSTALLATION	8

Introduction

DESCRIPTION:

The VHA Dementia Warning Signs Template should be used at any site for documenting observed or reported dementia warning signs and documenting the initiation of, or referral for, a dementia workup.

The template is intended for use by VA clinical staff primarily in primary care settings, i.e., Patient-Aligned Care Teams (PACT), Community-Based Outpatient Clinics (CBOC), Home-Based Primary Care (HBPC), and Geriatric Patient Aligned Care Teams (Geri-PACT).

When primary care team members observe dementia warning signs or have dementia warning signs reported to them for one of their patients, they should open the national VHA Dementia Warning Signs Template to document the observed or reported dementia warning signs and to document their initiation of a dementia workup, including referral(s) to specialty care providers as needed.

This template is to be used ONLY when warning signs have been identified. Please note: The presence of dementia warning signs does not mean the patient has dementia. However, an evaluation is recommended.

Thank you to the following test sites who installed and tested this reminder update:

East Orange Cleveland San Diego Black Hills

UPDATE_2_0_235 contains 1 Reminder Exchange entry: UPDATE_2_0_235 VA-DEMENTIA WARNING SIGNS with the following components:

TIU TEMPLATE FIELD

- 1 VA-BLANK IN PN TEXT
- 2 VA-DWS FIREARM SAFETY CLINICIAN LINK
- 3 VA-DWS FIREARM SAFETY LINK
- 4 VA-WORD PROCESS
- 5 VA-BLANK SP
- 6 VA-CLINICIAN FACT SHEET DCI
- 7 VA-NCP PREVENTATIVE SERVICES GUIDANCE LINK
- 8 VA-DEMENTIA STEERING COMMITTEE LINK

HEALTH FACTORS

- 9 VA-REMINDER UPDATES [C]
- 10 VA-UPDATE_2_0_235
- 11 VA-DEMENTIA REPORTED BY [C]
- 12 VA-DEMENTIA REPORTED BY

- 13 VA-DEMENTIA OTHER CONCERNS [C]
- 14 VA-DWS GUN UTD
- 15 VA-DWS GUN NO
- 16 VA-DWS GUN
- 17 VA-DWS DRIVE UTD
- 18 VA-DWS DRIVE NO
- 19 VA-DWS DRIVE
- 20 VA-DEMENTIA WARNING SIGNS [C]
- 21 VA-DWS OTHER WARNING SIGNS
- 22 VA-DWS SELF CARE
- 23 VA-DWS CONFUSED
- 24 VA-DWS DIRECTIONS
- 25 VA-DWS LOST
- 26 VA-DWS ASK SAME QUESTIONS
- 27 VA-DWS FAILS TO FOLLOW
- 28 VA-DWS FORGETFUL
- 29 VA-DWS UNKEMPT
- **30 VA-DWS APPOINTMENT**
- 31 VA-DWS REPORT [C]
- 32 VA-DWS COMPLETED
- 33 VA-DWS FOLLOW-UP COMPLETED

REMINDER TERM

34 VA-REMINDER UPDATE 2_0_235

HEALTH SUMMARY COMPONENT

PCE HEALTH FACTORS SELECTED

HEALTH SUMMARY TYPE

35 VA-DWS REPORTED

HEALTH SUMMARY OBJECTS

36 VA-DWS REPORTED (TIU)

TIU DOCUMENT DEFINITION

37 VA-DWS REPORTED

REMINDER DIALOG

38 VA-DEMENTIA WARNING SIGNS

Pre-Install

None

Install Details

This update is being distributed as a web host file. The address for the host file is: https://vaww.va.gov/reminders/docs/UPDATE_2_0_235.PRD

The file will be installed using Reminder Exchange, programmer access is not required. This update can be loaded with users on the system. Installation will take less than 5 minutes

Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

+	+ Next Screen - Prev	Screen ?	? More Actions	>>>
CFE	Create Exchange File Entry	LHF	Load Host File	
CHF	Create Host File	LMM	Load MailMan Messag	e
CMM	Create MailMan Message	LR	List Reminder Defin	itions
DFE	Delete Exchange File Entry	LWH	Load Web Host File	
IFE	Install Exchange File Entry	RI	Reminder Definition	Inquiry
IH	Installation History	RP	Repack	
Selec	t Action: Next Screen// LWH	Load Web	Host File	
Input	the URL for the .prd file:	https://		/UPDATE_2_0_235.PRD

At the <u>Select Action</u>: prompt, enter <u>LWH</u> for Load Web Host File At the <u>Input the url for the .prd file</u>: prompt, type the following web address:

https://Redacted/UPDATE_2_0_235.PRD

You should see a message that the file was successfully loaded.

https://

/UPDATE_2_0_235.PRD successfully loaded.

This will install one exchange entries in your exchange file

Added Reminder Exchange entry UPDATE_2_0_235 VA-DEMENTIA WARNING SIGNS

+Item	Entry	Source	Date Packed			
293	UPDATE_2_0_235 VA-DEMENTIA		03/03/2022@13:01			
	WARNING SIGNS					
294	UPDATE_2_0_235 VA-DEMENTIA		03/08/2022@06:20			
	WARNING SIGNS					
295	UPDATE_2_0_235 VA-DEMENTIA		03/23/2022@13:16			
	WARNING SIGNS					
296	UPDATE_2_0_236 VA-PROSTHETICS		11/04/2021@09:23			
297	UPDATE_2_0_237_VA-COVID-19		10/26/2021@04:22			
	TELEPHONE FOLLOW-UP UPDATE					
+	+ Next Screen - Prev Screen	n ?? More Actions	>>>			
CFE	Create Exchange File Entry	LHF Load Host File				
CHF	Create Host File	LMM Load MailMan Messag	je			
CMM	Create MailMan Message	LR List Reminder Defir	nitions			
DFE	Delete Exchange File Entry	LWH Load Web Host File				
IFE	Install Exchange File Entry	RI Reminder Definition	n Inquiry			
IH	Installation History	RP Repack				
Selec	Select Action: Next Screen// IFE Install Exchange File Entry					
Enter	a list or range of numbers (1-574	4): 295				

Search and locate an entry titled UPDATE_2_0_235 VA-DEMENTIA WARNING SIGNS

At the <u>Select Action</u> prompt, enter **IFE** for Install Exchange File Entry Enter the number that corresponds with your entry **UPDATE_2_0_235 VA-DEMENTIA WARNING SIGNS.** *(in this example it is entry 295 it will vary by site).* The date of the exchange file should be 03/23/2022.

Compone	ent			Category	Exists
Source:					
Date Packed:	03/23/2022@13:16:17				
Package Versi	on: 2.0P71				
Description:					
The following	g Clinical Reminder items	s were	selected 1	for packing:	
REMINDER DIAL	_0G				
VA-DEMENTI#	WARNING SIGNS				
REMINDER TERM	1				
VA-REMINDER	R UPDATE_2_0_235				
Keywords:					
ReyMondan					
Components:					
+ Ent	ter ?? for more actions				>>>
IA Install	all Components	IS	Install	Selected Component	
Select Action	n: Next Screen// IA				
At the Select	Action prompt, type IA for	Install	all Comp	onents and hit enter.	

Any items that are NEW in any categories, you should **INSTALL** There should not be any items that are pre-existing.

You will then be prompted to install the reminder dialog component, you will choose IA to install even though most items will have an X to indicate they already exist.

Packe	d reminder dialog: <mark>VA-DEMENTIA WARNING SIGNS [NATIONAL</mark>	DIALOG]	
		_	
Item	Seq. Dialog Findings	Туре	Exists
1	VA-DEMENTIA WARNING SIGNS	dialog	
2	1 VA-DEMENTIA WS Finding: *NONE*	element	
3	3 VA-GP DEMENTIA WS NOTE Finding: *NONE*	group	
4	4 VA-GP DEMENTIA WS SCREEN Finding: VA-DWS COMPLETED (HEALTH FACTOR)	group	x
5	4.5 VA-GP DEMENTIA WS REPORT	group	
6	4.5.4 VA-DEMENTIA WS QUESTIONS	element	
	Finding: VA-DWS ASK SAME QUESTIONS (HEALTH FACT	0R)	X
- 7	PXRM COMMENT	prompt	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DD	Dialog Details DT Dialog Text IS In	stall Selecte	d
DF	Dialog Findings DU Dialog Usage QU Qu	it	
DS	Dialog Summary IA Install All		
Selec	t Action: Next Screen// IA		

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – <u>VA-DEMENTIA WARNING</u> <u>SIGNS</u>

Pac	ke	ed reminder dialog: VA-	DEMENTIA	WARNING	SIGNS [NA	TIONAL	DIALOG]		
VA-	DE	MENTIA WARNING SIGNS (reminder	dialog)	installed	from	exchange	file.	
It	en	<u>ı Seq. Dialog Findings</u>					Тур	e	Exists
	1	VA-DEMENTIA WARNING	SIGNS				dial	og	Х
	2	1 VA-DEMENTIA WS Finding: *NONE*					eleme	nt	x
	3	3 VA-GP DEMENTIA WS Finding: *NONE*	NOTE				gro	up	x
	4	4 VA-GP DEMENTIA WS Finding: VA-DWS	SCREEN COMPLETI	ED (HEAL)	TH FACTOR)		gro	up	x x
	5	4.5 VA-GP DEMENTIA W Finding: *NONE*	s report	,	,		gro	up	х
	6	4.5.4 VA-DEMENTIA WS	QUESTIO				eleme	nt	X
	7	FINDING: VA-DWS	ASK SAM	E QUESTIC	MS (HEALT	H FACI	UK)	nt	×
	í	+ Next Screen -	Prev Sci	reen ?	More Act	ione	hi oiii	μι	^
nn		Dialog Details D	T Dial	na Text		S In	stall Sel	ected	
DF		Dialog Findings D	U Dialo	og Usade	0	U Qu	iit		
DS		Dialog Summary I	A Insta	all All					
Sel	.ec	t Action: Next Screen/	/ Q						

You will then be returned to this screen. At the Select Action prompt, type Q.

Component	Category	Exists
Source:		
Date Packed: 03/23/2022@13:16:17		
Package Version: 2.0P71		
Description: The following Clinical Reminder items were se REMINDER DIALOG VA-DEMENTIA WARNING SIGNS	lected for packing:	
REMINDER TERM VA-REMINDER UPDATE_2_0_235		
Keywords:		
Components:		
+ + Next Screen - Prev Screen ??	More Actions	>>>
IA Install all Components IS Select Action: Next Screen// Q	Install Selected Componer	nt

You will then be returned to this screen. At the Select Action prompt, type Q.

Install complete.

Post-Installation

 Make the Dialog template available to be attached to a Personal or Shared template or to be able to attach the dialog to a progress note title. Select Reminder Managers Menu <TEST ACCOUNT> Option: CP CPRS Reminder Configuration

CA Add/Edit Dominder C

- CA Add/Edit Reminder Categories
- CL CPRS Lookup Categories
- CS CPRS Cover Sheet Reminder List
- MH Mental Health Dialogs Active
- PN Progress Note Headers

RA Reminder GUI Resolution Active

TIU TIU Template Reminder Dialog Parameter

DEVL Evaluate Coversheet List on Dialog Finish

DL Default Outside Location

PT Position Reminder Text at Cursor

LINK Link Reminder Dialog to Template

TEST CPRS Coversheet Time Test

NP New Reminder Parameters

GEC GEC Status Check Active

WH WH Print Now Active

Select CPRS Reminder Configuration <TEST ACCOUNT> Option: **TIU** TIU Template Reminder Dialog Parameter

Reminder Dialogs allowed as Templates may be set for the following:

- 1 User USR [choose from NEW PERSON]
- 3 Service SRV [choose from SERVICE/SECTION]
- 4 Division DIV [SALT LAKE CITY]
- 5 System SYS [NATREM.FO-SLC.MED.VA.GOV]

Enter selection: **5** Choose the appropriate number for your site. Your site may do this by System or other levels. The example below uses SYSTEM level.

Setting Reminder Dialogs allowed as Templates for System: NATREM.FO-SLC.MED.VA.GOV

Select Display Sequence: ?

388 VA-COVID-19 TELEPHONE FOLLOW-UP

389 VA-SRA-CONSULT RESPONSE (D)

When you type a question mark above, you will see the list of #'s (with dialogs) that are already taken. Choose a number **NOT** on this list. For this example, 390 is not present so I will us<mark>e 390.</mark>

Select Display Sequence: 390

Are you adding 390 as a new Display Sequence? Yes// **YES** Display Sequence: 390// **<Enter>** 390

Clinical Reminder Dialog: VA-DEMENTIA WARNING SIGNS reminder dialog NATIONAL

...OK? Yes// <mark><Enter></mark> (Yes)

2. Add local consult orders if there are any. This group and all elements are local so if a site has more services that would be appropriate, more can be added. Any of the services listed above can be modified to match what is available at the local site. If you change the template display text make sure to also change the alternate text to match. If your site has no orders to add or decides not to add orders this group should be disabled.

26 4.25 Group: VAL-GP DEMENTIA WS REFERRAL Text: Optional Referrals:[BOX, SUPPRESS, SHOW, NO SELECTION]

27 4.25.5 Element: VAL-DEMENTIA WS REFERRAL GERIATRIC Text: A referral was made for a geriatric evaluation.

28 4.25.10 Element: VAL-DEMENTIA WS REFERRAL SOCIAL WORK Text: A referral was made for social work.

29 4.25.15 Element: VAL-DEMENTIA WS REFERRAL DEMENTIA CLINIC Text: A referral was made for the dementia clinic.

- 30 4.25.20 Element: VAL-DEMENTIA WS REFERRAL PCMHI
- 78 65.20.10 Group: VA-GP DWS PLAN REFERRED Text: [BOX, SUPPRESS, SHOW, NO SELECTION]
- 65.20.10.10 Element: VAL-DWS REFERRED INSTRUCTIONAL TEXT Text: Consider referral to (the following options will enter a consult to the corresponding service)
- 80 65.20.10.20 Element: VAL-DWS REFERRAL TO SOCIAL WORK Text: Social Work consult entered.
- 81 65.20.10.30 Element: VAL-DWS GERIATRIC CONSULT Text: Geriatric Clinic consult entered.
- 82 65.20.10.50 Element: VAL-DWS GERIATRIC PSYCH Text: Geriatric Psychiatric consult entered.
- 83 65.20.10.60 Element: VAL-DWS NEUROLOGY CONSULT Text: Neurology consult entered.

- 84 65.20.10.65 Element: VAL-DWS PCMHI Text: Referred to PCMHI.
- 85 65.20.10.70 Element: VAL-DWS DEMENTIA CLINIC Text: Dementia clinic consult entered.

3. An order set or individual quick orders can be added as finding items to this group. Make sure to change the text so the user will know what they are getting when they click this option. Local elements can also be added to this group if sites decide to have orders placed individually rather than all with one click. If your site decides not to add orders this group should be disabled.

 88 65.20.15 Group: VAL-GP DWS PLAN TESTING Text: Recommended routine labs/imaging for dementia
 *** ENTER LOCAL LAB ORDERS HERE ACCORDING TO LOCAL PREFERENCE SUGGESTED ROUTINE ORDERS INCLUDE: TSH, CBC, ELECTROLYTES, CA+, HEPATIC PANEL, DRUG AND ALCOHOL SCREEN, BUN, CREAT, GLUCOSE, B12, HIV

Other diagnostic testing:

- 4. An order set or individual quick orders can be added as finding items to this group. Make sure to change the text so the user will know what they are getting when they click this option. Local elements can also be added to this group if sites decide to have orders placed individually rather than all with one click. If your site has no orders to add or decides not to add orders this group should be disabled
 - 89 65.20.16 Group: VAL-GP DWS PLAN TESTING OTHER Text: Other diagnostic testing:
 - *** ENTER LOCAL ORDERS HERE ACCORDING TO LOCAL PREFERENCE SUGGESTED ORDERS INCLUDE: Neuropsych eval, brain imaging, consult to dementia specialist

Local elements can be added to this group to include any quick orders for the recommended testing for the items listed above. Work with your local specialists to determine what orders to include.

5. Elements can be added to this group to allow for finding items for some of these less common tests if your site has them available.

65.20.17 Group: VAL-GP DWS PLAN TESTING ADVANCED
Text:
Consider targeted testing for rare causes of cognitive
impairment if any of the following conditions are suspe
-Syphilis -Lyme Disease -Heavy Metal/Toxin exposure
-Endocrine disorders -Prion-related disorders such as
Creutzfeldt-Jacob Disease (CJD) -Meningeal carcinomato
-Central nervous system vasculitis
[NO BOX, SUPPRESS, SHOW, NO SELECTION]