



**Update \_2\_0\_264**

**Clinical Reminders**

**VA-SMART DIALOG AND MAMMOGRAM REMINDER  
UPDATE**

**Install Guide**

**June 2022**

Product Development  
Office of Information Technology  
Department of Veterans Affairs



**Table of Contents**

**INTRODUCTION ..... 1**

**INSTALL DETAILS ..... 11**

**INSTALL EXAMPLE..... 12**

**POST INSTALLATION..... 16**



# Introduction

Several revisions have been made to the VA-WH SMART BREAST IMAGING FOLLOW-UP reminder dialog and the VA-WH MAMMOGRAM SCREENING reminder definition.

1. The VA-WH NEXT BREAST PROCEDURE reminder term CONDITION has been updated. A Beginning Date/Time of 10.1.2020 has also been added to prevent the term from inappropriately evaluating procedures completed before the SMART functionality was deployed in CPRS 31b.
2. The following reminder terms have been moved from the reminder resolution and added to the reminder cohort logic as exclusions. These findings were not temporarily resolving the reminder correctly in conjunction with the Custom Due Date finding. Now, during the deferral and order period, the reminder will be N/A. It will come due again after the temporary resolution time has passed.
  - a. VA-MAMMOGRAM SCREEN DEFER
  - b. VA-WH MAMMOGRAM ORDER
  - c. VA-WH BREAST CARE ORDER HEALTH FACTOR
3. New finding FOUND TEXT has been added to display in the clinical maintenance to assist local staff understand temporary logic.
  - a. If the VA-BREAST CARE ORDER HEALTH FACTOR term is true, the following will display in clinical maintenance:

```
Cohort:
Reminder Term: VA-WH BREAST CARE ORDER HEALTH FACTOR
Health Factor: Women's Health Order Mammogram Screen HF
02/23/2022@12:49

* * * * *
*
* _
* _   A mammogram has been ordered in the past 90 days.
* _   This reminder will display as NOT APPLICABLE, even though
* _   it is applicable to this patient. If the patient
* _   does not have a mammogram within the 90 days, then the reminder
* _   will become DUE again.
* _
* * * * *
```

- b. If the VA-WH MAMMOGRAM SCREEN NOT INDICATED term is true, the following will display in clinical maintenance:

```
Cohort:
Reminder Term: VA-WH MAMMOGRAM SCREEN NOT INDICATED
Health Factor: Women's Health Mammogram Screen Not Indicated
02/23/2022@12:49
Comments: TEST

* * * * *
*
* _
* _      Mammogram screening is not indicated according to documentation.
* _
*
* * * * *
```

c. If the VA-WH MAMMOGRAM SCREEN DEFER term is true, the following will display in clinical maintenance:

```
Cohort:
Reminder Term: VA-WH MAMMOGRAM SCREEN DEFER
Health Factor: Women's Health Mammogram Declined
02/23/2022@12:49
Comments:

* * * * *
*
* _
* _      Mammogram Screening reminder has been deferred.
* _      During the 4 months after the deferral date, this reminder
* _      will display as NOT APPLICABLE, even though it is applicable
* _      to this patient. The reminder will become DUE again 4 months
* _      after the deferral was documented.
*
* * * * *
```

d. If the VA-WH MAMMOGRAM ORDER term is true, the following will display in clinical maintenance:

```
Cohort:
Reminder Term: VA-WH MAMMOGRAM ORDER
Orderable Item: MAMMOGRAM
03/21/2022@10:41 Status: active, Start date: 03/21/2022@10:41, Stop date:
04/20/2022 Duration: 0 D
*****
*
*- A mammogram has been ordered in the past 90 days.
*- While the patient has an active order, this reminder
*- will display as NOT APPLICABLE, even though it is applicable
*- to this patient. If the order is discontinued or the patient
*- does not have a mammogram within the 90 days, then the reminder
*- will become DUE again.
*
*****
```

- 4. Two incorrect URLs have been removed from the Mammogram Screening reference information.
- 5. Corrected an error that caused the reminder to be not applicable on patients <40 yo who had a mammogram and a new screening date documented in the VA-WH SMART BREAST IMAGING FOLLOW-UP note.
- 6. Branching logic has been added to the VA-WH MAMMOGRAM SCREENING reminder dialog to hide the “Change the frequency of mammograms for this patient” option if the VA-WH NEXT BREAST PROCEDURE reminder term is true. This indicates a future screening recommendation has been documented using the SMART dialog since CPRS 31b implementation. Once the mammogram reminder due date has been set using SMART, then the original frequency change options in the mammogram screening dialog do not impact the reminder logic.
- 7. The option to document “Patient breast images are normal” has been revised to “Patient breast images are benign/normal”.

Before:

```
Please summarize result:
 Patient breast images are normal
 Patient breast images are abnormal
 Patient breast images are incomplete
```

After:

Please summarize result:

- Patient breast images are benign/normal
- Patient breast images are abnormal
- Patient breast images are incomplete

UPDATE\_2\_0\_264 contains 1 Reminder Exchange entry:

### **VA-SMART DIALOG AND MAMMOGRAM REMINDER UPDATE**

8. Corrected an error occurring when a user selected 'no action required' under the next steps section.
9. Corrected a reminder resolution error occurring when a user selected "Other" in the Diagnostics options section.
10. A new general finding is included that will ensure the date a procedure report was received in SMART will be stored.
11. The "Written Report Available" section has been revised to ensure required fields are easier to complete.
12. A new section has been added to document breast density per the imaging report.
13. Options for documenting outside Breast Tomosynthesis procedures have been added.

Before:

Record prior or outside mammogram:

Written Report Available

NOTE: Obtain outside written report and forward for scanning into Vista Imaging and entry of BIRADS codes by Radiology

Outside report documentation:

Outside report

Received on: \*  /  / 2022

Location:

Select procedure to document (document next steps in next section):

Document follow-up actions:

Document how the patient was notified of the result

- Phone call
- In person
- Non-Registered Letter
- Registered Letter
- Secure Messaging
- Patient Not Notified: information forwarded to PACT/PC provider for patient notification
- Patient not contacted, chose date for tickler followup:

Additional Information:



Written Report Available

NOTE: Obtain outside written report and forward for scanning into Vista Imaging and entry of BIRADS codes by Radiology

Outside report documentation:

Outside report

Received on: \*  ...

Location:

Select procedure to document (document next steps in next section):

Mammogram Screening

Diagnostic Bilateral Mammogram:

Diagnostic Unilateral Mammogram:  
Side: \*

MRI

Ultrasound:  
Side:\*

Screening Breast Tomosynthesis

Diagnostic Bilateral Breast Tomosynthesis

Diagnostic Unilateral Breast Tomosynthesis:  
Side: \*

Biopsy:

Document breast density per the imaging report\*:

Dense (heterogeneously dense or extremely dense)

Not dense (predominately fatty or scattered fibroglandular tissue)

Unknown/not documented

After:

Document how the patient was notified of the result

- Phone call
- In person
- Non-Registered Letter
- Registered Letter
- Secure Messaging
- Information forwarded to PACT/PC provider for patient notification
- Patient not contacted, chose date for tickler followup:
- Patient already notified/not needed

Additional Information:

**The exchange file contains the following components:**

**TEMPLATE FIELDS**

- SMART WORD 40X4
- VA-WH SMART OUTSIDE OTHER
- VA-WH SMART ULTRASOUND TYPE
- VA-WH SMART UNILAT TYPE
- BLANK TEXT FIELD FOR DIALOGS
- VA-WH MAMMOGRAM GUIDANCE

**GENERAL FINDINGS**

TICKLER REMINDER COMMENT  
TICKLER REMINDER  
FOLLOW UP WITH PATIENT COMMENT  
FOLLOW UP WITH PATIENT  
WH PROCEDURE ACTION COMMENT  
RADIOLOGY REPORT RECEIVED  
NOTIFICATION NOT NEEDED  
EMAIL CONTACTED  
SECURE MESSAGING  
LETTER REGISTERED CONTACTED  
LETTER (CERTIFIED)  
LETTER CONTACTED  
LETTER  
WHO CONTACTED  
IN PERSON  
PHONE CALL  
BREAST CARE CASCADE NAME  
CLOSE CASCADE BREAST CARE  
BR CLOSE AND FOLLOW COMPLETE  
BR CLOSED PROCEDURE  
BREAST REFER FOR BIOPSY  
BIOPSY  
WV NOTIFICATION PATIENT  
OBTAIN PRIOR FILMS  
BREAST OBTAIN PRIOR FILMS  
BREAST BIOPSY ALREADY OBTAIN  
BIOPSY ALREADY OBTAINED  
BREAST CURRENTLY UNDER TREATMENT  
CURRENTLY UNDER TREATMENT  
CONSULT  
BREAST CONSULT  
BREAST REFER TO ONCOLOGIST  
REFER TO ONCOLOGIST  
BREAST REFER TO SURGEON  
REFER TO SURGEON  
IMAGE COMMENT  
BREAST NEED MRI  
RETURN FOR MRI  
BREAST NEED ULTRASOUND  
RETURN FOR ULTRASOUND  
BREAST NEED MAMMOGRAM  
RETURN FOR MAMMOGRAM  
BR NO NEXT TREATMENT DATE  
BR NO NEXT TREATMENT  
WH PATIENT  
BREAST next MAM 2Y  
AGE BASED SCREENING  
BREAST next MAM 1Y  
BREAST next MAM AT START AGE  
MAMMOGRAM SCREENING 2Y  
MAMMOGRAM SCREENING 1Y

BREAST UNSATISFACTORY RESULT  
BR UNSATISFACTORY RESULT  
BREAST ABNORMAL RESULT  
BR ABNORMAL RESULT  
BREAST NORMAL RESULT  
BR NEM RESULT  
BI-RAD 6 RESULT  
BI-RAD 6  
BI-RAD 5 RESULT  
BI-RAD 5  
BI-RAD 4 RESULT  
BI-RAD 4  
BI-RAD 3 RESULT  
BI-RAD 3  
BI-RAD 2 RESULT  
BI-RAD 2  
BI-RAD 1 RESULT  
BI-RAD 1  
BI-RAD 0 RESULT  
BI-RAD 0  
BREAST UNSATISFACTORY DENSE RESULT  
BREAST ABNORMAL DENSE RESULT  
BREAST NORMAL DENSE RESULT  
BI-RAD 6 DENSE RESULT  
BI-RAD 5 DENSE RESULT  
BI-RAD 4 DENSE RESULT  
BI-RAD 3 DENSE RESULT  
BI-RAD 2 DENSE RESULT  
BI-RAD 1 DENSE RESULT  
BI-RAD 0 DENSE RESULT  
BREAST TOMOSYNTHESIS UNILAT REPORT  
BR OUTSIDE REPORT PROCEDURE DATE  
BR OUTSIDE ACCESSION  
BREAST TOMOSYNTHESIS BILAT REPORT  
BREAST TOMOSYNTHESIS SCREENING REPORT  
OUTSIDE ULTRASOUND  
OUTSIDE MRI  
OUTSIDE MAMMOGRAM DX UNILAT  
OUTSIDE MAMMOGRAM DX BILAT  
MAMMOGRAM SCREENING REPORT  
OUTSIDE REPORT  
VIEW PROGRESS NOTE TEXT  
ALERT ID  
CLEAR ONE  
WV DIAGNOSIS  
WV PROCEDURE ACCESS  
MAMMOGRAM TEST  
VIEW LAST THREE BR TRMTS  
VIEW EPISODE  
SMART BI-RAD ALERT  
CLEAR OTHER ALERTS ID

CLEAR ALL FOR PROVIDER  
BR NEXT DUE AT AGE  
BR MAMMOGRAM SCREENING NEXT TRMT

## HEALTH FACTORS

VA-REMINDER UPDATES [C]  
VA-UPDATE\_2\_0\_264  
WH MAMMOGRAM [C]  
WH ORDER MAMMOGRAM SCREEN HF  
WH OUTSIDE INCOMPLETE MAMMOGRAM  
WH OUTSIDE ABNL MAMMOGRAM  
WH OUTSIDE NORMAL MAMMOGRAM  
VA-WH BREAST RESULT BIRAD 6  
VA-WH BREAST RESULT BIRAD 5  
VA-WH BREAST RESULT BIRAD 4  
VA-WH BREAST RESULT BIRAD 3  
VA-WH BREAST RESULT BIRAD 2  
VA-WH BREAST RESULT BIRAD 1  
VA-WH BREAST RESULT BIRAD 0  
VA-WH BREAST OTHER  
VA-WH MAMMOGRAM DX UNILATERAL  
VA-WH MAMMOGRAM DX BILATERAL  
VA-WH MAMMOGRAM SCREENING  
WH SMART OUTSIDE REPORT INFO  
WH MAMMOGRAM SCREEN NOT INDICATED  
WH BR CA SCREEN N/A 5 YRS-COMORBIDITIES  
WH BR CA SCREEN N/A 5 YRS-LE<5YRS  
VA LIFE EXPECTANCY [C]  
VA LIMITED LIFE EXPECTANCY  
WH BILATERAL MASTECTOMY  
WH MAMMOGRAM SCREEN FREQ - 1Y  
WH MAMMOGRAM SCREEN FREQ - 6M  
WH MAMMOGRAM DEFERRED  
VA-WH BR CA 40-44 WANTS SCREENING  
VA-WH BR CA 40-44 BEGIN AGE 45  
WH MAMMOGRAM SCREEN FREQ - 2Y  
WH MAMMOGRAM DECLINED  
WH BR CA 40-49 WANTS SCREENING  
WH MAMMOGRAM SCREEN FREQ - 4M  
VA-WH BR CA DEFER 1Y  
VA-WH BR CA DEFER 6M  
REMINDER FACTORS [C]  
INACTIVATE BREAST CANCER SCREEN  
WH ORDER MAMMOGRAM UNILAT HF  
WH ORDER MAMMOGRAM BILAT HF  
WH MAMMOGRAM OUTSIDE

## **REMINDER SPONSOR**

OFFICE OF PUBLIC HEALTH AND ENVIRONMENTAL HAZARDS AND THE WOMEN  
VETERANS HEALTH STRATEGIC HEALTH CARE GROUP  
VHA NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION (NCP)  
AND VHA OFFICE OF WOMEN'S HEALTH  
NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION

## **REMINDER TAXONOMY**

VA-MASTECTOMY  
VA-BREAST TUMOR  
VA-TERMINAL CANCER PATIENTS  
VA-MAMMOGRAM/SCREEN  
VA-WH BILATERAL MASTECTOMY

## **REMINDER TERM**

VA-REMINDER UPDATE\_2\_0\_264  
VA-WH BL ACTIVE BREAST IMAGE ALERT DATA  
VA-WH OPEN PROCEDURES PER CASCADE  
VA-WH BL BR ACTIVE ALERT ID  
VA-PATIENT IS MALE  
VA-WH MAMMOGRAM ALERT FOR MALE  
VA-WH BL BR BI-RAD 6  
VA-WH BL BR BI-RAD 5  
VA-WH BL BR BI-RAD 4  
VA-WH BL BR BI-RAD 3  
VA-WH BL BR BI-RAD 0  
VA-WH BL BR BI-RAD 2  
VA-WH BL BR BI-RAD 1  
VA-AGE >74  
VA-BL AGE>74  
VA-WH AGE BETWEEN 55 AND 74  
VA-WH AGE BETWEEN 45 AND 54  
VA-WH MAMMOGRAM AGE<55  
VA-WH OPEN BR PROCEDURES WITH BIRAD CODES  
VA-WH OPEN BR PROCEDURES WITHOUT BIRAD CODES AND DENSE  
VA-WH OPEN BR PROCEDURES WITH BIRAD CODES DENSE  
VA-WH BL GET BREAST IMAGES  
VA-WH BL BREAST CASCADE EXIST  
VA-WH NEXT BREAST PROCEDURE  
VA-AGE >39&<45  
VA-WH CHAMP VA PATIENT  
VA-WH BR CA 40-44 WANTS SCREEN TERM  
VA-WH BR CA 40-44 BEGIN AGE 45  
VA-WH MAMMOGRAM ORDER  
VA-WH HX BREAST CANCER/ABNORMAL MAM

VA-WH MAMMOGRAM SCREEN FREQ - 2Y  
VA-WH MAMMOGRAM SCREEN FREQ - 1Y  
VA-WH MAMMOGRAM SCREEN FREQ - 6M  
VA-WH MAMMOGRAM SCREEN FREQ - 4M  
VA-WH MAMMOGRAM UNSATISFACTORY IN RAD/WH PKG  
VA-WH MAMMOGRAM SCREEN DEFER  
VA-WH MAMMOGRAM SCREEN NOT INDICATED  
VA-WH BREAST CARE ORDER HEALTH FACTOR  
VA-WH MAMMOGRAM SCREEN DONE  
VA-WH MAMMOGRAM SCREEN IN RAD PKG  
VA-WH MAMMOGRAM SCREEN IN WH PKG  
VA-TERMINAL CANCER PATIENT  
VA-WH BILATERAL MASTECTOMY

### **REMINDER DEFINITIONS**

VA-WH BL BR ALERT WITH NO OPEN PROCEDURE  
VA-WH OPEN BREAST PROCEDURES PER CASCADE  
VA-WH MAMMOGRAM SCREENING

### **REMINDER DIALOG**

VA-WH SMART BREAST IMAGING FOLLOW-UP  
VA-WH MAMMOGRAM SCREENING

## **Pre-Installation Steps**

---

1. Record the orders mapped to the VA-WH SMART BREAST IMAGING FOLLOW-UP reminder dialog. The site will need to re-map these orders to the reminder dialog during the install.
  - a. Identify the order dialog mapped to the group VA-WH BR ADDITIONAL IMAGING. This order dialog will be used to replace the Q.BREAST TREATMENT MENU during the installation.
  - b. Identify the order dialog mapped to the element VA-WH BR CONSULT. This order dialog will be used to replace the Q.GMRCOR CONSULT during the installation.
2. Record the order dialog mapped (if one is mapped) to the VA-WH GP ORDERS - MAMMOGRAM SCREEN reminder group.
3. Create a copy of the reminder inquiry output for the VA-WH MAMMOGRAM SCREENING reminder using either CPRS or Vista. Use this copy to confirm all local mapping for reminder terms are correct after post-installation.

## **Install Details**

---

This update is being distributed as a web host file. The address for the host file is:

**[https://REDACTED/UPDATE\\_2\\_0\\_264.PRD](https://REDACTED/UPDATE_2_0_264.PRD)**

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 5 minutes.

## Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
CFE Create Exchange File Entry      LHF Load Host File
CHF Create Host File                LMM Load MailMan Message
CMM Create MailMan Message          LR List Reminder Definitions
DFE Delete Exchange File Entry      LWH Load Web Host File
IFE Install Exchange File Entry      RI Reminder Definition Inquiry
IH Installation History              RP Repack
Select Action: Next Screen// lwh    Load Web Host File
Input the URL for the .prd file: https:// UPDATE_2_0_264.PRD
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File  
At the **Input the url for the .prd file:** prompt, type the following web address:  
**[https://vaww.va.gov/reminders/docs/UPDATE\\_2\\_0\\_264.PRD](https://vaww.va.gov/reminders/docs/UPDATE_2_0_264.PRD)**

You should see a message at the top of your screen that the file successfully loaded.

Search and locate an entry titled **UPDATE\_2\_0\_264 VA-SMART DIALOG AND MAMMOGRAM REMINDER UPDATE** in reminder exchange.

+Item	Entry	Source	Date Packed
280	UPDATE_2_0_264 VA-SMART DIALOG AND MAMMOGRAM REMINDER UPDATE		03/30/2022@11:54
281	UPDATE_2_0_265 VA-COVID-19 TEMPLATE UPDATES		03/03/2022@07:57
282	UPDATE_2_0_267 GEC GERIATRIC EVALUATION		03/28/2022@05:59
283	UPDATE_2_0_27 VA-CHOICE DIALOGS		06/14/2017@11:33
284	UPDATE_2_0_272 VA-ATHLETES		02/22/2022@08:37

  

+Item	Entry	Source	?? More Actions	>>>
CFE	Create Exchange File Entry	LHF	Load Host File	
CHF	Create Host File	LMM	Load MailMan Message	
CMM	Create MailMan Message	LR	List Reminder Definitions	
DFE	Delete Exchange File Entry	LWH	Load Web Host File	
IFE	Install Exchange File Entry	RI	Reminder Definition Inquiry	
IH	Installation History	RP	Repack	

Select Action: Next Screen// IFE Install Exchange File Entry  
Enter a list or range of numbers (1-488): 280

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry **UPDATE\_2\_0\_264 VA-SMART DIALOG AND MAMMOGRAM REMINDER UPDATE** (in this example it is entry 280. It will vary by site). The date of the exchange file should be 03/30/2022.



```
Date Packed: 03/30/2022@11:54:28
Package Version: 2.0P71

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-WH MAMMOGRAM SCREENING
  VA-WH SMART BREAST IMAGING FOLLOW-UP

REMINDER TERM
  VA-REMINDER UPDATE_2_0_264

REMINDER DEFINITION
  VA-WH MAMMOGRAM SCREENING
Non-exchangeable order dialog(s):
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components          IS  Install Selected Component
Select Action: Next Screen// IA  Install all Components .....
```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.  
Select Action: Next Screen// **IA Install all Components**

- For all new entries, choose **Install**
- For REMINDER DIALOG LINK TYPE entries that exist, but the packed component is different, **OVERWRITE**
- For TERMS and TAXONOMIES that exist, but the packed component is different **SKIP**, except the following term. **Overwrite** this term if asked:
- **VA-WH NEXT BREAST PROCEDURE**
- For REMINDER DEFINITIONS that exist, but the packed component is different **OVERWRITE**.

You will be prompted to install the reminder dialog component (**VA-WH SMART BREAST IMAGING FOLLOW-UP**). You will choose to install all even though some components may have an X to indicate that they exist.

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-WH SMART BREAST IMAGING FOLLOW-UP	dialog	X		
2	10	VA-WH GP BR OPEN PROCEDURES* Finding: *NONE*	group	X		
3	10.10	VA-WH BR NO OPEN PROCEDURES* Finding: *NONE*	element	X		
4	10.15	VA-WH BR EPISODE OF CARE Finding: VIEW EPISODE (REMINDER GENERAL FINDING)	element	X		
5		PXRM GF VIEW BUTTON	prompt	X		
6	10.20	VA-WH BR LAST 3 TRMTS Finding: VIEW LAST THREE BR TRMTS (REMINDER GENERAL FINDING)	element	X		
7		PXRM GF VIEW BUTTON	prompt	X		
8	10.25	VA-WH BLANK LINE Finding: *NONE*	element	X		
9	10.30	VA-WH GP BR BI-RAD SELECTOR*	group	X		
<b>+</b> <b>+</b> Next Screen <b>-</b> Prev Screen <b>??</b> More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA    Install All						

At the **Select Action** prompt, type **IA** to install the dialog: **VA-WH SMART BREAST IMAGING FOLLOW-UP**.

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: Y// **Yes**

FINDING entry Q.BREAST TREATMENT MENU does not exist.

Enter a code from the list.

Select one of the following:

- D    Delete
- P    Replace with an existing entry
- Q    Quit the install

Enter response: **P**

Select ORDER DIALOG NAME: Enter local breast treatment menu identified in Pre-Installation steps.

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-WH SMART BREAST IMAGING FOLLOW-UP	dialog	X		
2	10	VA-WH GP BR OPEN PROCEDURES*	group	X		
		Finding: *NONE*				
3	10.10	VA-WH BR NO OPEN PROCEDURES*	element	X		
		Finding: *NONE*				
4	10.15	VA-WH BR EPISODE OF CARE	element	X		
		Finding: VIEW EPISODE (REMINDER GENERAL FINDING)		X		
5		PXRM GF VIEW BUTTON	prompt	X		
6	10.20	VA-WH BR LAST 3 TRMTS	element	X		
		Finding: VIEW LAST THREE BR TRMTS (REMINDER GENERAL FINDING)		X		
7		PXRM GF VIEW BUTTON	prompt	X		
8	10.25	VA-WH BLANK LINE	element	X		
		Finding: *NONE*				
9	10.30	VA-WH GP BR BI-RAD SELECTOR*	group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// Q Quit						

After completing this dialog install, you will type **Q**.

You will be prompted to install the reminder dialog component (**VA-WH MAMMOGRAM SCREENING**) You will choose to install all even though some components may have an X to indicate that they exist.

Packed reminder dialog: VA-WH MAMMOGRAM SCREENING

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-WH MAMMOGRAM SCREENING	dialog	X		
2	5	VA-WH MAMMOGRAM SCREENING HEADER*	group	X		
		Finding: *NONE*				
3	5.5	VA-WH GP MAM SCREEN HEADER*	group	X		
		Finding: *NONE*				
4	5.5.5	VA-WH MAM HEADER1 - SCREEN	element	X		
		Finding: *NONE*				
5	5.5.10	VA-WH GP MAM SCREEN HEADER-EXPANDED	group	X		
		Finding: *NONE*				
6	5.5.10.5	VA-WH MAM HEADER2 - SCREEN	element	X		
		Finding: *NONE*				
7	5.5.15	VA-WH DOTTED LINE	element	X		
		Finding: *NONE*				
8	5.5.20	VA-WH GP MAMMOGRAM SCREENING	group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// IA						

At the **Select Action** prompt, type **IA** to install the dialog: **VA-WH MAMMOGRAM SCREENING**).

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: Y// **Yes**



FINDING ITEM: [Add order dialog here.](#)

- 2. Using the reminder inquiry output from the pre-installation step, compare all reminder terms in the VA-WH MAMMOGRAM SCREENING reminder definition to confirm all local mapping is correct.**
- 3. Confirm all local orders are mapped correctly in the VA-WH SMART BREAST IMAGING FOLLOW-UP reminder dialog.**