



Update _2_0_299

Clinical Reminders

VA-POST ANESTHESIA SEDATION SCORE

Install Guide

March 2023

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION 1
INSTALL DETAILS 2
INSTALL EXAMPLE 2
POST INSTALLATION 10

Introduction

This reminder update will install a new reminder dialog for documentation of Post Anesthesia Sedation Scores using newly released MH instruments named PASS-1 and PASS-2. These instruments will auto score the answers to questions for Post Anesthesia rating.

As a part of this update, we were notified that EPRP was failing sites due to the incorrect verbiage being placed in the progress note when a PTSD-5 instrument was completed. This was updated and included in reminder dialog VA-MH PC-PTSD+I9 SCREEN. This extra piece in this update has nothing to do with Post Anesthesia.

UPDATE_2_0_299 contains 1 Reminder Exchange entry:

UPDATE_2_0_299 VA-AES POST ANESTHEISA/SEDATION SCORE

The exchange file contains the following components:

MH TESTS AND SURVEYS

- PASS-2
- PASS-1
- PC-PTSD-5

REMINDER DIALOG LINK TYPE

- 1 VA-PTSD POSITIVE WARM HANDOFF CHECK

TIU TEMPLATE FIELD

- 2 ONS STANDARDIZED NURSING DOCUMENTATION
- 3 VHA DIRECTIVE 1073

HEALTH FACTORS

- 4 VA-REMINDER UPDATES [C]
- 5 VA-UPDATE_2_0_299
- 6 VA-VAAES VAPAS [C]
- 7 VA-VAAES VAPAS PH2 TIME
- 8 VA-VAAES VAPAS PH1 TIME
- 9 VA-VAAES VAPAS VERSION 1.4

REMINDER SPONSOR

- 10 OFFICE OF NURSING INFORMATICS
- 11 OFFICE OF NURSING SERVICE
- 12 OFFICE OF MH AND SUICIDE PREVENTION
- 13 OFFICE OF MENTAL HEALTH SERVICES

REMINDER TERM

- 14 VA-REMINDER UPDATE_2_0_299

REMINDER DIALOG

- 15 VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS -
VERSION 1.0

Install Details

This update is being distributed as a web host file. The address for the host file is:
https://REDACTED /UPDATE_2_0_299.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 5 minutes.

Install Example

Step 1 – (must be completed first) Load the Mental Health Instruments

You will first load the instruments using the INSTRUMENT EXCHANGE Menu (YTQ INSTRUMENT EXCHANGE) in Vista. If you do not have access to this menu, you will need to go through your facility process to request a new Vista menu. The two instruments you will load are PASS-1 and PASS-2

THIS IS NOT REMINDER EXCHANGE.

ADDED NOTE NOT PART OF THIS ORIGINAL DOCUMENT:

Select OPTION NAME: MHA3 UTILITIES YTQ MHA3 MENU MHA3 Utilities

*** Mental Health ***
MHA3 Utilities

- 1 Print Test Form
- 2 Detailed Definition
- 3 Delete Patient Data
- 4 Stop/Re-Start Progress Notes for an Instrument
- 5 Exempt Test
- 6 Test Usage
- 7 XML Output
- 8 Instrument Exchange
- 10 Display Administration Details

Theses are PASS-1.JSON and PASS-2.JSON

Use this url for install REDACTED but they may be found here outside of the VA:

```

+          Enter ?? for more actions
CE (Create New Entry)   CH (Create Host File)   DR Trial Install <Dry Run>
RB (Rebuild Entry)     LH Load Host File       IE Install Exchange Entry
DL Delete Entry         LU Load from URL        BR Browse Specification
Select Action: Next Screen// lh Load Host File
Enter file name: _____/PASS-1.JSON

```

This will load: PASS1 for Reminder
 created on: 01/04/2023@13:34:30
 source: XXXXX@REDACTED

Description ---

(no description)

Do you want to continue? Yes//

PASS1 for Reminder loaded.

Press return to continue...

Once the PASS-1 is loaded, you will need to install

Entry	Source	Created
1	CCICM20 FOR REMINDER	
2	CR-PASS-1-B8	
3	CR-PASS-2-B4	
4	VFO20-READY-FOR-REMINDER	
5	YS*5.01*121 NEW INSTRUMENTS	
6	YS*5.01*121 UPDATES	
7	YS*5.01*123 INSTRUMENT UPDATE	
8	YS*5.01*134 T10C4	
9	YS*5.01*136	
10	YS*5.01*137	

```

+          Enter ?? for more actions
CE (Create New Entry)   CH (Create Host File)   DR Trial Install <Dry Run>
RB (Rebuild Entry)     LH Load Host File       IE Install Exchange Entry
DL Delete Entry         LU Load from URL        BR Browse Specification
Select Action: Next Screen// IE Install Exchange Entry
Enter a number (1-36): 2

```

Find entry named "CR-PASS-1-B8" and then select IE. In this instance, you would select number 2 to install. When prompted to Use verbose mode, you will choose YES

You will then repeat this process for PASS-2 instrument install

Use this url for install /srv/vista/patches/SOFTWARE/PASS-2.JSON

```

+          Enter ?? for more actions
CE (Create New Entry)   CH (Create Host File)   DR Trial Install <Dry Run>
RB (Rebuild Entry)     LH Load Host File       IE Install Exchange Entry
DL Delete Entry         LU Load from URL        BR Browse Specification
Select Action: Next Screen// LH Load Host File
Enter file name: Use this url for install _____/PASS-2.JSON

```


This will load: CR-PASS-2-B4
 created on: 11/17/2022@23:12:33
 source: XXXXX@REDACTED

Description ---

(no description)

Do you want to continue? Yes//
PASS1 for Reminder loaded.
Press return to continue...

Once the PASS-2 is loaded, you will need to install



Entry	Source	Created
1	CCICM20 FOR REMINDER	
2	CR-PASS-1-B8	
3	CR-PASS-2-B4	
4	VFQ20-READY-FOR-REMINDER	
5	YS*5.01*121 NEW INSTRUMENTS	
6	YS*5.01*121 UPDATES	
7	YS*5.01*123 INSTRUMENT UPDATE	
8	YS*5.01*134 T10C4	
9	YS*5.01*136	
10	YS*5.01*137	

Enter ?? for more actions

CE (Create New Entry)	CH (Create Host File)	DR Trial Install <Dry Run>
RB (Rebuild Entry)	LH Load Host File	IE Install Exchange Entry
DL Delete Entry	LU Load from URL	BR Browse Specification

Select Action: Next Screen// IE Install Exchange Entry
Enter a number (1-36): 2

Find entry named “CR-PASS-2-B4” and then select IE. In this instance, you would select number 3 to install. When prompted to Use verbose mode, you will choose YES

This completes the Mental Health instrument exchange install piece. You will now move to reminder exchange to install the remaining parts.

Go to Reminder exchange in Vista.
To Load the exchange file, use LWH. The URL is
https://REDACTED/UPDATE_2_0_299.PRD

```

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE Create Exchange File Entry      LHF Load Host File
CHF Create Host File                 LMM Load MailMan Message
CMM Create MailMan Message           LR  List Reminder Definitions
DFE Delete Exchange File Entry       LWH Load Web Host File
IFE Install Exchange File Entry       RI  Reminder Definition Inquiry
IH  Installation History              RP  Repack
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: https://[REDACTED]/UPDATE_2_0_299.prd
You will see a message indicating the file was loaded successfully.

```

Search and locate an entry titled UPDATE_2_0_299 VA-AES POST ANESTHESIA/SEDATION SCORE in reminder exchange.

+Item	Entry	Source	Date Packed
412	UPDATE_2_0_299 VA-AES POST ANESTHESIA/SEDATION SCORE		
413	UPDATE_2_0_30 VA-IRAQ & AFGHAN POST-DEPLOY		
414	UPDATE_2_0_300 VA-TH VVC UPDATE		
415	UPDATE_2_0_301 VA-VCL UPDATE		
416	UPDATE_2_0_302 VA-SEXUAL HEALTH		
417	UPDATE_2_0_302 VA-SEXUAL		

```

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE Create Exchange File Entry      LHF Load Host File
CHF Create Host File                 LMM Load MailMan Message
CMM Create MailMan Message           LR  List Reminder Definitions
DFE Delete Exchange File Entry       LWH Load Web Host File
IFE Install Exchange File Entry       RI  Reminder Definition Inquiry
IH  Installation History              RP  Repack
Select Action: Next Screen// IFE  Install Exchange File Entry
Enter a list or range of numbers (1-675): 412

```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry. Enter the number that corresponds with your entry titled UPDATE_2_0_299 VA-AES POST ANESTHESIA/SEDATION SCORE (in this example it is entry 412, it will vary by site) that is dated 01/26/2023.

```

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS - VERSION 1.0

REMINDER TERM
  VA-REMINDER UPDATE_2_0_299

Keywords:

Components:

+      Enter ?? for more actions      >>>
IA  Install all Components             IS  Install Selected Component
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.
Install all new components.
Overwrite any components that exist but are different

You will then be prompted to install the first dialog component.

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS - VERSION 1.0	dialog	X		
2	20	VA-GP VAAES POST ANESTHESIA/SEDATION SCORE-VAPAS - VER	group	X		
		Finding: VA-VAAES VAPAS VERSION 1.3 (HEALTH FACTOR)		X		
3	20.5	VA-BLANK LINE ALT PN TEXT ONLY	element	X		
		Finding: *NONE*				
4	25	VA-GP PASS PHASE 1	group	X		
		Finding: *NONE*				
5		PXRM PASS1 RESULT GROUP	rs.group	X		
6		PXRM PASS1 RESULT ELEMENT 1	rs.element	X		
7	25.3	VA-PASS PHASE 1	element	X		
		Finding: PASS-1 (MH TESTS AND SURVEY)		X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// IA Install All						

At the **Select Action** prompt, type **IA** to install the dialog – **VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS – VERSION 1.0**

Install reminder dialog and all components with no further changes: Y//ES

VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS - VERSION 1.0 (reminder dialog) ins						
Item	Seq.	Dialog Findings	Type	Exists		
1		VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS - VERSION 1.0	dialog	X		
2	20	VA-GP VAAES POST ANESTHESIA/SEDATION SCORE-VAPAS - VER Finding: VA-VAAES VAPAS VERSION 1.3 (HEALTH FACTOR)	group	X		
3	20.5	VA-BLANK LINE ALT PN TEXT ONLY Finding: *NONE*	element	X		
4	25	VA-GP PASS PHASE 1 Finding: *NONE*	group	X		
5		PXRM PASS1 RESULT GROUP	rs.group	X		
6		PXRM PASS1 RESULT ELEMENT 1	rs.element	X		
7	25.3	VA-PASS PHASE 1 Finding: PASS-1 (MH TESTS AND SURVEY)	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD		Dialog Details	DT	Dialog Text	IS	Install Selected
DF		Dialog Findings	DU	Dialog Usage	QU	Quit
DS		Dialog Summary	IA	Install All		
Select Action: Next Screen// Q						

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will then be prompted to install the second dialog component

Packed reminder dialog: VA-MH PC-PTSD+I9 SCREEN [NATIONAL DIALOG]

Item	Seq.	Dialog Findings	Type	Exists
1		VA-MH PC-PTSD+I9 SCREEN	element	X
		Finding: PC-PTSD-5 (MH TESTS AND SURVEY)		X
2		PXRM PC-PTSD-5 RESULT GROUP	rs.group	X
3		PXRM PC-PTSD-5 RESULT ELEMENT 1	rs.element	X
4		PXRM PC-PTSD-5 RESULT ELEMENT 2	rs.element	X

+ Next Screen		- Prev Screen		?? More Actions	
DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		
Select Action: Quit// IA Install All					

At the **Select Action** prompt, type **IA** to install the dialog – **VA-MH PC-PTSD+I9 SCREEN**

Install reminder dialog and all components with no further changes: Y//ES

Packed reminder dialog: VA-MH PC-PTSD+I9 SCREEN [NATIONAL DIALOG]

Item	Seq. Dialog Findings	Type	Exists
1	VA-MH PC-PTSD+I9 SCREEN Finding: PC-PTSD-5 (MH TESTS AND SURVEY)	element	X X
2	PXRM PC-PTSD-5 RESULT GROUP	rs.group	X
3	PXRM PC-PTSD-5 RESULT ELEMENT 1	rs.element	X
4	PXRM PC-PTSD-5 RESULT ELEMENT 2	rs.element	X

```
+ Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details  DT Dialog Text      IS Install Selected
DF Dialog Findings DU Dialog Usage     QU Quit
DS Dialog Summary  IA Install All
Select Action: Quit// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```
Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS - VERSION 1.0

REMINDER TERM
  VA-REMINDER UPDATE_2_0_299

Keywords:

Components:

+ + Next Screen  - Prev Screen  ?? More Actions  >>>
IA Install all Components      IS Install Selected Component
Select Action: Next Screen// Q  Q
```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.

Post Installation

1. **Make the Dialog templates available to be attached to a Personal or Shared templates or to be able to attach the dialog to a progress note title.**

Select Reminder Managers Menu <TEST ACCOUNT> Option: **CP CPRS**

Reminder Configuration

CA Add/Edit Reminder Categories
CL CPRS Lookup Categories
CS CPRS Cover Sheet Reminder List
MH Mental Health Dialogs Active
PN Progress Note Headers
RA Reminder GUI Resolution Active
TIU TIU Template Reminder Dialog Parameter
DEVL Evaluate Coversheet List on Dialog Finish
DL Default Outside Location
PT Position Reminder Text at Cursor
LINK Link Reminder Dialog to Template
TEST CPRS Coversheet Time Test
NP New Reminder Parameters
GEC GEC Status Check Active
WH WH Print Now Active

Select CPRS Reminder Configuration <TEST ACCOUNT> Option: **TIU TIU**

Template Reminder Dialog Parameter

Reminder Dialogs allowed as Templates may be set for the following:

- 1 User USR [choose from NEW PERSON]
- 3 Service SRV [choose from SERVICE/SECTION]
- 4 Division DIV [SALT LAKE CITY]
- 5 System SYS [NATREM.FO-SLC.MED.VA.GOV]

Enter selection: **5** *Choose the appropriate number for your site. Your site may do this by System or other levels. The example below uses SYSTEM level.*

Setting Reminder Dialogs allowed as Templates for System: NATREM.FO-SLC.MED.VA.GOV

Select Display Sequence: **?**

388 VA-COVID-19 TELEPHONE FOLLOW-UP
389 VA-SRA-CONSULT RESPONSE (D)

*When you type a question mark above, you will see the list of #'s (with dialogs) that are already taken. Choose a number **NOT** on this list. For this example, 390 is not present so I will use 390.*


Select Display Sequence: **390**

Are you adding 390 as a new Display Sequence? Yes// **YES**

Display Sequence: 390// **<Enter>** 390

Clinical Reminder Dialog: **VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS – VERSION 1.0** reminder dialog NATIONAL
...OK? Yes// **<Enter>** (Yes)

2. **Add the template to** your shared template folder so it can be used with any note title.
 - a In CPRS, open the Template Editor, expand Shared Templates, select the appropriate oncology folder (if applicable), then click New Template (top right-hand corner).



Type in the following information:

Name: **<Give your template a reasonable name>**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA-AES POST ANESTHESIA/SEDATION SCORE-VAPAS – VERSION 1.0**

— then click Apply.