



**Update \_2\_0\_367**

**Clinical Reminders**

**VA-ONC LUNG CANCER FOLLOW UP AND CONSULT**

**Install Guide**

**August 2023**

Product Development  
Office of Information Technology  
Department of Veterans Affairs

Table of Contents

<b>INTRODUCTION .....</b>	<b>1</b>
<b>PRE-INSTALL.....</b>	<b>1</b>
<b>INSTALL DETAILS .....</b>	<b>20</b>
<b>INSTALL EXAMPLE .....</b>	<b>20</b>
<b>POST INSTALLATION .....</b>	<b>27</b>

## Introduction

---

The purpose of these lung cancer follow-up and consult templates is to standardize follow-up care, regulate data capture, and to acquire the required documentation as determined by the National Oncology Program (NOP). Every VA with an Oncology program should be utilizing these notes to ensure standardized data capture for NOP dashboards.

In addition, the lung cancer clinical pathways are included the follow-up notes/consult reports to facilitate ease of use. Oncology Clinical Pathways (CPs) are decision-support tools that translate published evidence into standardized treatment process maps to help guide cancer treatment decisions. Oncology CPs are different from treatment guidelines, as CPs recommend a standard treatment based on efficacy, toxicity, and cost. VA Oncology CPs are developed by an interdisciplinary subject matter expert (SME) group and approved by national program offices.

The lung cancer follow-up reminder dialogs were built with the Cerner APSO (assessment, plan, subjective, objective) format. Furthermore, there are only 2 required mandatory fields. The 5 reminder dialog templates will be utilized by all VA medical centers with an Oncology service. Additionally, the 5 reminder dialog templates will also be used as consult reports and utilized by all Oncology providers:

- Nurse Practitioners
- Physician Assistants
- Residents
- Fellows
- Physicians

VHA Directive 1415 requires that all sites use the National Oncology templates. If your site does not provide oncology services, the site must request a waiver to exempt it from the requirement to use the National Oncology templates. The waiver can be requested by emailing REDACTED

## Pre-Install

---

Using the VistA menu Reminder Dialog Manager, Reminder Dialog, change your view to elements and do an inquire/print on the following reminder dialog elements, you will need to use this as a reference to use when adding the following orderable items during post-install. You may not have anything mapped to these elements. This was an optional step in the previous lung pathways install.

- VA-ONC RADIATION ONCOLOGY CONSULT
- VA-ONC LUNG REFER TO PALLIATIVE CARE
- VAL-ONC LUNG IVB POS 1ST OTHER
- VA-ONC HOSPICE REFERRAL
- VA-ONC LUNG SURGICAL CONSULT
- VA-EL ONC LUNG ON PDL1

- VA-EL ONC LUNG PTHWY NGS ON
- VA-EL ONC LUNG PTHWY LIQUID BIOPSY

UPDATE\_2\_0\_367 contains 1 Reminder Exchange entry:

**UPDATE\_2\_0\_367 VA-ONC LUNG CANCER FOLLOW UP AND CONSULT REPORT**

The exchange file contains the following components:

**REMINDER GENERAL FINDINGS**

VIEW PROGRESS NOTE TEXT

**TEMPLATE FIELD**

VA-ONC CN PREVIOUS NOTE VIEW  
 VA-ONC PRESUMPTIVE CONDITIONS  
 VA-ONC PATHWAYS SHAREPOINT  
 VA-ONC NPOP SHAREPOINT  
 VA-ONC CLINICAL TRIAL  
 VA-ONC OFF REASON TXT  
 VA-ONC TREATMENT REASON  
 VA-ONC WP 2 LINES  
 EDIT 55+  
 TELEHEALTH EMERGENCIES VA LOCATIONS  
 TELEHEALTH EMERGENCIES NON-VA LOCATIONS  
 VA-ONC VAMC SITE NAME  
 TEXT (1-50 CHAR)  
 VA-ONC CN EMOTIONAL WB  
 DATE  
 VA-ONC CN LUNG CHIEF COMPLAINT

**HEALTH FACTORS**

VA-REMINDER UPDATES [C]  
 VA-UPDATE\_2\_0\_367  
 VA-ONC CN [C]  
 VA-ONC CN LUNG  
 VA-ONC LUNG SCLC LIMITED STAGE 1ST LINE PATHWAY [C]  
 VA-ONC LUNG SCLC LIMITED SUB PS3-4 OFF PATH RSN  
 VA-ONC LUNG SCLC LIMITED SUB PS3-4 OFF PATH TX  
 VA-ONC LUNG SCLC LIMITED SUB PS3-4 MRI SURVEILLANCE  
 VA-ONC LUNG SCLC LIMITED SUB PS3-4 CLINICAL TRIAL  
 VA-ONC LUNG SCLC LIMITED SUB NO PROG OFF PATH RSN  
 VA-ONC LUNG SCLC LIMITED SUB NO PROG OFF PATH TX  
 VA-ONC LUNG SCLC LIMITED SUB NO PROG REFER RAD ONC  
 VA-ONC LUNG SCLC LIMITED SUB NO PROG CLINICAL TRIAL  
 VA-ONC LUNG SCLC LIMITED INITIAL OFF PATH RSN  
 VA-ONC LUNG SCLC LIMITED INITIAL OFF PATH TX  
 VA-ONC LUNG SCLC LIMITED INITIAL CIS+ETO+CHEST RAD  
 VA-ONC LUNG SCLC LIMITED INITIAL CLINICAL TRIAL

VA-ONC LUNG SCLC LIMITED INITIAL CONTRA OFF PATH RSN  
VA-ONC LUNG SCLC LIMITED INITIAL CONTRA OFF PATH TX  
VA-ONC LUNG SCLC LIMITED INITIAL CONTRA CARBO+ETO+CHEST RAD  
VA-ONC LUNG SCLC LIMITED INITIAL CONTRA CLINICAL TRIAL  
VA-ONC LUNG SQUAMOUS RELAPSE PATHWAY [C]  
VA-ONC LUNG SQR 2ND IMMUNO AMEN RAD OFF PATH REASON  
VA-ONC LUNG SQR 3RD LINE NOT AMEN RAD OFF PATH REASON  
VA-ONC LUNG SQR 3RD LINE NOT AMEN RAD OFF PATH TX  
VA-ONC LUNG SQR 3RD LINE NOT AMEN RAD REFER PALL CARE  
VA-ONC LUNG SQR 3RD LINE NOT AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG SQR 3RD LINE AMEN RAD OFF PATH REASON  
VA-ONC LUNG SQR 3RD LINE AMEN RAD OFF PATH TX  
VA-ONC LUNG SQR 3RD LINE AMEN RAD REFER RAD ONC  
VA-ONC LUNG SQR 3RD LINE AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG SQR 2ND IMMUNO NOT AMEN RAD OFF PATH REASON  
VA-ONC LUNG SQR 2ND IMMUNO NOT AMEN RAD OFF PATH TX  
VA-ONC LUNG SQR 2ND IMMUNO NOT AMEN RAD CARBO+PACLITAXEL  
VA-ONC LUNG SQR 2ND IMMUNO NOT AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG SQR 2ND IMMUNO AMEN RAD OFF PATH TX  
VA-ONC LUNG SQR 2ND IMMUNO AMEN RAD REFER RAD ONC  
VA-ONC LUNG SQR 2ND IMMUNO AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO NOT AMEN RAD OFF PATH RSN  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO NOT AMEN RAD OFF PATH TX  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO NOT AMEN RAD DOCETAXEL  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO NOT AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO AMEN RAD OFF PATH REASON  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO AMEN RAD OFF PATH TX  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO AMEN RAD REFER RAD ONC  
VA-ONC LUNG SQR 2ND CHEMO+/-IMMUNO AMEN RAD CLINICAL TRIAL  
ONCOLOGY (ONC) LUNG CANCER TREATMENTS [C]  
ONC LUNG PALLIATIVE REFERRAL  
VA-ONC LUNG SCLC RELAPSE PATHWAY [C]  
VA-ONC LUNG SCLC RELAPSE >6M OFF PATH RSN  
VA-ONC LUNG SCLC RELAPSE >6M OFF PATH TX  
VA-ONC LUNG SCLC RELAPSE >6M CARBO+ETO  
VA-ONC LUNG SCLC RELAPSE >6M CLINICAL TRIAL  
VA-ONC LUNG SCLC RELAPSE PS 3-4 OFF PATH RSN  
VA-ONC LUNG SCLC RELAPSE PS 3-4 OFF PATH TX  
VA-ONC LUNG SCLC RELAPSE PS 3-4 REFER HOSPICE  
VA-ONC LUNG SCLC RELAPSE PS 3-4 CLINICAL TRIAL  
VA-ONC LUNG SCLC RELAPSE PS 0-1-2 OFF PATHWAY RSN  
VA-ONC LUNG SCLC RELAPSE PS 0-1-2 OFF PATHWAY TX  
VA-ONC LUNG SCLC RELAPSE PS 0-1-2 LURBINECTEDIN  
VA-ONC LUNG SCLC RELAPSE PS 0-1-2 CLINICAL TRIAL  
VA-ONC LUNG SCLC RELAPSE >6M CARBO+ETO+ATEZ  
VA-ONC LUNG SCLC EXTENSIVE STAGE CLINICAL PATHWAY [C]  
VA-ONC LUNG SCLC EXT SUB NO PROG OFF PATH RSN  
VA-ONC LUNG SCLC EXT SUB NO PROG OFF PATH TX  
VA-ONC LUNG SCLC EXT SUB NO PROG MAIN ATEZO  
VA-ONC LUNG SCLC EXT SUB NO PROG CLINICAL TRIAL  
VA-ONC LUNG SCLC EXT INI IMMUNO OFF PATH RSN

VA-ONC LUNG SCLC EXT INI IMMUNO OFF PATH TX  
VA-ONC LUNG SCLC EXT INI IMMUNO REFER HOSPICE  
VA-ONC LUNG SCLC EXT INI IMMUNO CLINICAL TRIAL  
VA-ONC LUNG SCLC EXT INI CHEMO OFF PATH RSN  
VA-ONC LUNG SCLC EXT INI CHEMO OFF PATH TX  
VA-ONC LUNG SCLC EXT INI CHEMO CARBO+ETO  
VA-ONC LUNG SCLC EXT INI CHEMO CLINICAL TRIAL  
VA-ONC LUNG SCLC EXT INI CHEMO/IMMUNO OFF PATH RSN  
VA-ONC LUNG SCLC EXT INI CHEMO/IMMUNO OFF PATH TX  
VA-ONC LUNG SCLC EXT INI CHEMO/IMMUNO CARBO+ETO+ATEZO  
VA-ONC LUNG SCLC EXT INI CHEMO/IMMUNO CLINICAL TRIAL  
VA-ONC LUNG SCLC EXT SYMPTOM METS OFF PATH RSN  
VA-ONC LUNG SCLC EXT SYMPTOM METS OFF PATH TX  
VA-ONC LUNG SCLC EXT SYMPTOM METS CLINICAL TRIAL  
VA-ONC LUNG SCLC EXT SYMPTOM METS REFER RAD ONC  
VA-ONC LUNG SCLC INCIDENTAL DISCOVERY PATHWAY [C]  
VA-ONC LUNG SCLC INCIDENTAL INCOMPLETE RESECT OFF PATH RSN  
VA-ONC LUNG SCLC INCIDENTAL INCOMPLETE RESECT OFF PATH TX  
VA-ONC LUNG SCLC INCIDENTAL INCOMPLETE RESECT CIS+ETO+RAD  
VA-ONC LUNG SCLC INCIDENTAL INCOMPLETE RESECT CLINICAL TRIAL  
VA-ONC LUNG SCLC INCIDENTAL COMPLETE RESECT OFF PATH RSN  
VA-ONC LUNG SCLC INCIDENTAL COMPLETE RESECT OFF PATH TX  
VA-ONC LUNG SCLC INCIDENTAL COMPLETE RESECT CIS+ETO  
VA-ONC LUNG SCLC INCIDENTAL COMPLETE RESECT CLINICAL TRIAL  
VA-ONC LUNG NSCLC STAGE IVA M1B/M1A PATHWAY [C]  
VA-ONC LUNG IVA M1B/M1A AMEN DEF TX OFF PATHWAY REASON  
VA-ONC LUNG IVA M1B/M1A AMEN DEF TX OFF PATHWAY TX  
VA-ONC LUNG IVA M1B/M1A AMEN DEF TX REFER RAD ONC/ONC/SURG  
VA-ONC LUNG IVA M1B/M1A AMEN DEF TX CLINICAL TRIAL  
VA-ONC LUNG NON-SQUAMOUS RELAPSE PATHWAY [C]  
VA-ONC LUNG NSR 3RD LINE AMEN RAD OFF PATH RSN  
VA-ONC LUNG NSR 3RD LINE AMEN RAD OFF PATH TX  
VA-ONC LUNG NSR 3RD LINE AMEN RAD ONC  
VA-ONC LUNG NSR 3RD LINE AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG NSR 2ND IMMUNO AMEN RAD OFF PATH RSN  
VA-ONC LUNG NSR 2ND IMMUNO AMEN RAD OFF PATH TX  
VA-ONC LUNG NSR 2ND IMMUNO AMEN RAD ONC  
VA-ONC LUNG NSR 2ND IMMUNO AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO AMEN RAD OFF PATH RSN  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO AMEN RAD OFF PATH TX  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO AMEN RAD ONC  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO AMEN RAD CLINICAL TRIAL  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO EGFR EXON 20 INS MOBOCERTINIB  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG EGFREX20 MOBOCERTINIB  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO HER2 OFF PATH RSN  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO HER2 OFF PATH TX  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO HER2 TRASTUZUMAB  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO HER2 CLINICAL TRIAL  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG HER2 OFF PATH RSN  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG HER2 OFF PATH TX  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG HER2 TRASTUZUMAB

VA-ONC LUNG NSR 2ND S/P IMMUNO PROG HER2 CLINICAL TRIAL  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO KRAS G12 SOTORASIB  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG NO KRASG12/EGFREX20 OFF RSN  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG NO KRASG12/EGFRX20 OFF TX  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG NO KRASG12/EGFREX20 DOC  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG NO KRASG12/EGFREX20 CLIN TR  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG EGFREX20 OFF PATH REASON  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG EGFREX20 OFF PATH TREATMENT  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG EGFREX20 AMIVANTAMAB  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG EGFREX20 CLIN TRIAL  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO NO KRASG12/EGFREX20 OFF RSN  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO NO KRASG12/EGFREX20 OFF TX  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO NO KRASG12/EGFREX20  
DOCETAXEL  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO NO KRASG12/EGFREX20 CLIN TR  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG KRAS G12C OFF PATH RSN  
VA-ONC LUNG NSR 2ND IMMUNO ALONE OFF PATHWAY REASON  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO EGFR EXON 20 INS OFF PATH RSN  
VA-ONC LUNG NSR 3RD LINE PREV DOCETAXEL OFF PATHWAY REASON  
VA-ONC LUNG NSR 3RD LINE NO PREV DOCE TX OFF PATH REASON  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO KRAS G12 OFF PATH REASON  
VA-ONC LUNG NSR 3RD LINE PREV DOCETAXEL OFF PATHWAY TX  
VA-ONC LUNG NSR 3RD LINE PREV DOCETAXEL REFER TO HOSPICE  
VA-ONC LUNG NSR 3RD LINE PREV DOCETAXEL CLINICAL TRIAL  
VA-ONC LUNG NSR 3RD LINE NO PREV DOCE TX OFF PATH TX  
VA-ONC LUNG NSR 3RD LINE NO PREV DOCE TX DOCETAXEL  
VA-ONC LUNG NSR 3RD LINE NO PREV DOCE TX CLINICAL TRIAL  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG KRAS G12C OFF PATH TX  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG KRAS G12C SOTORASIB  
VA-ONC LUNG NSR 2ND S/P IMMUNO PROG KRAS G12C CLIN TRIAL  
VA-ONC LUNG NSR 2ND IMMUNO ALONE OFF PATHWAY TX  
VA-ONC LUNG NSR 2ND IMMUNO ALONE CARBOPLATIN+PEMET REXED  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO EGFR EXON 20 INS OFF PATH TX  
VA-ONC LUNG NSR 2ND IMMUNO ALONE CLINICAL TRIAL  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO EGFR EXON 20 INS CLIN TRIAL  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO KRAS G12 OFF PATH TX  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO EGFR EXON 20 INS AMIVANTAMAB  
VA-ONC LUNG NSR 2ND CHEMO+/-IMMUNO KRAS G12 CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB MUTATION NEGATIVE PATHWAY [C]  
VA-ONC LUNG IVA/IVB NEG MUT NO /CHEMO/IMMUNO OFF PATHWAY REASON  
VA-ONC LUNG IVA/IVB NEG MUT NO CHEMO/IMMUNO OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB NEG MUT NO CHEMO/IMMUNO HOSPICE  
VA-ONC LUNG IVA/IVB NEG MUT NO CHEMO/IMMUNO CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>=50 CHEMO PEMBROLIZUMAB ALONE  
VA-ONC LUNG IVA/IVB NEG MUT NO IMMUNO OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB NEG MUT NO IMMUNO CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB NEG MUT NO IMMUNO OFF PATHWAY REASON  
VA-ONC LUNG IVA/IVB NEG MUT NO IMMUNO CARBO+PEM  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>=50 CHEMO OFF PATHWAY REASON  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>1 NO CHEMO OFF PATHWAY REASON  
VA-ONC LUNG IVA/IVB NEG MUT PDL1<50 CHEMO OFF PATHWAY REASON

VA-ONC LUNG IVA/IVB NEG MUT PDL1>1 NO CHEMO OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>1 NO CHEMO PEMBRO ALONE  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>=50 CHEMO OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>=50 CHEMO CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>1 NO CHEMO CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB NEG MUT PDL1>=50 CHEMO EMBRO+CARBO+PEMETREX  
VA-ONC LUNG IVA/IVB NEG MUT PDL1<50 CHEMO OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB NEG MUT PDL1<50 CHEMO PEMBRO+CARBO+PEMETREX  
VA-ONC LUNG IVA/IVB NEG MUT PDL1<50 CHEMO CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB MUTATION POSITIVE PATHWAY [C]  
VA-ONC LUNG IVA/IVB SUB POS MUT EGFR 19/21 L858R OFF PATH REASON  
VA-ONC LUNG IVA/IVB SUB POS MUT EGFR 19/21 L858R FF PATH TX  
VA-ONC LUNG IVA/IVB SUB POS MUT EGFR 19/21 L858R DOCETAXEL  
VA-ONC LUNG IVA/IVB SUB POS MUT EGFR 19/21 L858R CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB SUB POS MUT ALL OTHER MUT OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB SUB POS MUT ALK GENE FUSION OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALL OTHER MUT OFF PATHWAY REASON  
VA-ONC LUNG IVA/IVB 2ND POS MUT EXON19/21 L858R OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALK GENE FUSION OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 1ST POS MUT ROS1 GENE FUSION OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 1ST POS MUT RET REARRANGE OFF PATHWAY REASON  
VA-ONC LUNG IVA/IVB 1ST POS MUT NTRK GENE FUSION OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 1ST POS MUT MET EXON 14 SKIP OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 1ST POS MUT EGFR 19/21 L858R OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 1ST POS MUT BRAF V600E OFF PATHWAY REASON  
VA-ONC LUNG IVA/IVB 1ST POS MUT ALK GENE FUSION OFF PATHWAY RSN  
VA-ONC LUNG IVA/IVB 1ST POS MUT NTRK GENE FUSION ENTRECTINIB  
VA-ONC LUNG IVA/IVB SUB POS MUT ALL OTHER MUT OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB SUB POS MUT ALL OTHER MUT DOCETAXEL  
VA-ONC LUNG IVA/IVB SUB POS MUT ALL OTHER MUT CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB SUB POS MUT ALK GENE FUSION OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB SUB POS MUT ALK GENE FUSION DOCETAXEL  
VA-ONC LUNG IVA/IVB SUB POS MUT ALK GENE FUSION CARBO+PEM+/-CI  
VA-ONC LUNG IVA/IVB SUB POS MUT ALK GENE FUSION CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALL OTHER MUT OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALL OTHER MUT CARBO+PEM+/-CI  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALL OTHER MUT CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALK GENE FUSION OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALK GENE FUSION LORLATINIB  
VA-ONC LUNG IVA/IVB 2ND POS MUT ALK GENE FUSION CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 2ND POS MUT EXON19/21 L858R OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 2ND POS MUT EXON19/21 L858R CARBO+PEM+/-CI  
VA-ONC LUNG IVA/IVB 2ND POS MUT EXON19/21 L858R CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 1ST POS MUT NON STANDARD NPOP CONSULT  
VA-ONC LUNG IVA/IVB 1ST POS MUT NTRK GENE FUSION OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 1ST POS MUT NTRK GENE FUSION CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 1ST POS MUT BRAF V600E OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 1ST POS MUT BRAF V600E DABRAFENIB+TRAMETINIB  
VA-ONC LUNG IVA/IVB 1ST POS MUT BRAF V600E CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 1ST POS MUT ROS1 GENE FUSION OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 1ST POS MUT ROS1 GENE FUSION ENTRECTINIB



VA-ONC LUNG IVA/IVB 1ST POS MUT ROS1 GENE FUSION CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 1ST POS MUT RET REARRANGE OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 1ST POS MUT RET REARRANGE SELPERCATNIB  
VA-ONC LUNG IVA/IVB 1ST POS MUT RET REARRANGE CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 1ST POS MUT MET EXON 14 SKIP OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 1ST POS MUT MET EXON 14 SKIP CAPMATINIB  
VA-ONC LUNG IVA/IVB 1ST POS MUT MET EXON 14 SKIP CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 1ST POS MUT ALK GENE FUSION OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 1ST POS MUT ALK GENE FUSION ALECTINIB  
VA-ONC LUNG IVA/IVB 1ST POS MUT ALK GENE FUSION CLINICAL TRIAL  
VA-ONC LUNG IVA/IVB 1ST POS MUT EGFR 19/21 L858R OFF PATHWAY TX  
VA-ONC LUNG IVA/IVB 1ST POS MUT EGFR 19/21 L858R OSIMERTINIB  
VA-ONC LUNG IVA/IVB 1ST POS MUT EGFR 19/21 L858R CLINICAL TRIAL  
VA-ONC LUNG SQUAMOUS STAGE IVB PATHWAY [C]  
VA-ONC LUNG SQUIVB PDL1 >=50% PEM ALONE  
VA-ONC LUNG SQUIVB NO IMMUNO\NO CHEMO OFF PATH REASON  
VA-ONC LUNG SQUIVB NO IMMUNO\NO CHEMO OFF PATH TX  
VA-ONC LUNG SQUIVB NO IMMUNO\NO CHEMO REFER HOSPICE  
VA-ONC LUNG SQUIVB NO IMMUNO\NO CHEMO CLINICAL TRIAL  
VA-ONC LUNG SQUIVB PD-L1>=1% IMMUNO ONLY OFF PATH REASON  
VA-ONC LUNG SQUIVB PD-L1>=1% IMMUNO ONLY OFF PATH TX  
VA-ONC LUNG SQUIVB PD-L1>=1% IMMUNO ONLY PEMBRO  
VA-ONC LUNG SQUIVB PD-L1>=1% IMMUNO ONLY CLINICAL TRIAL  
VA-ONC LUNG SQUIVB CHEMO ONLY OFF PATH REASON  
VA-ONC LUNG SQUIVB CHEMO ONLY OFF PATH TX  
VA-ONC LUNG SQUIVB CHEMO ONLY CARBO+PAC  
VA-ONC LUNG SQUIVB CHEMO ONLY CLINICAL TRIAL  
VA-ONC LUNG SQUIVB PDL1 >=50% OFF PATHWAY REASON  
VA-ONC LUNG SQUIVB PDL1 >=50% OFF PATHWAY TX  
VA-ONC LUNG SQUIVB PDL1 >=50% PEM+CARBO+PAC +PEM MAIN  
VA-ONC LUNG SQUIVB PDL1 >=50% CLINICAL TRIAL  
VA-ONC LUNG SQUIVB PDL1 0-49% OFF PATHWAY REASON  
VA-ONC LUNG SQUIVB PDL1 0-49% OFF PATHWAY TX  
VA-ONC LUNG SQUIVB PDL1 0-49% PEM+CARBO+PAC  
VA-ONC LUNG SQUIVB PDL1 0-49% CLINICAL TRIAL  
VA-ONC LUNG SQUIVB SYMPTOMATIC REFER RAD ONC  
VA-ONC LUNG METASTATIC NSCLC MOLECULAR TESTING [C]  
VA-ONC LUNG TISSUE SUFFICIENCY UNKNOWN  
VA-ONC LUNG CGP ORDER FUTURE DATE  
VA-ONC LUNG PD-L1 ORDER FUTURE DATE  
VA-ONC LUNG LIQUID BIOPSY: ON PATH  
VA-ONC LUNG TISSUE SUFFICIENT: OFF PATH  
VA-ONC LUNG REASSESS: ON PATH  
VA-ONC LUNG REPEAT BIOPSY: ON PATH  
VA-ONC LUNG TISSUE SUFFICIENT: ON PATH  
VA-ONC LUNG CGP ORDER: OFF PATH  
VA-ONC LUNG CGP ORDER: ON PATH  
VA-ONC LUNG PD-L1 ORDER: OFF PATH  
VA-ONC LUNG PD-L1 ORDER: ON PATH  
VA-ONC LUNG NSCLC STAGE III UNRESECTABLE CLINICAL PATHWAY [C]  
VA-ONC LUNG NSCLC III UNRESECT S/P CHEMO OFF PATH RSN

VA-ONC LUNG NSCLC III UNRESECT S/P CHEMO OFF PATH TX  
 VA-ONC LUNG NSCLC III UNRESECT S/P CHEMO SURVEILLANCE  
 VA-ONC LUNG NSCLC III UNRESECT S/P CHEMO CLINICAL TRIAL  
 VA-ONC LUNG NSCLC III UNRESECT SUB NO IMMUNO OFF PATH RSN  
 VA-ONC LUNG NSCLC III UNRESECT SUB NO IMMUNO OFF PATH TX  
 VA-ONC LUNG NSCLC III UNRESECT SUB NO IMMUNO SURVEILLANCE  
 VA-ONC LUNG NSCLC III UNRESECT SUB NO IMMUNO CLIN TRIAL  
 VA-ONC LUNG NSCLC III UNRESECT SUB IMMUNO OFF PATH RSN  
 VA-ONC LUNG NSCLC III UNRESECT SUB IMMUNO OFF PATH TX  
 VA-ONC LUNG NSCLC III UNRESECT SUB IMMUNO DURVALUMAB  
 VA-ONC LUNG NSCLC III UNRESECT SUB IMMUNO CLIN TRIAL  
 VA-ONC LUNG NSCLC III UNRESECT INI SQU OFF PATH RSN  
 VA-ONC LUNG NSCLC III UNRESECT INI SQU OFF PATH TX  
 VA-ONC LUNG NSCLC III UNRESECT INI SQU CARBO+PAC+RAD  
 VA-ONC LUNG NSCLC III UNRESECT INI SQU CLINICAL TRIAL  
 VA-ONC LUNG NSCLC III UNRESECT INI NSQ OFF PATH RSN  
 VA-ONC LUNG NSCLC III UNRESECT INI NSQ OFF PATH TX  
 VA-ONC LUNG NSCLC III UNRESECT INI NSQ CARBO+PEM+RAD  
 VA-ONC LUNG NSCLC III UNRESECT INI NSQ CLINICAL TRIAL  
 VA-ONC LUNG NSCLC III UNRESECT INI CON RAD OFF PATH RSN  
 VA-ONC LUNG NSCLC III UNRESECT INI CON RAD OFF PATH TX  
 VA-ONC LUNG NSCLC III UNRESECT INI CON RAD CARBO+PAC+RAD  
 VA-ONC LUNG NSCLC III UNRESECT INI CON RAD CLINICAL TRIAL  
 VA-ONC LUNG PANCOAST CLINICAL PATHWAY [C]  
 VA-ONC LUNG PAN NO SURG/YES CON CIS/ETO OFF PATH RSN  
 VA-ONC LUNG PAN NO SURG/YES CON CIS/ETO OFF PATH TX  
 VA-ONC LUNG PAN NO SURG/YES CON CIS/ETO ARBO+PAC+RAD+DURVALUMAB  
 VA-ONC LUNG PAN NO SURG/YES CON CIS/ETO CLINICAL TRIAL  
 VA-ONC LUNG PAN NO SURG/NO CON CIS/ETO OFF PATH RSN  
 VA-ONC LUNG PAN NO SURG/NO CON CIS/ETO OFF PATH TX  
 VA-ONC LUNG PAN NO SURG/NO CON CIS/ETO +RAD+DURVALUMAB  
 VA-ONC LUNG PAN SURG/YES CON CIS/ETO OFF PATH RSN  
 VA-ONC LUNG PAN NO SURG/NO CON CIS/ETO CLINICAL TRIAL  
 VA-ONC LUNG PAN SURG/YES CON CIS/ETO OFF PATH TX  
 VA-ONC LUNG PAN SURG/YES CON CIS/ETO CARBO+PAC+RAD+SURG  
 VA-ONC LUNG PAN SURG/YES CON CIS/ETO CLINICAL TRIAL  
 VA-ONC LUNG PAN SURG/NO CON CIS/ETO OFF PATH RSN  
 VA-ONC LUNG PAN SURG/NO CON CIS/ETO OFF PATH TX  
 VA-ONC LUNG PAN SURG/NO CON CIS/ETO +RAD+SURG  
 VA-ONC LUNG PAN SURG/NO CON CIS/ETO CLINICAL TRIAL  
 VA-ONC LUNG NSCLC IIIA CLINICAL PATHWAY [C]  
 VA-ONC LUNG NSCLC IIIA R1/R2 NO RRS OFF PATH RSN  
 VA-ONC LUNG NSCLC IIIA R1/R2 NO RRS OFF PATH TX  
 VA-ONC LUNG NSCLC IIIA R1/R2 NO RRS REF RADONC  
 VA-ONC LUNG NSCLC IIIA R1/R2 NO RRS CLIN TRIAL  
 VA-ONC LUNG NSCLC IIIA R1/R2 RRS OFF PATH RSN  
 VA-ONC LUNG NSCLC IIIA R1/R2 RRS OFF PATH TX  
 VA-ONC LUNG NSCLC IIIA R1/R2 RRS RERESECT  
 VA-ONC LUNG NSCLC IIIA R1/R2 RRS CLIN TRIAL  
 VA-ONC LUNG NSCLC IIIA R0 OFF PATH RSN  
 VA-ONC LUNG NSCLC IIIA R0 OFF PATH TX

VA-ONC LUNG NSCLC IIIA R0 SURVEILLANCE  
VA-ONC LUNG NSCLC IIIA R0 CLIN TRIAL  
VA-ONC LUNG NSCLC IIIA NO PROG OFF PATH RSN  
VA-ONC LUNG NSCLC IIIA NO PROG OFF PATH TX  
VA-ONC LUNG NSCLC IIIA NO PROG SURGERY  
VA-ONC LUNG NSCLC IIIA NO PROG CLIN TRIAL  
VA-ONC LUNG NSCLC IIIA INITIAL TX SQU OFF PATH RSN  
VA-ONC LUNG NSCLC IIIA INITIAL TX SQU OFF PATH TX  
VA-ONC LUNG NSCLC IIIA INITIAL TX SQU CIS+DOCETAXEL  
VA-ONC LUNG NSCLC IIIA INITIAL TX SQU CLIN TRIAL  
VA-ONC LUNG NSCLC IIIA INITIAL TX NSQ OFF PATH RSN  
VA-ONC LUNG NSCLC IIIA INITIAL TX NSQ OFF PATH TX  
VA-ONC LUNG NSCLC IIIA INITIAL TX NSQ CIS+PEMETREXED  
VA-ONC LUNG NSCLC IIIA INITIAL TX NSQ CLIN TRIAL  
VA-ONC LUNG NSCLC IIA/IIB CLINICAL PATHWAY [C]  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS SQU OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT SQU OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT SQU OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT SQU CIS+DOC+REF RAD ONC  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT SQU CLINICAL TRIAL  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT NSQ OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT NSQ OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT NSQ CIS+PEME+REF RAD ONC  
VA-ONC LUNG NSCLC IIA/IIB POST RESECT NSQ CLINICAL TRIAL  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS SQU OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS SQU CIS+DOC+RAD ONC  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS SQU CLIN TRIAL  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS NSQ OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS NSQ OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS NSQ CIS+PEME+RAD ONC  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 NO RRS NSQ CLIN TRIAL  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 RRS OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 RRS OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 RRS RERESECT  
VA-ONC LUNG NSCLC IIA/IIB R1/R2 RRS CLIN TRIAL  
VA-ONC LUNG NSCLC IIA/IIB R0 SQU OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB R0 SQU OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB R0 SQU CLIN TRIAL  
VA-ONC LUNG NSCLC IIA/IIB R0 NSQ OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB R0 NSQ OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB R0 NSQ CIS+PEME  
VA-ONC LUNG NSCLC IIA/IIB R0 SQU CIS+DOCE  
VA-ONC LUNG NSCLC IIA/IIB R0 NSQ CLIN TRIAL  
VA-ONC LUNG NSCLC IIA/IIB NO SURG OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB NO SURG OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB NO SURG REF ONC/RAD ONC  
VA-ONC LUNG NSCLC IIA/IIB NO SURG CLIN TRIAL  
VA-ONC LUNG NSCLC IIA/IIB SURG OFF PATH RSN  
VA-ONC LUNG NSCLC IIA/IIB SURG OFF PATH TX  
VA-ONC LUNG NSCLC IIA/IIB SURG RESECTION  
VA-ONC LUNG NSCLC IIA/IIB SURG CLINICAL TRIAL

VA-ONC LUNG NSCLC STAGE IA/IB PATHWAY[C]  
VA-ONC LUNG NSCLC IA/IB SURG SUBLOBAR RESECTION  
VA-ONC LUNG NSCLC IA/IB SURG SEGMENTECTOMY  
VA-ONC LUNG NSCLC IA/IB PATH NEG MARGINS SURVEILLANCE  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG NO RERESECT OFF PATH RSN  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG NO RERESECT OFF PATH TX  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG NO RERESECT REF RAD ONC  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG NO RERESECT CLINICAL TRIAL  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG RERESECT OFF PATH RSN  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG RERESECT OFF PATH TX  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG RERESECT  
VA-ONC LUNG NSCLC IA/IB PATH POS MARG RERESECT CLIN TRIAL  
VA-ONC LUNG NSCLC IA/IB PATH NEG MARGINS OFF PATH RSN  
VA-ONC LUNG NSCLC IA/IB PATH NEG MARGINS OFF PATH TX  
VA-ONC LUNG NSCLC IA/IB PATH NEGATIVE MARGINS SURVEILLANCE  
VA-ONC LUNG NSCLC IA/IB PATH NEG MARGINS CLIN TRIAL  
VA-ONC LUNG NSCLC PATH IB NEG MARG CLIN TRIAL  
VA-ONC LUNG NSCLC PATH IB POS MARG NO RERESECT OFF PATH RSN  
VA-ONC LUNG NSCLC PATH IB POS MARG NO RERESECT OFF PATH TX  
VA-ONC LUNG NSCLC PATH IB POS MARG NO RERESECT REF RAD ONC  
VA-ONC LUNG NSCLC PATH IB POS MARG NO RERESECT CLINICAL TRIAL  
VA-ONC LUNG NSCLC PATH IB POS MARG RERESECT OFF PATH RSN  
VA-ONC LUNG NSCLC PATH IB POS MARG RERESECT OFF PATH TX  
VA-ONC LUNG NSCLC PATH IB POS MARG RERESECT  
VA-ONC LUNG NSCLC PATH IB POS MARG RERESECT CLIN TRIAL  
VA-ONC LUNG NSCLC PATH IB NEG MARG OFF PATHWAY RSN  
VA-ONC LUNG NSCLC PATH IB NEG MARG OFF PATHWAY TX  
VA-ONC LUNG NSCLC PATH IB NEG MARG SURVEILLANCE  
VA-ONC LUNG NSCLC PATH IB NEG MARG CLINICAL TRIAL  
VA-ONC LUNG NSCLC PATH IA POS MARG NO RERESECT OFF PATH RSN  
VA-ONC LUNG NSCLC PATH IA POS MARG NO RERESECT OFF PATH TX  
VA-ONC LUNG NSCLC PATH IA POS MARG NO RERESECT REF RAD ONC  
VA-ONC LUNG NSCLC PATH IA POS MARG NO RERESECT CLINICAL TRIAL  
VA-ONC LUNG NSCLC PATH IA POS MARG RERESECT OFF PATH RSN  
VA-ONC LUNG NSCLC PATH IA POS MARG RERESECT OFF PATH TX  
VA-ONC LUNG NSCLC PATH IA POS MARG RERESECT  
VA-ONC LUNG NSCLC PATH IA POS MARG RERESECT CLIN TRIAL  
VA-ONC LUNG NSCLC PATH IA NEG MARGINS OFF PATH RSN  
VA-ONC LUNG NSCLC PATH IA NEG MARGINS OFF PATH TX  
VA-ONC LUNG NSCLC PATH IA NEG MARGINS SURVEILLANCE  
VA-ONC LUNG NSCLC PATH IA NEG MARGINS CLIN TRIAL  
VA-ONC LUNG NSCLC IA/IB SURG NO OFF PATH RSN  
VA-ONC LUNG NSCLC IA/IB SURG NO OFF PATH TX  
VA-ONC LUNG NSCLC IA/IB SURG NO REF RAD ONC  
VA-ONC LUNG NSCLC IA/IB SURG NO CLIN TRIAL  
VA-ONC LUNG NSCLC IA/IB SURG OFF PATH RSN  
VA-ONC LUNG NSCLC IA/IB SURG OFF PATH TX  
VA-ONC LUNG NSCLC IA/IB SURG LOBECTOMY  
VA-ONC LUNG NSCLC IA/IB SURG CLIN TRIAL  
ONC LUNG LLL-WEDGE RESECTION  
ONC LUNG LUL-WEDGE RESECTION

ONC LUNG RLL-WEDGE RESECTION  
ONC LUNG RML-WEDGE RESECTION  
ONC LUNG RUL-WEDGE RESECTION  
ONC LUNG LLL-SEGMENTECTOMY  
ONC LUNG RML-SEGMENTECTOMY  
ONC LUNG LUL-SEGMENTECTOMY  
ONC LUNG RUL-SEGMENTECTOMY  
ONC LUNG RIGHT SIDE-PNEUMONECTOMY  
ONC LUNG LEFT SIDE-PNEUMONECTOMY  
ONC LUNG LLL-LOBECTOMY  
ONC LUNG LUL-LOBECTOMY  
ONC LUNG RLL-LOBECTOMY  
ONC LUNG RML-LOBECTOMY  
ONC LUNG RUL-LOBECTOMY  
ONC LUNG TRAMETINIB STOP DATE  
ONC LUNG TRAMETINIB START DATE  
ONC LUNG DABRAFENIB STOP DATE  
ONC LUNG DABRAFENIB START DATE  
ONC LUNG BRIGATINIB STOP DATE  
ONC LUNG BRIGATINIB START DATE  
ONC LUNG ATEZOLUMUMAB STOP DATE  
ONC LUNG ATEZOLUMUMAB START DATE  
ONC LUNG DURVALUMAB STOP DATE  
ONC LUNG DURVALUMAB START DATE  
ONC LUNG CONSENT FOR CHEMO NOT OBTAINED  
ONC LUNG CONSENT FOR CHEMO OBTAINED  
ONC LUNG TOPOTECAN STOP DATE  
ONC LUNG TOPOTECAN START DATE  
ONC LUNG STEREOTACTIC BODY RT-UNKNOWN  
ONC LUNG RT INTENT DISCUSSED-UNKNOWN  
ONC LUNG CHEMO INTENT DISCUSSED-UNKNOWN  
ONC LUNG CLINICAL TRIAL  
ONC LUNG CLINICAL TRIAL STOP DATE  
ONC LUNG CLINICAL TRIAL START DATE  
ONC LUNG RT INTENT-SINGLE MODALITY  
ONC LUNG CHEMO INTENT-SINGLE MODALITY  
ONC LUNG RAMUCIRUMAB STOP DATE  
ONC LUNG RAMUCIRUMAB START DATE  
ONC LUNG PEMBROLIZUMAB STOP DATE  
ONC LUNG PEMBROLIZUMAB START DATE  
ONC LUNG NIVOLUMAB STOP DATE  
ONC LUNG NIVOLUMAB START DATE  
ONC LUNG ALECTINIB STOP DATE  
ONC LUNG ALECTINIB START DATE  
ONC LUNG NECITUMUMAB STOP DATE  
ONC LUNG NECITUMUMAB START DATE  
ONC LUNG OSIMERTINIB STOP DATE  
ONC LUNG OSIMERTINIB START DATE  
ONC LUNG GEFITINIB STOP DATE  
ONC LUNG GEFITINIB START DATE  
ONC LUNG AFATINIB STOP DATE

ONC LUNG AFATINIB START DATE  
ONC LUNG OTHER RT STOP DATE  
ONC LUNG OTHER RT START DATE  
ONC LUNG BONE RT STOP DATE  
ONC LUNG BONE RT START DATE  
ONC LUNG BRAIN RT STOP DATE  
ONC LUNG BRAIN RT START DATE  
ONC LUNG CHEST RT STOP DATE  
ONC LUNG CHEST RT START DATE  
ONC LUNG CHEMO OTHER STOP DATE  
ONC LUNG CHEMO OTHER START DATE  
ONC LUNG BEVACIZUMAB STOP DATE  
ONC LUNG PEMETREXED STOP DATE  
ONC LUNG PEMETREXED START DATE  
ONC LUNG CHEMO STOP DATE  
ONC LUNG CHEMO START DATE  
ONC LUNG BEVACIZUMAB START DATE  
ONC LUNG CERITINIB STOP DATE  
ONC LUNG CERITINIB START DATE  
ONC LUNG CRIZOTINIB STOP DATE  
ONC LUNG CRIZOTINIB START DATE  
ONC LUNG ERLOTINIB STOP DATE  
ONC LUNG ERLOTINIB START DATE  
ONC LUNG IRINOTECAN STOP DATE  
ONC LUNG IRINOTECAN START DATE  
ONC LUNG ETOPOSIDE STOP DATE  
ONC LUNG ETOPOSIDE START DATE  
ONC LUNG VINOURELBINE STOP DATE  
ONC LUNG VINOURELBINE START DATE  
ONC LUNG GEMITABINE STOP DATE  
ONC LUNG GEMITABINE START DATE  
ONC LUNG PACLITAXEL STOP DATE  
ONC LUNG PACLITAXEL START DATE  
ONC LUNG DOCETAXEL STOP DATE  
ONC LUNG DOCETAXEL START DATE  
ONC LUNG CARBOPLATIN STOP DATE  
ONC LUNG CARBOPLATIN START DATE  
ONC LUNG CISPLATIN STOP DATE  
ONC LUNG CISPLATIN START DATE  
ONC LUNG CHEMO TYPE-CERITINIB  
ONC LUNG RT TYPE-CONCURRENT  
ONC LUNG CHEMO TYPE-CRIZOTINIB  
ONC LUNG RT STOPPED-SITE CHANGE  
ONC LUNG RT STOPPED-OTHER  
ONC LUNG RT STOPPED-TOXICITY  
ONC LUNG RT STOPPED-PATIENT REQUEST  
ONC LUNG RT STOPPED-PROGRESSION  
ONC LUNG RT STOPPED-PLAN COMPLETED  
ONC LUNG REFERRAL-YES  
ONC LUNG PCU REFERRAL-NO  
ONC LUNG HOSPICE REFERRAL

ONC LUNG REFERRAL-PATIENT DECLINED  
ONC LUNG REFERRAL-NOT APPROPRIATE  
ONC LUNG REFERRAL-SERVICE NOT AVAILABLE  
ONC LUNG REFERRAL-OTHER  
ONC LUNG PRESENTED TO TUMOR BOARD-YES  
ONC LUNG PRESENTED TO TUMOR BOARD-NO  
ONC LUNG PRESENTED TO TUMOR BOARD-N/A  
ONC LUNG SURGERY TYPE-LOBECTOMY  
ONC LUNG SURGERY TYPE-PNEUMONECTOMY  
ONC LUNG SURGERY TYPE-SEGMENTECTOMY  
ONC LUNG SURGERY TYPE-WEDGE RESECTION  
ONC LUNG SURGERY OUTCOME-CLEAR MARGINS  
ONC LUNG SURGERY OUTCOME-+ MARGINS  
ONC LUNG SURGERY OUTCOME-NOT RESECTABLE  
ONC LUNG SURGERY OUTCOME-RESULT PENDING  
ONC LUNG SURGERY OUTCOME-NOT DETERMINED  
ONC LUNG SURGERY-NOT RECOMMENDED  
ONC LUNG SURGERY NO-PATIENT DECLINED  
ONC LUNG SURGERY NO-OTHER REASON  
ONC LUNG LN SAMPLING-ULTRASOUND  
ONC LUNG LN SAMPLING-MEDIASTINOSCOPY  
ONC LUNG LN SAMPLING-SURGERY  
ONC LUNG LN SAMPLING DONE-NO  
ONC LUNG CHEMO INTENT-ADJUVANT  
ONC LUNG CHEMO INTENT-PALLIATIVE  
ONC LUNG CHEMO INTENT-NOT DETERMINED  
ONC LUNG RT TYPE-ADJUVANT  
ONC LUNG CHEMO INTENT DISCUSSED-YES  
ONC LUNG CHEMO INTENT DISCUSSED-NO  
ONC LUNG RT TYPE-NEOADJUVANT  
ONC LUNG CHEMO TYPE-CISPLATIN  
ONC LUNG CHEMO TYPE-CARBOPLATIN  
ONC LUNG CHEMO TYPE-VINORELBINE  
ONC LUNG CHEMO TYPE-DOCETAXEL  
ONC LUNG CHEMO TYPE-PACLITAXEL  
ONC LUNG CHEMO TYPE-PEMETREXED  
ONC LUNG CHEMO TYPE-GEMCITABINE  
ONC LUNG CHEMO TYPE-ETOPOSIDE  
ONC LUNG CHEMO TYPE-IRINOTECAN  
ONC LUNG CHEMO TYPE-ERLOTINIB  
ONC LUNG CHEMO TYPE-BEVACIZUMAB  
ONC LUNG CHEMO TYPE-OTHER  
ONC LUNG CHEMO STOPPED-PLAN COMPLETED  
ONC LUNG CHEMO STOPPED-TOXICITY  
ONC LUNG CHEMO STOPPED-PROGRESSION  
ONC LUNG CHEMO STOPPED-PATIENT REQUEST  
ONC LUNG CHEMO STOPPED-OTHER  
ONC LUNG CHEMO NOT ADMIN-PT DECLINED  
ONC LUNG CHEMO NOT ADMIN-OTHER  
ONC LUNG CHEMO INTENT-CONCURRENT  
ONC LUNG DOUBLET NO-CONTRAINDICATED

ONC LUNG DOUBLET NO-DECLINED  
ONC LUNG DOUBLET NO-N/A  
ONC LUNG DOUBLET NO-OTHER  
ONC LUNG CHEMO INTENT-NEOADJUVANT  
ONC LUNG CHEMO NOT ADMIN-NOT RECOMMENDED  
ONC LUNG RT INTENT-PALLIATIVE  
ONC LUNG RT INTENT-CURATIVE  
ONC LUNG RT INTENT-NOT DETERMINED  
ONC LUNG RT INTENT DISCUSSED-YES  
ONC LUNG RT INTENT DISCUSSED-NO  
ONC LUNG STEREOTACTIC BODY RT-YES  
ONC LUNG STEREOTACTIC BODY RT-NO  
ONC LUNG RT SITE-CHEST  
ONC LUNG RT SITE-BRAIN  
ONC LUNG RT SITE-BONE  
ONC LUNG RT SITE-OTHER  
ONC LUNG RT NOT ADMIN-CONTRAINDICATED  
ONC LUNG RT NOT ADMIN-PT DECLINED  
ONC LUNG RT NOT ADMIN-OTHER  
ONCOLOGY (ONC) LUNG CANCER DIAGNOSES [C]  
VA-ONC LUNG METASTASIS OTHER  
VA-ONC LUNG METASTASIS CONTRALATERAL LUNG  
VA-ONC LUNG METASTASIS IPSILATERAL LUNG  
VA-ONC LUNG METASTASIS ADRENAL GLAND  
VA-ONC LUNG METASTASIS BONES  
VA-ONC LUNG METASTASIS LIVER  
VA-ONC LUNG METASTASIS BRAIN  
VA-ONC LUNG METASTASIS LYMPH NODES  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T0  
ONC LUNG CLINICAL PRIMARY TUMOR-T0  
ONC LUNG CANCER-RML  
ONC LUNG CANCER-RLL  
ONC LUNG CANCER-LUL  
ONC LUNG CANCER-LLL  
ONC LUNG CANCER-RUL  
ONC LUNG RLL-SEGMENTECTOMY  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T2  
ONC LUNG CLINICAL PRIMARY TUMOR-T2  
ONC LUNG CLINICAL PRIMARY TUMOR-T1  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T1  
ONC LUNG CLINICAL PRIMARY TUMOR-T1A(MI)  
ONC LUNG PATH PRIMARY TUMOR-T1A(MI)  
ONC LUNG CLIN DISTANT METASTASIS-M1C  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-TIS  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T1C  
ONC LUNG PATH DISTANT METASTASIS-M1C  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IVB  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IVA  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IIIC  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IA3  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IA2



ONC LUNG PATHOLOGICAL-SUMMARY STAGE IA1  
ONC LUNG CLINICAL PRIMARY TUMOR-TIS  
ONC LUNG CLINICAL PRIMARY TUMOR-T1C  
ONC LUNG CLINICAL-SUMMARY STAGE IVB  
ONC LUNG CLINICAL-SUMMARY STAGE IVA  
ONC LUNG CLINICAL-SUMMARY STAGE IIIC  
ONC LUNG CLINICAL-SUMMARY STAGE IA3  
ONC LUNG CLINICAL-SUMMARY STAGE IA2  
ONC LUNG CLINICAL-SUMMARY STAGE IA1  
ONC LUNG PRECISION ONCOLOGY  
ONC LUNG PD-L1-RESULTS PENDING  
ONC LUNG PD-L1-NOT DONE  
ONC LUNG PD-L1-NEGATIVE  
ONC LUNG PD-L1-POSITIVE  
ONC LUNG TISSUE SAMPLE SENT TO POP  
ONC LUNG TISSUE SAMPLE SENT  
ONC LUNG CLIN DISTANT METASTASIS-M1B  
ONC LUNG CLIN DISTANT METASTASIS-M1A  
ONC LUNG CLIN DISTANT METASTASIS-M1  
ONC LUNG CLIN DISTANT METASTASIS-M0  
ONC LUNG CLIN DISTANT METASTASIS-MX  
ONC LUNG CLIN REGIONAL LYMPH NODES-N3  
ONC LUNG CLIN REGIONAL LYMPH NODES-N2  
ONC LUNG CLIN REGIONAL LYMPH NODES-N1  
ONC LUNG CLIN REGIONAL LYMPH NODES-N0  
ONC LUNG CLIN REGIONAL LYMPH NODES-NX  
ONC LUNG CLINICAL PRIMARY TUMOR-T4  
ONC LUNG CLINICAL PRIMARY TUMOR-T3  
ONC LUNG CLINICAL PRIMARY TUMOR-T2B  
ONC LUNG CLINICAL PRIMARY TUMOR-T2A  
ONC LUNG CLINICAL PRIMARY TUMOR-T1B  
ONC LUNG CLINICAL PRIMARY TUMOR-T1A  
ONC LUNG CLINICAL PRIMARY TUMOR-TX  
ONC LUNG CLINICAL TNM STAGE DATE-  
ONC LUNG PATH DISTANT METASTASIS-M1B  
ONC LUNG PATH DISTANT METASTASIS-M1A  
ONC LUNG PATH DISTANT METASTASIS-M1  
ONC LUNG PATH DISTANT METASTASIS-M0  
ONC LUNG PATH DISTANT METASTASIS-MX  
ONC LUNG PATH REGIONAL LYMPH NODES-N3  
ONC LUNG PATH REGIONAL LYMPH NODES-N2  
ONC LUNG PATH REGIONAL LYMPH NODES-N1  
ONC LUNG PATH REGIONAL LYMPH NODES-N0  
ONC LUNG PATH REGIONAL LYMPH NODES-NX  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T4  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T3  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T2B  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T2A  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T1B  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-T1A  
ONC LUNG PATHOLOGICAL PRIMARY TUMOR-TX

ONC LUNG PATHOLOGICAL TNM STAGE DATE-  
ONC LUNG HISTOLOGY-OTHER  
ONC LUNG CLINICAL-SUMMARY STAGE 0  
ONC LUNG CLINICAL-SUMMARY STAGE I  
ONC LUNG CLINICAL-SUMMARY STAGE IA  
ONC LUNG CLINICAL-SUMMARY STAGE IB  
ONC LUNG CLINICAL-SUMMARY STAGE II  
ONC LUNG CLINICAL-SUMMARY STAGE IIA  
ONC LUNG CLINICAL-SUMMARY STAGE IIB  
ONC LUNG CLINICAL-SUMMARY STAGE III  
ONC LUNG CLINICAL-SUMMARY STAGE IIIA  
ONC LUNG CLINICAL-SUMMARY STAGE IIIB  
ONC LUNG CLINICAL-SUMMARY STAGE IV  
ONC LUNG CLINICAL-SUMMARY STAGE UNKNOWN  
ONC LUNG CLINICAL-SUMMARY STAGE LIM  
ONC LUNG CLINICAL-SUMMARY STAGE EXTENS  
ONC LUNG EGFR TEST-NOT DONE  
ONC LUNG EGFR TEST-RESULTS PENDING  
ONC LUNG EGFR TEST-WILD TYPE  
ONC LUNG EGFR TEST-MUTANT  
ONC LUNG ALK TRANSLOCATION-NOT DONE  
ONC LUNG ALK TRANSLOCATION-RESULT PENDING  
ONC LUNG ALK TRANSLOCATION-NORMAL  
ONC LUNG ALK TRANSLOCATION-REARRANGED  
ONC LUNG FIRST ABNORMAL RADIOLOGY DATE-  
ONC LUNG INITIAL DIAGNOSIS DATE-  
ONC LUNG HISTOLOGY-ADENOCARCINOMA  
ONC LUNG HISTOLOGY-SQUAMOUS CELL  
ONC LUNG HISTOLOGY-LARGE CELL  
ONC LUNG HISTOLOGY-NSCLC, NOS  
ONC LUNG HISTOLOGY-SCLC  
ONC LUNG ECOG-0  
ONC LUNG ECOG-1  
ONC LUNG ECOG-2  
ONC LUNG ECOG-3  
ONC LUNG ECOG-4  
ONC LUNG ECOG-5  
ONC LUNG ECOG-UNKNOWN  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE 0  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE I  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IA  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IB  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE II  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IIA  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IIB  
ONC LUNG PATHOLOGICAL-SUMM STAGE III  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IIIA  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IIIB  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE IV  
ONC LUNG PATHOLOGICAL-SUMMARY STAGE UNK  
ONC LUNG PATHOLOGICAL-SUMM STAGE LIM

ONC LUNG PATHOLOGICAL-SUMM STAGE EXT  
ONC LUNG MOLECULAR TEST-OTHER  
VA-ONCOLOGY PATHWAYS [C]  
VA-ONC LUNG SQUAMOUS STAGE IVB PATHWAY  
VA-ONC LUNG IVA/IVB MUTATION NEGATIVE PATHWAY  
VA-ONC LUNG IVA/IVB MUTATION POSITIVE PATHWAY  
VA-ONC LUNG NSCLC STAGE IVA M1B/M1A PATHWAY  
VA-ONC LUNG MOLECULAR TESTING  
VA-ONC GOALS OF CARE [C]  
VA-ONC GOALS OF CARE PALLIATIVE  
VA-ONC GOALS OF CARE CURATIVE  
VA-ONC PAIN MANAGEMENT PLAN [C]  
VA-ONC PAIN PLAN INITIAL  
VA-ONC PAIN PLAN MGT BY HEMONC  
VA-ONC PAIN PLAN ADEQUATE  
VA-ONC PAIN PLAN NEEDS ADJUSTMENT  
VA-ONC PAIN MGT PLAN OPIOID PRESCRIBED  
VA-ONC PAIN MGT PLAN BOWEL REGIMEN  
VA-ONC PAIN MGT PLAN NON OPIOID PLAN  
VA-ONC PAIN PLAN MANAGED BY  
VA-ONC PAIN PLAN DECLINED  
VA-ONC PAIN DENIES  
VA-ONC LUNG AJCC EDITION [C]  
VA-ONC LUNG AJCC 8TH EDITION  
VA-ONC LUNG AJCC 7TH EDITION  
VA-ONC CN TIME SPENT CHART REVIEW  
VA-ONC CN TIME SPENT FACE TO FACE  
VA-ONC CN TIME SPENT VVC CHART REVIEW  
VA-ONC CN TIME SPENT VVC FACE TO FACE  
VA-ONC CN TIME SPENT CVT FACE TO FACE  
VA-ONC CN TIME SPENT CVT CHART REVIEW  
VA-ONC LUNG NOTE TITLES [C]  
VA-ONC LUNG CANCER METASTASIS IVA-IVB NOTE  
VA-ONC LUNG SQUAMOUS 2ND/3RD LINE PATHWAY  
VA-ONC LUNG NON-SQUAMOUS 2ND/3RD LINE PATHWAY  
VA-ONC LUNG CANCER 2ND/3RD LINE METASTATIC NOTE  
VA-ONC LUNG NSCLC STAGE III UNRESECTABLE CLINICAL PATHWAY  
VA-ONC LUNG PANCOAST CLINICAL PATHWAY  
VA-ONC LUNG CANCER ADVANCED STAGE IIIB-IIIC & PANCOAST NOTE  
VA-ONC LUNG SCLC RELAPSE PATHWAY  
VA-ONC LUNG SCLC LIMITED STAGE 1ST LINE PATHWAY  
VA-ONC LUNG SCLC EXTENSIVE STAGE CLINICAL PATHWAY  
VA-ONC LUNG SCLC INCIDENTAL DISCOVERY PATHWAY  
VA-ONC LUNG CANCER SMALL CELL (SCLC) NOTE  
VA-ONC LUNG NSCLC STAGE IA/IB PATHWAY  
VA-ONC LUNG NSCLC IIA/IIB CLINICAL PATHWAY  
VA-ONC LUNG NSCLC IIAB/IIIA RESECTABLE CLINICAL PATHWAY [C]  
VA-ONC LUNG NSCLC IIAB/IIIA NO PROG AFTER 1ST LINE CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA NO PROG AFTER 1ST LINE SURGERY  
VA-ONC LUNG NSCLC IIAB/IIIA NO PROG AFTER 1ST LINE OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA NO PROG AFTER 1ST LINE OFF PATH RSN

VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU NO CIS CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU NO CIS OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU NO CIS OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU NO CIS PAC+CAR+NIV  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU CAND CIS OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU CAND CIS OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU CAND CIS PAC+CIS+NIV  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT SQU CAND CIS CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ NO CIS OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ NO CIS OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ NO CIS PEM+CAR+NIV  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ NO CIS CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ CAND CIS OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ CAND CIS OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ CAND CIS PEM+CIS+NIV  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT NSQ CAND CIS CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA EGFR/ALK MUT CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA EGFR/ALK SURGERY  
VA-ONC LUNG NSCLC IIAB/IIIA EGFR/ALK OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA EGFR/ALK OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R0 CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R0 SURVEILLANCE  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R0 OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R0 OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 NO RRS CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 NO RRS REF RADONC  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 NO RRS OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 NO RRS OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 RRS RERESECT  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 RRS CLIN TRIAL  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 RRS OFF PATH RSN  
VA-ONC LUNG NSCLC IIAB/IIIA NO MUT R1/R2 RRS OFF PATH TX  
VA-ONC LUNG NSCLC IIAB/IIIA SURG/NEO ADJ PERFORM MOLE TEST  
VA-ONC LUNG NSCLC IIAB/IIIA RESECTABLE CLINICAL PATHWAY  
VA-ONC LUNG CANCER EARLY STAGE IA-III A NOTE

**REMINDER SPONSOR**

VHA ONCOLOGY  
Office of Nursing Service

**REMINDER TERM**

VA-REMINDER UPDATE\_2\_0\_367  
VA-ONC CN BRANCH LUNG  
VA-ONC GOALS OF CARE RT  
VA-ONC PAIN EFFECTIVENESS RT  
VA-ONC LUNG METASTASIS RT  
VA-ONC LUNG HISTOLOGY RT  
VA-ONC LUNG PERFORMANCE STATUS RT

VA-ONC LUNG MOLECULAR TESTING RT  
VA-ONC LUNG CLINICAL STAGE RT  
VA-ONC LUNG DATE OF DIAGNOSIS RT  
VA-ONC ASU ATTENDING

## **HEALTH SUMMARY COMPONENT**

PCE HEALTH FACTORS SELECTED

## **HEALTH SUMMARY TYPE**

VA-ONC CN PREVIOUS LUNG (HS)  
VA-ONCOLOGY PATHWAYS LUNG (HS)  
VA-ONC CN LUNG HISTORY (HS)  
VA-ONC GOALS OF CARE HS  
VA-ONC PCA PAIN  
VA-ONC LUNG METASTASIS  
VA-ONC LUNG HISTOLOGY  
VA-ONC LUNG ECOG PS  
VA-ONC LUNG MOLECULAR TESTS  
VA-ONC LUNG STAGING  
VA-ONC LUNG DIAGNOSIS DATE

## **HEALTH SUMMARY OBJECTS**

VA-ONC CN PREVIOUS LUNG (TIU)  
VA-ONCOLOGY PATHWAYS LUNG (TIU)  
VA-ONC CN LUNG HISTORY (TIU)  
VA-ONC GOALS OF CARE HS (TIU)  
VA-ONC PCA PAIN (TIU)  
VA-ONC LUNG METASTASIS (TIU)  
VA-ONC LUNG HISTOLOGY  
VA-ONC LUNG ECOG PS  
VA-ONC LUNG MOLECULAR TESTS  
VA-ONC LUNG STAGING  
VA-ONC LUNG DIAGNOSIS DATE

## **TIU DOCUMENT DEFINITION**

VA-ONC CN PREVIOUS LUNG  
VA-ONCOLOGY PATHWAYS LUNG  
VA-ONC CN LUNG HISTORY  
VA-ONC GOALS OF CARE HS  
VA-ONC PCA PAIN  
VA-ONC LUNG METASTASIS  
VA-ONC LUNG HISTOLOGY  
VA-ONC LUNG ECOG PS  
VA-ONC LUNG MOLECULAR TEST  
VA-ONC LUNG STAGING  
VA-ONC LUNG DIAGNOSIS DATES

## REMINDER DIALOG

VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP  
VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP  
VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLOW-UP  
VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP  
VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA FOLLOW-UP

## Install Details

This update is being distributed as a web host file. The address for the host file is:  
[https://REDACTED/UPDATE\\_2\\_0\\_367.PRD](https://REDACTED/UPDATE_2_0_367.PRD)

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 15 minutes, total setup time could take close to an hour.

## Install Example

Go to Reminder exchange in Vista.

To Load the exchange file, use LWH. The URL is

[https://REDACTED/UPDATE\\_2\\_0\\_367.PRD](https://REDACTED/UPDATE_2_0_367.PRD)

```
+ + Next Screen - Prev Screen ?? More Actions >>>
CFE Create Exchange File Entry LHF Load Host File
CHF Create Host File LMM Load MailMan Message
CMM Create MailMan Message LR List Reminder Definitions
DFE Delete Exchange File Entry LWH Load Web Host File
IFE Install Exchange File Entry RI Reminder Definition Inquiry
IH Installation History RP Repack
Select Action: Next Screen// lwh Load Web Host File
Input the URL for the .prd file: https://REDACTED/UPDATE_2_0_367.PRD
```

You will see a message indicating the file was loaded successfully.

Search and locate an entry titled UPDATE\_2\_0\_367 VA-ONC LUNG CANCER FOLLOW UP AND CONSULT REPORT in reminder exchange.

```
+Item Entry Source Date Packed
546 UPDATE_2_0_367 VA-ONC LUNG CANCER FOLLOW UP AND CONSULT REPORT
547 UPDATE_2_0_37 MH UPDATES
548 UPDATE_2_0_39 VA-OSP SUICIDE SAFETY PLAN
549 UPDATE_2_0_3_2 TBI SCREENING
550 UPDATE_2_0_4 VA-VETERANS CHOICE/MAMMOGRAM/TBI UPDATE
551 UPDATE_2_0_4 VA-VETERANS
+ + Next Screen - Prev Screen ?? More Actions >>>
CFE Create Exchange File Entry LHF Load Host File
CHF Create Host File LMM Load MailMan Message
CMM Create MailMan Message LR List Reminder Definitions
DFE Delete Exchange File Entry LWH Load Web Host File
IFE Install Exchange File Entry RI Reminder Definition Inquiry
IH Installation History RP Repack
Select Action: Next Screen// IFE Install Exchange File Entry
Enter a list or range of numbers (1-774): 546
```

At the Select Action prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry titled UPDATE\_2\_0\_367 VA-ONC LUNG CANCER FOLLOW UP AND CONSULT REPORT (in this example it is entry 546, it will vary by site) that is dated 07/27/2023.

```

Date Packed: 08/07/2023 08:11:00
Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA FOLLOW-UP
VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP
VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLOW-UP
VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP
VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP

REMINDER TERM
VA-REMINDER UPDATE_2_0_367
Non-exchangeable TIU object(s):
+ Enter ?? for more actions >>>
IA Install all Components IS Install Selected Component
Select Action: Next Screen// IA
  
```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.  
**Install all new components.**  
**Overwrite any components that exist but are different.**

You will then be prompted to install the first dialog component

```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP [NA]
  
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP	dialog	
2	20	VA-GP ONC CN LUNG 1ST LINE METS MAIN Finding: VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)	group	X
3	20.10	VA-ONC CN LUNG PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
4	20.15	VA-ONC CN LUNG OBJECT Finding: *NONE*	element	X
5	20.25	VA-ONC CN LUNG PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
6	20.30	VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS Finding: *NONE*	element	X
+ + Next Screen - Prev Screen ?? More Actions				
DD		Dialog Details	DT Dialog Text	IS Install Selected
DF		Dialog Findings	DU Dialog Usage	QU Quit
DS		Dialog Summary	IA Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog – **VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP**  
 Install reminder dialog and all components with no further changes: Y//ES

```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP [NA
VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP (reminder dialog) installed
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB             dialog  X
      FOLLOW-UP
  2    20 VA-GP ONC CN LUNG 1ST LINE METS MAIN                group   X
      Finding: VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB NOTE
      (HEALTH FACTOR)
      Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)
  3    20.10 VA-ONC CN LUNG PREVIOUS HF VIEW                  element  X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)
  4    20.15 VA-ONC CN LUNG OBJECT                             element  X
      Finding: *NONE*
  5    20.25 VA-ONC CN LUNG PREVIOUS NOTE VIEW                element  X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)
  6    20.30 VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS        element  X
      Finding: *NONE*
+      + Next Screen    - Prev Screen    ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will then be prompted to install the next dialog component.

```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-U
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC       dialog
      FOLLOW-UP
  2    20 VA-GP ONC CN LUNG CANCER 2ND/3RD LINE METASTATIC   group
      MAIN
      Finding: VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC
      NOTE (HEALTH FACTOR)
      Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)
  3    20.10 VA-ONC CN LUNG PREVIOUS HF VIEW                  element  X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)
  4    20.15 VA-ONC CN LUNG OBJECT                             element  X
      Finding: *NONE*
  5    20.25 VA-ONC CN LUNG PREVIOUS NOTE VIEW                element  X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)
  6    20.30 VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS        element  X
+      + Next Screen    - Prev Screen    ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP**

Install reminder dialog and all components with no further changes: Y//ES



```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-U
VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP (reminder dialog) inst
Item  Seq.  Dialog Findings                                     Type  Exists
  1    VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC
      FOLLOW-UP                                           dialog  X
  2    20 VA-GP ONC CN LUNG CANCER 2ND/3RD LINE METASTATIC
      MAIN                                               group   X
      Finding: VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC
      NOTE (HEALTH FACTOR)                               X
      Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)       X
  3    20.10 VA-ONC CN LUNG PREVIOUS HF VIEW              element X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  4    20.15 VA-ONC CN LUNG OBJECT                         element X
      Finding: *NONE*
  5    20.25 VA-ONC CN LUNG PREVIOUS NOTE VIEW            element X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  6    20.30 VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS   element X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will then be prompted to install the next dialog component.

```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLO
Item  Seq.  Dialog Findings                                     Type  Exists
  1    VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST
      FOLLOW-UP                                           dialog
  2    20 VA-GP ONC CN LUNG ADVANCED STAGE MAIN           group
      Finding: VA-ONCOLOGY LUNG CANCER ADVANCED STAGE IIIB-IIIC &
      PANCOAST NOTE (HEALTH FACTOR)                       X
      Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)       X
  3    20.10 VA-ONC CN LUNG PREVIOUS HF VIEW              element
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  4    PXRM GF VIEW BUTTON                                prompt  X
  5    20.15 VA-ONC CN LUNG OBJECT                         element
      Finding: *NONE*
  6    20.25 VA-ONC CN LUNG PREVIOUS NOTE VIEW            element
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  7    PXRM GF VIEW BUTTON                                prompt  X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLOW-UP**

Install reminder dialog and all components with no further changes: Y//ES

```
Packed reminder dialog: VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLO
VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLOW-UP (reminder dialog) i
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLOW-UP	dialog	X
2	20	VA-GP ONC CN LUNG ADVANCED STAGE MAIN Finding: VA-ONCOLOGY LUNG CANCER ADVANCED STAGE IIIB-IIIC & PANCOAST NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)	group	X
3	20.10	VA-ONC CN LUNG PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
4		PXRM GF VIEW BUTTON	prompt	X
5	20.15	VA-ONC CN LUNG OBJECT Finding: *NONE*	element	X
6	20.25	VA-ONC CN LUNG PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
7		PXRM GF VIEW BUTTON	prompt	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details	DT	Dialog Text	IS Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU Quit
DS	Dialog Summary	IA	Install All	
Select Action: Next Screen// Q				

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will then be prompted to install the next dialog component

```
Packed reminder dialog: VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP [NATIONAL D
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP	dialog	
2	20	VA-GP ONC CN LUNG SCLC MAIN Finding: VA-ONCOLOGY LUNG CANCER SMALL CELL (SCLC) NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)	group	X
3	20.10	VA-ONC CN LUNG PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
4	20.15	VA-ONC CN LUNG OBJECT Finding: *NONE*	element	X
5	20.25	VA-ONC CN LUNG PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
6	20.30	VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS Finding: *NONE*	element	X
7	20.40	VA-GP ONC CN LUNG REQUIRED	group	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details	DT	Dialog Text	IS Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU Quit
DS	Dialog Summary	IA	Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog – **VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP**

Install reminder dialog and all components with no further changes: Y//ES

```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP [NATIONAL D
VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP (reminder dialog) installed from ex
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP          dialog  X
  2    20 VA-GP ONC CN LUNG SCLC MAIN                          group  X
      Finding: VA-ONCOLOGY LUNG CANCER SMALL CELL (SCLC) NOTE
      (HEALTH FACTOR)                                       X
      Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)          X
  3    20.10 VA-ONC CN LUNG PREVIOUS HF VIEW                 element X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  4    20.15 VA-ONC CN LUNG OBJECT                           element X
      Finding: *NONE*
  5    20.25 VA-ONC CN LUNG PREVIOUS NOTE VIEW              element X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  6    20.30 VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS      element X
      Finding: *NONE*
  7    20.40 VA-GP ONC CN LUNG REQUIRED                        group  X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will then be prompted to install the next dialog component

```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP [N
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-III A          dialog  X
      FOLLOW-UP
  2    20 VA-GP ONC CN LUNG EARLY MAIN                          group  X
      Finding: VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-III A NOTE
      (HEALTH FACTOR)                                       X
      Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)          X
  3    20.10 VA-ONC CN LUNG PREVIOUS HF VIEW                 element X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  4    20.15 VA-ONC CN LUNG OBJECT                           element X
      Finding: *NONE*
  5    20.25 VA-ONC CN LUNG PREVIOUS NOTE VIEW              element X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) X
  6    20.30 VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS      element X
      Finding: *NONE*
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP**

Install reminder dialog and all components with no further changes: Y//ES

```

Packed reminder dialog: VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA FOLLOW-UP [N
VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA FOLLOW-UP (reminder dialog) install
Item  Seq.  Dialog Findings                                     Type  Exists
  1    VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA          dialog  X
      FOLLOW-UP
  2    20 VA-GP ONC CN LUNG EARLY MAIN                        group   X
      Finding: VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA NOTE
      (HEALTH FACTOR)
      Add. Finding: VA-ONC CN LUNG (HEALTH FACTOR)
  3    20.10 VA-ONC CN LUNG PREVIOUS HF VIEW                element  X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)
  4    20.15 VA-ONC CN LUNG OBJECT                          element  X
      Finding: *NONE*
  5    20.25 VA-ONC CN LUNG PREVIOUS NOTE VIEW              element  X
      Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)
  6    20.30 VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS      element  X
      Finding: *NONE*
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```

Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA FOLLOW-UP
  VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP
  VA-ONCOLOGY LUNG CANCER STAGE IIIIB-IIIC & PANCOAST FOLLOW-UP
  VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP
  VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP

REMINDER TERM
  VA-REMINDER UPDATE_2_0_367
Non-exchangeable TIU object(s):
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components      IS  Install Selected Component
Select Action: Next Screen// Q

```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.

# Post Installation

---

## 1. Setup of Note Titles.

The note titles to create to associate with this template are:

ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA FOLLOW-UP  
ONCOLOGY LUNG CANCER STAGE IIIB-IIIC AND PANCOAST FOLLOW-UP  
ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP  
ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP  
ONCOLOGY LUNG CANCER SMALL CELL (SCLC) FOLLOW-UP

VHA Enterprise Standard Title: HEMATOLOGY AND ONCOLOGY **NOTE**

You can add these to your document class where progress notes used by Oncology are located. They should not be in your consult document class.

In addition to those notes, you will also need to create 5 consult note titles.

ONCOLOGY LUNG CANCER EARLY STAGE IA-IIIA CONSULT  
ONCOLOGY LUNG CANCER STAGE IIIB-IIIC AND PANCOAST CONSULT  
ONCOLOGY LUNG CANCER METASTASIS IVA-IVB CONSULT  
ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC CONSULT  
ONCOLOGY LUNG CANCER SMALL CELL (SCLC) CONSULT

VHA Enterprise Standard Title: HEMATOLOGY AND ONCOLOGY **CONSULT**

These note titles should be in your consult document class.

**Note:** Comprehensive information on Creating new Document Classes and Titles can be found beginning on page 54 of the TIU/ASU Implementation Guide at the following link: [http://www.va.gov/vdl/documents/Clinical/CPRS-Text Integration Utility \(TIU\)/tiuim.pdf](http://www.va.gov/vdl/documents/Clinical/CPRS-Text%20Integration%20Utility%20(TIU)/tiuim.pdf)

## 2. Make the Dialog template available to be attached to a Personal or Shared template or to be able to attach the dialog to a progress note title. (You can use the LINK feature in this menu for steps 2 and 3 if you choose to)

Select Reminder Managers Menu <TEST ACCOUNT> Option: **CP CPRS**

### Reminder Configuration

CA Add/Edit Reminder Categories  
CL CPRS Lookup Categories  
CS CPRS Cover Sheet Reminder List  
MH Mental Health Dialogs Active  
PN Progress Note Headers  
RA Reminder GUI Resolution Active  
**TIU TIU Template Reminder Dialog Parameter**  
DEVL Evaluate Coversheet List on Dialog Finish  
DL Default Outside Location  
PT Position Reminder Text at Cursor  
LINK Link Reminder Dialog to Template

TEST CPRS Coversheet Time Test  
NP New Reminder Parameters  
GEC GEC Status Check Active  
WH WH Print Now Active

Select CPRS Reminder Configuration <TEST ACCOUNT> Option: **TIU TIU**  
**Template Reminder Dialog Parameter**

Reminder Dialogs allowed as Templates may be set for the following:

- 1 User USR [choose from NEW PERSON]
- 3 Service SRV [choose from SERVICE/SECTION]
- 4 Division DIV [SALT LAKE CITY]
- 5 System SYS [NATREM.FO-SLC.MED.VA.GOV]

Enter selection: **5** Choose *the appropriate number for your site. Your site may do this by System or other levels. The example below uses SYSTEM level.*

Setting Reminder Dialogs allowed as Templates for System: NATREM.FO-SLC.MED.VA.GOV

Select Display Sequence:

- 388 VA-COVID-19 TELEPHONE FOLLOW-UP
- 389 VA-SRA-CONSULT RESPONSE (D)

*When you type a question mark above, you will see the list of #'s (with dialogs) that are already taken. Choose a number **NOT** on this list. For this example, 391 is not present so I will use 391.*

Select Display Sequence: **391**

Are you adding 391 as a new Display Sequence? Yes// **YES**

Display Sequence: 391// <Enter> 391

Clinical Reminder Dialog: **VA-ONCOLOGY LUNG CANCER STAGE IIIB-IIIC & PANCOAST FOLLOW-UP**

reminder dialog NATIONAL  
...OK? Yes// <Enter> (Yes)

Select Display Sequence: **392**

Are you adding 392 as a new Display Sequence? Yes// **YES**

Display Sequence: 392// <Enter> 392

Clinical Reminder Dialog: **VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP**

reminder dialog NATIONAL  
...OK? Yes// <Enter> (Yes)

Select Display Sequence: **393**

Are you adding 393 as a new Display Sequence? Yes// **YES**

Display Sequence: 393// <Enter> 393

Clinical Reminder Dialog: **VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP**

reminder dialog NATIONAL  
...OK? Yes// <Enter> (Yes)

Select Display Sequence: **394**

Are you adding 394 as a new Display Sequence? Yes// **YES**

Display Sequence: 394// **<Enter>** 394

Clinical Reminder Dialog: **VA-ONCOLOGY LUNG CANCER METASTASIS IVA-  
IVB FOLLOW-UP**

reminder dialog NATIONAL

...OK? Yes// **<Enter>** (Yes)

Select Display Sequence: **395**

Are you adding 395 as a new Display Sequence? Yes// **YES**

Display Sequence: 395// **<Enter>** 395


Clinical Reminder Dialog: **VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-  
IIIA FOLLOW-UP**

reminder dialog NATIONAL

...OK? Yes// **<Enter>** (Yes)


3. Associate the reminder dialogs with the note title in CPRS. In addition to connecting to the note titles, you should also add these to a shared template folder for your oncology service.

In CPRS, open the Template Editor, expand Document Titles, select the appropriate folder (if applicable), then click New Template (top right-hand corner).



- a. Type in the following information:  
 Name: **ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP**  
 Template Type: **Reminder Dialog**  
 Reminder Dialog: **VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP**


Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP** - then click Apply.



- b. Type in the following information:  
 Name: **ONCOLOGY LUNG CANCER STAGE III B-III C AND PANCOAST FOLLOW-UP**  
 Template Type: **Reminder Dialog**  
 Reminder Dialog: **VA-ONCOLOGY LUNG CANCER STAGE III B-III C AND PANCOAST FOLLOW-UP**

Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER STAGE III B-III C AND PANCOAST FOLLOW-UP** - then click Apply.






c. Type in the following information:

Name: **ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP**

Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP** - then click Apply.




d. Type in the following information:

Name: **ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA- ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP**

Next, type in the Associated Title name – **ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP** – then click Apply.




e. Type in the following information:

Name: **ONCOLOGY LUNG CANCER SMALL CELL (SCLC) FOLLOW-UP**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA-ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP**

Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER SMALL CELL (SCLC) FOLLOW-UP** – then click Apply.




f. Type in the following information:

Name: **ONCOLOGY LUNG CANCER EARLY STAGE IA-III A CONSULT**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA-ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP**

Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER EARLY STAGE IA-III A CONSULT** - then click Apply.




g. Type in the following information:

Name: **ONCOLOGY LUNG CANCER STAGE IIB-IIIIC AND PANCOAST CONSULT**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA-ONCOLOGY LUNG CANCER STAGE IIB-IIIIC AND PANCOAST FOLLOW-UP**

Next, type in the Associated Title name – **ONCOLOGY LUNG CANCER STAGE IIB-IIIIC AND PANCOAST CONSULT** – then click Apply.




h. Type in the following information:

Name: **ONCOLOGY LUNG CANCER METASTASIS IVA-IVB CONSULT**

Template Type: **Reminder Dialog**


Reminder Dialog: **VA-ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP**

Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER METASTASIS IVA-IVB CONSULT** – then click Apply.



- i. Type in the following information:  
 Name: **ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC CONSULT**  
 Template Type: **Reminder Dialog**  
 Reminder Dialog: **VA-ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP**

Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC CONSULT** – then click Apply.



- j. Type in the following information:  
 Name: **ONCOLOGY LUNG CANCER SMALL CELL (SCLC) CONSULT**  
 Template Type: **Reminder Dialog**  
 Reminder Dialog: **VA- ONCOLOGY LUNG CANCER SMALL CELL FOLLOW-UP**

Next, type in the Associated Title name - **ONCOLOGY LUNG CANCER SMALL CELL (SCLC) CONSULT** – then click Apply.

#### 4. Add objects to the reminder dialog.

Using the VistA menu **Reminder Dialog Management (DM)**, **Reminder Dialogs (DI)**, and **Change View (CV)** to **Dialog Elements (E)** view.

a. **Locate and edit the reminder dialog element:**

VAL-ONC CN PATIENT DEMOGRAPHICS

Edit the DIALOG/PROGRESS NOTE TEXT and **Replace the following data objects** with your local data objects if you do not have the included object.

- |PATIENT SEX|
- |PATIENT RACE|
- |PATIENT AGE|
- |PATIENT NAME

b. **Locate and edit the reminder dialog element:**

VAL-ONC CN VITALS TXT

Edit the DIALOG/PROGRESS NOTE TEXT and **replace the text** (add local object here) **with your local data object for the following:**

- Temperature
- Pulse
- Respiration
- Pulse Oximetry
- Blood Pressure
- Pain
- Height
- Weight

c. **Locate and edit the reminder dialog element:**

VAL-ONC CN LABS & IMAGING

EDIT the DIALOG/PROGRESS NOTE TEXT and **replace the text:** (please reference official report in record) with your local data objects for the most recent LABS & IMAGING.

**\*Or your site can choose to leave the verbiage “(please reference official report in record)”** without inserting data objects.

d. **Locate and edit the reminder dialog element:**

VAL-ONC CN VISIT VVC CONTACT HOME ADDRESS

EDIT the DIALOG/PROGRESS NOTE TEXT and **Replace the text** Home address: ADD local patient address object here **with your local data object.**

e. **Locate and edit the reminder dialog element:**

VAL-ONC CN VISIT VVC CONTACT PHONE NUMBER HOME/CELL.

EDIT the DIALOG/PROGRESS NOTE TEXT and **replace** (ADD patient xxx phone object here) **with your local data object the following:**

- Cell phone
- Home phone
- Work phone

5. **Suppress the resident supervision statement (if applicable):**

The template utilized branching logic via the element VAL-ONC CN ATTENDING OF RECORD to suppress documenting resident supervision in the dialog if needed.

\*This is a local element and can be disabled if not applicable to your site.

You should add user classes that do not need to document resident supervision (i.e. Attending Physician, etc.). To add user classes, edit the reminder term VA-ONC ASU ATTENDING. To add an additional user class, use the quote "" characters "CF.VA-ASU USER CLASS" and set the computed finding parameter of your new finding to the appropriate user class.

The reminder term is being released with the ASU class of ATTENDING PHYSICIAN and NURSE PRACTITIONER. Edit according to your local practices. If a user does not have one of the included user classes, the following will appear in the note.

The Attending of record for this patient is Dr. \* . This case has been reviewed by the Attending and they agree with the recommended plan.

6. **Enter your site name in the template field:** {FLD:VA-ONC VAMC SITE NAME}.

7. **Add consult quick orders as an additional finding to the elements identified in the pre-install:**

- VA-ONC RADIATION ONCOLOGY CONSULT – map radiation oncology consult
- VA-ONC LUNG REFER TO PALLIATIVE CARE - map palliative care consult
- VAL-ONC LUNG IVB POS 1ST OTHER – map IFC PRECISION ONCOLOGY ELECTRONIC CONSULT (element can also be inactivated if you don't have the Precision Oncology consult).
- VA-ONC HOSPICE REFERRAL - map hospice consult
- VA-ONC LUNG SURGICAL CONSULT – map thoracic surgery consult

**Add lab quick orders to these elements. Check with your local LIM. Your site may not have any lab orders for these tests; some sites use a consult for molecular testing. If your site does not have either, nothing else is needed.**

- VA-EL ONC LUNG ON PDL1 - map PD-L1 order
- VA-EL ONC LUNG PTHWY NGS ON – map Tumor Genomic testing order
- VA-EL ONC LUNG PTHWY LIQUID BIOPSY - map order for liquid biopsy

**8. Map the ten note titles to the HS Object [VA-ONC CN PREVIOUS LUNG]:**

NOTE TITLES:

- ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP
- ONCOLOGY LUNG CANCER STAGE IIIB-IIIC AND PANCOAST FOLLOW-UP
- ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP
- ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP
- ONCOLOGY LUNG CANCER SMALL CELL (SCLC) FOLLOW-UP
- ONCOLOGY LUNG CANCER EARLY STAGE IA-III A CONSULT
- ONCOLOGY LUNG CANCER STAGE IIIB-IIIC AND PANCOAST CONSULT
- ONCOLOGY LUNG CANCER METASTASIS IVA-IVB CONSULT
- ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC CONSULT
- ONCOLOGY LUNG CANCER SMALL CELL (SCLC) CONSULT

HS OBJECT:

[VA-ONC CN PREVIOUS LUNG]

Example:

Start with the TIU Maintenance Menu → Document Definitions (Manager) → Create TIU/Health Summary Objects → SL “VA-ONC CN PREVIOUS “ → Detailed Display/Edit HS Object

Map the note title the 10 note titles to the following object.

HS Object: VA-ONC CN PREVIOUS LUNG (TIU)  
Health Summary Type: VA-ONC CN PREVIOUS LUNG (HS)  
Report Period:  
Creator:

HS Object

Print Label: NO	Print Report Date and Time: NO
Print Blank Line after Label: NO	Print Confidentiality Banner: NO
Customized Header: YES	Print Report Date and Time: NO
Suppress Components w/o Data: YES	Print Component Header: NO

NO

Print Deceased Information: NO

Print Time-Occurrence Limits:

National Object: NO

Underline Component

Header: NO

Blank Line After Header: NO

Overwrite No Data: No Teleoncology note found

Select Item(s): Quit// **Edit HS Type**

Editing Health Summary Type VA-ONC CN PREVIOUS LUNG (HS)'

NAME: VA-ONC CN PREVIOUS LUNG (HS) Replace

TITLE:

SUPPRESS PRINT OF COMPONENTS WITHOUT DATA: yes//

LOCK:

OWNER:

Select COMPONENT: **SPN**

**SPN PROGRESS NOTES SELECTED SPN**

SUMMARY ORDER: 5// 5

OCCURRENCE LIMIT: **1**

TIME LIMIT:

HEADER NAME: Selected Prog Notes//

No selection items chosen.

Select new items one at a time in the sequence you want them displayed.  
You may select any number of items.

Select SELECTION ITEM: **ONCOLOGY LUNG CANCER EARLY  
STAGE IA-III A FOLLOW-UP**

Continue to add all the note titles listed above. When you are finished, you  
should have the following:

PROGRESS NOTES SELECTED (SPN)

- 1 ONCOLOGY LUNG CANCER EARLY STAGE IA-III A FOLLOW-UP
- 2 ONCOLOGY LUNG CANCER STAGE IIIB-IIIC AND PANCOAST FOLLOW-UP
- 3 ONCOLOGY LUNG CANCER METASTASIS IVA-IVB FOLLOW-UP
- 4 ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC FOLLOW-UP
- 5 ONCOLOGY LUNG CANCER SMALL CELL (SCLC) FOLLOW-UP
- 6 ONCOLOGY LUNG CANCER EARLY STAGE IA-III A CONSULT
- 7 ONCOLOGY LUNG CANCER STAGE IIIB-IIIC AND PANCOAST CONSULT
- 8 ONCOLOGY LUNG CANCER METASTASIS IVA-IVB CONSULT
- 9 ONCOLOGY LUNG CANCER 2ND/3RD LINE METASTATIC CONSULT
- 10 ONCOLOGY LUNG CANCER SMALL CELL (SCLC) CONSULT



9. **\*\*Optional\*\*** - The following note titles *may* be inactivated since the pathways now exist in the follow-up notes. Please ensure if you inactivate the note titles, that the templates are available in the Oncology shared templates folder and users at your site are aware they are being inactivated.

- ONCOLOGY-EARLY STAGE LUNG CANCER PATHWAYS
- ONCOLOGY-LOCALLY ADVANCED LUNG CANCER PATHWAYS
- ONCOLOGY-FIRST LINE STAGE IV METASTATIC LUNG CANCER PATHWAYS
- ONCOLOGY-LUNG CANCER RELAPSE PATHWAYS
- ONCOLOGY-SMALL CELL LUNG CANCER PATHWAYS

10. **\*\*Optional\*\*** - Add local reminder dialog groups and elements if needed. This template was built with the standard items needed for complete documentation. If a site would like to add guidance or information specific to that site, you can add text to the local reminder dialog group VAL-GP ONC CN LUNG LOCAL GROUP using the reminder dialog editor. This group is at the end of the dialog. This dialog group will install disabled, and you will need to remove the disable flag for it to be used.