



Update _2_0_380

Clinical Reminders

**VA-ONC PROSTATE CANCER FOLLOW UP AND
CONSULT TEMPLATES**

Install Guide

November 2023

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION 1
PRE-INSTALL..... 1
INSTALL DETAILS 16
INSTALL EXAMPLE 16
POST INSTALLATION 25

Introduction

The purpose of these prostate cancer follow-up and consult notes are to standardize follow-up care, regulate data capture, and to acquire the required documentation as determined by the National Oncology Program Office (NOPO). Every VA with an Oncology program should be utilizing these notes to ensure standardized data capture for NOPO dashboards.

In addition, the prostate cancer clinical pathways are included in the follow-up/consult templates to facilitate ease of use. Oncology Clinical Pathways (CPs) are decision-support tools that translate published evidence into standardized treatment process maps to help guide cancer treatment decisions. Oncology CPs are different from treatment guidelines, as CPs recommend a standard treatment based on efficacy, toxicity, and cost. VA Oncology CPs are developed by an interdisciplinary subject matter expert (SME) group and approved by national program offices.

The prostate cancer follow-up reminder dialogs were built with the Cerner APSO (assessment, plan, subjective, objective) format. Furthermore, there are only 2 required mandatory fields. The 3 reminder dialog templates will be utilized by all VA medical centers with an Oncology service. Additionally, the 3 reminder dialog templates will also be used as consult notes and utilized by all Oncology providers:

- Nurse Practitioners
- Physician Assistants
- Residents
- Fellows
- Physicians

VHA Directive 1415 requires that all sites use the National Oncology templates. If your site does not provide oncology services, the site must request a waiver to exempt it from the requirement to use the National Oncology templates. The waiver can be requested by emailing REDACTED

REDACTED URL for EDUCATION

Pre-Install

Check to see if you have installed and completed the post installation for Clinical Reminder Update_2_0_367

If not, STOP and install Clinical Reminder Update_2_0_367 before installing this update. If installing 367, do not add objects to the elements, you will do this in the post install of this update.

Using the VistA menu Reminder Dialog Manager, Reminder Dialog, change your view to elements and do an inquire/print on the following reminder dialog elements, you will need to use this as a reference for the quick orders/objects during post-install. You may not have any quick orders mapped to these elements.

VA-EL ONC PRCA ON PATHWAY
VA-EL ONC PRCA PTHWY NOT SUFF LIQUID BIOPSY
VA-ONC PALLIATIVE CARE CONSULT
VA-ONC RADIATION ONCOLOGY CONSULT
VA-ONC UROLOGY CONSULT
VAL-ONC CN PATIENT DEMOGRAPHICS
VAL-ONC CN ALLERGIES
VAL-ONC CN MED REC OBJECT
VAL-ONC CN VITALS TXT
VAL-ONC CN LABS & IMAGING
VAL-ONC CN VISIT VVC CONTACT HOME ADDRESS
VAL-ONC CN VISIT VVC CONTACT PHONE NUMBER HOME/CELL
VA-ONC HOSPICE REFERRAL

UPDATE_2_0_380 contains 1 Reminder Exchange entry named **UPDATE_2_0_380 VA-ONC PROSTATE CANCER FOLLOW UP**

This contains 3 updated reminder dialogs.

1. VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP
2. VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP
3. VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP

The exchange file contains the following components:

REMINDER GENERAL FINDINGS

- 1 TICKLER REMINDER COMMENT
- 2 TICKLER REMINDER
- 3 VIEW PROGRESS NOTE TEXT

TIU TEMPLATE FIELD

- 4 TEXT (1-50 CHAR)
- 5 VA-ONC CN EMOTIONAL WB
- 6 VA-ONC VAMC SITE NAME
- 7 DATE
- 8 VA-ONC WP 2 LINES
- 9 VA-ONC DAY/WEEK/MONTH/YEAR
- 10 EDIT 3*
- 11 VA-ONC PRESUMPTIVE CONDITIONS
- 12 VA-ONC OFF REASON TXT
- 13 VA-ONC NEITHER ON/OFF
- 14 VA-ONC PATHWAYS SHAREPOINT
- 15 VA-ONC INSTRUCTIONS ON/OFF PATHWAY
- 16 VA-ONC TREATMENT REASON
- 17 VA-ONC CLINICAL TRIAL

- 18 VA-ONC CN PREVIOUS NOTE VIEW
- 19 EDIT 55+
- 20 TELEHEALTH EMERGENCIES VA LOCATIONS
- 21 TELEHEALTH EMERGENCIES NON-VA LOCATIONS
- 22 VA-ONC CN PRCA CHIEF COMPLAINT
- 23 VA-ONC NPOP SHAREPOINT
- 24 VA-ONC CONFIRMATION TEXT

HEALTH FACTORS

- 25 VA-REMINDER UPDATES [C]
- 26 VA-UPDATE_2_0_380
- 27 VA-ONC CN [C]
- 28 VA-ONC CN PROSTATE
- 29 VA-ONC PRCA MOLECULAR TESTING PATHWAY [C]
- 30 VA-ONC PRCA GERM CGP ORDER: ON PATH
- 31 VA-CCGS CLINICAL CANCER GENETICS SERVICE [C]
- 32 VA-CCGS GERMLINE TESTING DECLINED
- 33 VA-ONC PRCA GERM CGP PT DECLINE: ON PATH
- 34 VA-ONC PRCA GERM CGP ORDER: OFF PATH
- 35 VA-ONC PRCA PRETEST CONSENT: OFF PATH
- 36 VA-ONC PRCA GERMLINE TEST CONSENT DISCUSSION IN FUTURE
- 37 VA-CCGS COMMON HEREDITARY PANEL ORDERED
- 38 VA-CCGS GERMLINE EXPEDITED TESTING
- 39 VA-CCGS GERMLINE PRETEST COUNSELING CONSULT OTHER
- 40 VA-CCGS GERMLINE PRETEST COUNSELING CONSULT CCGS
- 41 VA-CCGS GERMLINE PRETEST CONSENT DISCUSSION
- 42 VA-ONC PRCA PRETEST CONSENT: ON PATH
- 43 VA-CCGS GERMLINE PRETEST DOCUMENTATION
- 44 VA-ONC PRCA RADICAL PROSTATECTOMY PATHWAY [C]
- 45 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS >=5Y EBRT
- 46 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS <5Y OFF PATH RSN
- 47 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS <5Y OFF PATH TX
- 48 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS <5Y ADT
- 49 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS <5Y OBS
- 50 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS <5Y CLIN TRIAL
- 51 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS >=5Y OFF PATH RSN
- 52 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS >=5Y OFF PATH TX
- 53 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS >=5Y OBS
- 54 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS >=5Y EBRT+ADT
- 55 VA-ONC PRCA RAD PROSTATECTOMY NEG DISTANT METS >=5Y CLIN TRIAL
- 56 VA-ONC PRCA RADIATION THERAPY RECURRENCE PATHWAY [C]
- 57 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER SBRT
- 58 VA-ONC PRCA RAD TX RECUR S/P BIOPSY – CANCER OFF PATH RSN
- 59 VA-ONC PRCA RAD TX RECUR S/P BIOPSY – CANCER OFF PATH TX
- 60 VA-ONC PRCA RAD TX RECUR S/P BIOPSY – CANCER ADT
- 61 VA-ONC PRCA RAD TX RECUR S/P BIOPSY – CANCER OBSERVATION

62 VA-ONC PRCA RAD TX RECUR S/P BIOPSY – CANCER CLIN TRIAL
63 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER OFF PATH RSN
64 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER OFF PATH TX
65 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER HIFU
66 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER ADT
67 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER CRYOTHERAPY
68 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER BRACHYOTHERAPY
69 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER RAD PROSTATE
70 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER OBSERVATION
71 VA-ONC PRCA RAD TX RECUR S/P BIOPSY + CANCER CLIN TRIAL
72 VA-ONC PRCA RAD TX RECUR LOCAL NO METS NO OFF PATH RSN
73 VA-ONC PRCA RAD TX RECUR LOCAL NO METS NO OFF PATH TX
74 VA-ONC PRCA RAD TX RECUR LOCAL NO METS NO OBSERVATION
75 VA-ONC PRCA RAD TX RECUR LOCAL NO METS NO CLIN TRIAL
76 VA-ONC PRCA RAD TX RECUR LOCAL YES METS NO OFF PATH RSN
77 VA-ONC PRCA RAD TX RECUR LOCAL YES METS NO OFF PATH TX
78 VA-ONC PRCA RAD TX RECUR LOCAL YES METS NO BIOPSY
79 VA-ONC PRCA RAD TX RECUR LOCAL YES METS NO CLIN TRIAL
80 VA-ONC PRCA TISSUE SUFFICIENCY UNKNOWN
81 VA-ONC PRCA TUMOR GENOMIC TEST ORDER FUTURE DATE
82 VA-ONC PRCA TISSUE SUFFICIENT: OFF PATH
83 VA-ONC PRCA REASSESS: ON PATH
84 VA-ONC PRCA LIQUID BIOPSY: ON PATH
85 VA-ONC PRCA REPEAT BIOPSY: ON PATH
86 VA-ONC PRCA TISSUE SUFFICIENT: ON PATH
87 VA-ONC PRCA CGP ORDER: ON PATH
88 VA-ONC PRCA CGP ORDER: OFF PATH
89 VA-ONC PRCA CRPC M1 PATHWAY [C]
90 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE LU177+ADT
91 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA – OFF PATH RSN
92 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA – OFF PATH TX
93 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA – ABIRATERONE
94 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA – CLIN TRIAL
95 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA + OFF PATH RSN
96 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA + OFF PATH TX
97 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA + OFF PATHWAY
98 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA + OLAPARIB/ABIRATERONE
99 VA-ONC PRCA CRPC M1 1ST ECOG <=2 BRCA + CLIN TRIAL
100 VA-ONC PRCA CRPC M1 3RD ECOG <=2 LU177 OFF PATH RSN
101 VA-ONC PRCA CRPC M1 3RD ECOG <=2 LU177 OFF PATH TX
102 VA-ONC PRCA CRPC M1 3RD ECOG <=2 LU177+ADT
103 VA-ONC PRCA CRPC M1 3RD ECOG <=2 LU177 CLIN TRIAL
104 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA OFF PATH RSN
105 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA OFF PATH TX
106 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA LU177 OFF PATH RSN
107 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA LU177 CLIN TRIAL

108 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE OFF PATH RSN
 109 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE OFF PATH TX
 110 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE CLINICAL TRIAL
 111 VA-ONC PRCA CRPC M1 2ND ECOG <=2 OFF PATH RSN
 112 VA-ONC PRCA CRPC M1 1ST NEUROENDOCRINE OFF PATH RSN
 113 VA-ONC PRCA CRPC M1 1ST ECOG >2 OFF PATH RSN
 114 VA-ONC PRCA CRPC M1 1ST ECOG >2 CLINICAL TRIAL
 115 VA-ONC PRCA CRPC M1 1ST ECOG <=2 CONTRA ABI OFF PATH RSN
 116 VA-ONC PRCA CRPC M1 1ST ECOG <=2 CONTRA ABI OFF PATH TX
 117 VA-ONC PRCA CRPC M1 1ST ECOG <=2 CONTRA ABI CLINICAL TRIAL
 118 VA-ONC PRCA CRPC M1 1ST ECOG <=2 OFF PATH RSN
 119 VA-ONC PRCA CRPC M1 1ST NEUROENDOCRINE CLINICAL TRIAL
 120 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA CLINICAL TRIAL
 121 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CLIN TRIAL
 122 VA-ONC PRCA CRPC M1 1ST ECOG <=2 CLINICAL TRIAL
 123 VA-ONC PRCA CRPC M1 1ST NEUROENDOCRINE OFF PATH TX
 124 VA-ONC PRCA CRPC M1 1ST NEUROENDOCRINE PLATINUM
 125 VA-ONC PRCA CRPC M1 1ST ECOG >2 OFF PATH TX
 126 VA-ONC PRCA CRPC M1 1ST ECOG >2 REFER PALLIATIVE CARE
 127 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA LU177 OFF PATH TX
 128 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA OLAPARIB+ADT
 129 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA RADIUM 223+ADT
 130 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA ABI+PRED+ADT
 131 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA ENZALUT+ADT
 132 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA CABA PEMBRO+ADT
 133 VA-ONC PRCA CRPC M1 3RD ECOG <=2 CONTRA LU177 CABA+ADT
 134 VA-ONC PRCA CRPC M1 2ND ECOG <=2 OFF PATH TX
 135 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE OLAPARIB+ADT
 136 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE RADIUM 223+ADT
 137 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE CABAZIT+ADT
 138 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE ABI+PRED+ADT
 139 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE ENZALUT+ADT
 140 VA-ONC PRCA CRPC M1 2ND ECOG <=2 CONTRA DOCE PEMBRO+ADT
 141 VA-ONC PRCA CRPC M1 2ND ECOG <=2 DOCETAXEL+ADT
 142 VA-ONC PRCA CRPC M1 1ST ECOG <=2 OFF PATH TX
 143 VA-ONC PRCA CRPC M1 1ST ECOG <=2 CONTRA ABI DOCE+ADT
 144 VA-ONC PRCA CRPC M1 1ST ECOG <=2 CONTRA ABI ENZALUT+ADT
 145 VA-ONC PRCA CRPC M1 1ST ECOG <=2 ABI+PRED+ADT
 146 VA-ONC PRCA CRPC M0 PATHWAY [C]
 147 VA-ONC PRCA CRPC M0 <=10M APALUTAMIDE+ADT
 148 VA-ONC PRCA CRPC M0 >10M OFF PATH RSN
 149 VA-ONC PRCA CRPC M0 <=10M OFF PATH RSN
 150 VA-ONC PRCA CRPC M0 >10M CLIN TRIAL
 151 VA-ONC PRCA CRPC M0 <=10M CLIN TRIAL
 152 VA-ONC PRCA CRPC M0 >10M OFF PATH TX
 153 VA-ONC PRCA CRPC M0 >10M CONTINUE ADT

154 VA-ONC PRCA CRPC M0 <=10M OFF PATH TX
155 VA-ONC PRCA CRPC M0 <=10M DAROLUTAMIDE+ADT
156 VA-ONC PRCA CSPC M1 PATHWAY [C]
157 VA-ONC PRCA CSPC M1 HV NO TRIPLE YES ABI OFF PATH RSN
158 VA-ONC PRCA CSPC M1 HV NO TRIPLE YES ABI OFF PATH TX
159 VA-ONC PRCA CSPC M1 HV NO TRIPLE YES ABI CLINICAL TRIAL
160 VA-ONC PRCA CSPC M1 HV NO TRIPLE/ABI OFF PATH RSN
161 VA-ONC PRCA CSPC M1 HV NO TRIPLE/ABI OFF PATH TX
162 VA-ONC PRCA CSPC M1 HV NO TRIPLE/ABI CLINICAL TRIAL
163 VA-ONC PRCA CSPC M1 LV NO ABI OFF PATH RSN
164 VA-ONC PRCA CSPC M1 LV CAND ABI OFF PATH RSN
165 VA-ONC PRCA CSPC M1 HV TRIPLE OFF PATH RSN
166 VA-ONC PRCA CSPC M1 HV TRIPLE CLINICAL TRIAL
167 VA-ONC PRCA CSPC M1 LV CAND ABI CLINICAL TRIAL
168 VA-ONC PRCA CSPC M1 LV NO ABI CLINICAL TRIAL
169 VA-ONC PRCA CSPC M1 HV TRIPLE OFF PATH TX
170 VA-ONC PRCA CSPC M1 HV NO TRIPLE YES ABI ADT+ABI+PRED
171 VA-ONC PRCA CSPC M1 HV NO TRIPLE/ABI ADT+ENZALUTAMIDE
172 VA-ONC PRCA CSPC M1 HV TRIPLE DAROLUTAMIDE+DOCE+ADT
173 VA-ONC PRCA CSPC M1 LV NO ABI OFF PATH TX
174 VA-ONC PRCA CSPC M1 LV NO ABI RAD ONC+ADT+ENZALUTAMIDE
175 VA-ONC PRCA CSPC M1 LV CAND ABI RAD ONC+ADT+ABI+PRED
176 VA-ONC PRCA CSPC M1 LV CAND ABI OFF PATH TX
177 VA-ONC PRCA REGIONAL RISK PATHWAY [C]
178 VA-ONC PRCA REGION LE>=5Y ADT
179 VA-ONC PRCA REGION LE>=5Y RAD TX+ADT
180 VA-ONC PRCA REGION LE<5Y OFF PATH RSN
181 VA-ONC PRCA REGION LE<5Y OFF PATH TX
182 VA-ONC PRCA REGION LE<5Y ADT
183 VA-ONC PRCA REGION LE<5Y OBSERVATION
184 VA-ONC PRCA REGION LE<5Y CLIN TRIAL
185 VA-ONC PRCA REGION LE>=5Y OFF PATH RSN
186 VA-ONC PRCA REGION LE>=5Y ADT+ABI
187 VA-ONC PRCA REGION LE>=5Y RAD TX+ADT+ABI
188 VA-ONC PRCA REGION LE>=5Y OFF PATH TX
189 VA-ONC PRCA REGION LE>=5Y CLIN TRIAL
190 VA-ONC PRCA HIGH/VERY HIGH RISK PATHWAY [C]
191 VA-ONC PRCA HI/VHI LE >=5Y RAD+ADT
192 VA-ONC PRCA HI/VHI NODE METS OFF PATH RSN
193 VA-ONC PRCA HI/VHI NODE METS OFF PATH TX
194 VA-ONC PRCA HI/VHI NODE METS ADT+/-EBRT
195 VA-ONC PRCA HI/VHI NODE METS CLINICAL TRIAL
196 VA-ONC PRCA HI/VHI T3/PSA/+MAR OFF PATH RSN
197 VA-ONC PRCA HI/VHI T3/PSA/+MAR OFF PATH TX
198 VA-ONC PRCA HI/VHI T3/PSA/+MAR REFER RAD ONC
199 VA-ONC PRCA HI/VHI T3/PSA/+MAR CLINICAL TRIAL

200 VA-ONC PRCA HI/VHI LE <5Y OFF PATH RSN
 201 VA-ONC PRCA HI/VHI LE <5Y OFF PATH TX
 202 VA-ONC PRCA HI/VHI LE <5Y ADT
 203 VA-ONC PRCA HI/VHI LE <5Y RADIATION TX
 204 VA-ONC PRCA HI/VHI LE <5Y OBSERVATION
 205 VA-ONC PRCA HI/VHI LE <5Y CLINICAL TRIAL
 206 VA-ONC PRCA HI/VHI LE >=5Y OFF PATH RSN
 207 VA-ONC PRCA HI/VHI LE >=5Y OFF PATH TX
 208 VA-ONC PRCA HI/VHI LE >=5Y PROSTECTOMY+NODES
 209 VA-ONC PRCA HI/VHI LE >=5Y RAD+ADT+ABI
 210 VA-ONC PRCA HI/VHI LE >=5Y CLINICAL TRIAL
 211 VA-ONC PRCA UNFAVORABLE INTERMEDIATE RISK PATHWAY [C]
 212 VA-ONC PRCA UNFAVOR INTER LE >=10Y AFTER SURG T3/PSA/+MAR/METS
 213 VA-ONC PRCA UNFAVOR INTER LE >=10Y OFF PATH RSN
 214 VA-ONC PRCA UNFAVOR INTER LE >=10Y OFF PATH TX
 215 VA-ONC PRCA UNFAVOR INTER LE <10Y OFF PATH RSN
 216 VA-ONC PRCA UNFAVOR INTER LE <10Y OFF PATH TX
 217 VA-ONC PRCA UNFAVOR INTER LE <10Y RAD TX + ADT
 218 VA-ONC PRCA UNFAVOR INTER LE <10Y OBSERVATION
 219 VA-ONC PRCA UNFAVOR INTER LE <10Y CLINICAL TRIAL
 220 VA-ONC PRCA UNFAVOR INTER LE >=10Y RAD PROSTATECTOMY
 221 VA-ONC PRCA UNFAVOR INTER LE >=10Y RAD TX+ADT
 222 VA-ONC PRCA UNFAVOR INTER LE >=10Y CLINICAL TRIAL
 223 VA-ONC PRCA FAVORABLE INTERMEDIATE RISK PATHWAY [C]
 224 VA-ONC PRCA FI LE >=10Y ACTIVE SURV MRI Q 1-2Y
 225 VA-ONC PRCA FI LE >=10Y ACTIVE SURV PROG MRI Q 1-2Y
 226 VA-ONC PRCA FI LE >=10Y ACTIVE SURV PROG BIOPSY Q 1-2Y
 227 VA-ONC PRCA FI LE >=10Y ACTIVE SURV REPEAT BIOPSY Q 1-2Y
 228 VA-ONC PRCA FI LE >=10Y ACTIVE SURV PSA Q 6-12M
 229 VA-ONC PRCA FI LE >=10Y ACTIVE SURV DRE Q 6M
 230 VA-ONC PRCA FI LE >=10Y ACTIVE SURV CONCERN PROGRESS
 231 VA-ONC PRCA FI LE >=10Y AFTER SURG T3/PSA/+MAR/METS
 232 VA-ONC PRCA FI LE <10Y OFF PATH RSN
 233 VA-ONC PRCA FI LE <10Y OFF PATH TX
 234 VA-ONC PRCA FI LE <10Y RADIATION THERAPY
 235 VA-ONC PRCA FI LE <10Y OBSERVATION
 236 VA-ONC PRCA FI LE <10Y CLINICAL TRIAL
 237 VA-ONC PRCA FI LE >=10Y OFF PATH RSN
 238 VA-ONC PRCA FI LE >=10Y OFF PATH TX
 239 VA-ONC PRCA FI LE >=10Y RADICAL PROSTATECTOMY
 240 VA-ONC PRCA FI LE >=10Y RADIATION THERAPY
 241 VA-ONC PRCA FI LE >=10Y ACTIVE SURVEILLANCE
 242 VA-ONC PRCA FI LE >=10Y CLINICAL TRIAL
 243 VA-ONC PRCA LOW RISK PATHWAY [C]
 244 VA-ONC PRCA LR LE >=10Y ACTIVE SURV PROG MRI
 245 VA-ONC PRCA LR LE >=10Y ACTIVE SURV PROG BIOPSY

246 VA-ONC PRCA LR LE >=10Y ACTIVE SURV MRI Q 1-2Y
247 VA-ONC PRCA LR LE >=10Y ACTIVE SURV REPEAT BIOPSY Q 1-2Y
248 VA-ONC PRCA LR LE >=10Y ACTIVE SURV PSA Q 6-12M
249 VA-ONC PRCA LR LE >=10Y ACTIVE SURV DRE Q 6M
250 VA-ONC PRCA LR LE >=10Y ACTIVE SURV CONCERN PROGRESS
251 VA-ONC PRCA LR LE <10Y OFF PATHWAY RSN
252 VA-ONC PRCA LR LE <10Y OFF PATHWAY TX
253 VA-ONC PRCA LR LE <10Y OBSERVATION
254 VA-ONC PRCA LR LE <10Y CLINICAL TRIAL
255 VA-ONC PRCA LR LE >=10Y OFF PATHWAY RSN
256 VA-ONC PRCA LR LE >=10Y OFF PATHWAY TX
257 VA-ONC PRCA LR LE >=10Y RADICAL PROSTATECTOMY
258 VA-ONC PRCA LR LE >=10Y RADIATION THERAPY
259 VA-ONC PRCA LR LE >=10Y ACTIVE SURVEILLANCE
260 VA-ONC PRCA LR LE >=10Y CLINICAL TRIAL
261 VA-ONC PRCA VERY LOW RISK PATHWAY [C]
262 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURV PROG MRI
263 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURV PROG BIOPSY
264 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURV MRI Q 1-2Y
265 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURV REPEAT BIOPSY Q 1-2Y
266 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURV PSA Q 6-12M
267 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURV DRE Q 6M
268 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURV CONCERN PROGRESS
269 VA-ONC PRCA VLR LE >20Y ACTIVE SURV CONCERN PROGRESS
270 VA-ONC PRCA VLR LE >20Y ACTIVE SURV PROG MRI
271 VA-ONC PRCA VLR LE >20Y ACTIVE SURV PROG BIOPSY
272 VA-ONC PRCA VLR LE >20Y ACTIVE SURV MRI Q 1-2Y
273 VA-ONC PRCA VLR LE >20Y ACTIVE SURV REPEAT BIOPSY Q 1-2Y
274 VA-ONC PRCA VLR LE >20Y ACTIVE SURV PSA Q 6-12M
275 VA-ONC PRCA VLR LE >20Y ACTIVE SURV DRE Q 6M
276 VA-ONC PRCA VLR LE <10Y OFF PATHWAY RSN
277 VA-ONC PRCA VLR LE <10Y OFF PATHWAY TX
278 VA-ONC PRCA VLR LE <10Y OBSERVATION
279 VA-ONC PRCA VLR LE <10Y CLINICAL TRIAL
280 VA-ONC PRCA VLR LE 10-20Y OFF PATHWAY RSN
281 VA-ONC PRCA VLR LE 10-20Y OFF PATHWAY TX
282 VA-ONC PRCA VLR LE 10-20Y ACTIVE SURVEILLANCE
283 VA-ONC PRCA VLR LE 10-20Y CLINICAL TRIAL
284 VA-ONC PRCA VLR LE >20Y OFF PATHWAY RSN
285 VA-ONC PRCA VLR LE >20Y OFF PATHWAY TX
286 VA-ONC PRCA VLR LE >20Y RADICAL PROSTATECTOMY
287 VA-ONC PRCA VLR LE >20Y RADIATION THERAPY
288 VA-ONC PRCA VLR LE >20Y ACTIVE SURVEILLANCE
289 VA-ONC PRCA VLR LE >20Y CLINICAL TRIAL
290 VA-ONCOLOGY PATHWAYS [C]
291 VA-ONC PRCA HIGH/VERY HIGH RISK PATHWAY

292 VA-ONC PRCA VERY LOW RISK PATHWAY
293 VA-ONC PRCA LOW RISK PATHWAY
294 VA-ONC PRCA FAVORABLE INTERMEDIATE RISK PATHWAY
295 VA-ONC PRCA UNFAVORABLE INTERMEDIATE RISK PATHWAY
296 PROSTATE CANCER (ONC PCA) DIAGNOSES [C]
297 VA-ONC PCA METASTASIS BONES
298 VA-ONC PCA SURGICAL MARGINS NOT ASSESSED
299 VA-ONC PCA HISTOLOGY INTRADUCTAL
300 VA-ONC PCA METASTASIS LYMPH NODES DISTANT
301 VA-ONC PCA METASTASIS LYMPH NODES REGIONAL
302 VA-ONC PROSTATE PATHOLOGY DATE
303 VA-ONC PCA SURGICAL MARGINS LIMITED
304 VA-ONC PCA SURGICAL MARGINS NEGATIVE
305 VA-ONC PCA SURGICAL MARGINS POSITIVE
306 VA-ONC PCA LOCATION OF SYSTEMATIC POSITIVE CORE SAMPLES
307 VA-ONC PCA NUMBER OF SYSTEMATIC POSITIVE SAMPLES
308 VA-ONC PCA NUMBER OF SYSTEMATIC CORE SAMPLES
309 VA-ONC PCA LOCATION OF TARGETED POSITIVE CORE SAMPLES
310 VA-ONC PCA NUMBER OF TARGETED POSITIVE SAMPLES
311 VA-ONC PCA NUMBER OF TARGETED CORE SAMPLES
312 VA-ONC PCA BIOPSY SYSTEMATIC CORE
313 VA-ONC PCA BIOPSY TARGETED CORE
314 VA-ONC PCA NUMBER OF TRANSPERINEAL POSITIVE SAMPLES
315 VA-ONC PCA NUMBER OF TRANSRECTAL POSITIVE SAMPLES
316 VA-ONC PCA NUMBER OF TRANSPERINEAL CORE SAMPLES
317 VA-ONC PCA NUMBER OF TRANSRECTAL CORE SAMPLES
318 VA-ONC PCA INITIAL PSA
319 VA-ONC PCA INITIAL TESTOSTERONE
320 VA-ONC PCA PATHOLOGICAL PRIMARY TUMOR-T2C
321 VA-ONC PCA PATHOLOGICAL PRIMARY TUMOR-T2B
322 VA-ONC PCA PATHOLOGICAL PRIMARY TUMOR-T2A
323 VA-ONC PCA PATHOLOGICAL PRIMARY TUMOR-T1C
324 VA-ONC PCA PATHOLOGICAL PRIMARY TUMOR-T1B
325 VA-ONC PCA PATHOLOGICAL PRIMARY TUMOR-T1A
326 VA-ONC PCA PATHOLOGICAL PRIMARY TUMOR-T1
327 VA-ONC PCA GRADE GROUP 5
328 VA-ONC PCA GRADE GROUP 4
329 VA-ONC PCA GRADE GROUP 3
330 VA-ONC PCA GRADE GROUP 2
331 VA-ONC PCA GRADE GROUP 1
332 VA-ONC PCA METASTASIS OTHER
333 VA-ONC PCA METASTASIS BRAIN
334 VA-ONC PCA METASTASIS ADRENAL GLAND
335 VA-ONC PCA METASTASIS LIVER
336 VA-ONC PCA METASTASIS LUNGS
337 VA-ONC PCA METASTASIS LYMPH NODES

338 VA-PROSTATE CA RISK: UNFAV INTERMED
339 VA-PROSTATE CA RISK: METASTATIC
340 VA-PROSTATE CA RISK: REGIONAL
341 VA-PROSTATE CA RISK: VERY HIGH
342 VA-PROSTATE CA RISK: HIGH
343 VA-PROSTATE CA RISK: FAV INTERMED
344 VA-PROSTATE CA RISK: LOW
345 VA-PROSTATE CA RISK: VERY LOW
346 VA-PROSTATE CA STAGE
347 VA-PROSTATE CA GLEASON SCORE
348 ONC PCA PATHOLOGICAL-SUMMARY STAGE IIIA
349 ONC PCA CLINICAL-SUMMARY STAGE IIIA
350 ONC PCA PATHOLOGICAL-SUMMARY STAGE IVB
351 ONC PCA PATHOLOGICAL-SUMMARY STAGE IVA
352 ONC PCA PATHOLOGICAL-SUMMARY STAGE IIIC
353 ONC PCA PATHOLOGICAL-SUMMARY STAGE IIIB
354 ONC PCA PATHOLOGICAL-SUMMARY STAGE IIC
355 ONC PCA CLINICAL-SUMMARY STAGE IVB
356 ONC PCA CLINICAL-SUMMARY STAGE IVA
357 ONC PCA CLINICAL-SUMMARY STAGE IIIC
358 ONC PCA CLINICAL-SUMMARY STAGE IIIB
359 ONC PCA CLINICAL-SUMMARY STAGE IIC
360 ONC PCA PATHOLOGICAL PRIMARY TUMOR-T3
361 ONC PCA PATHOLOGICAL PRIMARY TUMOR-T2
362 ONC PCA CLINICAL PRIMARY TUMOR-T3
363 ONC PCA CLINICAL PRIMARY TUMOR-T2
364 ONC PCA CLINICAL PRIMARY TUMOR-T1
365 ONC PCA PRECISION ONCOLOGY
366 ONC PCA TISSUE SAMPLE SENT TO POP
367 ONC PCA TISSUE SAMPLE SENT
368 ONC PCA PATHOLOGICAL PRIMARY TUMOR-T0
369 ONC PCA PATHOLOGICAL PRIMARY TUMOR-TX
370 ONC PCA PATHOLOGICAL STAGE DATE
371 ONC PCA CLINICAL METASTASIS-M1C
372 ONC PCA CLINICAL METASTASIS-M1B
373 ONC PCA CLINICAL METASTASIS-M1A
374 ONC PCA CLINICAL METASTASIS-M0
375 ONC PCA CLINICAL LYMPH NODES-N1
376 ONC PCA CLINICAL LYMPH NODES-N0
377 ONC PCA CLINICAL LYMPH NODES-NX
378 ONC PCA PATHOLOGICAL METASTASIS-M1C
379 ONC PCA PATHOLOGICAL METASTASIS-M1B
380 ONC PCA PATHOLOGICAL METASTASIS-M1A
381 ONC PCA PATHOLOGICAL METASTASIS-M0
382 ONC PCA PATHOLOGICAL LYMPH NODES-N1
383 ONC PCA PATHOLOGICAL LYMPH NODES-N0

384 ONC PCA PATHOLOGICAL LYMPH NODES-NX
385 ONC PCA PATHOLOGICAL PRIMARY TUMOR-T4
386 ONC PCA PATHOLOGICAL PRIMARY TUMOR-T3B
387 ONC PCA PATHOLOGICAL PRIMARY TUMOR-T3A
388 ONC PCA CLINICAL PRIMARY TUMOR-T4
389 ONC PCA CLINICAL PRIMARY TUMOR-T3B
390 ONC PCA CLINICAL PRIMARY TUMOR-T3A
391 ONC PCA CLINICAL PRIMARY TUMOR-T2C
392 ONC PCA CLINICAL PRIMARY TUMOR-T2B
393 ONC PCA CLINICAL PRIMARY TUMOR-T2A
394 ONC PCA CLINICAL PRIMARY TUMOR-T1C
395 ONC PCA CLINICAL PRIMARY TUMOR-T1B
396 ONC PCA CLINICAL PRIMARY TUMOR-T1A
397 ONC PCA CLINICAL PRIMARY TUMOR-T0
398 ONC PCA CLINICAL PRIMARY TUMOR-TX
399 ONC PCA CLINICAL STAGE DATE-
400 ONC PCA >12 POSITIVE-TRANSPERINEAL
401 ONC PCA >12 POSITIVE-TRANSRECTAL
402 ONC PCA >12 CORE-TRANSRECTAL
403 ONC PCA >12 CORE-TRANSPERINEAL
404 ONC PCA ADDITIONAL TESTOSTERONE
405 ONC PCA CORE SAMPLES-TRANSPERINEAL
406 ONC PCA POSITIVE SAMPLES-TRANSPERINEAL
407 ONC PCA BIOPSY-TRANSRECTAL
408 ONC PCA BIOPSY-TRANSPERINEAL
409 ONC PCA LAB TEST-OTHER:
410 ONC PCA ECOG-1
411 ONC PCA ECOG-0
412 ONC PCA ECOG-2
413 ONC PCA ECOG-3
414 ONC PCA ECOG-4
415 ONC PCA ECOG-5
416 ONC PCA 1 CORE-TRANSPERINEAL
417 ONC PCA 2 CORE-TRANSPERINEAL
418 ONC PCA 3 CORE-TRANSPERINEAL
419 ONC PCA 4 CORE-TRANSPERINEAL
420 ONC PCA 5 CORE-TRANSPERINEAL
421 ONC PCA 6 CORE-TRANSPERINEAL
422 ONC PCA 7 CORE-TRANSPERINEAL
423 ONC PCA 8 CORE-TRANSPERINEAL
424 ONC PCA 9 CORE-TRANSPERINEAL
425 ONC PCA 10 CORE-TRANSPERINEAL
426 ONC PCA 11 CORE-TRANSPERINEAL
427 ONC PCA 12 CORE-TRANSPERINEAL
428 ONC PCA 1 POS SAMPLE-TRANSPERINEAL
429 ONC PCA 2 POS SAMPLE-TRANSPERINEAL

430 ONC PCA 3 POS SAMPLE-TRANSPERINEAL
431 ONC PCA 4 POS SAMPLE-TRANSPERINEAL
432 ONC PCA 5 POS SAMPLE-TRANSPERINEAL
433 ONC PCA 6 POS SAMPLE-TRANSPERINEAL
434 ONC PCA 7 POS SAMPLE-TRANSPERINEAL
435 ONC PCA 8 POS SAMPLE-TRANSPERINEAL
436 ONC PCA 9 POS SAMPLE-TRANSPERINEAL
437 ONC PCA 10 POS SAMPLE-TRANSPERINEAL
438 ONC PCA 11 POS SAMPLE-TRANSPERINEAL
439 ONC PCA 12 POS SAMPLE-TRANSPERINEAL
440 ONC PCA ADDITIONAL PSA
441 ONC PCA 1 CORE-TRANSRECTAL
442 ONC PCA 2 CORE-TRANSRECTAL
443 ONC PCA 3 CORE-TRANSRECTAL
444 ONC PCA 4 CORE-TRANSRECTAL
445 ONC PCA 5 POS SAMPLE-TRANSRECTAL
446 ONC PCA 6 CORE-TRANSRECTAL
447 ONC PCA 7 CORE-TRANSRECTAL
448 ONC PCA 8 CORE-TRANSRECTAL
449 ONC PCA 9 CORE-TRANSRECTAL
450 ONC PCA 10 CORE-TRANSRECTAL
451 ONC PCA 11 CORE-TRANSRECTAL
452 ONC PCA 12 CORE-TRANSRECTAL
453 ONC PCA 1 POS SAMPLE-TRANSRECTAL
454 ONC PCA 2 POS SAMPLE-TRANSRECTAL
455 ONC PCA 3 POS SAMPLE-TRANSRECTAL
456 ONC PCA 4 POS SAMPLE-TRANSRECTAL
457 ONC PCA 5 CORE-TRANSRECTAL
458 ONC PCA 6 POS SAMPLE-TRANSRECTAL
459 ONC PCA 7 POS SAMPLE-TRANSRECTAL
460 ONC PCA 8 POS SAMPLE-TRANSRECTAL
461 ONC PCA 9 POS SAMPLE-TRANSRECTAL
462 ONC PCA 10 POS SAMPLE-TRANSRECTAL
463 ONC PCA 11 POS SAMPLE-TRANSRECTAL
464 ONC PCA 12 POS SAMPLE-TRANSRECTAL
465 ONC PCA DIAGNOSIS DATE-
466 ONC PCA CLINICAL-SUMMARY STAGE I
467 ONC PCA CLINICAL-SUMMARY STAGE IIA
468 ONC PCA CLINICAL-SUMMARY STAGE IIB
469 ONC PCA CLINICAL-SUMMARY STAGE III
470 ONC PCA CLINICAL-SUMMARY STAGE IV
471 ONC PCA PATHOLOGICAL-SUMMARY STAGE I
472 ONC PCA PATHOLOGICAL-SUMMARY STAGE IIA
473 ONC PCA PATHOLOGICAL-SUMMARY STAGE IIB
474 ONC PCA PATHOLOGICAL-SUMMARY STAGE III
475 ONC PCA PATHOLOGICAL-SUMMARY STAGE IV

476 ONC PCA HISTOLOGY-ADENOCARCINOMA
477 ONC PCA HISTOLOGY-SMALL CELL
478 ONC PCA HISTOLOGY-OTHER:
479 ONC PCA GLEASON SCORE ONE-1
480 ONC PCA GLEASON SCORE ONE-2
481 ONC PCA GLEASON SCORE ONE-3
482 ONC PCA GLEASON SCORE ONE-4
483 ONC PCA GLEASON SCORE ONE-5
484 ONC PCA GLEASON SCORE TWO-1
485 ONC PCA GLEASON SCORE TWO-2
486 ONC PCA GLEASON SCORE TWO-3
487 ONC PCA GLEASON SCORE TWO-4
488 ONC PCA GLEASON SCORE TWO-5
489 ONC PCA GLEASON SCORE COMBINED-2
490 ONC PCA GLEASON SCORE COMBINED-3
491 ONC PCA GLEASON SCORE COMBINED-4
492 ONC PCA GLEASON SCORE COMBINED-5
493 ONC PCA GLEASON SCORE COMBINED-6
494 ONC PCA GLEASON SCORE COMBINED-7
495 ONC PCA GLEASON SCORE COMBINED-8
496 ONC PCA GLEASON SCORE COMBINED-9
497 ONC PCA GLEASON SCORE COMBINED-10
498 ONC PCA RISK STATUS-LOW
499 ONC PCA RISK STATUS-INTERMEDIATE
500 ONC PCA RISK STATUS-HIGH
501 ONC PCA ECOG-UNKNOWN
502 PROSTATE CANCER (ONC PCA) TREATMENTS [C]
503 VA-ONC PCA SURGICAL RESECTION OTHER
504 VA-ONC PCA SURGICAL RESECTION TURP
505 VA-ONC PCA SURGICAL RESECTION CYSTOPROSTATECTOMY
506 VA-ONC PCA SURGICAL RESECTION OPEN
507 VA-ONC PCA SURGICAL RESECTION ROBOTIC
508 VA-ONC GOALS OF CARE [C]
509 VA-ONC GOALS OF CARE PALLIATIVE
510 VA-ONC GOALS OF CARE CURATIVE
511 VA-ONC PAIN MANAGEMENT PLAN [C]
512 VA-ONC PAIN PLAN INITIAL
513 VA-ONC PAIN PLAN MGT BY HEMONC
514 VA-ONC PAIN PLAN ADEQUATE
515 VA-ONC PAIN PLAN NEEDS ADJUSTMENT
516 VA-ONC PAIN MGT PLAN OPIOID PRESCRIBED
517 VA-ONC PAIN MGT PLAN BOWEL REGIMEN
518 VA-ONC PAIN MGT PLAN NON OPIOID PLAN
519 VA-ONC PAIN PLAN MANAGED BY
520 VA-ONC PAIN PLAN DECLINED
521 VA-ONC PAIN DENIES

522 VA-ONC CN TIME SPENT CHART REVIEW
523 VA-ONC CN TIME SPENT FACE TO FACE
524 VA-ONC CN TIME SPENT VVC CHART REVIEW
525 VA-ONC CN TIME SPENT VVC FACE TO FACE
526 VA-ONC CN TIME SPENT CVT FACE TO FACE
527 VA-ONC CN TIME SPENT CVT CHART REVIEW
528 VA-ONC PCA NOTE TITLES [C]
529 VA-ONC PCA CANCER LOCALIZED NOTE
530 VA-ONC PRCA RADIATION THERAPY RECURRENCE PATHWAY
531 VA-ONC PRCA RADICAL PROSTATECTOMY PATHWAY
532 VA-ONC PCA CANCER RECURRENCE/PROSTATECTOMY NOTE
533 VA-ONC PRCA MOLECULAR TESTING
534 VA-ONC PRCA CSPC M1 PATHWAY
535 VA-ONC PRCA REGIONAL RISK PATHWAY
536 VA-ONC PRCA CRPC M1 PATHWAY
537 VA-ONC PRCA CRPC M0 PATHWAY
538 VA-ONC PCA CANCER REGIONAL/METASTATIC NOTE

REMINDER SPONSOR

539 VHA ONCOLOGY
540 OFFICE OF PATIENT CARE SERVICES
541 Office of Nursing Service

REMINDER COMPUTED FINDINGS

VA-ASU USER CLASS

REMINDER TERM

542 VA-REMINDER UPDATE_2_0_380
543 VA-ONC CN BRANCH PROSTATE
544 VA-ONC PROSTATE METASTASIS RT
545 VA-ONC PCA PATHOLOGY RT
546 VA-ONC PCA ECOG RT
547 VA-ONC PCA STAGING RT
548 VA-ONC PCA DX DATE RT
549 VA-ONC GOALS OF CARE RT
550 VA-ONC PAIN EFFECTIVENESS RT
551 VA-ONC ASU ATTENDING

HEALTH SUMMARY COMPONENT

PCE HEALTH FACTORS SELECTED
LAB SURGICAL PATHOLOGY
LAB TESTS SELECTED

HEALTH SUMMARY TYPE

552 VA-ONCOLOGY PRCA PATHWAYS (HS)
553 VA-ONC CN PREVIOUS PROSTATE

554 VA-ONC CN PRCA HISTORY (HS)
555 VA-ONC PROSTATE METASTASIS
556 VA-ONC PCA SURG PATH REPORT
557 VA-ONC PCA PATHOLOGY
558 VA-ONC PCA TESTOSTERONE
559 VA-ONC PCA LAST PSA
560 VA-ONC PCA ECOG
561 VA-ONC PCA STAGING
562 VA-ONC PCA DX DATE
563 VA-ONC GOALS OF CARE HS
564 VA-ONC PCA PAIN

HEALTH SUMMARY OBJECTS

565 VA-ONCOLOGY PRCA PATHWAYS (TIU)
566 VA-ONC CN PREVIOUS PROSTATE (TIU)
567 VA-ONC CN PRCA HISTORY (TIU)
568 VA-ONC PROSTATE METASTASIS (TIU)
569 VA-ONC PCA SURG PATH REPORT
570 VA-ONC PCA PATHOLOGY
571 VA-ONC PCA TESTOSTERONE
572 VA-ONC PCA LAST PSA
573 VA-ONC PCA ECOG
574 VA-ONC PCA STAGING
575 VA-ONC PCA TS DX DATE
576 VA-ONC GOALS OF CARE HS (TIU)
577 VA-ONC PCA PAIN (TIU)

TIU DOCUMENT DEFINITION

578 VA-ONCOLOGY PRCA PATHWAYS
579 VA-ONC CN PREVIOUS PROSTATE
580 VA-ONC CN PRCA HISTORY
581 VA-ONC PROSTATE METASTASIS
582 VA-ONC PCA SURG PATH REPORT
583 VA-ONC PCA PATHOLOGY
584 VA-ONC PCA TESTOSTERONE
585 VA-ONC PCA LAST PSA
586 VA-ONC PCA ECOG
587 VA-ONC PCA STAGING
588 VA-ONC PCA DX DATE
589 VA-ONC GOALS OF CARE HS
590 VA-ONC PCA PAIN

ACTIVE MEDICATIONS
ALLERGIES/ADR
PATIENT SEX
PATIENT RACE
PATIENT AGE

PATIENT NAME

REMINDER DIALOG

- 591 VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP
- 592 VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP
- 593 VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP

Install Details

This update is being distributed as a web host file. The address for the host file is:
https://REDACTED/UPDATE_2_0_380.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====


This update can be loaded with users on the system. Installation will take less than 15 minutes.
Post install will take about an hour.

Install Example

Go to Reminder exchange in Vista.

To Load the exchange file, use LWH. The URL is:

https://REDACTED/UPDATE_2_0_380.PRD



```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History              RP   Repack
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file
https://REDACTED/UPDATE_2_0_380.PRD
```

You should receive a message that the file was successfully loaded into your reminder exchange file.

Search and locate exchange file name UPDATE_2_0_380 VA-ONC PROSTATE CANCER FOLLOW UP

+Item	Entry	Source	Date Packed
522	UPDATE_2_0_380 VA-ONC PROSTATE CANCER FOLLOW UP		
523	UPDATE_2_0_39 VA-OSP SUICIDE SAFETY PLAN		
524	UPDATE_2_0_3_2 TBI SCREENING		
525	UPDATE_2_0_4 VA-VETERANS CHOICE/MAMMOGRAM/TBI UPDATE		
526	UPDATE_2_0_40 VA-HERPES ZOSTER IMMUNIZATION UPDATE		
527	UPDATE_2_0_40 VA-HERPES		
+ + Next Screen - Prev Screen ?? More Actions >>>			
CFE	Create Exchange File Entry	LHF	Load Host File
CHF	Create Host File	LMM	Load MailMan Message
CMM	Create MailMan Message	LR	List Reminder Definitions
DFE	Delete Exchange File Entry	LWH	Load Web Host File
IFE	Install Exchange File Entry	RI	Reminder Definition Inquiry
IH	Installation History	RP	Repack
Select Action: Next Screen// IFE Install Exchange File Entry			
Enter a list or range of numbers (1-758): 522			

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry
Enter the number that corresponds with your entry titled UPDATE_2_0_380 VA-ONC
PROSTATE CANCER FOLLOW UP (in this example it is entry 522, it will vary by site) that is
dated 11/06/2023.

```

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP
  VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP
  VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP

REMINDER TERM
  VA-REMINDER UPDATE_2_0_380
Non-exchangeable order dialog(s):
-----
                                Name: SD RTC
+    Enter ?? for more actions    >>>
IA   Install all Components        IS   Install Selected Component
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.

For all new entries in any category, you will choose **I to Install**

All items below that are **SKIPS are because you completed Update 2_0_367 and Update _2_0_20 (March 2017 release)** **Do not skip if you have not completed these updates.**

For Template Fields that exist but are different, **OVERWRITE**
EXCEPT: VA-ONC VAMC SITE NAME, **SKIP** this item

For Reminder Terms that exist, but are different, **OVERWRITE**
EXCEPT: VA-ONC ASU ATTENDING, **MERGE** this item

For Health Factors that exist, but are different, **OVERWRITE**

For Health Summary Type that exist, but are different, **OVERWRITE**
EXCEPT: VA-ONC PCA TESTOSTERONE **SKIP** this item
VA-ONC PCA LAST PSA **SKIP** this item.

For Health Summary Objects that exist, but are different, **OVERWRITE**
EXCEPT: VA-ONC PCA TESTOSTERONE **SKIP** this item
VA-ONC PCA LAST PSA **SKIP** this item.

TIU Document Definitions that exist, but are different: **OVERWRITE**
EXCEPT: VA-ONC PCA TESTOSTERONE **SKIP** this item
VA-ONC PCA LAST PSA **SKIP** this item.

You will then be prompted to install the first dialog component.

Packed reminder dialog: VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-U

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP	dialog			
2	20	VA-GP ONC CN PRCA REGIONAL/METASTATIC Finding: VA-ONC PCA CANCER REGIONAL/METASTATIC NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	group	X		
3	20.10	VA-ONC CN PRCA PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X		
4		PXRM GF VIEW BUTTON	prompt	X		
5	20.15	VA-ONC CN PRCA OBJECT Finding: *NONE*	element			
6	20.20	VA-ONC CN PRCA PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA						

Select IA to install all.

Install reminder dialog and all components with no further changes: Y//Y

VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP (reminder dialog) inst

Item	Seq.	Dialog Findings	Type	Exists
1		VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP	dialog	X
2	20	VA-GP ONC CN PRCA REGIONAL/METASTATIC Finding: VA-ONC PCA CANCER REGIONAL/METASTATIC NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	group	X
3	20.10	VA-ONC CN PRCA PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
4		PXRM GF VIEW BUTTON	prompt	X
5	20.15	VA-ONC CN PRCA OBJECT Finding: *NONE*	element	X
6	20.20	VA-ONC CN PRCA PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	element	X
+ + Next Screen - Prev Screen ?? More Actions				
DD		Dialog Details	DT Dialog Text	IS Install Selected
DF		Dialog Findings	DU Dialog Usage	QU Quit
DS		Dialog Summary	IA Install All	
Select Action: Next Screen// Q				

After install, select Q for quit to move to the next dialog.

You will then be prompted to install the second dialog component.

Packed reminder dialog: VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOL

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP	dialog			
2	20	VA-GP ONC CN PRCA RECURRENCE/PROSTATECTOMY Finding: VA-ONC PCA CANCER RECURRENCE/PROSTATECTOMY NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	group	X		
3	20.10	VA-ONC CN PRCA PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X		
4	20.15	VA-ONC CN PRCA OBJECT Finding: *NONE*	element	X		
5	20.20	VA-ONC CN PRCA PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	element	X		
6	20.30	VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS	element	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA						

Select IA to install the dialog

Install reminder dialog and all components with no further changes: Y//Y

VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP (reminder dialog)						
Item	Seq.	Dialog	Findings	Type	Exists	
1		VA-ONCOLOGY PROSTATE CANCER		dialog	X	
		RECURRENCE/PROSTATECTOMY FOLLOW-UP				
2	20	VA-GP ONC CN PRCA RECURRENCE/PROSTATECTOMY		group	X	
		Finding: VA-ONC PCA CANCER RECURRENCE/PROSTATECTOMY NOTE			X	
		(HEALTH FACTOR)				
		Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)			X	
3	20.10	VA-ONC CN PRCA PREVIOUS HF VIEW		element	X	
		Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)			X	
4	20.15	VA-ONC CN PRCA OBJECT		element	X	
		Finding: *NONE*				
5	20.20	VA-ONC CN PRCA PREVIOUS NOTE VIEW		element	X	
		Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)			X	
		Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)			X	
6	20.30	VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS		element	X	
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// Q						

Select Q to quit and move to the next dialog.

You will then be prompted to install the last dialog component.

Packed reminder dialog: VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP [NATIONA

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP	dialog			
2	20	VA-GP ONC CN PRCA LOCALIZED Finding: VA-ONC PCA CANCER LOCALIZED NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	group	X X X		
3	20.10	VA-ONC CN PRCA PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X X		
4	20.15	VA-ONC CN PRCA OBJECT Finding: *NONE*	element	X		
5	20.20	VA-ONC CN PRCA PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	element	X X X		
6	20.30	VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS Finding: *NONE*	element	X		
7	20.40	VA-GP ONC CN PRCA REQUIRED	group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA						

Select IA to install the dialog

Install reminder dialog and all components with no further changes: Y//Y

VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP (reminder dialog) installed from

Item	Seq.	Dialog Findings	Type	Exists
1		VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP	dialog	X
2	20	VA-GP ONC CN PRCA LOCALIZED Finding: VA-ONC PCA CANCER LOCALIZED NOTE (HEALTH FACTOR) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	group	X X X
3	20.10	VA-ONC CN PRCA PREVIOUS HF VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X X
4	20.15	VA-ONC CN PRCA OBJECT Finding: *NONE*	element	X
5	20.20	VA-ONC CN PRCA PREVIOUS NOTE VIEW Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING) Add. Finding: VA-ONC CN PROSTATE (HEALTH FACTOR)	element	X X X
6	20.30	VA-ONC CN ONCOLOGICAL DOCUMENTATION TOOLS Finding: *NONE*	element	X
7	20.40	VA-GP ONC CN PRCA REQUIRED	group	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details		DT Dialog Text	IS Install Selected
DF	Dialog Findings		DU Dialog Usage	QU Quit
DS	Dialog Summary		IA Install All	
Select Action: Next Screen// Q				

Select Q to quit and you will be returned to the following screen

```

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP
  VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP
  VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP

REMINDER TERM
  VA-REMINDER UPDATE_2_0_380
Non-exchangeable order dialog(s):
-----
                        Name: SD RTC
+      + Next Screen      - Prev Screen      ?? More Actions
IA  Install all Components                      IS  Install Selected Component
Select Action: Next Screen// Q
  
```

Select Q for quit to exit the install

Post Installation

1. Create Note Titles:

The note titles to create to associate with these templates are:

ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP
ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP
ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP

VHA Enterprise Standard Title: HEMATOLOGY AND ONCOLOGY NOTE

You can add these to your document class where progress notes used by Oncology are located. **These titles should not be in your consult document class.**

In addition, you will also need to create 3 consult note titles:

ONCOLOGY PROSTATE CANCER LOCALIZED CONSULT
ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY CONSULT
ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC CONSULT

VHA Enterprise Standard Title: HEMATOLOGY AND ONCOLOGY CONSULT

These note titles should be in your consult document class.

***Note:** Comprehensive information on Creating new Document Classes and Titles can be found beginning on page 54 of the TIU/ASU Implementation Guide at the following link:
[http://www.va.gov/vdl/documents/Clinical/CPRS-Text_Integration_UTILITY_\(TIU\)/tiuim.pdf](http://www.va.gov/vdl/documents/Clinical/CPRS-Text_Integration_UTILITY_(TIU)/tiuim.pdf)*

2. Make the Dialog template available to be attached to a Personal or Shared template or to be able to attach the dialog to a progress note title. (You can use the LINK feature in this menu for steps 2 and 3 if you choose to)

Select Reminder Managers Menu <TEST ACCOUNT> Option:

CP CPRS Reminder Configuration
CA Add/Edit Reminder Categories
CL CPRS Lookup Categories
CS CPRS Cover Sheet Reminder List
MH Mental Health Dialogs Active
PN Progress Note Headers
RA Reminder GUI Resolution Active
TIU TIU Template Reminder Dialog Parameter
DEVL Evaluate Coversheet List on Dialog Finish
DL Default Outside Location
PT Position Reminder Text at Cursor
LINK Link Reminder Dialog to Template
TEST CPRS Coversheet Time Test
NP New Reminder Parameters
GEC GEC Status Check Active
WH WH Print Now Active

Select CPRS Reminder Configuration <TEST ACCOUNT> Option:

TIU TIU Template Reminder Dialog Parameter

Reminder Dialogs allowed as Templates may be set for the following:

- 1 User USR [choose from NEW PERSON]
- 3 Service SRV [choose from SERVICE/SECTION]
- 4 Division DIV [SALT LAKE CITY]
- 5 System SYS [NATREM.FO-SLC.MED.VA.GOV]

Enter selection: **5** Choose the appropriate number for your site. Your site may do this by System or other levels. The example below uses SYSTEM level.

Setting Reminder Dialogs allowed as Templates for System:

NATREM.FO-SLC.MED.VA.GOV

Select Display Sequence:

388 VA-COVID-19 TELEPHONE FOLLOW-UP

389 VA-SRA-CONSULT RESPONSE (D)

When you type a question mark above, you will see the list of #'s (with dialogs) that are already taken. Choose a number NOT on this list. For this example, 391 is not present so I will use 391.

Select Display Sequence: **391**

Are you adding 391 as a new Display Sequence? Yes// **YES**

Display Sequence: 391// <Enter> 391

Clinical Reminder Dialog: **VA- ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP** reminder dialog NATIONAL ...OK? Yes// <Enter> (Yes)

Select Display Sequence: **392**

Are you adding 392 as a new Display Sequence? Yes// **YES**

Display Sequence: 392// <Enter> 392

Clinical Reminder Dialog: **VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP** reminder dialog NATIONAL ...OK? Yes// <Enter> (Yes)

Select Display Sequence: **393**


Are you adding 393 as a new Display Sequence? Yes// **YES**

Display Sequence: 393// <Enter> 393

Clinical Reminder Dialog: **VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP** reminder dialog NATIONAL ...OK? Yes// <Enter> (Yes)

3. **Associate the reminder dialogs with the note title in CPRS.**

In CPRS, open the Template Editor, expand Document Titles, select the appropriate folder (if applicable), then click New Template (top right-hand corner).



a. Type in the following information:


Name: **ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA- ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP**

Next, type in the Associated Title name –

ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP - then click Apply.



b. Type in the following information:


Name: **ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP**

Template Type: **Reminder Dialog**

Reminder Dialog **VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP**

Next, type in the Associated Title name –

ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP - then click Apply.



c. Type in the following information:


Name: **ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP**
 Template Type: **Reminder Dialog**
 Reminder Dialog: **VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP**

Next, type in the Associated Title name – **ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP** - then click Apply.

d. Type in the following information:

Name: **ONCOLOGY PROSTATE CANCER LOCALIZED CONSULT**
 Template Type: **Reminder Dialog**
 Reminder Dialog: **VA-ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP**


Next, type in the Associated Title name – **ONCOLOGY PROSTATE CANCER LOCALIZED CONSULT** – then click Apply.



e. Type in the following information:

Name: **ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY CONSULT**
 Template Type: **Reminder Dialog**
 Reminder Dialog: **VA-ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP**

Next, type in the Associated Title name –
ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY CONSULT –
 then click Apply.



- f. Type in the following information:
 Name: **ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC CONSULT**
 Template Type: **Reminder Dialog**
 Reminder Dialog: **VA-ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP**

Next, type in the Associated Title name –
ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC CONSULT - then click Apply.

4. Add objects to the reminder dialog.
 Using the VistA menu Reminder Dialog Management (DM), Reminder Dialogs (DI), and Change View (CV) to Dialog Elements (E) view.

- a. Locate and edit the reminder dialog element:
 VAL-ONC CN PATIENT DEMOGRAPHICS

Replace the following data objects with your local data objects if you do not have the included objects.

- |PATIENT SEX|
- |PATIENT RACE|
- |PATIENT AGE|
- |PATIENT NAME|

- b. Locate and edit the reminder dialog element:
 VAL-ONC CN VITALS TXT

Replace the text (add local object here) with your local data objects for the following:

- Temperature
- Pulse
- Respiration
- Pulse Oximetry
- Blood Pressure
- Pain
- Height
- Weight

c. Locate and edit the reminder dialog element:

VAL-ONC CN LABS & IMAGING

You can replace the text “(please reference official report in record)” with your local data objects for the most recent LABS & IMAGING. Check with your local Oncologists to determine what labs and/or imaging that they would like in the note.

***Or your site can choose to leave the verbiage “(please reference official report in record)” without inserting data objects.**

d. Locate and edit the reminder dialog element:

VAL-ONC CN VISIT VVC CONTACT HOME ADDRESS

Replace the text: (ADD local patient address object here) **with your local data object.**

e. Locate and edit the reminder dialog element:

VAL-ONC CN VISIT VVC CONTACT PHONE NUMBER HOME/CELL.

Replace the text (ADD patient xxx phone object here) **with your local data objects** for the following:

- Cell phone
- Home phone
- Work phone

f. Locate and edit the reminder dialog element VAL-ONC CN MED REC OBJECT

The reminder dialog element above is being sent out with object [VAL-ONC CN MED REC OBJECT]. If your site has their own object that is being use, you may want to replace this with your local object.

5. Suppress the resident supervision statement (if applicable):

The template utilizes branching logic via the element VAL-ONC CN ATTENDING OF RECORD to suppress documenting resident supervision in the dialog if needed.

*If **MERGED** in install, you do not need to complete this step.

*This is a local element and can be disabled if not applicable to your site.

You should add user classes that **do not** need to document resident supervision (i.e., Attending Physician, etc.). To add user classes, edit the reminder term VA-ONC ASU ATTENDING.

**To add an additional user class, use the quote “characters
“CF.VA-ASU USER CLASS”

and set the computed finding parameter of your new finding to the appropriate user class.

The reminder term is being released with the ASU class of ATTENDING PHYSICIAN and NURSE PRACTITIONER. Edit according to your local practices.

If a user does not have one of the included user classes, the following will appear in the note.

The Attending of record for this patient is Dr. * . This case has been reviewed by the Attending and they agree with the recommended plan.

6. **Enter your site name in the template field: {FLD:VA-ONC VAMC SITE NAME}.**

SKIP adding your site name in the template field if previously done in Clinical Reminder Update_2_0_367.

7. **Add orders as an additional finding to the elements identified in the pre-install (optional):**

VA-ONC RADIATION ONCOLOGY CONSULT – radiation oncology consult.

VA-ONC HOSPICE REFERRAL - hospice consult.

VA-ONC UROLOGY CONSULT–urology consult.

8. **Add consult quick order as an additional finding to the element:**

VA-ONC PRCA GENETIC PRETEST OPTION CONSULT CCGS

-if **NTO site**, map IFC CVT-TELEONCOLOGY TO DURHAM VAMC (if you have this consult, you are an NTO site)

-if **not NTO site**, map IFC PRECISION ONCOLOGY ELECTRONIC CONSULT


9. **Add orders as an additional finding to the element:**

VA-ONC PRCA GENETIC PRETEST OPTION OTHER

-map community care genetic consult, or any other genetic consult your site uses; may need to create order menu if have more than one genetics service (i.e. Community care and Genomic Medicine Service)

**If you don't have genetic counseling consults, then inactivate this element.


Ex:



10. Add values to Additional Prompt:

Navigate to reminder dialogs. Change view to Additional Prompts. Search for prompt:
VA-ONC CCGS GENETIC SERVICE

Add community care and any other genetic counseling that you have at your facility. If you don't have any, then don't add anything.



Ex:

11. Map HS Objects:

**You may already have lab tests in these objects. If you do, and skipped them on install, go to #12.

- | VA-ONC PCA LAST PSA|
- | VA-ONC PCA TESTOSTERONE|.

Start with the TIU Maintenance Menu → Document Definitions (Manager) → Create TIU/Health Summary Objects → Detailed Display/Edit HS Object → Edit HS Type.

```

HS Object: VA-ONC PCA LAST PSA
Object Label: Prostate Specific Antigen (PSA)
Health Summary Type: VA-ONC PCA LAST PSA
Report Period:
Creator:
HS Object
Print Label: YES
Print Report Date and Time: NO
Print Blank Line after Label: NO
Print Confidentiality Banner: NO
Customized Header: YES
Print Report Date and Time: NO
Suppress Components w/o Data: NO
Print Component Header: NO
Print Deceased Information: NO
Print Time-Occurrence Limits: NO
National Object: NO
Underline Component Header: NO
Blank Line After Header: NO
Overwrite No Data: No PSA data available
Enter ?? for more actions
Edit HS Object
Change HS Type
Inquire about a HS Type
Edit HS Type
Select Item(s): Quit//

```

Map your local PSA lab test(s) to the following object | VA-ONC PCA LAST PSA |

Select Item(s): Quit// **Edit HS Type**

Editing Health Summary Type VA-ONC PCA LAST PSA

NAME: VA-ONC PCA LAST PSA Replace
TITLE:
SUPPRESS PRINT OF COMPONENTS WITHOUT DATA: no//
LOCK:
OWNER:
Select COMPONENT: **SLT**
SLT LAB TESTS SELECTED SLT
SUMMARY ORDER: 5// 5
OCCURRENCE LIMIT: **1**
TIME LIMIT:
HEADER NAME: Selected Lab Tests//

No selection items chosen.

Select new items one at a time in the sequence you want them displayed.
You may select any number of items.

Select SELECTION ITEM: **PROSTATE SPECIFIC ANTIGEN**

Continue to add **PSA, Free**; **PSA, % Free**; & **PSA, Total**. Check with your local LIM for the correct lab tests if these tests do not exist in your system. When you are finished, you should have the following:

LAB TESTS SELECTED (SLT)

Type Name: PSA results

Abb	Ord	Component Name	Max Occ	Hos Loc	ICD Text	Pro Nar	CPT Mod	Selection
SLT	5	Lab Tests Selected						PROSTATE SPECIFIC ANTIGEN PSA, FREE [spl] PSA, % FREE [spl] PSA, TOTAL [spl]

```

HS Object: VA-ONC PCA TESTOSTERONE
Object Label: Testosterone
Health Summary Type: VA-ONC PCA TESTOSTERONE
Report Period:
Creator:
HS Object
Print Label: YES
Print Blank Line after Label: NO
Customized Header: YES
Suppress Components w/o Data: NO
Print Deceased Information: NO
National Object: NO
Print Report Date and Time: NO
Print Confidentiality Banner: NO
Print Report Date and Time: NO
Print Component Header: NO
Print Time-Occurrence Limits: NO
Underline Component Header: NO
Blank Line After Header: NO
Overwrite No Data: No TESTOSTERONE data available
Enter ?? for more actions
Edit HS Object
Change HS Type
Select Item(s): Quit//
Inquire about a HS Type
Edit HS Type

```

Map your local Testosterone lab tests to the following object | VA-ONC PCA TESTOSTERONE |

Select Item(s): Quit// **Edit HS Type**

Editing Health Summary Type VA-ONC PCA TESTOSTERONE

NAME: VA-ONC PCA TESTOSTERONE Replace
 TITLE:
 SUPPRESS PRINT OF COMPONENTS WITHOUT DATA: no//
 LOCK:
 OWNER:
 Select COMPONENT: **SLT**
SLT LAB TESTS SELECTED SLT
 SUMMARY ORDER: 5// 5
 OCCURRENCE LIMIT: **1**
 TIME LIMIT:
 HEADER NAME: Selected Lab Tests//

No selection items chosen.
 Select new items one at a time in the sequence you want them displayed.
 You may select any number of items.

Select SELECTION ITEM: **FREE TESTOSTERONE**

Continue to add **TESTOSTERONE, TOTAL**. Check with your local LIM for the correct lab tests if these tests do not exist in your system. When you are finished, you should have the following:

LAB TESTS SELECTED (SLT)

Type Name: Testosterone Results

Abb	Ord	Component Name	Max Occ	Time	Hos Loc	ICD Text	Pro Nar	CPT Mod	Selection
SLT	5	Lab Tests Selected							FREE TESTOSTERONE [QST] TESTOSTERONE, TOTAL

12. Map the six note titles to HS Object:

|VA-ONC CN PREVIOUS PROSTATE|:

NOTE TITLES:

- ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP
- ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP
- ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP
- ONCOLOGY PROSTATE CANCER LOCALIZED CONSULT
- ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY CONSULT
- ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC CONSULT

HS OBJECT:

|VA-ONC CN PREVIOUS PROSTATE|

Example:

Start with the TIU Maintenance Menu → Document Definitions (Manager) → Create TIU/Health Summary Objects → SL “VA-ONC CN PREVIOUS PROSTATE “ → Detailed Display/Edit HS Object.

Map the note title the 6 note titles to the following object.

HS Object: VA-ONC CN PREVIOUS PROSTATE (TIU)
 Health Summary Type: VA-ONC CN PREVIOUS PROSTATE (HS)
 Report Period:
 Creator:

HS Object

Print Label: NO	Print Report Date and Time: NO
Print Blank Line after Label: NO	Print Confidentiality Banner: NO
Customized Header: YES	Print Report Date and Time: NO
Suppress Components w/o Data: YES	Print Component Header: NO
Print Deceased Information: NO	Print Time-Occurrence Limits: NO

National Object: NO
Blank Line After Header: NO
Overwrite No Data: No Teleoncology note found

Underline Component Header: NO

Select Item(s): Quit// **Edit HS Type**

Editing Health Summary Type VA-ONC CN PREVIOUS PROSTATE (HS)'

NAME: VA-ONC CN PREVIOUS PROSTATE (HS) Replace

TITLE:

SUPPRESS PRINT OF COMPONENTS WITHOUT DATA: yes//

LOCK:

OWNER:

Select COMPONENT: **SPN**

SPN PROGRESS NOTES SELECTED SPN

SUMMARY ORDER: 5// 5

OCCURRENCE LIMIT: **1**

TIME LIMIT:

HEADER NAME: Selected Prog Notes//

No selection items chosen.

Select new items one at a time in the sequence you want them displayed.
You may select any number of items.

Select SELECTION ITEM: **ONCOLOGY PROSTATE CANCER
LOCALIZED FOLLOW-UP**

Continue to add all the note titles listed above. When you are finished, you should have the following:

PROGRESS NOTES SELECTED (SPN)

- 1 ONCOLOGY PROSTATE CANCER LOCALIZED FOLLOW-UP
- 2 ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY FOLLOW-UP
- 3 ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC FOLLOW-UP
- 4 ONCOLOGY PROSTATE CANCER LOCALIZED CONSULT
- 5 ONCOLOGY PROSTATE CANCER RECURRENCE/PROSTATECTOMY CONSULT
- 6 ONCOLOGY PROSTATE CANCER REGIONAL/METASTATIC CONSULT

13. ****Optional**** - Add local reminder dialog groups and elements if needed. This template was built with the standard items needed for complete documentation. If a site would like to add guidance or information specific to that site, you can add text to the local reminder dialog group **VAL-GP ONC CN PROSTATE LOCAL GROUP** using the reminder dialog editor.

This group is at the end of the dialog. This dialog group will install disabled, and you will need to remove the disable flag for it to be used.