



**Update\_2\_0\_24**

**Clinical Reminders**

**VA-MISC UPDATE 3**

**INSTALLATION and SETUP GUIDE**

**May 2017**

Product Development  
Department of Veterans Affairs

# Contents

---

<b>Introduction.....</b>	<b>3</b>
<b>Install Details .....</b>	<b>13</b>
<b>Pre-Install .....</b>	<b>21</b>
<b>Install Example .....</b>	<b>21</b>
<b>Post-Install Set-up Instructions.....</b>	<b>38</b>

## Introduction

---

UPDATE\_2\_0\_24 includes multiple reminder dialog updates. It contains Reminder Exchange entry:  
**UPDATE\_2\_0\_24 MISC UPDATE 3**

This update will update the following reminder content:

### **REMINDER DEFINITION**

VA-BL AUDC 1 TO 4 NO AUD DX  
VA-BL ALCOHOL WITHIN SAFE LIMIT  
VA-BL AUD-C>7/NO AUD DX  
VA-BL AUD-C 5-7/NO AUD DX  
VA-BL AUD-C>0 & AUD DX  
VA-FOLLOW-UP POS PTSD/DEP COMBINATION  
VA-BL SHOW COMBINED PTSD/DEPRESSION DISPO  
VA-BL SHOW PTSD DISPO  
VA-BL SHOW DEPRESSION DISPO  
VA-PTSD SCREENING  
VA-DEPRESSION SCREENING

### **REMINDER DIALOG**

VA-GP AAA PRIOR DIAGNOSIS  
VA-PTSD SCREEN  
VA-DEPRESSION SCREEN  
VA-GP ALC OFFER TREATMENT  
VA-AUD-C POS BRIEF ALCOHOL COUNSELING  
VA-ALCOHOL USE SCREENING (AUDIT-C)  
VA-ALCOHOL USE SCREENING (AUDIT-C) RESPONSE  
REQUIRED  
VA-HPV SERIES 2 GIVEN PREVIOUSLY  
VA-HPV SERIES 1 GIVEN PREVIOUSLY  
VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT  
VA-AH-BPR PURPOSE STATEMENT ADDL INFO  
VA-AH-BPR PURPOSE STATEMENT  
VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER  
VA-GP WH PCC EXPOSURE SUB ABUSE HX  
VA-GP WH PCC ASSESSMENT AND PLAN  
VA-TDR-DG READER DX LESION ACT KER (A)  
VA-FOLLOW-UP POS PTSD/DEP COMBINATION  
VA-PTSD SCREENING

### **TIU TEMPLATE FIELD**

BLANK TEXT FIELD FOR DIALOGS  
BLANK SPACE1  
VA-URL PHQ9  
VA-DEP CPG

BLANK TEXT FOR DIALOGS  
YES/NO  
NUMB1-16  
VA-SUD CPG  
VA-URL AUDIT-C  
VA-BENE TRAV INTER-FACILITY TRANSFER  
AH-BPR BURN PIT REGISTRY  
AH-BPR SELF-ASSESSMENT PT RESPONSES  
VA-BENE TRAV MEDICAL JUSTIFICATION STATEMENT  
DISPLAY  
VA-BENE TRAV WORD PROCESSING  
VA-BENE TRAV SPECIAL MODE OTHER NEEDS O2  
VA-BENE TRAV SPECIAL MODE OTHER NEEDS IV  
VA-BENE TRAV TYPE WC  
VA-BENE TRAV WORD PROCESS RQD  
VA-BENE TRAV WORD INDENT REQ  
VA-BENE TRAV SM ALS DISPLAY  
VA-BENE TRAV SPECIAL MODE EST TIME FRAME  
STATEMENT DISPLAY  
VA-BENEFIT TRAVEL  
VA-BENE DATE REQUIRED  
VA-BENE TRAV CLINICAL CONDITION  
VA-BENE TRAV YES NO  
VA-BENE TRAV TYPE OF CARE PROVIDED  
VA-BENE 45 BOX  
VA-BENE TRAV TRANSPORTATION CC MODES  
VA-BENE TRAV POV NOT ACCESSABLE  
VA-BENE TRAV MEDICALLY NECESSARY  
VA-MH FREQ  
VA-MH WHEN  
VA-PTSD CPG  
VA-PTSD PCL INFO

## **HEALTH FACTORS**

VA-REMINDER UPDATES  
VA-UPDATE\_2\_0\_24  
AAA SCREENING AND F/U  
ABD AORTIC ANEURYSM >5.4 CM  
ABD AORTIC ANEURYSM 4.0-5.4 CM  
ABD AORTIC ANEURYSM 3.0-3.9 CM  
MENTAL HEALTH  
REFUSED PTSD SCREEN  
SEVERE CHRONIC COGNITIVE IMPAIRMENT  
UNABLE TO SCREEN-ACUTE MED CONDITION  
UNABLE TO SCREEN - ACUTE ILLNESS  
REFUSED DEPRESSION SCREENING

ALCOHOL USE  
REFUSED REFERRAL TO SUBSTANCE ABUSE  
VA-ALCOHOL OTHER  
VA-ALCOHOL REFER TO SUD PROGRAM  
VA-ALCOHOL REFER TO PCMH  
VA-ALCOHOL TREAT IN PRIMARY CARE  
VA-ALCOHOL SELF HELP/AA MEETINGS  
VA-ALCOHOL USE AND MEDICAL PROBLEMS  
ALCOHOL - ADVISE TO ABSTAIN  
ALCOHOL - ADVISE ON SAFE LIMITS  
ALCOHOL COUNSELING - PT. RESPONSE  
ALCOHOL - RESPONSE TO COUNSELING  
ALCOHOL USE WITHIN SAFE LIMITS  
VA-ALCOHOL ABSTAIN OTHER ADVICE  
REFUSED ALCOHOL USE SCREENING  
REMINDER FACTORS  
LIFE EXPECTANCY < 1 YEAR  
VA-BENE SM TYPE  
VA-BENE SM TYPE WC  
VA-BENE SM BLS REASON  
VA-BENE SM BLS RIGID NECK BRACE  
VA-BENE SM BLS VRE  
VA-BENE SM BLS OTHER  
VA-BENE SM BLS 02  
VA-BENE SM BLS NG TUBE  
VA-BENE SM BLS MRSA  
VA-BENE SM BLS WEARING HALO  
VA-BENE SM BLS FOLEY  
VA-BENE SM BLS DNR  
VA-BENE SM BLS DNI  
VA-BENE SM BLS BARIATRIC STRETCHER  
VA-BENE SM TYPE BLS TYPE  
VA-BENE SM TYPE BLS TYPE RESTRAINT  
VA-BENE SM TYPE BLS TYPE CARE  
VA-BENE SM TYPE BLS TYPE STRETCHER  
VA-BENE SM TYPE BLS  
VA-BENE ALS REASON  
VA-BENE SM ALS VRE+  
VA-BENE SM ALS RIGID NECK BRACE  
VA-BENE SM ALS RESTRAINTS/HIGH-RISK  
VA-BENE SM ALS OTHER  
VA-BENE SM ALS OPEN IV  
VA-BENE SM ALS OPEN CENTRAL LINE  
VA-BENE SM ALS OXYGEN  
VA-BENE SM ALS NG TUBE  
VA-BENE SM ALS MONITOR

VA-BENE SM ALS INTRAOUSSEOUS LINE  
VA-BENE SM ALS FOLEY  
VA-BENE SM ALS ENDOTRACHEAL INTUBATION  
VA-BENE SM ALS EMERGENT TRANSPORT  
VA-BENE SM ALS DNR  
VA-BENE SM ALS DNI  
VA-BENE SM ALS HIGH CPR  
VA-BENE SM ALS CHEST TUBE  
VA-BENE SM ALS BARIATRIC STRETCHER  
VA-BENE SM ALS AIR AMBULANCE  
VA-BENE SM TYPE ALS  
VA-BENE SM FREQUENCY  
VA-BENE SM FREQUENCY ROUND TRIP  
VA-BENE SM FREQUENCY 1 WAY  
VA-BENE SM DURATION  
VA-BENE SM 1 YEAR  
VA-BENE 6 MONTHS  
VA-BENE SM 3 MONTHS  
VA-BENE SM ONE TIME  
VA-BENE SM MEDICAL JUSTIFICATION  
VA-BENE SM MEDICAL JUSTIFICATION NO  
VA-BENE SM MEDICAL JUSTIFICATION YES  
VA-BENE SM FACILITY TRANSFER  
VA-BENE SM FACILITY TRANSFER NO  
VA-BENE SM FACILITY TRANSFER YES  
VA-WH PCC ASSESSMENT CAT  
VA-WH PCC ASSESSMENT PLAN  
VA-TDR CATEGORY-PROBLEM A  
VA-TDR DX LESION A-ACTINIC KERATOSIS  
SUICIDE RISK SCREEN  
PRIOR SUICIDE ATTEMPT  
NO PRIOR SUICIDE ATTEMPT  
SUICIDE PLAN  
NO SUICIDE PLAN  
SUICIDAL THOUGHTS  
NO THOUGHTS OF SUICIDE  
HOPELESSNESS PRESENT/FUTURE  
NO HOPELESSNESS  
PTSD ASSESSMENT COMPLETE  
DEPRESSION ASSESS COMPLETE  
VA-PTSD/DEP OTHER DISPOSITIONS  
VA-REFUSE INTERVENTION PTSD/DEP  
VA-MANAGE PTSD/DEP IN PRIMARY CARE  
REFERRAL TO MHC NON-URGENT  
REFERRAL TO MHC IMMEDIATE (EMERGENT)  
VA-OUTSIDE EVAL/TREATMENT FOR PTSD/DEP

VA-CURRENT F/U OR RX FOR PTSD/DEP SCR  
VA-NO INTERVENTION FOR PTSD/DEP SCR  
PTSD SCREEN POSITIVE OTHER  
REFUSES MH REFERRAL FOR PTSD  
PTSD TO BE MANAGED IN PC  
OUTSIDE EVAL/TREATMENT FOR PTSD  
CURRENT F/U OR RX FOR PTSD  
NO INTERVENTION NEEDED FOR PTSD SCR  
DEPRESSION SCREEN POSITIVE PC MANAGED  
DEPRESSION SCREEN POSITIVE OTHER  
REFUSE DEPRESSION INTERVENTION  
OUTSIDE EVAL/TREATMENT FOR DEPRESSION  
CURRENT F/U OR RX FOR DEPRESSION  
NO INTERVENTION NEEDED FOR DEP SCR  
LIFE EXPECTANCY < 6 MONTHS  
REFER FOR ALCOHOL TREATMENT  
NON-DRINKER (NO ALCOHOL FOR >1 YR)  
ALCOHOL - PRIOR TREATMENT  
ALCOHOL - ADDL ASSESSMENT OF USE  
ALCOHOL - ASSESS USE  
ALCOHOL - NO PRIOR TREATMENT

#### **EDUCATION TOPICS**

ALCOHOL USE AND MEDICAL PROBLEMS  
MEDICAL PROBLEMS OF ALCOHOL (SCREENING)

#### **REMINDER SPONSOR**

VA National Center for Health Promotion and Disease Prevention (NCP)  
National Center for Health Promotion and Disease Prevention (NCP)  
Mental Health Services  
National Clinical Practice Guideline Council  
Office of Patient Care Services  
Office of Public Health and Environmental Hazards  
Office of Nursing Service  
Office of Telehealth Services Store and Forward Telehealth  
Office of Quality & Performance  
Mental Health and Behavioral Science Strategic Group  
VHA National Center for Health Promotion and Disease Prevention (NCP) and VHA  
Office of Women's Health

#### **REMINDER LOCATION LIST**

VA-SUD CLINICS

#### **REMINDER TAXONOMY**

VA-PNEUMOC DZ RISK - IMMUNOCOMPROMISED  
VA-MAMMOGRAM/SCREEN  
VA-ABD AORTIC ANEURYSM  
VA-AUD DIAGNOSIS  
VA-DEMENTIA DIAGNOSIS  
VA-PTSD DX OUTPT VISIT  
VA-TERMINAL CANCER PATIENTS  
VA-DEPRESSION DX OUTPT VISIT

**REMINDER TERM**

VA-REMINDER UPDATE\_2\_0\_24  
VA-BL AUDC SCORE 1 TO 4 NO AUD DX  
VA-SEX  
VA-BL AUDIT-C >7  
VA-BL AUDC 5 SAFE LIMITS  
VA-AUDC SCORE=0  
VA-BL AUD-C>7/NO AUD DX  
VA-BL AUD-C 5-7/NO AUD DX  
VA-BL AUD DIAGNOSIS  
VA-PTSD SCREEN POSITIVE  
VA-PRIOR SUICIDE ATTEMPT  
VA-NO PRIOR SUICIDE ATTEMPT  
VA-SUICIDE PLAN  
VA-NO SUICIDE PLAN  
VA-SUICIDAL THOUGHTS  
VA-NO THOUGHTS OF SUICIDE  
VA-HOPELESSNESS  
VA-NO HOPELESSNESS  
VA-PTSD/DEPRESSION ASSESSMENT COMPLETE  
VA-DEPRESSION SCREEN POSITIVE  
VA-BL FOLLOW-UP POS PTSD/DEP COMBINATION  
VA-BL SHOW COMBINED PTSD/DEPRESSION DISPO  
VA-BL SHOW PTSD DISPO  
VA-BL SHOW DEPRESSION DISPO  
VA-VETERAN PATIENT  
VA-MH ACUTE ILLNESS  
VA-COGNITIVE IMPAIRMENT  
VA-PTSD DIAGNOSIS  
VA-LIFE EXPECTANCY < 6 MONTHS  
VA-REFUSED PTSD SCREEN  
VA-PTSD SCREEN  
VA-REFUSED DEPRESSION SCREENING  
VA-DEP PHQ2  
VA-DEPRESSION DIAGNOSIS  
VA-POS ALCOHOL SCREEN  
VA-REFUSED ALCOHOL SCREENING  
VA-ALCOHOL USE SCREEN

VA-ALC USE INTERVENTIONS  
VA-AUD DIAGNOSIS  
VA-SUD CLINIC VISIT  
VA-ALCOHOL NONE PAST 1YR  
VA-ALC MEDICAL PROBLEMS  
VA-ALCOHOL USE INFORMATION

### 1. Update to Reminder Term: VA-AAA SCREENING NOT APPLICABLE

Finding item 1, Health Factor IMAGING FOR AAA N/A CURRENTLY beginning date and time changed to T-1Y

#### Before:

Findings:  
Finding Item: IMAGING FOR AAA N/A CURRENTLY (FI(1)=HF(866))  
Finding Type: HEALTH FACTOR  
Health Factor Category: IMAGING  
Beginning Date/Time: T-6M

#### After:

Findings:  
Finding Item: IMAGING FOR AAA N/A CURRENTLY (FI(1)=HF(493))  
Finding Type: HEALTH FACTOR  
Health Factor Category: IMAGING  
Beginning Date/Time: T-1Y

### 2. Update of Reminder Dialog Group: VA-TDR-DG READER DX LESION ACT KER (A)

HF changed from VA-TDI DX LESION A-ACTINIC KERATOSIS to VA-TDR DX LESION A-ACTINIC KERATOSIS as finding item for reminder dialog element VA-TDR-DG READER DX LESION ACT KER (A).

#### Before:

Name: VA-TDR-DG READER DX LESION ACT KER  
Finding Item: HF VA-TDI DX LESION A-ACTINIC KERATOSIS

#### After:

Name: VA-TDR-DG READER DX LESION ACT KER  
Finding Item: HF VA-TDR DX LESION A-ACTINIC KERATOSIS

### 3. Update to Reminder Dialog Group: VA-GP WH PCC ASSESSMENT AND PLAN

Dialog group VA-GP WH PCC ASSESSMENT AND PLAN updated to add the alternate progress note text and remove the Roman Numeral XII to be consistent with other sections of the dialog.

#### Before:

Name: VA-GP WH PCC ASSESSMENT AND PLAN  
Dialog text: XII. PRECONCEPTION CARE ASSESSMENT AND PLAN  
Alternate P/N text:

#### After:


Name: VA-GP WH PCC ASSESSMENT AND PLAN  
Dialog text: XII. PRECONCEPTION CARE ASSESSMENT AND PLAN  
Alternate P/N text: PRECONCEPTION CARE ASSESSMENT AND PLAN

#### 4. Update to Reminder Dialog Elements: VA-AH-BPR PURPOSE STATEMENT and VA-AH-BPR PURPOSE STATEMENT ADDL INFO


The two TXML fields containing URLs within the elements : VA-AH-BPR PURPOSE STATEMENT and VA-AH-BPR PURPOSE STATEMENT ADDL INFO were updated.

##### Before:

1.




2.




##### After:

1.



2.



## 5. Update to Reminder Dialog Group: VA-GP WH PCC EXPOSURE SUB ABUSE HX

Dialog group VA-GP WH PCC EXPOSURE SUB ABUSE HX updated to correct a spelling error. The words illicit and personal were spelled incorrectly.

### Before:

Name: VA-GP WH PCC EXPOSURE SUB ABUSE HX  
Dialog text: Substance abuse: If patient has history of  
illicit drug use or abuse of legal drugs, counsel about  
presonal and pregnancy related risks and importance  
of cessation before and during pregnancy.


### After:

Name: VA-GP WH PCC EXPOSURE SUB ABUSE HX  
Dialog text: Substance abuse: If patient has history of  
illicit drug use or abuse of legal drugs,  
counsel about personal and pregnancy related risks  
and importance of cessation before and during  
pregnancy.


## 6. Update to Reminder Dialog Element: VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT

The URL in the txml field contained in the element VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT was updated to reflect the updated directive name and web link

### Before:



### After:



## 7. Update to Reminder Dialog Group: VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSER

Dialog group VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSER updated to edit the dialog text and reference the new VA Directive.

Before:

Resources  
[VA Directive 2007-015 - Inter-facility Transfer Policy](#) VHA Directive has been update  
Dialog needs to reflect this update.

INTER-FACILITY TRANSFER:  
Is this Special Mode Transport in relation to an Inter-facility Transfer in accordance with VA  
Directive 2007-015?

[ Yes ]  
 [ No ]

After:

Resources  
[VHA Directive 1094 - Inter-facility Transfer Policy](#) UPDATED DIALOG  
Displays New VHA Directive Text  
and has correct Hyperlink

INTER-FACILITY TRANSFER:  
Is this Special Mode Transport in relation to an Inter-facility Transfer in accordance with VHA  
Directive 1094 Inter-Facility Transfer Policy?

[ Yes ]  
 [ No ]

8. Update to Reminder Taxonomy VA-MAMMOGRAM/SCREEN to show the updated CPT codes of 77065, 77066 and 77067. Codes 77055, 77056 and 77057 were inactivated 1/1/2017.
9. Update of Reminder Dialog Elements . Progress note test read "...dose **if** the HPV vaccine has been documented previously....". Text now reads "...dose **of** the HPV vaccine has been documented previously..."
  - a. VA-HPV SERIES 1 GIVEN PREVIOUSLY
  - b. VA-HPV SERIES 2 GIVEN PREVIOUSLY
10. A new Reminder Dialog was added to allow sites that do the Alcohol Screening and the Follow-Up in the same dialog instead of using 2 different dialogs. The new dialog is VA-ALCOHOL USE SCREENING (AUDIT-C) WITH FOLLOW-UP. When the new web host file is installed, your reminder definition VA-ALCOHOL USE SCREEN (AUDIT-C) will not have a dialog associated with it. If your site wants to configure it so user can perform both the screening and the follow up in the same dialog, you will link VA-ALCOHOL USE SCREENING (AUDIT-C) WITH FOLLOW-UP dialog to your Alcohol Screening reminder. If your site will use 2 dialog system, then you will need to link VA-ALCOHOL USE SCREENING (AUDIT-C) to your Alcohol Screening reminder.
11. Reminder Dialog Group VA-GP ALC OFFER TREATMENT was made LOCAL instead of NATIONAL to all facilities to customize. This group is part of Reminder Dialogs VA-AUD-C POS BRIEF ALCOHOL COUNSELING and VA-ALCOHOL USE SCREENING (AUDIT-C) WITH FOLLOW-UP
12. It was found that it was possible for the Alcohol Screening and Alcohol Follow-up reminder to both be due at the same time if an old positive AUDC existed and the follow-up had never been completed. To alleviate this, a beginning date/time entry of T-14D was added to FI(1) and FI(10) of the reminder definition VA-ALCOHOL AUDIT-C POSITIVE F/U EVAL so that if the Alcohol Screen is positive and a year old (without the follow-up being done), the Follow-up Reminder will no longer be DUE. When this scenario occurs, the VA-ALCOHOL USE

SCREEN (AUDIT-C) will then become DUE NOW Custom cohort logic was changed to CUSTOMIZED PATIENT COHORT LOGIC (OPTIONAL):  
 (FI(1)&FF(9)&'(FF(7)!FF(3))!((FI(1)&(FF(7)!FF(3))!((FI(9)&FI(10))))&'(FI(1)&FI(7)&FF(6)) For CLIN2 us only →created an exchange entry in NATREM prior to V3 VA-ALCOHOL USE SCREEN (AUDIT-C) prior to update\_2\_0\_24\_3

13. Check box for Dispositions when a Depression/PTSD screen was positive has been removed and user will now be required to make at least one selection for Dispositions.

**Before:**

DEPRESSION SCREEN DISPOSITIONS

I have reviewed the results of the Depression screen and have evaluated the patient including inquiry about feelings of hopelessness, suicidal thoughts, suicide plan if thoughts are present, and prior suicide attempts. Based on the evaluation, the following disposition plan will be implemented:

**After:**

PTSD/DEPRESSION SCREEN DISPOSITIONS

I have reviewed the results of the PTSD and DEPRESSION screens and have evaluated the patient including inquiry about feelings of hopelessness, suicidal plan if thoughts are present, and prior suicide attempts. Based on the evaluation, the following disposition plan will be implemented:

14. Radio button for “feedback” section on Alcohol Follow-up reminder removed. Now required to select a problem addressed. This was done by adding element VA-ALC SCREEN POS MEDICAL PROBLEMS

**Before:**

Brief Alcohol Counseling:

FEEDBACK (Required)

Feedback regarding relationship of alcohol to patient's health issues and risks were reviewed with patient.  
 Check specific problem(s) addressed:

**After:**

FEEDBACK (Required)

Feedback regarding relationship of alcohol to patient's health issues and risks were reviewed with patient

Check specific problems addressed. (One or more required)

\*

high blood pressure  heart disease  liver disease  seizures  injury  medication interactions  
 depression and anxiety  insomnia  bleeding from the stomach  stroke  dementia  cancers  Other:

15. Alcohol follow up.....radio button for advised patient to abstain from drinking.....now requires a selection

**Before:**

ADVISE (One is required)

Advised patient to drink within recommended limits: No more than 2 drinks/day on average and no more than 4 drinks on any one day AND no more than 14 drinks per week.

Advised patient to abstain from drinking alcohol due to contraindications.  
 Patient should abstain due to:

Medication interaction  Medical condition  Operating heavy machinery  Pregnant/attempting pregnancy  Other:

**After:**

ADVISE (One is required)

Advised patient to drink within recommended limits: No more than 2 drinks/day on average and no more than 4 drinks on any one day AND no more than 14 drinks per week.

Advised patient to abstain from drinking alcohol due to contraindications.  
 Patient should abstain due to:

\*

Medication interaction  Medical condition  Operating heavy machinery  Pregnant/attempting pregnancy  
 Other:

16. The following two taxonomies were updated:
  - a. VA-IMAGING FOR AAA (NON-SPECIFIC)
  - b. VA-PNEUMOC DZ RISK – IMMUNOCOMPROMISED

The update includes an addition to the taxonomy VA-IMAGING FOR AAA (NON-SPECIFIC) which is included in the term that is finding 9 below. This may not completely resolve the AAA screening reminder since anything that is a “non-specific” CPT code requires either that the provider evaluate the study and check off that AAA screening was completed or that the radiology procedure is included in either finding 8 or finding 12.

If you have a new radiology procedure that is always used for AAA screening and associated with this new CPT code, please be sure to put the radiology procedure in the reminder term that is finding 8 below.


RT VA-IMAGING FOR AAA (CPT NON-SPECIFIC PROC)	Finding # 9
RT VA-IMAGING FOR AAA (RAD NON-SPECIFIC PROC)	Finding # 12
RT VA-IMAGING FOR AAA (RAD SPECIFIC PROC)	Finding # 8

VA-PNEUMOC DZ RISK – IMMUNOCOMPROMISED had SNOMED code 707147002 added.  
Asplenia (disorder)


17. Reminder Dialog group VA-PTSD SCREEN POS REVIEW EVAL GP had incorrect grammar in dialog text. Changed from “I have reviewed the results of the PTSD screen and have evaluated **the the** patient” to remove the extra “the”.
18. An issue was brought up that when the AUDC score is 1 to 4 and patient doesn’t have an AUD diagnosis, there are times when the provider would still want to do a follow-up. The current dialog won’t allow that to happen.




This was changed to the below screen shot when using the dialog for screening and follow-up



This was changed to the below screen shot when using the stand alone follow-up for positive alcohol. You will have to navigate to the Not Applicable folder to be able to invoke the dialog since it is not currently due.



19. Removed the heading 'SCREEN FOR ALCOHOL' when you click 'Declined alcohol screening' checkbox



20. Other grammatical fixes made for alcohol follow up

- 21. Added option for staff to be able to choose that AUDC was positive and the Primary Care provider was notified. This dialog element (VAL-POS AUDC NOTIFY PC PROVIDER )was made local so that sites can configure whether or not they want to require a primary provider be documented using the PXR COMMENT box.

Annual screening for problem drinking, using the AUDIT-C is required for all patients.  
 -Screening questions should be asked verbatim -

AUD-C

Unable to Screen

[VA/DOD CPG for Mgmt of Substance Use Disorders](#)  
[AUDIT-C Questionnaire](#)

---

Follow up Alcohol Screen (Physician, Mental Health, PA, APN, SW, RN only)  
 Options when a patient does not have an AUD diagnosis

Patient has AUDC score of 5-7 without a diagnosis of Alcohol Use Dependence (AUD)  
 Patient has AUDC score >7 without a diagnosis of Alcohol Use Dependence (AUD)  
 Patient had an AUDC score of 0 and does not require a follow-up  
 Patient had an AUDC score of 1 to 4 but still needs alcohol counseling  
 Patient unable to be screened for alcohol use  
 Patient had an AUDC score that was positive. Primary Care provider notified  
 PC Provider notified:

- 22. Pop up in CPRS was removed when a score of 3 or 4 for females or a score of 4 for men on AUDC
- 23. Added text to meet EPRP reviewer requirements for showing concern for patients alcohol use. Concern expressed to patient about their alcohol use. Patient to be followed up with at next ambulatory visit.
- 24. Text was changed for radio button when patient was being screened by a non-credentialed provider and the screen was positive. This is a change from item 21 above.

Follow up Alcohol Screen (Physician, Mental Health, PA, APN, SW, RN only)  
 Patient has AUD diagnosis within the past year. Follow up is required for any AUDC score >0

Patient has AUDC score >0 and has diagnosis of Alcohol Use Dependence (AUD)  
 Patient had a negative AUDC score and does not require a follow-up  
 Patient unable to be screened for alcohol use  
 Patient had an AUDC score that was positive. Licensed provider notified of need for follow-up

- 25. It was decided to give sites flexibility on whether they want to force user to make a choice on the Alcohol Screening dialog. To do this, we made this group VA-GP ALC SCREEN AUDC local. The name was changed to VAL-GP ALC SCREEN AUDC. To force the user to make a choice on the dialog, you will edit the group and set the EXCLUDE FROM PROGRESS NOTE entry to NO or no entry. This was sent out to NOT force users to make a choice.

26. It was found that in the VA-AAA SCREENING dialog, when user clicked box for “Known Abdominal Aortic Aneurysm”, all 3 ICD10 Diagnosis were being entered into PCE. This was because dialog group VA-GP AAA PRIOR DIAGNOSIS was set to NO PICK LIST in the Taxonomy Pick List field. This was changed to be ICD Diagnosis Only. Below are screenshots before and after the change.

**Before:**

Reminder Resolution: Screen for Abd Aortic Aneurysm

VHA National Center for Health Promotion and Disease Prevention (NCP) recommends that men between the ages of 65 and 75 years who are current or former smokers undergo a one-time ultrasound screening for abdominal aortic aneurysm.

[NCP Recommendations on Screening for AAA](#)

Order Ultrasound

Record Prior or Outside Ultrasound, CT or MRI Results for AAA screening  
(Be sure that the testing done was adequate to fully evaluate for AAA)

Known Abdominal Aortic Aneurysm  
choose one if known

Abdominal aortic aneurysm present - 3.0-3.9 cm.

Abdominal aortic aneurysm present - 4.0-5.4 cm.

Abdominal aortic aneurysm present - >5.4 cm.

Repair of AAA previously

Not candidate for AAA repair at this time: review in 1 year.

EXCLUSIONS

Refuses Ultrasound for Screening for AAA

Clear Clinical Maint Visit Info < Back Next > Finish Cancel

Screen for Abd Aortic Aneurysm:  
The patient is known to have an abdominal aortic aneurysm.

Diagnoses: Abdominal Aortic Aneurysm, Ruptured (ICD-10-CM I71.3), Abdominal Aortic Aneurysm, without Rupture (ICD-10-CM I71.4), Dissection of Abdominal Aorta (ICD-10-CM I71.02)

\* Indicates a Required Field

## After:

Reminder Resolution: Screen for Abd Aortic Aneurysm

VHA National Center for Health Promotion and Disease Prevention (NCP) recommends that men between the ages of 65 and 75 years who are current or former smokers undergo a one-time ultrasound screening for abdominal aortic aneurysm.

[NCP Recommendations on Screening for AAA](#)

Order Ultrasound

Record Prior or Outside Ultrasound, CT or MRI Results for AAA screening  
(Be sure that the testing done was adequate to fully evaluate for AAA)

Known Abdominal Aortic Aneurysm

Selectable Diagnoses: \* ((None Selected))

choose one if known

Abdominal aortic aneurysm present - >5.4 cm.

Dissection Of Abdominal Aorta (I71.02)

Abdominal Aortic Aneurysm, Ruptured (I71.3)

Abdominal Aortic Aneurysm, Without Rupture (I71.4)

OK Cancel

Repair of AAA previously

Not candidate for AAA repair at this time: review in 1 year.

EXCLUSIONS

Clear Clinical Maint Visit Info < Back Next > Finish Cancel

Screen for Abd Aortic Aneurysm:  
The patient is known to have an abdominal aortic aneurysm.

<No encounter information entered>

\* Indicates a Required Field

27. Added Reminder Term VA-UPDATE\_2\_0\_24 that will add health factor named VA-UPDATE\_2\_0\_24. This health factor is used to track installation of the exchange file.

## Install Details

This update is being distributed as a web host file. The address for the host file is:

**[http://vista.med.va.gov/reminders/UPDATE\\_2\\_0\\_24.PRD](http://vista.med.va.gov/reminders/UPDATE_2_0_24.PRD)**

The file will be installed using Reminder Exchange, programmer access is not required.

**Installation:**

=====

This update can be loaded with users on the system. Installation will take approximately an hour.

## Pre-Install

---


1. Using the Vista menu, Reminder Dialog Management, Reminder Dialogs, change view to Group, inquire into group VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER YES and identify the local orderable item used as a finding item. During the installation, when this prompt “FINDING entry Q. ORZ BENE TRAV OPTIONAL INTER-FACILITY does not exist. Select one of the following:” appears during installation, you will replace with your local QO for beneficiary travel inter-facility transfer. Your facility may not have an order as a finding item for this group. If you do not have an order as a finding item, when this prompt “FINDING entry Q. ORZ BENE TRAV OPTIONAL INTER-FACILITY does not exist. Select one of the following:” appears during installation, you will delete the finding item.
2. Using the Vista menu, Reminder Dialog Management, Reminder Dialogs, change view to Group, inquire into group VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER NO and identify the local orderable item used as a finding item. During the installation, when this prompt “FINDING entry Q. ORZ VA BENEFICIARY TRAVEL does not exist. Select one of the following:” appears during installation, you will replace with your local QO for beneficiary travel
3. Using the Vista menu, Reminder Term Management, Inquire about Reminder Term, inquire into the reminder term VA-AAA SCREENING NOT APPLICABLE. Save the inquiry for mapping post install if needed. You will OVERWRITE this term during the install and may need to see what you previously had mapped.
4. Here are the general rules to use when doing the install
  - a. If a reminder component does not exist, use the INSTALL action
  - b. If a reminder component exists, but the packed component is different, use the OVERWRITE action except when the component is a Reminder Term. If the component is a Reminder Term, you should use the MERGE option to make sure you don't lose your local mapped items. With this install, there is an exception to that rule with the VA-AAA SCREENING NOT APPLICABLE. You will use the OVERWRITE action for this one.

## Install Example

---

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      IH  Installation History
CHF  Create Host File                 LHF  Load Host File
CMM  Create MailMan Message           LMM  Load MailMan Message
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI  Reminder Definition Inquiry
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: http://vista.med.va.gov/reminders/UPDATE_2_0_24.PRD
```



At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type or copy and paste the following web address:

[http://vista.med.va.gov/reminders/UPDATE\\_2\\_0\\_24.PRD](http://vista.med.va.gov/reminders/UPDATE_2_0_24.PRD)

`http://vista.med.va.gov/reminders/UPDATE_2_0_24.PRD` successfully loaded.

You should see a message at the top of your screen that the file successfully loaded.

Search and locate the entry that begins with UPDATE\_2\_0\_24 in reminder exchange. If you possibly have installed more than once, you want to use the entry with date of 5/2/2017@12:34:26

+Item	Entry	Source	Date Packed
285	UPDATE_2_0_24 MAMMOGRAM TAXONOMY UPDATE	MONTGOMERY@ZZ ALBANY	03/22/2017@07:49
286	UPDATE_2_0_24 MISC UPDATE 3	MONTGOMERY@SALT LAKE CI	03/22/2017@07:43
287	UPDATE_2_0_24 MISC UPDATE 3	MONTGOMERY@SALT LAKE CI	03/28/2017@10:13
288	UPDATE_2_0_24 MISC UPDATE 3	MONTGOMERY@SALT LAKE CI	04/18/2017@14:58
289	UPDATE_2_0_24 MISC UPDATE 3	MONTGOMERY@SALT LAKE CI	04/19/2017@08:47
290	UPDATE_2_0_24 MISC UPDATE 3	MONTGOMERY@SALT LAKE CI	05/02/2017@12:34
291	UPDATE_2_0_4 VA-VETERANS CHOICE/MAMMOGRAM/TBI UPDATE	GRIFFITH@SALT LAKE CI	12/07/2015@06:03
292	UPDATE_2_0_5 MISC UPDATE 1	GRIFFITH@SALT LAKE CI	05/03/2016@09:42

+      + Next Screen   - Prev Screen   ?? More Actions      >>>

CFE Create Exchange File Entry      IH    Installation History  
CHF Create Host File                    LHF   Load Host File  
CMM Create MailMan Message            LMM   Load MailMan Message  
DFE Delete Exchange File Entry        LWH   Load Web Host File  
IFE Install Exchange File Entry        RI    Reminder Definition Inquiry  
Select Action: Next Screen// ife    Install Exchange File Entry  
Enter a list or range of numbers (1-462): 290

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entries titled UPDATE\_2\_0\_24 MISC UPDATE 3 *(in this example it is entry 290, it will vary by site).*

Component	Category	Exists
Source: MONTGOMERY,ALAN at SALT LAKE CITY		
Date Packed: 05/02/2017@12:34:26		
Package Version: 2.0P35		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-GP AAA PRIOR DIAGNOSIS		
VA-DEPRESSION SCREEN		
VA-AH-BPR PURPOSE STATEMENT ADDL INFO		
VA-PTSD SCREEN		
VA-AH-BPR PURPOSE STATEMENT		
VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT		
VA-GP ALC OFFER TREATMENT		
VA-GP WH PCC ASSESSMENT AND PLAN		
VA-HPV SERIES 1 GIVEN PREVIOUSLY		
+ Enter ?? for more actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// IA		

At the **Select Action** prompt, type **IA** for Install All Component and hit enter.

During the install, you will see multiple prompts. Please follow the guidance below.

1. If a reminder component does not exist, use the INSTALL action
2. If a reminder component exists, but the packed component is different, use the OVERWRITE action except when the component is a Reminder Term. If the component is a Reminder Term, you should use the MERGE option to make sure you don't lose your local mapped items. With this install, there is an exception to that rule with the **VA-AAA SCREENING NOT APPLICABLE**. You will use the OVERWRITE action for this one Reminder Term.

**There are multiple dialogs included in this update, you will be prompted several times to install the dialog components**

Item	Seq.	Dialog Findings	Type	Exists		
Dialog Components Mar 22, 2017@10:47:02 Page: 1 of 2						
Packed reminder dialog: VA-GP ALC OFFER TREATMENT						
1		VA-GP ALC OFFER TREATMENT Finding: *NONE*	group			
2	5	VA-ALC SELF-HELP/AA MEETINGS Finding: VA-ALCOHOL SELF HELP/AA MEETINGS (HEALTH FACTOR)	element	X		
3		PXRM COMMENT	prompt	X		
4	10	VA-ALC TREAT IN PRIMARY CARE Finding: VA-ALCOHOL TREAT IN PRIMARY CARE (HEALTH FACTOR)	element	X		
5		PXRM COMMENT	prompt	X		
6	15	VA-ALC REFER TO PCMHI Finding: VA-ALCOHOL REFER TO PCMHI (HEALTH FACTOR)	element	X		
7		PXRM COMMENT	prompt	X		
+      + Next Screen    - Prev Screen    ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA						

At the **Select Action** prompt, type **IA** to install the dialog – **VA-GP ALC OFFER TREATMENT**

All dialog components for VA-GP ALC OFFER TREATMENT are new.

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Item	Seq. Dialog Findings	Type	Exists		
1	VA-ALCOHOL USE SCREENING (AUDIT-C) WITH FOLLOW-UP	dialog			
2	5 VA-TEXT ALCOHOL MISUSE SCREEN INSTRUCTIONS Finding: *NONE*	element			
3	10 VA-GP ALC SCREEN AUDC Finding: *NONE*	group	X		
4	10.5 VA-MH AUDC Finding: AUDC (MH TESTS AND SURVEY)	element	X		
5	PXRM AUDC RESULT GROUP	rs.group	X		
6	PXRM AUDC RESULT ELEMENT 1	rs.element	X		
7	PXRM AUDC RESULT ELEMENT 2	rs.element	X		
8	PXRM AUDC RESULT ELEMENT 3	rs.element			
9	PXRM AUDC RESULT ELEMENT 4	rs.element			
10	10.10 VA-TEXT BLANK LINE WITH TEMPLATE FIELD	element	X		
+      + Next Screen    - Prev Screen    ?? More Actions					
DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		
Select Action: Next Screen// IA    Install All					
Install reminder dialog and all components with no further changes: Y//					

At the **Select Action** prompt, type **IA** to install the dialog – **VA-ALCOHOL USE SCREENING (AUDIT-C) WITH FOLLOW-UP**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-ALCOHOL USE SCREENING (AUDIT-C)

Item	Seq.	Dialog Findings	Type	Exists
1		VA-ALCOHOL USE SCREENING (AUDIT-C)	dialog	X
2	5	VA-TEXT ALCOHOL MISUSE SCREEN INSTRUCTIONS Finding: *NONE*	element	X
3	10	VA-GP ALC SCREEN AUDC Finding: *NONE*	group	X
4	10.5	VA-MH AUDC Finding: AUDC (MH TESTS AND SURVEY)	element	X
5		PXRM AUDC RESULT GROUP	rs.group	X
6		PXRM AUDC RESULT ELEMENT 1	rs.element	X
7		PXRM AUDC RESULT ELEMENT 2	rs.element	X
8		PXRM AUDC RESULT ELEMENT 3	rs.element	X
9		PXRM AUDC RESULT ELEMENT 4	rs.element	X
10	10.10	VA-TEXT BLANK LINE WITH TEMPLATE FIELD	element	X

```

+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
    
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-ALCOHOL USE SCREENING (AUDIT-C)**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-HPV SERIES 2 GIVEN PREVIOUSLY

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-HPV SERIES 2 GIVEN PREVIOUSLY Finding: *NONE*	element			
----- REPLACEMENT ITEMS -----						
2	1	VA-GP PATIENT HAS NO AUD DX OPTIONS Finding: *NONE*	group	X		
3	1.10	VA-GP ALCOHOL SCREEN/FOLLOW-UP 1 Finding: *NONE*	group	X		
4	1.10.5	VA-ALC FOLLOW-UP 1 INTRO TEXT Finding: *NONE*	element	X		
5	1.10.10	BLANK SPACE1 Finding: *NONE*	element	X		
6	1.10.20	VA-GP ALC ADVICE1	group	X		
+      + Next Screen    - Prev Screen    ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA    Install All						

At the **Select Action** prompt, type **IA** to install the dialog – **VA-HPV SERIES 2 GIVEN PREVIOUSLY**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-HPV SERIES 1 GIVEN PREVIOUSLY

Item	Seq.	Dialog Findings	Type	Exists
1		VA-HPV SERIES 1 GIVEN PREVIOUSLY Finding: *NONE*	element	
----- REPLACEMENT ITEMS -----				
2	1	VA-GP PATIENT HAS NO AUD DX OPTIONS Finding: *NONE*	group	X
3	1.10	VA-GP ALCOHOL SCREEN/FOLLOW-UP 1 Finding: *NONE*	group	X
4	1.10.5	VA-ALC FOLLOW-UP 1 INTRO TEXT Finding: *NONE*	element	X
5	1.10.10	BLANK SPACE1 Finding: *NONE*	element	X
6	1.10.20	VA-GP ALC ADVICE1	group	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details		DT Dialog Text	IS Install Selected
DF	Dialog Findings		DU Dialog Usage	QU Quit
DS	Dialog Summary		IA Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog – **VA-HPV SERIES 1 GIVEN PREVIOUSLY**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

```

Packed reminder dialog: VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT
+-----+-----+-----+-----+-----+
Item  Seq.  Dialog Findings                                     Type  Exists
+-----+-----+-----+-----+-----+
  1    VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT
       Finding: *NONE*
+-----+-----+-----+-----+-----+
+ Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details  DT  Dialog Text      IS  Install Selected
DF  Dialog Findings  DU  Dialog Usage    OU  Quit
DS  Dialog Summary  IA  Install All
Select Action: Quit// IA

```

At the **Select Action** prompt, type **IA** to install the dialog **VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

```

Packed reminder dialog: VA-AH-BPR PURPOSE STATEMENT ADDL INFO
+-----+-----+-----+-----+-----+
Item  Seq.  Dialog Findings                                     Type  Exists
+-----+-----+-----+-----+-----+
  1    VA-AH-BPR PURPOSE STATEMENT ADDL INFO
       Finding: *NONE*
+-----+-----+-----+-----+-----+
+ Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details  DT  Dialog Text      IS  Install Selected
DF  Dialog Findings  DU  Dialog Usage    OU  Quit
DS  Dialog Summary  IA  Install All
Select Action: Quit// IA

```

At the **Select Action** prompt, type **IA** to install the dialog **VA-AH-BPR PURPOSE STATEMENT ADDL INFO**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

```
Packed reminder dialog: VA-AH-BPR PURPOSE STATEMENT
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-AH-BPR PURPOSE STATEMENT Finding: *NONE*	group	X
2	5	VA-AH-BPR PURPOSE STATEMENT ADDL INFO CHECKBOX Finding: *NONE*	group	X
3	5.5	VA-AH-BPR PURPOSE STATEMENT ADDL INFO Finding: *NONE*	element	X

```

+ Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details  DT Dialog Text      IS Install Selected
DF Dialog Findings DU Dialog Usage     OU Quit
DS Dialog Summary IA Install All
Select Action: Quit// Q

```

At the **Select Action** prompt, type **IA** to install the dialog **VA-AH-BPR PURPOSE STATEMENT**  
Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

```
Packed reminder dialog: VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER Finding: *NONE*	group	X
2	10	VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER YES Finding: ORZ VA BENE TRAV OPTIONAL INTER-FACILITY (ORDER DIALOG) Add. Finding: VA-BENE SM FACILITY TRANSFER YES (HEALTH FACTOR)	group	X
3		PXRM VISIT DATE FORCED TODAY	forced	X
4	10.10	VA-BENE TRAV VA-BENE TRAV TEXT Finding: *NONE*	element	X
5	20	VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER NO	group	X

```

+ + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details  DT Dialog Text      IS Install Selected
DF Dialog Findings DU Dialog Usage     OU Quit
DS Dialog Summary IA Install All
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog **VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER**

Install reminder dialog and all components with no further changes: Y// ES  
 After installing, you will need to address the following prompts.

FINDING entry Q.ORZ VA BENEFICIARY TRAVEL does not exist.

Select one of the following:

- D Delete
- P Replace with an existing entry
- Q Quit the install

Enter response: **p Replace with an existing entry**

Select ORDER DIALOG NAME:

**Choose P and replace with the order identified in the pre install for your beneficiary travel quick order**

FINDING entry Q.ORZ VA BENE TRAV OPTIONAL INTER-FACILITY does not exist.

Select one of the following:

- D Delete
- P Replace with an existing entry
- Q Quit the install

**Choose P and replace with the order identified in the pre install for your INTER FACILITY beneficiary travel quick order. If you did not have one then choose D to delete.**

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-GP WH PCC EXPOSURE SUB ABUSE HX

Item	Seq.	Dialog	Findings	Type	Exists
1		VA-GP WH PCC EXPOSURE SUB ABUSE HX	Finding: *NONE*	group	X
2	5	VA-EL WH PCC EXPOSURE SUB ABUSE TREATMENT	Finding: *NONE*	element	X
3	10	VA-EL WH PCC EXPOSURE SUB ABUSE TREATMENT NO	Finding: *NONE*	element	X

+ Next Screen    - Prev Screen    ?? More Actions

DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	OU	Quit
DS	Dialog Summary	IA	Install All		

Select Action: Quit// IA

At the **Select Action** prompt, type **IA** to install the dialog **VA-GP WH PCC EXPOSURE SUB ABUSE HX**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-GP WH PCC ASSESSMENT AND PLAN

Item	Seq.	Dialog Findings	Type	Exists
1		VA-GP WH PCC ASSESSMENT AND PLAN Finding: *NONE*	group	X
2	5	VA-EL WH PCC ASSESSMENT 2 Finding: *NONE*	element	X
3	10	VA-EL WH PCC ASSESSMENT 3 Finding: *NONE*	element	X
4		PXRM COMMENT	prompt	X
5	15	VA-EL WH PCC ASSESSMENT 4 Finding: VA-WH PCC ASSESSMENT PLAN (HEALTH FACTOR)	element	X
6		PXRM COMMENT	prompt	X
7	20	VA-BLANK LINE IN PM	element	X

+ Next Screen - Prev Screen ?? More Actions

DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		

Select Action: Next Screen// IA

At the **Select Action** prompt, type **IA** to install the dialog **VA-GP WH PCC ASSESSMENT AND PLAN**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-TDR-DG READER DX LESION ACT KER (A)

Item	Seq.	Dialog Findings	Type	Exists
1		VA-TDR-DG READER DX LESION ACT KER (A) Finding: VA-TDR DX LESION A-ACTINIC KERATOSIS (HEALTH FACTOR)	group	X
2		VA-TDR PXRM DX LESION ACTINIC KERATOSIS	prompt	X

+ Next Screen - Prev Screen ?? More Actions

DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		

Select Action: Quit// IA

At the **Select Action** prompt, type **IA** to install the dialog **VA-TDR-DG READER DX LESION ACT KER (A)**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-FOLLOW-UP POS PTSD/DEP COMBINATION	dialog			
2	5	VA-PTSD/DEPRESSION COMBINATION MASTER GP* Finding: *NONE*	group			
3	5.5	VA-POS PTSD/DEP HEADER Finding: *NONE*	element			
4	5.10	VA-SUICIDE RISK EVALUATION GP Finding: *NONE*	group			
5	5.10.5	VA-SUICIDE RISK EVAL QUESTION 1 GP Finding: *NONE*	group			
6	5.10.5.5	VA-SUICIDE RISK EVAL QUESTION 1 NO Finding: NO HOPELESSNESS (HEALTH FACTOR)	element	X		
7	5.10.5.10	VA-SUICIDE RISK EVAL QUESTION 1 YES Finding: HOPELESSNESS PRESENT/FUTURE (HEALTH FACTOR)	element	X		
8	5.10.10	VA-SUICIDE RISK EVAL QUESTION 2 GP	group			
+      + Next Screen    - Prev Screen    ?? More Actions						
DD	Dialog Details		DT	Dialog Text	IS	Install Selected
DF	Dialog Findings		DU	Dialog Usage	QU	Quit
DS	Dialog Summary		IA	Install All		
Select Action: Next Screen// IA						

At the **Select Action** prompt, type **IA** to install the dialog **VA-FOLLOW-UP POS PTSD/DEP COMBINATION**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-PTSD SCREENING [NATIONAL DIALOG]

Item	Seq.	Dialog Findings	Type	Exists
1		VA-PTSD SCREENING	dialog	X
2	10	VA-PTSD SCREEN Finding: *NONE*	group	X
3	10.10	VA-MH PC PTSD SCREEN Finding: PC PTSD (MH TESTS AND SURVEY)	element	X
4		PXRM PTSD RESULT GROUP	rs.group	X
5		PXRM PTSD RESULT ELEMENT 1	rs.element	X
6		PXRM PTSD RESULT ELEMENT 2	rs.element	X
7	10.15	VA-TEXT BLANK LINE WITH TEMPLATE FIELD Finding: *NONE*	element	X
8	10.20	VA-GP PTSD UNABLE TO SCREEN Finding: *NONE*	group	X
9	10.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details		DT Dialog Text	IS Install Selected
DF	Dialog Findings		DU Dialog Usage	QU Quit
DS	Dialog Summary		IA Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog **VA-PTSD SCREENING**  
 Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-DEPRESSION SCREEN [NATIONAL DIALOG]

Item	Seq.	Dialog Findings	Type	Exists
1		VA-DEPRESSION SCREEN	dialog	X
2	5	VA-DEP PHQ2 Finding: *NONE*	group	X
3	5.10	VA-MH PHQ2 Finding: PHQ-2 (MH TESTS AND SURVEY)	element	X
4		PXRM PHQ2 RESULT GROUP	rs.group	X
5		PXRM PHQ2 RESULT ELEMENT 1	rs.element	X
6		PXRM PHQ2 RESULT ELEMENT 2	rs.element	X
7	5.15	VA-TEXT BLANK LINE WITH TEMPLATE FIELD Finding: *NONE*	element	X
8	5.20	VA-GP DEP UNABLE TO SCREEN Finding: *NONE*	group	X
9	5.20.5	VA-HF ACUTE ILLNESS Finding: UNABLE TO SCREEN - ACUTE ILLNESS (HEALTH FACTOR)	element	X

+ + Next Screen - Prev Screen ?? More Actions

DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		

Select Action: Next Screen// IA

At the **Select Action** prompt, type **IA** to install the dialog **VA-DEPRESSION SCREEN**  
Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q to move to the next dialog to install

Packed reminder dialog: VA-AUD-C POS BRIEF ALCOHOL COUNSELING [NATIONAL DIALOG]

Item	Seq.	Dialog Findings	Type	Exists
1		VA-AUD-C POS BRIEF ALCOHOL COUNSELING	dialog	
2	3	VA-AUDC SCORE=0 Finding: *NONE*	element	
3	4	VA-GP ALCOHOL FOLLOW-UP 0 Finding: *NONE*	group	X
4	4.5	TEXT AUDC PRIOR RESULTS Finding: *NONE*	element	X
5	4.10	VA-GP ALCOHOL AUDIT C POSITIVE Finding: *NONE*	group	X
6	4.10.10	VA-GP ALCOHOL AUDIT C 5-7* Finding: *NONE*	group	X
7	4.10.10.10	VA-GP ALC AUDIT C 5-7 Finding: *NONE*	group	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details		DT Dialog Text	IS Install Selected
DF	Dialog Findings		DU Dialog Usage	QU Quit
DS	Dialog Summary		IA Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog **VA-AUD-C POS BRIEF ALCOHOL COUNSELING**

Install reminder dialog and all components with no further changes: Y// ES

After completing this dialog install, you will type Q and be prompted for the items below.

HEALTH SUMMARY TYPE entry named OB AUDC already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// overwrite the current entry  
Are you sure you want to overwrite? N// y YES  
HEALTH SUMMARY TYPE entry named OB PC-PTSD PCL already exists but the  
packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// overwrite the current entry  
Are you sure you want to overwrite? N// y YES

HEALTH SUMMARY TYPE entry OB PHQ2 PHQ9 is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall  
HEALTH SUMMARY OBJECTS entry named AUDC (TIU) already exists and the  
packed component is identical, skipping.

HEALTH SUMMARY OBJECTS entry named PC-PTSD PCL (TIU) already exists  
and the packed component is identical, skipping.

HEALTH SUMMARY OBJECTS entry PHQ2 PHQ9 (TIU) is NEW,  
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall  
TIU DOCUMENT DEFINITION entry named AUDIT-C already exists and the  
packed component is identical, skipping.

TIU DOCUMENT DEFINITION entry named PC-PTSD PCL already exists and  
the packed component is identical, skipping.

TIU DOCUMENT DEFINITION entry PHQ2 PHQ9 is NEW,

what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall

```
Exchange File Components      May 03, 2017@11:13:57      Page: 1 of 25
-----
Component                      Category      Exists
-----
Source:      MONTGOMERY,ALAN at SALT LAKE CITY
Date Packed: 05/02/2017@12:34:26
Package Version: 2.0P35

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
VA-GP AAA PRIOR DIAGNOSIS
VA-DEPRESSION SCREEN
VA-AH-BPR PURPOSE STATEMENT ADDL INFO
VA-PTSD SCREEN
VA-AH-BPR PURPOSE STATEMENT
VA-BENE TRAV SPECIAL MODE RESOURCES ELEMENT
VA-GP ALC OFFER TREATMENT
VA-GP WH PCC ASSESSMENT AND PLAN
VA-HPV SERIES 1 GIVEN PREVIOUSLY
+      + Next Screen  - Prev Screen  ?? More Actions      >>>
IA  Install all Components      IS  Install Selected Component
Select Action: Next Screen// Q
```

After all the dialogs and components have been installed, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

## Post-Install Set-up Instructions

---

1. VA-AAA reminder dialog

- a. Check that any local components mapped to the reminder term VA-AAA SCREENING NOT APPLICABLE prior to install are present post install. **Review Pre-Install inquiry to assist in identifying the correct order.**

2. VA-BENE TRAVEL reminder dialog

- a. Review reminder dialog group GROUP VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER YES and confirm that your local facility order is included as a finding item. **Review Pre-Install inquiry to assist in identifying the correct order.**  
**Note: If you did not have an order mapped pre-install to this group, you do not need to map one now. This inclusion of an order on this group depends on local policy.**
- b. Review reminder dialog group GROUP VA-BENE TRAV GP OPTIONAL INTER-FACILITY TRANSFER NO and confirm that your local facility order is included as a finding item. **Review Pre-Install inquiry to assist in identifying the correct order.**

3. Link reminder dialog to reminder definition VA-ALCOHOL SCREEN (AUDIT-C). You will choose the dialog that best suits your facility’s practice on whether you will have the Alcohol Follow-up as part of the screening dialog.

Using Reminders Manager Menu, Reminder Dialog Management (DM), Reminder dialogs (DI), in reminder view, search and locate VA-ALCOHOL USE SCREEN. Choose the corresponding number (in this example it is 164 but it will vary by site).

```

REMINDER VIEW (ALL REMINDERS BY NAME)
+Item Reminder Name                               Linked Dialog Name & Dialog Status
164 VA-ALCOHOL USE SCREEN (AUDIT-C)
165 VA-ANTIPSYCHOTIC MED SIDE EFF EVAL VA-AIMS
166 VA-BENZODIAZEPINE (CN302) BUT NOT
167 VA-BL ALCOHOL SCREEN
168 VA-BL ALCOHOL WITHIN SAFE LIMIT VA-AUD-C POS BRIEF ALCOHOL
169 VA-BL AUD-C 5-7/NO AUD DX
170 VA-BL AUD-C>0 & AUD DX
171 VA-BL AUD-C>7/NO AUD DX
172 VA-BL AUDC 1 TO 4 NO AUD DX
173 VA-BL DEPRESSION SCREEN
174 VA-BL LDL CV DISEASE ONLY
175 VA-BL LDL DIABETES DX ONLY
176 VA-BL LDL DM AND CV DISEASE
177 VA-BL MULTIPLE SSD AND MST YES
+ + Next Screen - Prev Screen ?? More Actions >>>
AR All reminders LR Linked Reminders QU Quit
CV Change View RN Name/Print Name
Select Item: Next Screen// 164

```

At the Select Item prompt, Choose LR to link reminder, at the LINKED REMINDER DIALOG prompt, type VA-ALCOHOL USE:

Select Item: Quit// LR Link Reminder

REMINDER NAME: VA-ALCOHOL USE SCREEN (AUDIT-C)

LINKED REMINDER DIALOG: VA-ALCOHOL USE

- 1 VA-ALCOHOL USE SCREENING (AUDIT-C) reminder dialog NATIONAL
- 2 VA-ALCOHOL USE SCREENING (AUDIT-C) WITH FOLLOW-UP reminder Dialog

You will have two choices. Choose item number 1 if you only want the alcohol use screen dialog to display. Choose item number 2 if you want the combined alcohol use screen with follow-up dialog to display.

```

REMINDER NAME: VA-ALCOHOL USE SCREEN (AUDIT-C)
Item Dialog Name          Latest Update  Other Linked Reminders
This reminder is linked to dialog:
  1  VA-ALCOHOL USE SCREENING (AUDIT-
+ Next Screen  - Prev Screen  ?? More Actions  >>>
AD  Autogenerate Dialog      QU  Quit
LR  Link Reminder
Select Item: Quit// 0

```

Type Q to quit.

- 4. It is possible that your reminder dialogs may not be linked to the appropriate reminder definitions. Please ensure that they are linked appropriately

**NOTE:** If you require further technical assistance, please notify your local IT support to log a national CA Service Desk Manager (SDM) ticket (previously a Remedy™ ticket) or contact the VA Service Desk at 1-888-596-4357 and have them submit a national CA ticket to the Incident Area: **NTL.APP.VISTA.CLINICAL REMINDERS 2\_0** and we will contact you.