

## Care Management Setup Checklist

| ACTIVITY  | STATUS | RESPONSIBLE | COMMENTS |
|---|--------|-------------|----------|
| 1. After the package is installed, there is no required set up unless it is intended to implement this package immediately.   |        |             |          |
| <p>A. <b>Assign a Care Management (CM) option appropriate for the user. Following is a list of CM options:</b><br/> ORRCMM CLINICIAN<br/> ORRCMM CLINICIAN QUERY<br/> ORRCMM NURSING<br/> ORRCMM NURSE QUERY<br/> ORRCMM CM CLINICIAN NURSE<br/> ORRCMM CM ALL</p>  |        |             |          |
| <p>B. <b>Assign Care Management Perspective to User:</b><br/> Menu option: CPRS CONFIGURATION (IRM) then<br/> HealtheVet Desktop Configuration menu.<br/> Following is a list of perspectives:<br/> ORRCMM CLINICIAN<br/> ORRCMM CLINICIAN QUERY<br/> ORRCMM NURSING<br/> ORRCMM NURSE QUERY<br/> ORRCMM CM CLINICIAN NURSE<br/> ORRCMM CM ALL</p>                        |        |             |          |
| <p>C. <b>Set the Activation Date:</b><br/> Menu option: CPRS CONFIGURATION MENU –General<br/> Parameter Tools Type EP to select Edit Parameter Values.<br/> Type ORRC ACTIVATION DATE. You can set this<br/> parameter at a user, Service, Division or System level. Set<br/> this value to today’s date, or to the date you expect new<br/> users to begin using CM.</p> |        |             |          |

| ACTIVITY  | STATUS | RESPONSIBLE | COMMENTS |
|---|--------|-------------|----------|
| <p><b>D. Security Keys:</b><br/>ORRC QUERY RESULT EXPORT key should be assigned to users who need to export or print reports generated through the CM Query Tool perspective</p>  |        |             |          |
| <p><b>E. Use the Care Management Installation Guide</b> and follow the steps outlined for installing server side and client side application.</p> <ol style="list-style-type: none"> <li>1) JRE 1.4.2_12 needs to be on each workstation running Care Management.</li> <li>2) Ensure that every workstation has the Sentillion Desktop Components installed (this is usually done when cloning).</li> </ol> |        |             |          |