

## Clinical Reminders Setup Checklist

ACTIVITY	STATUS	RESPONSIBLE	COMMENTS
Required packages to be set up:			
<p><b>A. File Numbers</b></p> <p>811.7 REMINDER CATEGORY            811.4 REMINDER COMPUTED FINDINGS            810.8 REMINDER COUNTING GROUP            811.9 REMINDER DEFINITION            801.41 REMINDER DIALOG            801.5 REMINDER DIALOG PATIENT ASSOCIATION            811.8 REMINDER EXCHANGE            810.7 REMINDER EXTRACT COUNTING RULE            810.2 REMINDER EXTRACT DEFINITION            810.3 REMINDER EXTRACT SUMMARY            801.43 REMINDER FINDING ITEM PARAMETER            801.45 REMINDER FINDING TYPE PARAMETER            802.4 REMINDER FUNCTION FINDING FUNCTIONS            801.55 REMINDER GEC DIALOG ASSOCIATION</p> <p><b>HISTORY</b></p> <p>801.42 REMINDER GUI PROCESS            810.4 REMINDER LIST RULE            810.9 REMINDER LOCATION LIST            810.5 REMINDER PATIENT LIST            810.1 REMINDER REPORT TEMPLATE            801.9 REMINDER RESOLUTION STATUS            811.6 REMINDER SPONSOR            811.2 REMINDER TAXONOMY            811.5 REMINDER TERM            800 CLINICAL REMINDER PARAMETERS</p>			

ACTIVITY	STATUS	RESPONSIBLE	COMMENTS
<p><b>B. Client installation:</b>                      System must be fully patched. Reminders interfaces with many other VistA packages through the use of findings. For reminders to be able to evaluate findings properly, the clinical reminder index global (^PXRMINDX) must be present. This is installed with patch PXR*1.5*12</p>			
<p><b>C. Security Key:</b> PXRMANAGER - This key is assigned to users who are responsible for managing clinical reminders.</p>			
<p><b>D. Menu Option:</b> PXRMANAGERS MENU</p>			
<p><b>E. Editing System Parameters</b></p>			